

US EPA ARCHIVE DOCUMENT

REVISED

NOT TO BE COPIED

Procedure No. 1047 Copy No. 1 Page No. 1 of 7 Pages

Title GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SBP-1382) RESIDUES

IN FOOD PRODUCTS - REVISED

Date 9/9/82

Signed A. J. Furmanec *RF*

CONFIDENTIAL

SUMMARY

A general analytical method for the determination of Resmethrin residues in diverse types of food products is described. The procedure has been successfully applied to the following substances; sliced american cheese, sliced ham, cream pie, white flour, polished rice, white bread, popcorn kernels, saltine crackers, sliced roast beef, chicken meat, skin, fat, hearts, liver, gizzard and eggs.

Whole milk and sugar require a different sample extraction, details are in notes 1 and 2. Recoveries of Resmethrin added to these substances at the 1 ppm level range from 85 to 101% with a detection limit of 0.1 ppm.

Sample preparation involves extraction in a blender with acetonitrile, partitioning between acetonitrile and petroleum ether and clean-up on Florisil^R and Alumin AR columns.

Three chromatographic columns are mandated. The 3% OV-1 column is the primary column. If interfering peaks are encountered chromatography on 3% OV-210 or on "mixed column" (50:50 mixture of 3% OV-17 on 80-100 mesh Varaport 30 and 3% OV-1 on 100-120 mesh Gas Chrom Q) will separate the interferences or vice-versa.

NOT TO BE COPIED

REVISED
Procedure No. 1047 Copy No. 1 Page No. 2 of 7 Pages

Title GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SBP-1382) RESIDUES
IN FOOD PRODUCTS - REVISED

Date 9/9/82

Signed *D. Furmanec RP*

EQUIPMENT

CONFIDENTIAL

Hewlett Packard Gas Chromatograph 5840A
Automatic sampler 7671A Hewlett-Packard
GC Glass Column 6 ft. x 2 mm I.D. 1/4" O.D.
Blender (Waring)
Flash-evaporator (Searle)
Separatory Funnel 2 liter
Round Bottom Flask, 500 ml. 250 ml.
Pear-shaped Flask, 50 ml.
Clean-up column glass, 450 mm long x 15 mm I.D.
Pipet, 1.0 ml., 5.0 ml., 25.0 ml.
Buchner Funnel 70 mm Dia.
Suction Flask, 1000 ml.
Whatman #1 filter paper 70 mm.

MATERIALS:

Acetonitrile, HPLC grade
Hexane, HPLC grade
Petroleum Ether, Pesticide grade
Methylene Chloride, Pesticide grade
Ether Anhydrous, Reagent A.C.S.
Acetone, HPLC grade
Sodium Sulfate, Anhydrous - Reagent grade
Sodium Chloride, Reagent grade
3% OV-1 on 100-120 Mesh on Gas Chrom Q (Applied Science Lab., Inc.)
3% OV-210 on 100-120 Mesh on Gas Chrom Q (Applied Science Lab., Inc.)
3% OV-17 on 80-100 Mesh Varaport 30 (Varian Aerograph)
Florasil 60-100 PR Mesh (Floridin Company)
Alumin AR CC-10 100-200 Mesh Activity Grade 4 (Mallinckrodt Chemical Works)
Resmethrin Standard (Penick Corporation)

NOT TO BE COPIED

REVISÉD
Procedure No. 1047 Copy No. 1 Page No. 3 of 7 Pages

Title GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SBP-1382) RESIDUES
IN FOOD PRODUCTS - REVISED

Date 9/9/82

CONFIDENTIAL

Signed D. Furmanec

PROCEDURE:

Weigh accurately 25.0 gm. of a representative food sample into blender (If milk or sugar see notes 1 and 2). Add 100 ml. of acetonitrile and activate the blender for 5 minutes. Decant the solvent through Whatman #1 filter paper in a Buchner funnel using gentle suction. Repeat the extraction in the blender two additional times using 100 ml and 50 ml portions of acetonitrile and filter each time through the same Whatman #1 filter paper. Transfer combined extracts to a 2 liter separatory funnel.

Acetonitrile - Petroleum Ether Partition

Add 100 ml. petroleum ether to the acetonitrile filtrate and shake well. Then add 75 ml. saturated salt solution, 700 ml. of distilled water and shake again. Transfer the petroleum ether layer to 250 ml. round bottom flask. Extract two additional times with 100 ml. and 50 ml. portions of petroleum ether. Each time remove petroleum ether in the same round bottom flask on a flash evaporator under vacuum at 50°C to an oily residue. Quantitatively transfer the oily residue to 25 ml. volumetric flask with small portions of hexane and dilute to volume.

COLUMN CLEAN-UP:

Florisil^R Column

Pack a chromatographic column (450 mm l x 15 mm I.D.) with 10 gm. of Florisil^R (activated at 600°C). Add a 10 mm layer of anhydrous sodium sulfate, to the top of the column, pre-wet column with hexane and drain to the top of the bed. Transfer 5.0 ml. sample aliquot to the top of the column and drain into the bed of the column.

Elute the column with 100 ml. of solvent mixture 50% ether - 50% hexane, and collect the total eluate. Concentrate the eluate on flash-evaporator under vacuum at 50°C to an oily residue.

NOT TO BE COPIED

REVISÉD
Procedure No. 1047 Copy No. 1 Page No. 4 of 7 Pages

Title GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SBP-1382) RESIDUES

IN FOOD PRODUCTS - REVISED

Date 9/9/82

~~CONFIDENTIAL~~

Signed

D. J. Furmanec RP

Alumin AR Column:

Pack a chromatographic column (450 mm l x 15 mm I.D.) with 20 gm of Alumin AR. Add a 10 mm. layer of anhydrous sodium sulfate to the top of the column. Pre-wet column with hexane and drain to the top of the bed. Transfer quantitatively the residue with three 5 ml. portions of hexane to the top of the column and drain each time into the bed of the column. Elute the column with 100 ml. hexane and discard the eluate. Continue eluation with 100 ml. of solvent mixture 25 ml. methylene chloride and 75 ml. hexane and collect the eluate in 250 ml. round bottom flask. Evaporate the eluate on a flash-evaporator under vacuum at 50°C to about 5 ml. Transfer quantitatively the residue with 2 x 5 ml. hexane to 50 ml. pear-shaped flask and evaporate the hexane on a flash eváporator again.

Sample Preparation:

Add 1.0 ml. hexane to the pear shaped flask and dissolve the residue.

Standard Preparation:

Dissolve 25 mg. of standard SBP-1382 accurately weighed in a 250 ml. volumetric flask and dilute to volume with hexane. Dilute 5.0 ml. aliquot to 100 ml. with hexane. Final concentration 0.005 mg/ml.

GLC Analysis:

Equipment and Conditions:

Instrument:

Hewlett Packard Gas Chromatograph 5840A

Column:

- a) 6 ft. x 2 mm. I.D. 1/4" O.D. glass packed with 3% OV-1 on 100-120 mesh Gas Chrom Q
- b) 6 ft. x 2 mm. I.D. 1/4" O.D. glass packed with 3% OV-210 on 100-120 Mesh Gas Chrom Q
- c) "Mixed Column" - 50:50 mixture of 3% OV-17 on 80-100 mesh Varaport 30 and 3% OV-1 on 80-120 mesh Gas Chrom Q.

Temperature:

Injector: 250°C

Detector (F.I.D.) 250°C

NOT TO BE COPIED

Procedure No. REVISIED 1047 Copy No. 1 Page No. 5 of 7 Pages

Title GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SBP-1382) RESIDUES
IN FOOD PRODUCTS - REVISED

Date 9/9/82

Signed D. Fulmanee RP

CONFIDENTIAL

Column Oven Temperature:

- I. Primary Column: 3% OV-1 180-250°C at the rate of 5°C per minute and held at 250°C for 15 min.
- II. First Alternate Column: 3% OV-210 100-250°C at the rate of 10°C per minute and held at 250°C for 15 min
- III. Second Alternate Column: "Mixed Column" - 50:50 mixture of 3% OV-17 on 80-100 mesh Varaport 30 and 3% OV-1 on 100-120 mesh Gas Chrom Q.
(A) 180-250°C at the rate of 5°C per minute and held at 250°C for 15 min.
(B) Isothermally at 230°C.

Carrier Gas:

Helium at 30ml per minute
Hydrogen pressure - 20 p.s.i.g.
Air pressure - 40 p.s.i.g.

Attenuation:

2 (4)

Determination:

Condition the column before use by injecting the standard several times. Inject the sample and check for interfering peaks. If interfering peaks are encountered use the alternate column (3% OV-210). Inject 5.0 ul of the standard and sample preparation until sufficient chromatograms of each are obtained.

NOT TO BE COPIED

REVISID
Procedure No. 1047 Copy No. 1 Page No. 6 of 7 Pages

Title GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SBP-1382) RESIDUES

IN FOOD PRODUCTS - REVISED

Date 9/9/82

Signed R. J. Furmanec *RF*

~~CONFIDENTIAL~~

Calculation:

Calculation are based on equal volume injection using peak areas. From the chromatograms of the standard and sample preparation calculate ppm of SBP-1382 in the sample by the following formula:

$$\frac{A(s)}{A(std)} \times \frac{W(std)}{W(s)} \times 10^6 = \text{ppm of SBP-1382}$$

Where:

$A(s)$ = peak area of SBP-1382 in the sample

$A(std)$ = peak area of SBP-1382 in the standard

$W(s)$ = peak of the sample in mgs.

$W(std)$ = weight of the standard in mgs./ml. in final dilution.

REVISED

NOT TO BE COPIED

Procedure No. 1047 Copy No. 1 Page No. 7 of 7 Pages

Title GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SBP-1382) RESIDUES
IN FOOD PRODUCTS - REVISED

Date 9/9/82

Signed R. J. Furmanee 

CONFIDENTIAL

NOTE #1 (MILK)

Weigh 25 gm. of milk sample in 125 ml. Erlenmeyer flask. Add 50 ml. 5% acetone in acetonitrile and shake the mixture vigorously for 2 min. Decant the organic phase through Whatman #1 filter paper into a 250 ml. round bottom flask. Repeat three additional extractions with 50 ml. 5% acetone in acetonitrile each time shaking for 2 min. and filtering through the same Whatman #1 filter paper. Strip about 25 ml. of the organic phase to remove most of the acetone on flash-evaporator transfer quantitatively to 2 liter separatory funnel and proceed with acetonitrile petroleum ether partitioning as described on page 3.

NOTE #2 (SUGAR)

Dissolve 25 gm. of sugar in 500 ml. distilled water in 2 l. separatory funnel and proceed with acetonitrile petroleum ether partitioning as described on page 3.

PENICK

CORPORATION

1050 Wall Street West, Lyndhurst, New Jersey 07071 • Telephone: (201) 935-6600 • Telex: 133525

December 21, 1982

Mr. Franklin D. R. Gee (PM#17)
Insecticide and Rodenticide Branch
Registration Division (TS-767)
Environmental Protection Agency
401 M Street S.W.
Washington, D.C. 20460

Subject: Analytical Method
Penick's Procedure No.
1047, Revised


Dear Mr. Gee:

Authorization is hereby given to US EPA for trial of the subject method and for publication in the Pesticide Analytical Method Book.

Please inform Ms. Martha Bradley, Chemical Branch (703-557-7377).

Should you need further assistance, call us.

Sincerely,



Emily F. Hedal
Supervisor Regulatory Affairs
Pesticides Technical Support Group

EFH:PH

cc: M. deVries
S. Kutch

NOT TO BE COPIED

Procedure No. 1047 Copy No. 1 Page No. 1 of 7 Pages

Title GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SBP-1322)
RESIDUES IN FOOD PRODUCTS

Date 11/14/80

Signed *[Signature]* RP

SUMMARY:

A general analytical method for determination of Resmethrin residues in diverse types of food products is described. The procedure has been successfully applied to the following substances; sliced american cheese, sliced ham, cream pie, white flour, polished rice, white bread, popcorn kernels, saltine crackers, sliced roast beef, chicken meat, skin, fat, hearts, liver, gizzard and eggs.

Whole milk and sugar require a different sample extraction, details are in notes 1 and 2. Recoveries of Resmethrin added to these substances at the 1 ppm level range from 85 to 101% with a detection limit of 0.1 ppm.

Sample preparation involves extraction in a blender with acetonitrile, partitioning between acetonitrile and petroleum ether and clean-up on Florisil^R and Alumin AR columns.

Two gas chromatographic columns are mandated. The 3% OV-1 column is the primary column, if interfering peaks are encountered chromatography on 3% OV-210 will separate the interferences or vice-versa.

**PENICK
CONFIDENTIAL**

PENICK CONFIDENTIAL

NOT TO BE COPIED

Procedure No. 1047 Copy No. 1 Page No. 2 of 7 PagesTitle GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SBP-1322)RESIDUES IN FOOD PRODUCTSDate 11/14/80Signed R. Furman RPEquipment:

Hewlett Packard Gas Chromatograph 5840A
 Automatic sampler 7671A Hewlett-Packard
 GC Glass Column 6 ft. x 2 mm I.D. 1/4" O.D.
 Blender (Waring)
 Flash-evaporator (Searle)
 Separatory Funnel 2 liter
 Round Bottom Flask, 500 ml. 250 ml
 Pear-shaped Flask, 50 ml
 Clean-up column glass, 450 mm long x 15 mm I.D.
 Pipet, 1.0 ml., 5.0 ml., 25.0 ml.
 Buchner Funnel 70 mm Dia.
 Suction Flask, 1000 ml.
 Whatman #1 filter paper 70 mm.

Materials:

Acetonitrile, HPLC grade
 Hexane, HPLC grade
 Petroleum Ether, Pesticide grade
 Methylene Chloride, Pesticide grade
 Ether Anhydrous, Reagent A.C.S.
 Acetone, HPLC grade
 Sodium Sulfate, Anhydrous, Reagent grade
 Sodium Chloride, Reagent grade
 3% OV-1 on 100-120 Mesh on Gas Chrom Q (Applied Science Lab., Inc)
 3% OV-210 on 100-120 Mesh on Gas Chrom Q (Applied Science Lab., Inc)
 Florasil 60-100 PR Mesh (Floridin Company)
 Alumin AR CC-10 100-200 Mesh Activity Grade 4 (Mallinckrodt
 Chemical Works)
 Resmethrin Standard (Penick Corporation)

PENICK CONFIDENTIAL

NOT TO BE COPIED

Procedure No. 1047 Copy No. 1 Page No. 3 of 7 PagesTitle GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (S2P-1322)
RESIDUES IN FOOD PRODUCTSDate 11/14/80Signed R. J. Fulmer RPProcedure:

Weigh accurately 25.0 gm. of a representative food sample into blender (If milk or sugar see notes 1 and 2). Add 100 ml of acetonitrile and activate the blender for 5 minutes. Decant the solvent through Whatman #1 filter paper in a Buchner funnel using gentle suction. Repeat the extraction in the blender two additional times using 100 ml and 50 ml portions of acetonitrile and filter each time through the same Whatman #1 filter paper. Transfer combined extracts to a 2 liter seapratory funnel.

Acetonitrile - Petroleum Ether Partition

Add 100 ml petroleum ether to the acetonitrile filtrate and shake well. Then add 75 ml saturated salt solution, 700 ml of distilled water and shake again. Transfer the petroleum ether layer to 250 ml round bottom flask. Extract two additional times with 100 ml and 50 ml portions of petroleum ether. Each time remove petroleum ether in the same round bottom flask on a flash evaporator under vacuum at 50°C to an oily residue. Quantitatively transfer the oily residue to 25 ml volumetric flask with small portions of hexane and dilute to volume.

PENICK CONFIDENTIAL

NOT TO BE COPIED

Procedure No. 1047 Copy No. 1 Page No. 4 of 7 PagesTitle GLC. PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SBP-1322)
RESIDUES IN FOOD PRODUCTSDate 11/14/80Signed R. J. Furman RPColumn Clean-upFlorisil^R Column:

Pack a chromatographic column (450 mm 1 x 15 mm I.D.) with 10 gm. of Florisil^R (activated at 600°C). Add a 10 mm layer of anhydrous sodium sulfate, to the top of the column, prewet column with hexane and drain to the top of the bed. Transfer 5.0 ml sample aliquot to the top of the column and drain into the bed of the column.

Elute the column with 100 ml of solvent mixture 50% ether - 50% hexane, and collect the total eluate. Concentrate the eluate on flash-evaporator under vacuum at 50°C to an oily residue.

Alumin AR Column:

Pack a chromatographic column (450 mm 1 x 15 mm I.D.) with 20 gm of Alumin AR.

Add a 10 mm layer of anhydrous sodium sulfate to the top of the column. Prewet column with hexane and drain to the top of the bed. Transfer quantitatively the residue with three 5 ml portions of hexane to the top of the column and drain each time into the bed of the column. Elute the column with 100 ml hexane and discard the eluate. Continue eluation with 100 ml of solvent mixture 25 ml methylene chloride and 75 ml hexane and collect the eluate in 250 ml. round bottom flask. Evaporate the eluate on a flash-evaporator under vacuum at 50°C to about 5 ml. Transfer quantitatively the residue with 2 x 5 ml hexane to 50 ml pear-shaped flask and evaporate the hexane on a flash evaporator again.

PENICK CONFIDENTIAL

Procedure No. 1047 Copy No. 1 Page No. 5 of 7 Pages

Title GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SBP-1382)
RESIDUES IN FOOD PRODUCTS

Date 11/14/80

Signed G. Furman RP

Sample Preparation:

Add 1.0 ml hexane to the pear shaped flask and dissolve the residue.

PENICK CONFIDENTIAL

Standard Preparation:

Dissolve 25 mg of standard SBP-1382, accurately weighed in a 250 ml volumetric flask and dilute to volume with hexane. Dilute 5.0 ml aliquot to 100 ml with hexane. Final concentration 0.005 mg/ml.

GLC Analysis:

Equipment and Conditions:

Instrument: Hewlett Packard Gas Chromatograph 5840A

- Column:
- a) 6 ft. x 2 mm I.D. 1/4" O.D. glass packed with 3% OV-1 on 100-120 mesh Gas Chrom Q
 - b) 6 ft. x 2 mm. I.D. 1/4" O.D. glass packed with 3% OV-210 on 100-120 Mesh Gas Chrom Q

Temperature: Injector 250°C
Detector (F.I.D.) 250°C

Column Oven Temperature:

Primary Column: 3% OV-1 180-250°C at the rate of 5°C per minute and held at 250°C for 15 min.

Alternate Column: 3% OV-210 100-250°C at the rate of 10°C per minute and held at 250°C for 15 min.

Carrier Gas: Helium at 30 ml per min.
Hydrogen pressure - 20 p.s.i.g.
Air pressure - 40 p.s.i.g.

Attenuation: 2 (3)

Procedure No. 1047 Copy No. 1 Page No. 6 of 7 PagesTitle GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SBP-1382)RESIDUES IN FOOD PRODUCTSDate 11/14/80Signed R. J. [Signature] RPDetermination:

Condition the column before use by injecting the standard several times. Inject the sample and check for interfering peaks. If interfering peaks are encountered use the alternate column (3% OV-210). Inject 5.0 ul of the standard and sample preparation until sufficient chromatograms of each are obtained.

Calculation:

Calculations are based on equal volume injection using peak areas. From the chromatograms of the standard and sample preparation calculate ppm of SBP-1382 in the sample by the following formula:

$$\frac{A(s)}{A(std)} \times \frac{W(std)}{W(s)} \times 10^6 = \text{ppm of SBP-1382}$$

Where:

- $A(s)$ = peak area of SBP-1382 in the sample
 $A(std)$ = peak area of SBP-1382 in the standard
 $W(s)$ = weight of the sample in mgs
 $W(std)$ = weight of standard in mgs/ml in final dilution

PENICK CONFIDENTIAL


Procedure No. 1047 Copy No. 1 Page No. 7 of 7 PagesTitle GLC PROCEDURE FOR THE DETERMINATION OF RESMETHRIN (SPP-1322)
RESIDUES IN FOOD PRODUCTSDate 11/14/80Signed H. Furmanee R.P.NOTE #1 (MILK)

Weigh 25 gm of milk sample in 125 ml Erlenmeyer flask. Add 50 ml 5% acetone in acetonitrile and shake the mixture vigorously for 2 min. Decant the organic phase through Whatman #1 filter paper into a 250 ml round bottom flask. Repeat three additional extractions with 50 ml 5% acetone in acetonitrile each time shaking for 2 min. and filtering through the same Whatman #1 filter paper. Strip about 25 ml of the organic phase to remove most of the acetone on flash-evaporator, transfer quantitatively to 2 liter separatory funnel and procede with acetonitrile petroleum ether partitioning as described on page 3.

NOTE #2 (SUGAR)


Dissolve 25 gm of sugar in 500 ml distilled water in 2 l. separatory funnel and procede with acetonitrile petroleum ether partitioning as described on page 3.

PENICK CONFIDENTIAL



3.3
3.70

Retention Time	Area	Height	Width
15.57	1000	1000	1000
16.58	1000	1000	1000
17.59	1000	1000	1000




DATE: 11/11/77
 TIME: 11:00 AM
 OPERATOR:

NO. 111111

NO. 111111

Retention Time	Area	Height	Width
6.86	1000	1000	1000
6.98	1000	1000	1000
8.32	1000	1000	1000
8.38	1000	1000	1000
10.33	1000	1000	1000
11.36	1000	1000	1000

TC 13(1) X




SAMPLE NAME Sliced white Bread
 SAMPLE NO. Treated (or) wrapped
 INSTRUMENT HP 5840A GC
 COLUMN 37, OV-1 I

CONCENTRATION

HP 5840A GC
 HP 5840A GC
 HP 5840A GC


HP 5840A GC



HP 5840A GC
 HP 5840A GC
 HP 5840A GC

START


0.25



TO-M (2) 11

START


0.25



SOURCE: Milk
 SAMPLE NO: Treated (1K) open
 INSTRUMENT: HP 5840A GC
 COLUMN: 3700V-210 II

5 HP RUN # 4
 BOTTLE 24
 AREA %

RT	AREA	AREA %
2.35	830600000	99.642
4.14	2206	0.000
4.94	521	0.000
6.19	799	0.000
6.35	131	0.000
7.31	664	0.000
7.50	873	0.000
7.66	363	0.000
7.87	127	0.000
8.30	324	0.000
8.53	1138	0.000
8.82	506	0.000
9.10	795	0.000
9.27	1194	0.000
9.63	2000	0.000
9.76	1974	0.000
9.92	1876	0.000
11.06	1910	0.000
11.32	335	0.000
11.77	3029	0.000



STATION NAME Midi
 SAMPLE NO. Unstrat'd wrapped
 INSTRUMENT HP 5840A GC
 COLUMN 37.00-1 I

CONFIDENTIAL

A 10 MINUTE 1.0000 BT &
 10.00 15.00 20.00 25.00 30.00
 10.00 15.00 20.00 25.00 30.00

MIN	AREA	CONC
10.00	1.000000	0.0000
15.00	1.000000	0.0000
17.00	1.000000	0.0000
20.00	1.000000	0.0000
22.00	1.000000	0.0000
25.00	1.000000	0.0000

10.00 15.00 20.00 25.00 30.00
 10.00 15.00 20.00 25.00 30.00

SD *[Signature]*

