

US EPA ARCHIVE DOCUMENT

MOBAY CHEMICAL CORPORATION
AGRICULTURAL DIVISION

Research and Development Department

TITLE: A Method for the Determination of ¹⁴C-GUTHION and GUTHION Oxygen Analog in Bovine Tissues and Milk Utilizing Gas Chromatography and High Pressure Liquid Chromatography


AUTHORS: J. P. Wargo, R. J. Pollock¹ and R. R. Gronberg

ABSTRACT: A method is described for the analysis of GUTHION and GUTHION oxygen analog in bovine tissues and milk. The method for tissues and milk involves an initial extraction using acetone and dichloromethane (DCM) followed by an acetonitrile:hexane partition clean-up. The initial extraction for fat utilized hexane and acetonitrile followed by a partition clean-up. An additional solvent partition clean-up followed by a silica gel column clean-up and another acetonitrile:hexane partition clean-up was performed prior to gas liquid chromatographic analysis for GUTHION and high pressure liquid chromatographic analysis for GUTHION oxygen analog. The sensitivity of the method is 0.01 ppm for tissues and 0.001 ppm for milk.

DATE: April 10, 1978

NOTEBOOK
REFERENCE: 78-R-183

APPROVED BY


T. Bill Waggover

¹Analytical Development Corporation
Monument, Colorado

CONFIDENTIAL

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

GUTHION is a Reg. TM of the Parent Company of Farbenfabriken Bayer GmbH, Leverkusen.

66439

RESEARCH REPORT

To: Chemagro Agricultural Division
Mobay Chemical Corporation
P. O. Box 4913 - Hawthorn Road
Kansas City, Missouri 64120

ADC Project: #378

Date: April 10, 1978

ANALYTICAL METHOD FOR THE DETERMINATION OF
RESIDUES OF GUTHION AND GUTHION OXYGEN ANALOG IN
BOVINE TISSUES (LIVER, MUSCLE, FAT, KIDNEY) AND MILK

I. INTRODUCTION

This report contains an analytical method to determine residues of Guthion and Guthion oxygen analog in bovine tissues and milk. The sensitivity of detection is 0.01 ppm in tissue and 0.001 ppm in milk. The parent Guthion is quantitated on a gas chromatograph equipped with a flame photometric detector, while the oxygen analog is quantitated on an HPLC at a wavelength of 280 nm.

Recovery data and chromatograms obtained from the validation analyses of each tissue and milk are contained in this report.* The data were obtained by analyzing the samples according to the method contained in this report.

CONFIDENTIAL

II. MATERIALS AND METHOD

A. Apparatus Required

1. Assorted laboratory glassware.
2. Waring (explosion proof) Laboratory Blendor.
3. Gas Chromatograph, Tracor 222, equipped with a flame photometric detector in the phosphorous mode.
4. High Pressure Liquid Chromatograph.
5. Rotary vacuum evaporator.

B. Reagents Required

1. Acetone - Nanograde[®] - Registered trademark of Mallinckrodt Chemical Works, St. Louis, Missouri.
2. Acetonitrile - Nanograde[®].
3. Benzene - Nanograde[®].

*See Chemagro Report No. 66440.

PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

II. MATERIALS AND METHODB. Reagents Required (cont.)

4. Dichloromethane - Nanograde[®].
5. Ethyl Acetate - Nanograde[®].
6. Filter paper, Whatman No. 31.
7. Filter paper, Whatman 2V fluted.
8. Glass wool.
9. Hexane - Nanograde[®].
10. Methanol - Nanograde[®].
11. Silica Gel - 10% water; heat the silica gel (80/100 mesh) 24 hours at 120°C; deactivate 10% with water (90g silica gel and 10 ml distilled water); mix well and allow the deactivated silica gel to stand 24 hours before use.
12. Sodium Chloride (10% w/v in distilled water).
13. Sodium Sulfate, - granular anhydrous - benzene wash before use.
14. Super Cel, Hyflow.

CONFIDENTIAL

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATIONC. Special Reagents1. PS Guthion Spiking Solutions

- a. Solution "A" - Dissolve 100.0 mg of recrystallized PS Guthion in benzene and dilute to volume in a 50 ml volumetric flask. This solution contains 2 mg PS Guthion per ml.
- b. Solution "B" - Pipet 5.0 ml of Solution "A" into a 250 ml volumetric flask and dilute to volume with benzene. This solution contains 40 µg of PS Guthion per ml.
- c. Solution "C" - Pipet 1.0 ml of Solution "A" into a 200 ml volumetric flask and dilute to volume with benzene. This solution contains 10 µg of PS Guthion per ml.

II. MATERIALS AND METHOD

PROPERTY OF

C. Special ReagentsCHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION1. PS Guthion Spiking Solutions (cont.)

d. Solution "D" - Pipet 5.0 ml of Solution "B" into a 100 ml volumetric flask and dilute to volume with benzene. This solution contains 2 μ g of PS Guthion per ml.

2. PO Guthion Spiking Solutions

a. Solution "A" - Dissolve 95.0 mg of recrystallized PO Guthion in benzene and dilute to volume in a 50 ml volumetric flask. This solution contains 2 mg PO Guthion equivalents per ml.

b. Solution "B" - Pipet 5.0 ml of Solution "A" into a 250 ml volumetric flask and dilute to volume with benzene. This solution contains 40 μ g of PO Guthion per ml.

c. Solution "C" - Pipet 1.0 ml of Solution "A" into a 200 ml volumetric flask and dilute to volume with benzene. This solution contains 10 μ g of PO Guthion per ml.

d. Solution "D" - Pipet 5.0 ml of Solution "B" into a 100 ml volumetric flask and dilute to volume with benzene. This solution contains 2 μ g of PO Guthion per ml.

3. PS Guthion Standard Solutions (For GC Analysis)

a. Solution #1 - Dissolve 10.0 mg of recrystallized Guthion in benzene and dilute to volume in a 500 ml volumetric flask. This solution contains 20 μ g of PS Guthion per ml. Store in a freezer.

b. Solution #2 - Pipet 5.0 ml of Solution #1 into a 50 ml volumetric flask and dilute to volume with ethyl acetate. This solution contains 2.0 μ g of PS Guthion per ml. (The stability of Guthion of ethyl acetate for long duration is not known.)

CONFIDENTIAL

II. MATERIALS AND METHOD

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATIONC. Special Reagents (cont.)4. PO Guthion (Equivalent) Standard Solution (For LC Analysis)

- a. Solution #1 - Dissolve 9.50 mg of recrystallized PO Guthion in 50% acetonitrile/dichloromethane and dilute to volume in a 500 ml volumetric flask. This solution contains 20 µg PO Guthion equivalents per ml. Store in a freezer.
- b. Solution #2 - Pipet 5.0 ml of Solution #1 into a 50 ml volumetric flask and dilute to volume with 50% acetonitrile/dichloromethane. This solution contains 2.0 µg of PO Guthion per ml.

D. Detailed Procedure

1. Initial Extraction (except for fat; see Procedure Notes, pages 9 & 10)
 - a. Weigh 200 grams of milk or 100 grams of tissue into a blender jar and add 5 grams Hyflow Super Cel (10 grams for milk).
 - b. Add 200 ml of acetone and blend 2 minutes (400 ml for milk).
 - c. Filter under vacuum through Whatman No. 31 filter paper.
 - d. Return the filter cake and paper to the blender jar, add 200 ml of dichloromethane (DCM) and blend for 5 minutes.
 - e. Filter under vacuum through Whatman No. 31 filter paper and wash the filter cake with 100 ml DCM.
 - f. Transfer the combined filtrates to a 1,000 ml separatory funnel. Rinse the suction flask with 100 ml DCM and add to the combined filtrates.
 - g. Shake the separatory funnel vigorously for 30 seconds. Allow the layers to separate and drain lower organic phase through Whatman 2V fluted filter paper into a 1,000 ml flat bottom (FB) flask.

CONFIDENTIAL

II. MATERIALS AND METHOD

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATIOND. Detailed Procedure1. Initial Extraction (cont.)

h. Evaporate the extract just to dryness (oil).

2. Acetonitrile Hexane Partition

a. Dissolve the residue in 200 ml of hexane pre-saturated with acetonitrile (ACN), and transfer to a 500 ml separatory funnel.

b. Extract the hexane solution with 100 ml of ACN presaturated with hexane.

c. Drain lower ACN phase into a second 500 ml separatory funnel containing 200 ml hexane pre-saturated with ACN. Shake for one minute.

d. Draw off the lower ACN phase from the second separatory funnel into a 500 ml ~~RB~~ flask.e. Repeat Steps b through d twice, with 50 ml portions of ACN presaturated with hexane. ^{RB}

f. Evaporate the combined ACN fractions just to dryness.

3. Dichloromethane/Aqueous Partition

a. Dissolve the residue in 25 ml methanol and add 100 ml of 10% NaCl in distilled water. Transfer to a 250 ml separatory funnel.

b. Extract with 3 x 50 ml portions of DCM.

c. Drain lower layers into a 250 ml ~~RB~~ flask. ^{RB}

d. Evaporate the combined DCM fractions just to dryness.

e. Dissolve residue in 10 ml benzene.

66439

Research Report (cont.) - Chemagro - (ADC#378) - April 10, 1978

- 6 -

CONFIDENTIAL

II. MATERIALS AND METHOD

D. Detailed Procedure (cont.)

PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

4. Chromatographic Column Cleanup

- a. Tamp a plug of glass wool into the bottom of a 46 cm x 2 cm bell column having a 300 ml reservoir.
- b. Fill the column to the bell with benzene and begin a slow drip. Slowly pour 15 grams of 10% deactivated silica gel into the benzene. Add 10 grams of Na_2SO_4 to the top of the silica gel and allow the benzene to drain to the top of the Na_2SO_4 .
- c. Transfer the residue dissolved in benzene from above (II.D.3.e.) to the column. Start column drip rate at one per second into a beaker.
- d. Just as the top of the solvent reaches the silica gel bed, add a 15 ml benzene rinse of sample flask to column.
- e. Just as the top of the solvent reaches the silica gel bed, add 25 ml benzene.
- f. Just as the top of the solvent reaches the silica gel bed, change receivers to a 250 ml FB flask and add 150 ml DCM and collect. This fraction contains PS Guthion. (Label the flasks "A" Fraction.)
- g. Evaporate the DCM just to dryness and proceed with the ACN Hexane Partition #2 in Section D.5. below.
- h. Change receivers to a second 500 ml FB flask and add 250 ml 25% ACN in DCM and collect. This fraction contains the oxygen analog, PO Guthion. (Label the flasks "B" Fraction.)
- i. Evaporate the ACN/DCM just to dryness and proceed with ACN Hexane Partition #2 in Section D.5. below.

CONFIDENTIAL

II. MATERIALS AND METHOD

PROPERTY OF

D. Detailed Procedure (cont.)CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION5. Acetonitrile Hexane Partition #2

- a. Dissolve the residue in 100 ml of hexane presaturated with ACN and transfer to a 250 ml separatory funnel.
- b. Extract the hexane solution with 50 ml of ACN presaturated with hexane.
- c. Drain lower ACN phase into a second 250 ml separatory funnel containing 100 ml hexane presaturated with ACN. Shake.
- d. Draw off the lower ACN phase from the second separatory funnel into a 250 ml FB flask.
- e. Repeat Steps b through d twice, with 25 ml portions of ACN presaturated with hexane.
- f. Evaporate the combined ACN fractions just to dryness.
- g. Proceed to Section D.6. with the "A" Fractions (PS Guthion Fraction).
- h. Proceed to Section D.7. with the "B" Fractions (PO Guthion Fraction).

6. Gas Liquid Chromatography

Column: 24" x 2 mm i.d. borosilicate glass column packed with 5% OV-210 on 80/100 mesh Chromosorb W (HP).

Column Conditioning: Attach column to injection port fitting. Purge the column for 15 minutes, 40 cc/min helium. Cap end of column and heat to 250°C for 1 hour. Cool, remove cap, purge 2 hours, 40 cc/min helium at 249°C. Attach column to detector and proceed with analysis.

- 8 - CONFIDENTIAL

II. MATERIALS AND METHOD

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATIOND. Detailed Procedure6. Gas Liquid Chromatography (cont.)

Gases: Rotometer Settings

N₂ = 4.0O₂ = 12.0

Air = 100

H₂ = 150Temperature °C: Column: 225°C
Injection Port: 250°C

Sensitivity: Instrumental response for the standard (2 ng/μl equivalent to 0.1 ppm in tissues or 0.01 ppm in milk) must be sufficient to obtain greater than 50% full scale response.

- a. Dissolve the sample from Section D.5.g. above in ethyl acetate, 5 ml volume for tissue samples and 1 ml volume for milk. [If the sample contains a sizeable residue (volume), transfer quantitatively to a graduated tube and dilute to volume]. The sample or standard are generally injected in a volume of 5. μl. A greater volume may be used to obtain the desired sensitivity.
- b. Identify Guthion by its retention time and measure the peak produced.

Calculations:

$$\text{ppm} = \frac{\text{sample peak}}{\text{standard peak}} \times \text{Std}^1 \times \text{sample dilution}$$

¹Use 0.1 for tissues and 0.01 for milk calculations.7. High Pressure Liquid ChromatographyColumn: Partisil PXS 10/25 Silica Gel Whatman, 25 cm
x 4.6 mm i.d.

Mobile Phase: 50% ACN/50% DCM

CONFIDENTIAL

PROPERTY OF

II. MATERIALS AND METHODCHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATIOND. Detailed Procedure7. High Pressure Liquid Chromatography (cont.)

Pressure: approximately 500 psi

Sensitivity Setting: 0.50 AUFS for milk
0.25 AUFS for tissue

Wavelength: 280 nm

Pre-Column: HC Pellosil, 6 cm x 4 mm i.d.

Instrument: Glenco HPLC System I, Glenco Scientific,
Inc., Houston, Texas

- a. Dissolve the sample from Section D.5.h. above in 50% ACN/50% DCM, 5 ml volume for tissue samples and 1 ml volume for milk. Filter through a 0.5 micron Milli-Pore filter using a syringe and Swinny adaptor. Inject 100 μ l samples interspersed with 100 μ l samples of a 2 ng/ μ l standard.
- b. Identify PO Guthion (oxygen analog) by its retention time and measure the peak produced.

Calculations:

$$\text{ppm} = \frac{\text{sample peak}}{\text{standard peak}} \times \text{Std}^1 \times \text{sample dilution}$$

¹Use 0.1 for tissues and 0.01 for milk calculations.8. Procedure Notesa. Initial Extraction for Fat (Only)

- (1) Weigh 100 grams of fat into a blender jar and add 10 grams of Hyflow Super Cel.
- (2) Add 250 ml hexane and blend for 2 minutes.
- (3) Filter under vacuum through Whatman No. 31 filter paper.

CONFIDENTIAL

II. MATERIALS AND METHODD. Detailed Procedure8. Procedure Notesa. Initial Extraction for Fat (Only) (cont.)

- (4) Return the filter cake to the blender jar, add 250 ml hexane and blend for 5 minutes.
- (5) Filter under vacuum through Whatman No. 31 filter paper and wash the filter cake with 150 ml of acetonitrile.
- (6) Transfer the combined filtrates to a 1,000 ml separatory funnel. Rinse the suction flask with 50 ml acetonitrile and add to the combined filtrate.
- (7) Shake the separatory funnel vigorously for one minute and allow the layers to separate.
- (8) Proceed with Section D.2.c. and follow the procedure to the end (D.7.b.).

9. Recoveries Through Method

The following tables (I through V) (Appendix I) and enclosed chromatograms (Appendix II)* show recovery values and chromatograms obtained for the analysis of each matrix for PS Guthion and PO Guthion. Typical control, recovery and standard chromatograms are included in this report.

ANALYTICAL DEVELOPMENT CORPORATION

By: Durene V. Gilsdorf
 Durene V. Gilsdorf, M.T. (ASCP)
 Research Assistant

& Sharon V. Stone
 Sharon V. Stone
 Research Technician

66439

Research Report (cont.) - Chemagro - (ADC#378) - April 10, 1978

- 11 -

Approved by: Joseph P. Wargo, Jr.
Joseph P. Wargo, Jr.
Technical Manager, R&D

DVG:as

CONFIDENTIAL

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

66439

C O N F I D E N T I A L

**PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION**

A P P E N D I X I

66439

TABLE I
MUSCLE TISSUE (BOVINE) RECOVERIES

A. PS Guthion - Analyzed by GC/FPD

<u>Sample No.</u>	<u>Sample Fortification</u>	<u>ppm Found</u>	<u>Corrected ppm</u>	<u>% Recovery</u>
1a	Muscle Control (100 grams)	0.0101)	0.0079 (Avg)	
2a	Muscle Control	0.0057)		
3a	Muscle Control + 0.05 ppm PS Guthion	0.0344	0.0265	53.0
4a	Muscle Control + 0.05 ppm PS Guthion	0.0476	0.0397	79.4
5a	Muscle Control + 0.10 ppm PS Guthion	0.1117	0.1038	103.8
6a	Muscle Control + 0.10 ppm PS Guthion	0.0987	0.0908	90.8

B. PO Guthion - Analyzed by HPLC

<u>Sample No.</u>	<u>Sample Fortification</u>	<u>ppm Found</u>	<u>Corrected ppm</u>	<u>% Recovery</u>
1b	Muscle Control (100 grams)	0.0000	0.0000	---
2b	Muscle Control	0.0000	0.0000	---
3b	Muscle Control + 0.05 ppm PO Guthion	0.0422	0.0422	84.4
4b	Muscle Control + 0.05 ppm PO Guthion	0.0430	0.0430	86.0
5b	Muscle Control + 0.10 ppm PO Guthion	0.0734	0.0734	73.4
6b	Muscle Control + 0.10 ppm PO Guthion	0.0796	0.0796	79.6

REFERENCE: NBR #A-31

CONFIDENTIAL

PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

TABLE II
FAT TISSUE (BOVINE) RECOVERIES

A. PS Guthion - Analyzed by GC/FPD

<u>Sample No.</u>	<u>Sample Fortification</u>	<u>ppm Found</u>	<u>% Recovery</u>
1a	Fat Control (100 grams)	0.0000	---
2a	Fat Control	0.0000	---
3a	Fat Control + 0.05 ppm PS Guthion	0.0437	87.4
4a	Fat Control + 0.05 ppm PS Guthion	0.0459	91.8
5a	Fat Control + 0.1 ppm PS Guthion	0.0974	97.4
6a	Fat Control + 0.1 ppm PS Guthion	0.0980	98.0

CONFIDENTIAL

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION.
MOBAY CHEMICAL CORPORATION

B. PO Guthion - Analyzed by HPLC

<u>Sample No.</u>	<u>Sample Fortification</u>	<u>ppm Found</u>	<u>% Recovery</u>
1b	Fat Control (100 grams)	0.0000	---
2b	Fat Control	0.0000	---
3b	Fat Control + 0.05 ppm PO Guthion	0.0474	94.8
4b	Fat Control + 0.05 ppm PO Guthion	0.0520	104.0
5b	Fat Control + 0.1 ppm PO Guthion	0.0820	82.0
6b	Fat Control + 0.1 ppm PO Guthion	0.0890	89.0

REFERENCE: NBR #A-31

TABLE III
MILK (BOVINE) RECOVERIES

CONFIDENTIAL

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

A. PS Guthion - Analyzed by GC/FPD

Sample No.	Sample Fortification	ppm Found	% Recovery
1a	Milk Control (200 grams)	0.0000	---
2a	Milk Control	0.0000	---
3a	Milk Control + 0.005 ppm PS Guthion	0.0039	78.0
4a	Milk Control + 0.005 ppm PS Guthion	0.0041	82.0
5a	Milk Control + 0.01 ppm PS Guthion	0.0083	83.0
6a	Milk Control + 0.01 ppm PS Guthion	0.0074	74.0
7a	Milk Control + 0.02 ppm PS Guthion	0.0189	95.0
8a	Milk Control + 0.02 ppm PS Guthion	0.0162	81.0

B. PO Guthion - Analyzed by HPLC

Sample No.	Sample Fortification	ppm Found	Corrected ppm	% Recovery
1b	Milk Control (200 grams)	0.0024)	0.0022 (Avg)	
2b	Milk Control	0.0020)		
3b	Milk Control + 0.005 ppm PO Guthion	0.0053	0.0031	62.0
4b	Milk Control + 0.005 ppm PO Guthion	0.0057	0.0035	70.0
5b	Milk Control + 0.01 ppm PO Guthion	0.0100	0.0078	78.0
6b	Milk Control + 0.01 ppm PO Guthion	0.0106	0.0084	84.0
7b	Milk Control + 0.02 ppm PO Guthion	0.0192	0.0170	85.0
8b	Milk Control + 0.02 ppm PO Guthion	0.0182	0.0160	80.0

REFERENCE: NBR #A-31

CONFIDENTIAL

66439

PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION
TABLE IV

KIDNEY TISSUE (BOVINE) RECOVERIES

A. PS Guthion - Analyzed by GC/FPD

<u>Sample No.</u>	<u>Sample Fortification</u>	<u>ppm Found</u>	<u>% Recovery</u>
1a	Kidney Control (100 grams)	0.0000	---
2a	Kidney Control	0.0000	---
3a	Kidney Control + 0.05 ppm PS Guthion	0.0557	111.4
4a	Kidney Control + 0.05 ppm PS Guthion	0.0559	111.8
5a	Kidney Control + 0.1 ppm PS Guthion	0.1049	104.9
6a	Kidney Control + 0.1 ppm PS Guthion	0.0980	98.0

B. PO Guthion - Analyzed by HPLC

<u>Sample No.</u>	<u>Sample Fortification</u>	<u>ppm Found</u>	<u>% Recovery</u>
1b	Kidney Control (100 grams)	0.0000	---
2b	Kidney Control	0.0000	---
3b	Kidney Control + 0.05 ppm PO Guthion	0.0411	82.2
4b	Kidney Control + 0.05 ppm PO Guthion	0.0453	90.6
5b	Kidney Control + 0.1 ppm PO Guthion	0.0818	81.8
6b	Kidney Control + 0.1 ppm PO Guthion	0.0901	90.1

REFERENCE: NBR #A-31

66439

TABLE V
LIVER TISSUE (BOVINE) RECOVERIES

CONFIDENTIAL

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATIONA. PS Guthion - Analyzed by GC/FPD

<u>Sample No.</u>	<u>Sample Fortification</u>	<u>ppm Found</u>	<u>% Recovery</u>
1a	Liver Control (100 grams)	0.0000	---
2a	Liver Control	0.0000	---
3a	Liver Control + 0.05 ppm PS Guthion	0.0421	84.2
4a	Liver Control + 0.05 ppm PS Guthion	0.0447	89.2
5a	Liver Control + 0.1 ppm PS Guthion	0.0944	94.4
6a	Liver Control + 0.1 ppm PS Guthion	0.0806	80.6

B. PO Guthion - Analyzed by HPLC

<u>Sample No.</u>	<u>Sample Fortification</u>	<u>ppm Found</u>	<u>% Recovery</u>
1b	Liver Control (100 grams)	0.0000	---
2b	Liver Control	0.0000	---
3b	Liver Control + 0.05 ppm PO Guthion	0.0450	90.0
4b	Liver Control + 0.05 ppm PO Guthion	0.0449	89.8
5b	Liver Control + 0.1 ppm PO Guthion	0.0844	84.4
6b	Liver Control + 0.1 ppm PO Guthion	0.0866	86.6

REFERENCE: NBR #A-31

66439

Sample No. 3-24-78
PS BUTYRON LIVER


CONTROL

Wt. A-31 DA #19

PEAK HEIGHT IN MM. 0.0/197.0

CONFIDENTIAL

PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION



0.0mm

66439

Sample No. 3-24-78

PS GUTHION LIVER

FORTIFIED AT 0.1 ppm PS

WGR A-31 07/19

PEAK HEIGHT IN mm

150.8 / 187.0

CONFIDENTIAL

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

Baseline

$\frac{150.8}{187.0} = 0.8069 \text{ ppm}$
80.69 ppm

150.8 mm

66439

Sample & no. _____ Date: 4-4-78

PO GUTHON - LIVEE

CONTROL

NBR - A-31 ~~CO~~ # 23


PEAK HEIGHT IN MM

0.0 / 159.0

CONFIDENTIAL

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION



66439

Sample # Date 4-4-78
PO GUTTION LIVER FORTIFIED

Pt 0.1 ppm PO

PO #31 PO #23

PEAK HEIGHT IN CM

136.0 / 157.0

-131

0.866
0.612

CONFIDENTIAL


PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

66439

Sample & To DE 3-24-78
P.S. (LITHIUM) KIDNEY
CONTROL

MBR A-31 09 #18
PEAK HEIGHT IN MM

0.0/152.3



CONFIDENTIAL

PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

66439

00

66439

Sample & No. PS GUTHON Date 3-24-78
KIDNEY


FORTIFIED AT. 0.1ppm PS

NBR A31 PG #18

PEAK HEIGHT IN MM 155.2 W. 148.0

CONFIDENTIAL

PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION



66439

Sample No. 3-22-78

PO GUTTION KIDNEY

FORTIFIED AT 0.1 ppm PO

A 31 P 9 7 18

PEAK HEIGHT IN CM

155.0 / 172.0

H6

Kidney

1000

100

10

10

1

153

6.5901

90.1%

CONFIDENTIAL

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION



66439

Sample no. 3-7-78
PS (GUTHRIE) MUSCLE
CONTROL

HR A-31 9.0 / 158.0
PEAK HEIGHT IN mm

Muscle Control

CONFIDENTIAL

PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOSAY CHEMICAL CORPORATION



$\frac{158.0}{2} = \frac{9.0}{2/5} = 0.0051 \mu\mu$
27

66439

Sample No. 3-7-78
PS - GUTTION MUSCLE
FORTIFIED @ 0.1 ppm PS

#34 A-31
PEAK HEIGHT IN MM

184.5 / 187.0

CONFIDENTIAL

PROPERTY OF
CHEMACRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

10/19 + 10/20/78

10/19

AS/5mm

$$\frac{184.5}{2} = \frac{184.5}{4} = 0.09225$$

$$0.09225 = 9.225\%$$

66439

DATE: 3-2-78
P.O. - GUTHION MUSCLE
CONTROL

NO. A-31

PEAK HEIGHT IN MM

0.0/204.0

26.5
0.0
7.0
1.0
0.0
1.0
1.0

10
20
30
40
50
60
70
80
90
100

CONFIDENTIAL

PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION


7.0

66439

Sample & tin Date 3-2-78
PO GUTHON MUSCLE
FORTIFIED @ 0.1ppm PO
TBR A-31 Chrom. No.
PEAK HEIGHT IN mm Sample/Std
164.0/206.0

CONFIDENTIAL

PROPERTY OF
CHEMASRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION



66439

3-1-52
PS - GUTHERAL - F-F-T

CONFIDENTIAL

NO. P-31
PEAK HEIGHT 1.0

0.0 / 149.0

FM
Corr. 2 (Ind. P)

CONFIDENTIAL

PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION



5/10/52

0.000

0.0000 ppm

06439

Date: 3-6-78
 P. CUTTING FAT
 FORTIFIED PT. D. 1 gm P
 NDR: P-31 Gram No. _____
 PEAK HEIGHT IN MM Sample/Stg. _____
 120.2 / 123.8

CONFIDENTIAL

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

5 ml / 5 ml

120.2 mm

123.8	120.2
4	75

32 0.0974 → 97.47

66439


Sample & No. PD-GUTHRIE Date 3-6-78
CAUTION

NBR. F-31 Chrom. No. 1011
PEAK-HEIGHT IN MIN.
0.0-158.0

#1. C. H. P. a. 100.0

CONFIDENTIAL

PROPERTY OF
CHEMACRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION



66439

Sample & No. Date 3-6-72

PC - BUTYRAN FAT


FORTIFIED @ 0.1 ppm PC

NBR - A-31 Chrom. No.

PEAK HEIGHT IN mm
Sample/Std.

123.0/150.0

MS Sample No. 861-01
0.1 ppm
1000/0



CONFIDENTIAL

PROPERTY OF
CHEMACRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

66439

PS GUTTION - 3-16-78
MILK

CONTROL

A 31 pg. # 13

PEAK HEIGHT 10.0mm

0.0/169.0

CONFIDENTIAL

PROPERTY OF

CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

66439

PS GUTHRIE 3-16-78
MILK
FORTIFIED @ 0.01ppm PS
A-31 P# 13

156.0 / 187.0

CONFIDENTIAL

PROPERTY OF

CHEMAGRO-AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

66439

Sample & No. Date 3-10-78

PO-GUTHON MILK


CONTROL

Vol. A 31 PRSE 13

Area (sq in) Sample/Field

18.0 / 75.0

#1 Milk Contact road/pond



CONFIDENTIAL

PROPERTY OF
CHEMAGRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

66439

3-7-78
PS GUTHON STANDARD

2ug/ml

A-31
PEAK HEIGHT 10 mm

187.0

CONFIDENTIAL

PROPERTY OF
CHICAGO AGRICULTURAL DIVISION
MIDWAY CHEMICAL CORPORATION

12/1/78

RS

66439

Sample No. PO GUTHION Date 3-2-78
STANDARD

2.19 ml

CON. NO. A-31

PEAK HEIGHT IN mm

206.0

PO STD GUTHION 2.19 ml

206

CONFIDENTIAL

PROPERTY OF
CHEMACRO AGRICULTURAL DIVISION
MOBAY CHEMICAL CORPORATION

200

↑

40

