

US EPA ARCHIVE DOCUMENT

TITLE: "The Determination of Sulfamate Residues"

AUTHOR: Harlan L. Pease

MAY 8 1964

ABSTRACT


An analytical method for determining residues of ammonium sulfamate in apples and pears is described. The method is based on the sulfamate-nitrous acid reaction, liberating sulfate which is reduced to H_2S and measured spectrophotometrically after treating with zinc, p-aminodimethylaniline and ferric chloride to form methyleneblue. Blanks of 0.1 ppm were obtained on untreated apples and pears and recoveries of about 85% were realized when known quantities of ammonium sulfamate were added to untreated fruit.

THE DETERMINATION OF SULFAMATE RESIDUES

By: Harlan L. Pease

INTRODUCTION

A method has been developed which could have wide applicability for the determination of sulfamate residues. It was developed specifically for determining residues of ammonium sulfamate in fruit. The method has a sensitivity of 10 μ g of ammonium sulfamate and shows a blank of 0.1 ppm on untreated apple and pear samples. The method is based on the ability to measure microgram quantities of sulfate resulting from the well-known reaction of sulfamic acid with sodium nitrite in an acid medium.



Before effecting this reaction it is necessary to, first, separate the sulfamate residues from the fruit tissue by extraction with water; second, remove interfering colors from this extract with activated carbon; and, third, separate sulfates from the water extract by alumina chromatographic techniques (6). The sulfamate is then reacted with sodium nitrite to form sulfate, as shown above, and this sulfate reduced to hydrogen sulfide using a reducing acid mixture (1). The sulfide is distilled into alkaline zinc acetate solution forming zinc sulfide which, in turn, is treated with p-aminodimethylaniline and ferric chloride to form methyleneblue (2). Methyleneblue has an absorption maximum at 665 m μ .

REAGENTS AND APPARATUS

✓ Standard Potassium Sulfate Solution (Equivalent to 200 $\mu\text{g/ml}$.

of $\text{NH}_4\text{SO}_3\text{NH}_2$) - Dissolve 0.305 grams of reagent grade K_2SO_4 in deionized water, transfer to a 1000 ml. volumetric flask, dilute to the mark with deionized water and mix well.

✓ Standard Potassium Sulfate Solution (Equivalent to 20 $\mu\text{g/ml}$.

of $\text{NH}_4\text{SO}_3\text{NH}_2$) - Pipet 10 ml. of the above solution into a 100 ml. volumetric flask, dilute to volume with deionized water and mix thoroughly.

Water (deionized). Barnstead bantam Demineralizer, Fisher Scientific Co., Pittsburgh, Pa., No. 9-933 or equivalent.

✓ Sodium Hydroxide, 12% aqueous solution (deionized H_2O).

✓ Hydrochloric Acid, reagent grade; total S as SO_4
<0.0002%; 1 N (deionized H_2O).

✓ Ammonium Hydroxide, reagent grade; total S as SO_4
<0.0002%; 1 N; 0.1 N (deionized H_2O).

✓ Sodium Nitrite, granular, reagent grade; total S as SO_4
<0.01%.

✓ Zinc Acetate, 1% aqueous solution (deionized H_2O).

✓ p-Aminodimethylaniline Sulfate, 0.1% in 5.5M hydrochloric acid (5 + 6). Eastman Organic Chemical No. 1333

✓ N,N-dimethyl-p-phenylenediamine sulfate (deionized H_2O).

✓ Ferric Chloride, 0.023M in 1.2M hydrochloric acid (1 + 9)
(deionized H_2O).

Reducing Acid Mixture, use fresh, unopened bottles of both hydriodic acid and hypophosphorous acid. Both are somewhat unstable; hydriodic acid is easily oxidized to give free iodine and hypophosphorous acid, upon standing, will slowly oxidize to orthophosphoric acid. (60)

Add 320 ml. of 47% hydriodic acid (Fisher, Certified Reagent A-134), 320 ml. of conc. hydrochloric acid and 100 ml. of 30% hypophosphorous acid (J.T. Baker, purified 0178 (50%) diluted to 30%) to a 1000 ml. Erlenmeyer flask. Insert a nitrogen line, and with a gentle stream of nitrogen bubbling through the mixture (40-50 ml./min.), heat to boiling; boil gently for 20 minutes. Cool to room temperature and transfer to an automatic pipettor, similar to Scientific Glass Company No. J-2130 with a 20 ml. capacity pipet, previously flushed with nitrogen. It is important that the volumes and the concentrations of acids used conform with the specified directions to assure complete reduction efficiency when the reducing acid is applied to sulfate samples at about 100°C.

Hypophosphorous acid, if heated too strongly (above 130°C), will undergo autooxidation and give off phosphine gas. Phosphine gas, if present in the reducing mixture, will distill into the zinc acetate flask and interfere with the methyleneblue formation. It is recommended that each new lot of acid be tested by analyzing aliquots of the

standard K_2SO_4 solution, as described under Procedure, before analyzing unknowns. If low or erratic results are obtained with the standard K_2SO_4 solution, a new lot of reducing acid mixture should be prepared. The efficiency of the reducing acid mixture should be checked periodically.

Nitrogen - Cylinder supply, with necessary regulator, pressure reducing valve and needle control valve.

✓ "Darco" G-60 - Activated carbon: Darco Corporation, New York.

✓ Alumina - Acid-washed. With stirring, digest reagent grade Aluminum Oxide (Merck) 71707 at room temperature for 1 hour with 1N HCl. Wash with water, allowing the larger particles to settle and decant off the fines until the residue settles in less than 1 minute for a fall height of 10 cm.

Indicator paper - pHydrion: Micro Essential Lab, Brooklyn, N.Y.

Chromatographic Column, Acid-washed Alumina (2 columns are required and 4 are desired for routine use) - Insert a plug of glass wool into the lower end of a 50 ml. buret and add a slurry of the acid-washed alumina until a column of about 20-25 cm. is formed. Place the tip of the buret through a one-hole rubber stopper, insert into a filtering flask and apply suction, washing down any particles that may be on the side of the buret with distilled

water. Insert a second plug of glass wool at the top of the column and firmly press it into contact with the hard surface of the oxide column. Eliminate sulfate contamination by washing with 50 ml. of 1N HCl, followed by 100 ml. of deionized water, 50 ml. of 1N NH₄OH and 50 ml. of N/10 NH₄OH. Continue to wash with an additional 50 ml. of 1N NH₄OH and 50 ml. of N/10 NH₄OH, followed with 100 ml. of deionized water. Never allow the column to be left dry. Always fill the buret with deionized water and let stand until ready for use.

Before using the column for the first time, it should be checked to determine if all sulfate contamination has been eliminated. This is done by making an additional wash of the column with 20 ml. of 1N NH₄OH and 50 ml. of N/10 NH₄OH. Determine the presence of sulfate using the total volume of effluent and the procedure described on pages 9 and 10. If sulfate is detected, as indicated by a final absorbance reading greater than 0.020, repeat washing the column with the NH₄OH solutions as often as necessary. Secondly, each column should be checked to establish the volume of NH₄OH solutions required to elute the sulfate. It has been our recent experience that two washes, each consisting of 20 ml. of 1N NH₄OH and 50 ml. of N/10 NH₄OH, is satisfactory. To determine the necessary elution volumes, add 2.5 mL of the standard K₂SO₄ solution

equivalent to 200 $\mu\text{g/ml}$. of $\text{NH}_4\text{SO}_3\text{NH}_2$ to 100 ml. of 1N HCl and pour this solution through the chromatographic column as described under Procedure, page 9. Elute the sulfate as suggested and analyze the ammonium sulfate solution as described on page 9 and 10, using a 1/5 aliquot. If necessary, additional NH_4OH may be used to completely elute all sulfate.

Blendor, Waring, equipped with standard size (1000 ml.) glass jar.

Spectrophotometer, Beckman Model B, or its equivalent with 5 cm. cell.

Centrifuge, International Size 1, Model BE 50, or its equivalent, equipped with either 250 ml. or 500 ml. capacity centrifuge bottles.

Stirring Apparatus, magnetic.

Reducing Apparatus, see Figure.

PROCEDURE

Weigh a 100 g. sample of diced fruit into a Waring Blendor jar, add 125 ml. of distilled water and blend for 3 to 4 minutes. Transfer quantitatively to either a 250 ml. or 500 ml. centrifuge bottle, centrifuge for 10 minutes at 2000 RPM and carefully decant the water extract through glass wool into a 600 ml. beaker. Rinse the blendor jar with an additional 125 ml. portion of distilled water, collecting the rinse in the centrifuge

bottle. Stopper with a cork plug, shake thoroughly for two minutes, and centrifuge as before, combining the water extract with that in the 600 ml. beaker. Repeat the extraction one more time, using another 125 ml. portion of distilled water and combine the extract with those previously collected.

Add 2 g. of "Darco" to the combined extracts and heat to 60°C. on a hot plate in about 10 minutes. Using a vacuum, filter through Whatman #1 filter paper into a flask containing 5 ml. of 12% NaOH solution. Wash the beaker with two 5 ml. portions of distilled water. The direct filtration of the "Darco" slurry is very slow. Allow overnight for filtration. (The procedure should not be interrupted for any length of time before this point is reached.) Pour the clear filtrate into a 400 ml. beaker using several small portions of distilled water as wash. The pH of the solution should be 11.5 to 12; if not, add additional 12% NaOH solution and carefully concentrate on a hot plate to a volume of about 90 ml. Check the pH of the solution with pHydrion paper during the concentration step. The pH should remain above 7; if not, add additional 12% NaOH. The concentration should be completed in 1 to 2 hours if the solution is maintained at a boil. An appreciably longer heating time could result in hydrolysis of the $\text{NH}_4\text{SO}_3\text{NH}_2$, giving low results.

Cool in an ice bath, add 10 ml. of conc. HCl, measured by graduate, two teaspoonsful of acid washed (wet) Al_2O_3 , and stir continuously on a magnetic stirrer for 1 hour. Allow to settle and

filter with vacuum through #40 filter paper. Wash the Al_2O_3 remaining in the beaker by resuspending in 5 ml. of distilled water, allowing to settle, and decanting the supernatant onto the filter. Repeat this washing operation a second time. It should be kept in mind that the HCl acid concentration must remain between 2 and 5% for complete adsorption of SO_4 on Al_2O_3 , and therefore the wash volumes should not exceed those prescribed.

Remove the final traces of sulfate using the alumina column prepared according to the directions under Reagents and Apparatus. Drain the water from the column, wash with 20 ml. of 1N HCl, and elute the sample through the column. Wash the column with 20 ml. of 1N HCl, combine with the sample eluate and transfer to a 400 ml. beaker using a minimum amount of distilled water as wash. (As all undesirable SO_4 has now been removed this is a convenient stopping point.) Before reusing this column, it must be regenerated as described below. This column should be reserved for sulfate clean-up.

To regenerate the alumina column, continue to wash with 50 ml. of 1N HCl, followed in order with 100 ml. of deionized water, 20 ml. of 1N NH_4OH , and 50 ml. of N/10 NH_4OH . Repeat the NH_4OH wash using 20 ml. of 1N solution and 50 ml. of N/10 solution. Complete the wash with 100 ml. of deionized water. The column is now ready for reuse after a final washing with 20 ml. of 1N HCl just prior to the addition of the sample. No change in the adsorption power of the column has been observed after 100 analyses (6).

Add 0.1 gram of reagent grade NaNO_2 to the solution in the 400 ml. beaker and place on a magnetic stirrer for 30 minutes, thus allowing sufficient time for the nitrous acid to react.

Pass this solution, now containing sulfate resulting from the nitrous acid-sulfamate reaction, through a second alumina column, prepared as described under Reagents and Apparatus, drained, and washed with 20 ml. of 1N HCl. Wash the beaker and stirring bar with 50 ml. of 1N HCl and add to the column, followed by 100 ml. of deionized water. Discard the eluate.

Elute the sulfate from the column by adding 20 ml. of 1N NH_4OH solution followed by 50 ml. of N/10 NH_4OH solution. Repeat using a second 20 ml. portion of 1N NH_4OH and 50 ml. of N/10 NH_4OH . Measure the volume of the combined ammonium sulfate solution and transfer $1/5$ the volume to a 200 ml. round-bottom flask equipped with a 29/26 joint and 10/30 side arm. Ammonium sulfamate in the range of 0.5 ppm to 8 ppm may be determined using the above aliquot. If maximum sensitivity is desired, transfer the entire ammonium sulfate solution to the reaction flask. Add 3 to 4 glass beads and concentrate to about 5 ml. using an electric heating mantle. Continue to concentrate to a volume of less than $1/2$ ml. using a steam bath.

Assemble the reducing apparatus as illustrated (Figure) with the sulfate solution to be analyzed already in the reaction flask and the receiving flask containing 130 ml. of 1% zinc acetate solution plus 5 ml. of 12% sodium hydroxide solution.

With the cooling water running through the condenser, add 20 ml. of the reducing acid mixture to the reaction flask via the dropping funnel. Immediately insert the nitrogen line allowing a constant flow of nitrogen to bubble gently (40-50 ml./min.) through the system. Heat the contents of the flask to boiling (5-6 min.) and continue to boil under reflux for exactly 5 minutes. Continue the nitrogen flow for an additional 5-6 minutes before disconnecting the equipment.

Without delay, remove the dropping tube from the receiving flask, washing the end of the tube with small amounts of deionized water and add, by pipet, 25 ml. of 0.1% *p*-aminodimethylaniline solution; swirl to dissolve the solids. Add, by pipet, 5 ml. of 0.023 M ferric chloride solution, mix, and allow to stand for 10 minutes. Dilute to the mark with deionized water, mix well and allow to stand for at least 20 minutes. (The solution must be kept out of direct sunlight and standing time should not exceed 2 hours.)

Prepare a reagent blank solution by adding to a 250 ml. volumetric flask, 130 ml. of 1% zinc acetate, 5 ml. of 12% sodium hydroxide, 25 ml. of 0.1% *p*-aminodimethylaniline solution, swirl to dissolve the solids, and then add 5 ml. of 0.023 M ferric chloride solution. Mix and allow to stand for 10 minutes before diluting to the mark with deionized water. Mix well and allow to stand for at least 20 more minutes, as with the sample solutions.

Determine the absorbance of the sample solution at 665 m μ , using 5 cm. cells and the reagent blank solution in the reference cell. The amount of ammonium sulfamate in the sample is then determined from a calibration curve covering ammonium sulfamate concentrations over the range of 0 to 160 μ g. The calibration curve is prepared using aliquots of a standard potassium sulfate solution, calculated as equivalent amounts of ammonium sulfamate, that have been subjected to the reducing procedure described above. Aliquots representing 20 μ g, 60 μ g, 100 μ g, and 160 μ g of NH₄SO₃NH₂ are suggested for calibration purpose. The volume of an aliquot must be reduced to less than 1/2 ml., as with the sample, before reduction.

DISCUSSION

Ammonium sulfamate hydrolyzes relatively easily in acid medium but is more stable in alkaline systems. During those steps in the procedure in which the fruit extract must be acidified (while in contact with the alumina), it is necessary to proceed without undue delay. This is no longer a problem once the sulfamate has been converted to sulfate. It is also necessary to proceed immediately with the color development procedure at the conclusion of the reduction step and after complete absorption of the liberated hydrogen sulfide. Sulfide is readily oxidized to sulfite on contact with atmospheric oxygen.

To obtain low and consistent blanks it is necessary to (a) use very pure reagents; (b) use deionized water as specified;

and (c) pass the fruit extract through an alumina column to remove final traces of sulfate prior to the conversion of sulfamate to sulfate.

Several reducing acid mixtures were considered:

HI-H₃PO₂-HCl (1), (3), (5); SnCl₂·H₃PO₄ (4); and titanium-H₃PO₄ (7), the former being chosen because of its simplicity in preparation and handling. It is very important, however, that the mixture be prepared as directed under Reagents and Apparatus if optimum reducing conditions are to exist. It is equally important that the volume of the sulfate solution be reduced to 1/2 ml. or less, prior to addition of the reducing acid solution. The times and nitrogen flow rate suggested during the reduction steps are also critical.

RESULTS

Blanks of 0.1 ppm as ammonium sulfamate, inclusive of equipment and reagents, were obtained on several varieties of apples and one of pears (Table I). Such sulfur-containing pesticides as parathion, captan, thiram and "Guthion" do not interfere. Recovery of ammonium sulfamate over the range of 0.1 ppm to 6 ppm when added to untreated apples averaged 87%. Recoveries obtained on pears over the range of 0.2 ppm to 1 ppm averaged 83% (Table II). Equally satisfactory results should be obtained on a wide variety of materials.

E. I. du Pont de Nemours & Company
Experimental Station
Industrial and Biochemicals Department
Wilmington, Delaware
May 8, 1964

TABLE I

Blanks of Untreated Fruit

Apparent Ammonium Sulfamate Level

<u>Fruit</u> <u>(Apples)</u>	<u>Sample Weight</u> <u>Analyzed</u>	<u>Apparent</u> <u>$\mu\text{g NH}_4\text{SO}_3\text{NH}_2$</u>	<u>Blank</u> <u>ppm Residue</u>
Rome	100 g.	12	0.12
Rome	100 g.	10	0.10
Rome	100 g.	14	0.14
Rome	100 g.	9	0.09
Rome	100 g.	11	0.11
Rome	100 g.	10	0.10
Stayman	100 g.	15	0.15
Stayman	100 g.	10	0.10
Stayman	100 g.	11	0.11
Stayman	100 g.	12	0.12
Stayman	100 g.	10	0.10
Red Delicious	100 g.	10	0.10
	Average	11 μg	0.11 ppm
(Pears)			
D'Anjou	100 g.	14	0.14
D'Anjou	100 g.	11	0.11
	Average	13 μg	0.13 ppm

TABLE II

Recovery of Known Amounts of Ammonium Sulfamate
Added to Untreated Fruit Samples

<u>Fruit</u>	<u>Sample Weight</u>	<u>μg NH₄SO₃NH₂ Added</u>	<u>μg NH₄SO₃NH₂* Found</u>	<u>% Recovery</u>	<u>ppm Level</u>
Apples	100 g.	10	9	90	0.1
Apples	100 g.	20	18	90	0.2
Apples	100 g.	30	27	90	0.3
Apples	100 g.	40	32	80	0.4
Apples	100 g.	50	44	88	0.5
Apples	100 g.	60	55	92	0.6
Apples	100 g.	100	92	92	1.0
Apples	100 g.	100	84	84	1.0
Apples	100 g.	200	169	85	2.0
Apples	100 g.	300	249	83	3.0
Apples	100 g.	400	354	89	4.0
Apples	100 g.	500	404	81	5.0
Apples	100 g.	600	539	<u>90</u>	6.0
			Average	-	87%
			Average Deviation	-	3.6%
Pears	100 g.	20	16	80	0.2
Pears	100 g.	40	33	83	0.4
Pears	100 g.	60	49	82	0.6
Pears	100 g.	100	87	<u>87</u>	1.0
			Average	-	83%
			Average Deviation	-	2.0%

*Corrected for average blank.

REFERENCES

1. Archer, E. E., *Analyst*, 81 (1956) 181.
2. Fogo, J. K. and Popowsky, M., *Anal. Chem.*, 21 (1949) 732.
3. Horton, A. D. and Thomason, P.E., *Anal. Chem.*, 23 (1951) 1859.
4. Kiba, T.; Takagi, T.; Yoshimura, Y.; and Kishi, I.; *Anal. Abstr.*, 3 (1956) 2708.
5. Luke, L.L., *Anal. Chem.*, 21 (1949) 1369.
6. Nydahl, F., *Anal. Chem.*, 26 (1954) 580.
7. Suzuki, S.; Harimaya, K.; Tsuji, N.; and Yamaoka, N.; *Bull. Soc. Chem., Japan*, 30 (1957) 771.

SULFATE REDUCING APPARATUS

