

US EPA ARCHIVE DOCUMENT

445370-58

PBJ

VOLUME 60 OF 65 OF SUBMISSION

CGA-245704

TITLE

ANALYTICAL METHOD FOR THE DETERMINATION OF TOTAL RESIDUES OF CGA-245704
AS CGA-210007 IN TOBACCO, LEAFY VEGETABLES, AND FRUITING VEGETABLES
BY COLUMN SWITCHING HIGH PERFORMANCE LIQUID CHROMATOGRAPHY

DATA REQUIREMENT

EPA Guideline No 860-1340

AUTHOR

C G HAYWORTH, SCIENTIST I

COMPLETION DATE

March 24, 1998

PERFORMING LABORATORY

Human Safety Department
Novartis Crop Protection, Inc
Greensboro, NC 27419

LABORATORY PROJECT IDENTIFICATION

AG-671A Supercedes AG-671
Novartis Number 491-97

SUBMITTER/SPONSOR

Novartis Crop Protection, Inc.
Post Office Box 18300
Greensboro, NC 27419-8300

VOLUME 1 OF 1 OF STUDY

PAGE 1 OF 144

STATEMENT OF NO DATA CONFIDENTIALITY CLAIM

No claim of confidentiality is made for any information contained in this study on the basis of its falling within the scope of FIFRA Section 10 (d)(1)(A), (B) or (C)

Company Novartis Crop Protection, Inc


Company Representative H Lee Hubbard, Ph.D.

Title Senior Regulatory Manager

Signature *H Lee Hubbard* Date 3/27/99

These data are the property of Novartis Crop Protection, Inc and, as such, are considered confidential for all purposes other than compliance with FIFRA Section 10


Submission of these data in compliance with FIFRA does not constitute a waiver of any right to confidentiality that may exist under any other statute or in any other country



SUBMITTER/SPONSOR Novartis Crop Protection, Inc. P O Box 38300, Greensboro, NC 27419-8300

STATEMENT CONCERNING GOOD LABORATORY PRACTICES

The document in this submission, is an analytical method and is not considered a final report because validation of the method is not included. Therefore, certification of compliance with Environmental Protection Agency's Good Laboratory Practice Standards (40 CFR Part 160, October 16, 1989) is not required.



Robert K. Williams
Agent of Submitter/Sponsor
And Manager, Residue Chemistry
Human Safety Department

3/27/98
Date

SUBMITTER/SPONSOR Novartis Crop Protection, Inc. P. O. Box 18300, Greensboro, NC 27419

TABLE OF CONTENTS OF SUBMITTAL VOLUME

<u>TITLE</u>	<u>PAGE NO</u>
Analytical Method No AG-671A Analytical Method for the Determination of Total Residues of CGA-245704 as CGA-210007 in Tobacco, Leafy Vegetables, and Fruiting Vegetables by Column Switching High Performance Liquid Chromatography (SUPERCEDES Analytical Method No AG-671)	5
Appendix I Separate Documents Accompanying This Report	82
Appendix II Separate Documents Accompanying This Submission	143
Appendix III Previously Submitted Documents with EPA MRID Numbers	144

SUBMITTER/SPONSOR Novartis Crop Protection, Inc. Post Office Box 18300, Greensboro, NC 27419

Novartis Crop Protection, Inc
Human Safety Department

Analytical Method for the Determination of Total Residues of
CGA-245704 as CGA-210007 in Tobacco, Leafy Vegetables, and
Fruiting Vegetables by Column Switching High Performance
Liquid Chromatography

Analytical Method No. 671A

Project Number 436001

Submitted By: M. C. Grunenwald

Title Chemist IV

Signature *M.C. Grunenwald*

Date 3/24/98

Approved By Robert K. Williams

Title Manager, Residue Chemistry

Signature: *[Signature]*

Date 3/24/98

TABLE OF CONTENTS

	<u>Page No.</u>
I. SUMMARY AND INTRODUCTION	6
A. SCOPE	6
B. PRINCIPLE	6
II. MATERIALS AND METHODS	7
A. APPARATUS	7
B. REAGENTS	8
C. ANALYTICAL PROCEDURE	10
1.0 Sample Preparation	10
2.0 Extraction	10
3.0 Partition Cleanup	11
4.0 Analyte Concentration	13
D. INSTRUMENTATION	14
1.0 Description and Operating Conditions	14
2.0 Calibration and Standardization	15
E. INTERFERENCES	16
F. CONFIRMATORY TECHNIQUE	16
G. TIME REQUIRED FOR ANALYSIS	17
H. MODIFICATIONS AND POTENTIAL PROBLEMS	17
I. PREPARATION OF STANDARD SOLUTIONS AND SAMPLE FORTIFICATION PROCEDURES	19
1.0 Preparation of Analytical Standard	19
2.0 Preparation of Fortification Standard	19
3.0 Sample Fortification Procedures	20

TABLE OF CONTENTS
(Continued)

	<u>Page</u> <u>No.</u>
J. METHODS OF CALCULATIONS	21
1.0 Determination of Sample Residues..	21
2.0 Determination of Procedural Recoveries .	22
3.0 Determination of Corrected Sample Residues	22
III. RESULTS AND DISCUSSION	23
A. Accuracy	23
B. Precision	25
C. Extractability	25
D. Accountability	26
IV. CONCLUSION	27
V. CERTIFICATION	28
VI. TABLES AND FIGURES	29
TABLE I. LIQUID CHROMATOGRAPHIC OPERATING CONDITIONS FOR ANALYSIS OF CGA-210007	29
TABLE II. SOURCE OF ACCURACY AND PRECISION DATA	30
TABLE III. TOBACCO MULTICHROM WORKSHEET	31
TABLE IV. LETTUCE MULTICHROM WORKSHEET	32
TABLE V. TOMATO MULTICHROM WORKSHEET	33
TABLE VI. RECOVERY DATA FROM CONTROL SUBSTRATES FORTIFIED WITH CGA-245704 AND ANALYZED BY ANALYTICAL METHOD AG-671A	34

TABLE OF CONTENTS
(Continued)

	<u>Page</u> <u>No.</u>
TABLE VII	STATISTICAL RESULTS FROM THE ANALYSIS OF FORTIFIED CONTROLS BASED ON FORTIFICATION LEVEL. 38
TABLE VIII.	EXTRACTABILITY AND ACCOUNTABILITY DATA FROM ANALYSES OF RADIOLABELED TOBACCO, LETTUCE, AND TOMATO ANALYZED BY ANALYTICAL METHOD AG-671A.... 40
TABLE IX.	COMPARISON OF VALIDATION EXTRACTABILITY AND ACCOUNTABILITY WITH QUANTITATION OF METABOLITE FRACTIONS BY 2-DIMENSIONAL THIN LAYER CHROMATOGRAPHY... 41
TABLE X.	COMPARATIVE RESULTS OF LC/MS DETECTION VERSUS HPLC/UV DETECTION IN TOBACCO, LETTUCE, AND TOMATO SUBSTRATES..... 42
TABLE XI.	OPERATING CONDITIONS FOR THE CONFIRMATION OF CGA-245704 (AS CGA-210007) BY HIGH PERFORMANCE LIQUID CHROMATOGRAPHY WITH MASS SPECTROMETRIC DETECTION..... 43
FIGURE 1	STRUCTURES AND CHEMICAL NAMES OF CGA-245704 AND CGA-210007 44
FIGURE 2.	FLOW DIAGRAM FOR ANALYTICAL METHOD AG-671A 45
FIGURE 3.	SCHEMATIC DIAGRAM OF THE HPLC COLUMN SWITCHING SYSTEM AND REPRESENTATIVE CHROMATOGRAM OF STANDARD PROFILE FROM COLUMN 1..... 46
FIGURE 4.	REPRESENTATIVE HPLC CHROMATOGRAMS OF CGA-210007 STANDARDS 48

TABLE OF CONTENTS
(Continued)

		<u>Page</u> <u>No</u>
FIGURE 5	REPRESENTATIVE CALIBRATION CURVE USING THE CGA-210007 STANDARDS FROM FIGURE 4.....	52
FIGURE 6.	REAGENT BLANKS CHROMATOGRAMS	53
FIGURE 7.	HPLC CHROMATOGRAMS FROM THE ANALYSIS OF ¹⁴ C-TOBACCO	57
FIGURE 8.	HPLC CHROMATOGRAMS FROM THE ANALYSIS OF ¹⁴ C-LETTUCE	60
FIGURE 9.	HPLC CHROMATOGRAMS FROM THE ANALYSIS OF ¹⁴ C-TOMATO	63
FIGURE 10.	REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF FLUE-CURED TOBACCO.....	66
FIGURE 11.	REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF HEAD LETTUCE.....	68
FIGURE 12.	REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF UNWASHED TOMATOES....	70
FIGURE 13.	REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF TOMATO WET POMACE.....	72
FIGURE 14	REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF TOMATO PASTE	74
VII.	REFERENCES	76

I. SUMMARY AND INTRODUCTION

A. Scope

This method is used for the determination of total residues of CGA-245704 as CGA-210007 in tobacco, leafy vegetables, and fruiting vegetables (including processed fractions). The limit of detection (LOD), as determined by the smallest standard concentration injected, is 0.35 ng for CGA-210007. The limit of quantitation (LOQ), as demonstrated by the smallest acceptable fortification level used in the study, is 0.02 ppm CGA-245704 equivalents for all substrates. The chemical names and structures for CGA-245704 and CGA-210007 are presented in Figure 1.

Analytical Method AG-671 was a re-issue of REM 172.02² and REM 172.11² to incorporate additional substrates, update formatting, lower the LOQ, include significant purification, procedural, and instrument modifications, and to report procedural and radiovalidation data. Analytical Method AG-671A is a re-issue of AG-671 to incorporate instrumental conditions and supporting data for a confirmatory procedure utilizing HPLC with mass spectrometric detection (LC/MS). Analytical Method AG-671A supercedes Analytical Method AG-671, all method references herein will be cited as Analytical Method AG-671A.

B. Principle

CGA-245704 and metabolites containing the CGA-210007 moiety are hydrolyzed to CGA-210007 by heating a representative crop sample in the presence of sodium hydroxide. Methanol is added and the mixture shaken. The extract is purified on Chem Elut and C₁₈ columns followed by an acidic liquid-liquid partition. Additional purification and final determination are performed using column switching high performance liquid chromatography (HPLC) with UV-detection or conformational HPLC analyses with mass spectrometric detection (LC/MS). A flow diagram for the method is presented in Figure 2.

II. MATERIALS AND METHODS

Below is a compilation of the materials needed to run this method. Standard laboratory equipment is not listed.

A. Apparatus

- 1.0 Lab Shaker, IKA Labortechnik KS501 Digital, IKA Works, Cat. #2526401, or equivalent
- 2.0 Lab Oven, Fisher Scientific, Model #506G, Cat. #13-246-506G, or equivalent
- 3.0 Coors Porcelain with Fixed Perforated Plate (Büchner funnel), Fisher Scientific, Cat. #10-356C, or equivalent
- 4.0 Filter Flask, Pyrex Brand with Tubulation, Fisher Scientific, Cat. #10-180E, or equivalent
- 5.0 Glass Microfibre Filters, circles 70 mm, Fisher Scientific, Cat. #1827 070, or equivalent
- 6.0 Disposable Pasture pipettes with 2-mL bulb
- 7.0 Graduated Cylinders, 250, 100, 50, 20 mL
- 8.0 Kontes Recovery Flasks with outer joint, 1000 mL (teardrop flasks), Fisher Scientific, Cat. #K608675-4724, or equivalent
- 9.0 Volumetric pipettes, 0.5-10 mL
- 10.0 Graduated pipettes, 0.5-10 mL
- 11.0 Beakers, 100 mL
- 12.0 Rotary evaporator, Büchi, Fisher Scientific, Cat. #09-548-105F, or equivalent
- 13.0 C₁₈ Mega Bond Elute SPE column, 5 g/20 cc, Varian, Cat. #1225-6023

- 14.0 Vacuum Manifold, Waters, Cat. #WAT054760, or equivalent
- 15.0 Aspirator Pump, Cole-Parmer Instrument Company, Cat. #H07049-02, or equivalent
- 16.0 Chem Elut 20-mL capacity, Varian, Cat. #12198008
- 17.0 Erlenmeyer Flasks 250, 125 mL
- 18.0 Separatory Funnels, 500, 250, 125 mL with glass stopcocks
- 19.0 Round Bottom Flask, 250 mL
- 20.0 Pyrex Brand Crucible Holders (carbon filter, tubes), Fisher Scientific, Cat. #08-261A, or equivalent
- 21.0 Sanitary Cotton, Fisher Scientific, Cat #07900, used as a filter, inserted into the neck of a crucible holder, approximately 1" x 1" x 1" in size.
- 22.0 Vortex Maxi Mix Plus, Fisher Scientific, Cat. #12-815-18, or equivalent
- 23.0 Vials for LC injection, Wheaton 2-mL glass, or equivalent
- 24.0 pH Stick, Haake-Büchler, Fisher Scientific, Cat. #13-643-50, or equivalent
- 25.0 Ultrasonic Cleaner, Cole-Parmer Instrument Company, Cat. #H08895-05, or equivalent
- 26.0 Rotary Evaporator Trap with inner vapor tube at bottom, Kontes Cat. #K570200-0124, 100 mL, 24/40

3. Reagents

- 1.0 Sodium Hydroxide (NaOH) 1N Volumetric Solution, reagent grade, J.T. Baker, Cat. #5635-02, or equivalent
- 2.0 Methanol (MeOH), HPLC grade, Fisher Scientific, Cat. #A452-4, or equivalent

- 3.0 HPLC-Grade Water either purified in house with a HYDRO™ purification system or obtained from manufacturer
- 4.0 Sodium Chloride (NaCl), Certified ACS grade, Fisher Scientific Cat #S-421-3, or equivalent
- 5.0 Saturated Sodium Chloride (Sat. NaCl), add sodium chloride to reagent-grade water until maximum solubility is obtained
- 6.0 Hydrochloric Acid (HCl) concentrated, reagent grade, Fisher Scientific Cat. #A144-212, or equivalent
- 7.0 1M HCl, carefully mix 82.5 mL of conc. HCl with reagent-grade water and bring to a 1-liter final volume
- 8.0 Phosphoric Acid (H₃PO₄) 85%, Certified ACS grade, Fisher Scientific, Cat. #A242, or equivalent
- 9.0 0.02M H₃PO₄, mix 2.31 g of H₃PO₄ with reagent-grade water and bring to a 1-liter final volume
- 10.0 Dichloromethane (DCM), HPLC grade, Fisher Scientific, Cat. #D143-4, or equivalent
- 11.0 Ethyl Acetate (EtOAc), HPLC grade, Fisher Scientific, Cat. #E-195, or equivalent
- 12.0 9:1 DCM/EtOAc (v/v), combine 3600 mL of DCM with 400 mL of EtOAc
- 13.0 Acetonitrile (ACN), HPLC grade, Fisher Scientific, Cat. #998-4, or equivalent
- 14.0 Mobile Phase: 80:20 0.02M H₃PO₄/ACN (v/v), mix 14.76 g of H₃PO₄ (85%) with reagent-grade water and bring to a volume of 6400 mL. Solution pH range should fall between 2.8-3.0. Combine the 0.02M H₃PO₄ solution with 1600 mL of ACN for a total volume of 8.0 liters. Degas with helium purge.

- 15.0 Mobile Phase: 75:25 0.02M H₃PO₄/ACN (v/v), mix 13.84 g of H₃PO₄ (85%) with reagent-grade water and bring to a volume of 6000 mL. Solution pH range should fall between 2.8-3.0. Combine the 0.02M H₃PO₄ solution with 2000 mL of ACN for a total volume of 8 liters. Degas with helium purge.
- 16.0 CGA-245704 Analytical Standard, Novartis Crop Protection, Inc., P. O. Box 18300, Greensboro, NC 27419
- 17.0 CGA-210007 Analytical Standard, Novartis Crop Protection, Inc., P. O. Box 18300, Greensboro, NC 27419

C. Analytical Procedure

1.0 Sample Preparation

Samples are received and stored frozen at -20°C. Samples are prepared under the general guidelines of the U.S. Food and Drug Administration Pesticide Analytical Manual Volume I, Section 141 and Novartis SOP 7.21 Revision 2³ (or the current revision).

2.0 Extraction

- 2.1 Weigh a 4-g subsample from a well-homogenized chopped or ground crop sample into a 250-mL Erlenmeyer flask. Fortify samples with CGA-245704 if needed (see Section II.I.3.0 for fortification details).
- 2.2 Add 36 mL of distilled water and 4 mL of 1N sodium hydroxide to the substrate (the pH should be ≥ 12 , however, it should not be necessary to check this with every sample). Stopper the flask, swirl the mixture, and allow it to shake (mechanical shaker) at high speed for approximately 10 minutes. Place the stoppered flask in a preheated oven.

(65-70°C) for 30 minutes. Swirl the flask approximately every 10 minutes.

- 2.3 Allow the mixture to cool to room temperature (this step can be expedited by use of a cold tapwater bath), and add 80 mL of methanol to the flask. Seal the flask and shake (mechanical shaker) the sample at high speed for 30 minutes. Filter under vacuum by carefully pouring the sample through a glass fiber filter paper inserted into a Buchner funnel assembly. After filtering, rinse the sampling flask with approximately 3 x 20 mL of methanol. Use each rinse to wash the filter cake; stir and compress the filter cake during this rinsing to ensure quantitative transfer.
- 2.4 Transfer the filtrate to a 1000-mL teardrop flask and remove the methanol, rotating slowly to avoid bumping, with a water aspirated rotary evaporator while maintaining the water bath temperature at approximately 30°C. Reduce the sample volume to approximately 50 mL and then increase rotation to full speed until approximately 8-10 mL of sample remains.

3.0 Partition Cleanup

- 3.1 Chem Elut Column: Add 8 mL of 1M HCl to the sample from Section II.C.2.4, swirl, and load onto a Chem Elut Column. Wait approximately 10-20 minutes for sample absorption, rinse the 1000-mL teardrop flask with 100 mL 9:1 dichloromethane/ethyl acetate, elute the column with the rinse from the flask, and collect into a 250-mL Erlenmeyer flask. (A small amount of pressure may be applied to the top of the Chem Elut column to efficiently drain the column, if necessary.) Transfer the eluant to a 500-mL

separatory funnel. Rinse the 250-mL Erlenmeyer flask with 10 mL of 9:1 dichloromethane/ethyl acetate and add the rinse to the 500-mL separatory funnel.

- 3.2 Alkaline Partition: Add 20 mL of 1N sodium hydroxide to the column eluant, stopper the separatory funnel, and shake **gently** to minimize emulsions (use frequent venting) Allow 10-15 minutes for the phases to separate **completely**. (See Section II.H.1.0, Modifications and Potential Problems, if emulsions develop from difficult substrates, i.e. cured tobacco.) Drain off the lower, organic layer and discard appropriately
- 3.3 C₁₈ Solid Phase Extraction. Attach the C₁₈ Mega Bond Elut column onto a vacuum manifold. Precondition the cartridge by rinsing with 15 mL of methanol, followed by 20 mL of water, and discard appropriately. **Do not let column bed go dry at anytime.** Load the upper, aqueous fraction from the alkaline partition (Section II.C.3.2) and elute quickly through the cartridge under maximum vacuum collecting the filtrate in a 100-mL beaker. Rinse the cartridge with 50 mL of water and collect the rinse in the same 100-mL beaker. (See Section II.H.2.0, Modifications and Potential Problems.)
- 3.4 Acidic Partition: Transfer all of the filtrate from the beaker back into the 500-mL separatory funnel used for the alkaline partitioning (Section II.C.3.2). Rinse the beaker with 40 mL of saturated sodium chloride solution, and then add the rinse to the separatory funnel. Add 4 mL of phosphoric acid (85%) to the funnel, and shake well (use frequent venting) to mix thoroughly.

- 3.4.1 Partition the solution with 50 mL of ethyl acetate and wait approximately 10 minutes to allow the phases to separate completely. Drain the lower, aqueous layer into a 125-mL Erlenmeyer flask. Drain the upper, organic layer and any emulsion into a 125-mL separatory funnel. Transfer the aqueous layer back to the 500-mL separatory funnel. Rinse the 125-mL Erlenmeyer flask with 50 mL of ethyl acetate and add the rinse to the 500-mL separatory funnel containing the aqueous layer. Partition and allow the layers to separate, and then drain off the lower, aqueous layer and discard appropriately. Combine the upper, organic layer and any emulsion or residual aqueous material with the organic fraction in the 125-mL separatory funnel. (Drain and discard appropriately any aqueous layer that may have settled from the combined organic layers upon sitting.)
- 3.4.2 Insert a small amount of cotton (1" x 1" x 1") into the stem of a crucible holder (carbon filter tube). Filter the organic layer through the cotton and collect the filtrate in a 250-mL round bottom flask. Rinse the crucible holder with 10 mL of ethyl acetate draining it through the cotton filter

4.0 Analyte Concentration

- 4.1 Evaporate the contents of the round bottom flask from Section II.C.3.4.2

to complete dryness using a rotary evaporator (equipped with a trap adapter), rotating slowly at room temperature (~25°C). Add the appropriate final volume of 0.02 M H₃PO₄ to the sample and vortex or sonicate with an ultrasonic cleaner to ensure dissolution of the analyte.

- 4.2 Transfer an aliquot of the final fraction to an injection vial for HPLC analysis. Analyze the samples according to Section II.D., Instrumentation or II.F., Confirmatory Technique

D. Instrumentation

1.0 Description and Operating Conditions

An aliquot of the final sample fraction from Section II.C.4.2 is analyzed by column switching high performance liquid chromatography for residues of CGA-210007. The instrumental conditions are described in Table I and a diagram of the switching valve setup is presented in Figure 3. Control of the switching valve is accomplished via time-programmed contact closures of the detector, autoinjector, or other timing source.

- 1.1 Determine the retention time of CGA-210007 on column 1 (Spherisorb SAX/SB) by connecting the outlet of column 1 directly to the detector (by passing the switching valve) and injecting 25 ng of the analyte. (Inject 50 µl of the 0.5 ng/µl CGA-210007 standard solution prepared in Section II.I.1.1.). Record the starting and ending retention times of the eluting CGA-210007 analyte. Repeat the cuttime determination until the starting and ending retention times are consistent between analyses. The cuttime (ending retention time - starting

retention time) should be \leq 1 minute if possible.

- 1.2 Connect the HPLC system as shown in Figure 3 for **column switching**. Program the valve to switch to the INJECT position at the beginning of the CGA-210007 analyte peak and return to the LOAD position at the end of the analyte peak, as determined in Section II.D.1.1.
- 1.3 Inject 25 ng of CGA-210007 to determine its retention time through the two columns and to confirm that the valve time programming is working properly.

2.0 Calibration and Standardization

- 2.1 The HPLC system is calibrated with each analytical run by checking the cuttime and comparing it to the retention time and detector response relative to previous runs. Retention times should not vary by more than 5% and detector response should not vary more than 10% between runs. Otherwise the system should be inspected and proper maintenance or recalibration should be performed to minimize this variation. (See Section II.H.4.0, Modifications and Potential Problems.)
- 2.2 The HPLC system is standardized by injecting 50- μ L aliquots of the standard solutions of CGA-210007 in a working range of 0.35-25 ng/injection during residue analyses. A linear regression function is generated from the data comparing detector response and ng injected. Representative standardization data and chromatograms are presented in Figures 4. The corresponding calibration curve is presented in Figure 5.

E. Interferences

1.0 Analysis of control samples from the tobacco, leafy vegetable, and fruiting vegetable crop groups showed no significant (>0.01 ppm) interferences at the analyte retention time with a limit of quantitation of 0.02 ppm. No interferences have been observed in reagent blanks.

F. Confirmatory Technique

A confirmatory procedure for the detection and quantitation of CGA-245704 (as CGA-210007) has been developed that employs HPLC with mass spectrometric detection (LC/MS), in compliance with the confirmatory method requirement of US EPA Residue Chemistry Guidelines, OPPTS 860.1340, "Residue Analytical Methods." This technique utilizes the same extracts produced by the primary method, and therefore can be used in the event of chromatographic or other interferences with the HPLC/UV technique. The LC/MS technique can confirm detections of CGA-210007 to levels of <0.05 ppm in tobacco, and <0.02 ppm in both lettuce and tomatoes (fruit, wet pomace, and paste).

A validation of this technique was performed by analyzing identical sample sets on both detection systems simultaneously. In addition to calibration standards, a typical validation set consisted of one untreated control, two or more fortified controls, and two or more field samples containing incurred residues of CGA-245704.

After analysis by negative ion electrospray LC/MS (ES-), the peak area measurements of the 179 m/z ion (m-H⁻) were transcribed into the "Manual Data Entry" mode of the Multichrom Worksheet program for complete calculation of residues in the crop substrates.

Comparative results of the LC/MS detection technique versus HPLC/UV detection in the tobacco, lettuce, and tomato substrates are presented in Table X.

* Samples of flue-cured tobacco were obtained from an in-progress storage stability study and were fortified with CGA-210007 rather than CGA-245704.

Representative chromatograms of CGA-210007 in flue-dried tobacco, head lettuce, unwashed tomatoes, wet pomace, and tomato paste are presented in Figures 10-14, respectively. Each set of representative chromatograms includes a calibration standard, an untreated control, a fortified control (procedural recovery), and a treated sample containing field-incurred residues.

Complete HPLC and mass spectrometric conditions are presented in Table XI.

G. Time Required For Analysis

The extraction and cleanup of a set of 6 samples may be completed by a skilled analyst in approximately 6-8 hours. Automated injection of the samples can be performed overnight.

H. Modifications and Potential Problems

1.0 Any difficult emulsions formed during the sample work up can be minimized by swirling the emulsion with a glass stirring rod and allowing time for complete separation. If the emulsion is persistent (i.e. Section II.C.3 2), the complete sample should be transferred to a 250-ml polyurethane centrifuge bottle. Cap the bottle and proceed with centrifuging of sample(s) by slowly ramping to approximately 9,000 rpm for approximately 10 minutes and then slowly reduce speed. Decant the sample back into the separatory funnel and rinse the centrifuge bottle with an appropriate solvent; add the rinse to the separatory funnel containing the sample. Keep any emulsions with the analyte fraction until the final separation is conducted.

2.0 Use of the C₁₈ Mega Bond Elut column was necessary to minimize persistent emulsions encountered in substrates. The procedure is intended to remove unwanted nonpolar and basic components. The elution should be

conducted at maximum vacuum to minimize interaction of the ionized analyte with the SPE stationary phase.

- 3.0 A previous report noted that CGA-245704 is very unstable in acetonitrile solutions stored in light at room temperature. Therefore, the CGA-245704 standard solutions were stored in amber bottles during the validation of this method. The a.1 of a stock solution containing 100 µg/mL of CGA-245704 in acetonitrile (stored in the dark at -20°C), should be stable for at least 6 months
- 4.0 If minimal variability is noted in retention time and detector response between daily cuttime calibrations, then it is not necessary to calibrate on a daily basis. However if this condition changes, then daily calibrations should be begun until stability is achieved.
- 5.0 If an interference peak is present at the retention time of the CGA-210007 analyte during method validation or field trial sample analyses and this peak prevents accurate quantitation of CGA-210007, then reduce the starting and ending "cuttime" in portions of about 15 seconds each. Proceed until the interference disappears or is reduced to an insignificant height. If this does not remove the interference or analyses are being done for tolerance enforcement requirements, then proceed with quantitation by use of the confirmatory method.
- 6.0 The mobile phases for the HPLC systems if prepared in large quantity, can be recycled during the analyses. However, if problems develop such as increased baseline noise, control interference peak, or shifting retention times, this practice should be discontinued and new mobile phase prepared. The mobile phases should not be recycled during the determination of the analyte cuttime on column 1, as this would result in cross-contamination of the mobile phases.

- 7.0 LC/MS: Decreased procedural recoveries were observed when a sample injection volume of 100 μ L was used during method development. It is suspected that ionization of CGA-210007 is suppressed when injecting more than 50 μ L of a crop substrate matrix. The use of a column effluent splitter installed between the analytical HPLC column and the MS detector, with a split ratio optimized for the size of the sample injection, is recommended. Refer to Table XI for a description of the split ratio and injection volume used in this study.
- 8.0 LC/MS Since this method monitors the CGA-210007 [m-H] ion at 179 m/z, acetic acid (60 m/z) should not be used as a component of the LC/MS mobile phase, as it creates a spectral interference by reassociating at 179 m/z, representing the acetic acid [m-H]-2m negative ion.

I. Preparation of Standard Solutions and Sample Fortification Procedures

1.0 Preparation of Analytical Standard

Weigh 10 mg of CGA-210007 into a 100-mL volumetric flask, and dilute to the 100-mL mark with acetonitrile. The concentration of this stock solution is 100 ng/ μ L. Store the stock injection solution in the refrigerator when not in use.

- 1.1 Make serial dilutions of the 100 ng/ μ L stock solution with 0.02M H₃PO₄ until working solutions containing 0.007 ng/ μ L to 0.50 ng/ μ L are obtained. Store these injection standards in the refrigerator when not in use.

2.0 Preparation of Fortification Standard

Weigh 10 mg of CGA-245704 analytical standard into a 100-mL volumetric flask and dilute to the 100-mL mark with acetonitrile. The concentration of this

stock solution is 100 ng/ μ L. Store the stock fortification standard in the freezer when not in use.

- 2 1 Make serial dilutions of the 100 μ g/mL standard solution with acetonitrile to give a series of fortification standards in a range of 0.02 ng/ μ L to 1.0 ng/ μ L of CGA-245704. Store these fortification standards in the refrigerator when not in use.

3.0 Sample Fortification Procedures

The method is validated for each set of samples analyzed by including an untreated control sample and one or more control samples fortified prior to extraction with CGA-245704.

- 3.1 Add between 0.50 mL to 2.0 mL of the appropriate standard solution of CGA-245704 to 4 g of the crop substrate prior to the addition of the hydrolysis solution (Section II.C 2.1) Allow the sample to sit for a few minutes before adding the aqueous NaOH required for hydrolysis. If a higher concentration spike is desired, adjust the concentration of the fortification solution so that no more than 2 mL of solution is added to the substrate.
- 3.2 Analyze the fortified samples (procedural recoveries) through the same method procedures as for the treated samples.
- 3.3 See Section II.J.2.0 for the procedure to calculate the ppm of CGA-245704 in the samples.

J. Methods of Calculations

1.0 Determination of Sample Residues
(Tolerance Enforcement Data)

Inject 50- μ L aliquots of the final fractions prepared in Section II.C.4 2 into the HPLC under the same conditions as the standards. Make appropriate dilutions of the samples (if necessary) with 0.02M H_3PO_4 to bring the sample peak heights within the range of the standard curve. Compare the peak heights observed in the treated (unknown analyte quantity) samples to the standard curve or enter the peak height into the least squares program to determine the nanograms of CGA-210007 in the injected aliquot. Representative chromatograms for reagent blanks, control, and procedural recovery crop samples are shown in Figures 4-9.

To calculate the residue results, the mg injected must first be calculated according to Equation 1:

$$(1) \text{ Sample weight injected (mg)} = \frac{(G)(V_2)(V_1)}{(V_1 V_e + G(M/100))} \times CF_1 \times CF_2$$

- G = sample weight extracted (g)
- V_a = aliquot volume (mL)
- V_1 = injection volume (μ L)
- V_2 = volume of final fraction (mL)
- V_e = extraction volume (mL)
- M = % moisture of sample
- CF_1 = conversion factor to convert grams of sample to milligrams ($CF_1 = 1000/1$)
- CF_2 = conversion factor to convert microliters to milliliters ($CF_2 = 1/1000$)

In this method, the total extraction aliquot is used for analyses therefore, $V_a = [V_e + G(M/100)]$, and these values cancel out of Equation 1.

Determine the concentration (ppm) of the analyte in residue samples (tolerance enforcement) according to Equation 2

$$(2) \text{ Analyte (ppm)} = \frac{(\text{ng analyte found})}{(\text{mg sample injected})}$$

To convert the ppm found for CGA-210007 to parent CGA-245704 equivalents, multiply the results from Equation 2 by 1.17. This ratio (F) is determined from the molecular weights as defined below.

$$(3) \text{ Ratio of molecular weights (F)} = \frac{\text{MW parent (210.28 g/mol)}}{\text{MW of the (180.18 g/mol) CGA-210007 metabolite}} = 1.17$$

2.0 Determination of Procedural Recoveries
(Method Validation and Residue Field Trial Data)

- 2.1 Calculate the ng of CGA-210007 and mg sample weight injected using the procedures described in Section II.J.1.0. Calculate the final ppm values of the control and fortified samples according to Equation 2 (Section II.J.1.0).

Determine the recovery factor by first subtracting the background detector response, if any, in the control sample from the analyte response in the recovery sample. Calculate the recovery factor as a percentage (R) by Equation 4:

$$(4) R = \frac{(\text{ppm analyte found} - \text{ppm control})}{\text{ppm analyte added}} \times 100$$

Convert the ppm found of CGA-210007 to parent CGA-245704 equivalents by multiplying by F, the ratio of the molecular weights.

3.0 Determination of Corrected Sample Residues
(Method Validation and Residue Field Trial Data)

- 3.1 Determine the ng of CGA-210007 injected and calculate the mg injected using the procedures described in Section II.J.1.0.

Calculate the final ppm values of the residue samples according to Equation 2 (Section II.J 1 C). Calculate the ppm found of CGA-210007 to parent CGA-245705 equivalents by multiplying by F, the ratio of the molecular weights. Calculate the method recovery according to Equation 4 (Section II.J 2.1). Use the recovery, or mean recovery if two or more fortified samples are analyzed, to correct the concentration (ppm) of each analyte for the treated sample according to the following equation.

$$(5) \text{ Corrected Residue (ppm)} = \frac{\text{ng analyte found}}{\text{mg sample}} \times \frac{1}{R}$$

Note: Recovery adjustments for controls, recovery corrections for small control values (≤ 0.02 ppm), and recoveries for residue concentrations are not used in tolerance enforcement calculations. The recovery correction is used to correct the amount of analyte found in the fortified (spiked) samples. The recovery correction is utilized for method validation and for the correction of residue data.

III. RESULTS AND DISCUSSION

Analytical Method AG-671A was evaluated for accuracy, precision, extractability, and accountability by analyses of residue field trial samples and metabolism radiolabeled samples. A comprehensive report of these details is presented in ABR-97117, Final Report Amendment Number 1⁴. A summary of the data and conclusions are outlined below.

A. Accuracy

Untreated tobacco, leafy vegetable and fruiting vegetable (including tomato processing fractions) crop group substrates were fortified with

CGA-245704 at levels ranging from 0.02-20 ppm. Recovery results from fortified control samples (Table VI) were used to calculate accuracy in terms of a mean, standard deviation (sd), and coefficient of variation (CV%) = $sd/mean \times 100$. Procedural recovery data from fortified control field trial samples, along with procedural recovery data from validation sets analyzed herein (Table II), were compiled to assess the overall accuracy of Analytical Method AG-671A.

Procedural recoveries were calculated by dividing the concentration of CGA-245704 measured in each fortified sample by the concentration added to the sample. Representative Multichrom analysis sheets (spreadsheets) showing data for all of the intermediate calculations are presented in Tables III-V. Substrate, fortification levels, and residues in the fortified controls are presented in Table VI and summarized below.

Commodity	CGA-245704	CGA-245704	
	Mean Recovered (%)	Mean	Standard Deviation
Tobacco (Green, Uncured)	84	9.1	(n=17)
Tobacco (Cured)	88	1.1	(n= 7)
Head Lettuce	94	10	(n=14)
Leaf Lettuce	90	9.3	(n=14)
Spinach	89	8.3	(n=15)
Celery	90	9.2	(n=17)
Tomato	87	6.1	(n=15)
Tomato Processing	100	12	(n=25)
Bell Peppers	84	4.7	(n= 9)
Hot Peppers	96	1.2	(n= 9)

No significant interferences (>0.01 ppm) were found at or near the retention time of CGA-210007 in the reagent blanks (Figure 6) or the control samples. An interference (>0.01 ppm) was present at the retention time of CGA-210007 in the tomato processing puree (0.01 ppm) and paste (0.02 ppm) sample sets. This apparent contamination was detected in the tomato fruit of one of the field trials (02-FR-030-96) and increased in magnitude with concentration of the processing fractions. Subsequent verification of the tomato puree residue by liquid chromatography/mass spectroscopy indicates the presence of a CGA-210007 contaminate in the control sample. Results from an additional field trial (0W-FR-412-96) showed no indication of any significant (>0.01 ppm) residue or interference in the tomato processing control samples. HPLC

chromatograms for the ¹⁴C-tobacco, ¹⁴C-lettuce and ¹⁴C-tomato samples (control, fortified, and radiolabeled samples) are shown in Figures 7-9. Representative chromatograms of fortified field trial and processing (control and procedural recovery) samples are presented in ABR-97117, Final Report Amendment Number 1.

These acceptable procedural recovery results confirm that the analytical procedure outlined in Analytical Method AG-671A is adequate for recovery of residues and is satisfactory for use as a tolerance enforcement method.

B. Precision

Precision of the method was determined by calculating the mean, standard deviation, and coefficient of variation of results from the analysis of fortified control (Table VI) and ¹⁴C-incurred samples. Table VII contains the statistical results for recovery values from fortification experiments. Table VIII contains the statistical results for the triplicate analyses of ¹⁴C-incurred tobacco, lettuce, and tomato samples. The precision of Analytical Method AG-671A is demonstrated by its CV values and by the close agreement of triplicate sample results.

C. Extractability

The extractability of the TRR by Analytical Method AG-671A was determined by comparing the total ppm of ¹⁴C-residue in the sample from combustion analysis to the ppm ¹⁴C-residue in the initial sample extract (Section II.C 2.2 3) of the method. The formula for the determination of percent extractability is:

$$\% \text{ extractability} = \frac{\text{Total } ^{14}\text{C-residue (ppm) in sample extract by LSC}}{\text{TRR (ppm)}} \times 100$$

Table VIII summarizes the TRR values, concentration of ¹⁴C-residue in sample extracts, and the total extractability. The TRR values were determined by combustion analysis in the metabolism studies and are reported in ¹⁴C-CGA-245704 equivalents. The concentration of the ¹⁴C-extractable residue was calculated by LSC analyses of the extracts and

reported in ppm ^{14}C -CGA-245704 equivalents. The total extractabilities were 108%, 87%, and 99% for the ^{14}C -tobacco, ^{14}C -lettuce, and ^{14}C -tomato, respectively. These extractability results are consistent with quantitation by 2-dimensional thin layer chromatography used in the metabolism studies. The validation extractability results are compared with the metabolism extractability results in Table IX.

D. Accountability

The accountability of Analytical Method AG-671A is determined by comparing the total ppm of ^{14}C -residue found in the sample by combustion analysis to the ppm of CGA-245704 equivalent residue found in the final extract by HPLC analysis. The formula for the determination of percent accountability is:

$$\% \text{ accountability} = \frac{\text{Total residue (ppm) in sample extract determined by HPLC analysis}}{\text{TRR (ppm)}} \times 100$$

The total mean accountabilities were 78%, 33%, and 71% for the ^{14}C -tobacco, ^{14}C -lettuce, and ^{14}C -tomato, respectively. These accountability results are consistent with the 2-dimensional thin layer chromatography results obtained in the respective metabolism studies after hydrolysis of the sample to free CGA-210007. In Analytical Method AG-671A, the hydrolysis step is based on treatment of the sample with 1N NaOH/heat and would most likely result in a greater conversion to total CGA-210007⁵⁻⁷ than the metabolism 0.1N NaOH/room temperature treatment. The accountability results for lettuce were less than the results obtained in the metabolism study after quantifying by 2-dimensional thin layer chromatography. The metabolism accountability for lettuce was based on the surface rinse radioactivity and the lettuce head penetrated radioactivity. The hydrolysis conditions for conversion of penetrated radioactivity were conducted in 0.1N NaOH (room temperature) and cellulase. The accountability results reported herein, were based on an unwashed lettuce sample. The difference in magnitude of the lettuce accountability results can be attributed to

the different metabolism profile for lettuce. The validation accountability results are compared with the metabolism accountability results in Table IX.

IV. CONCLUSION

Novartis Analytical Method AG-671A is shown to be a valid method for the determination and confirmation of CGA-245704 in RAC and processing substrates from the tobacco, leafy and fruiting vegetable crop groups. The validity is demonstrated by the acceptable accuracy and precision obtained on numerous procedural recovery samples (radiovalidation and field trial sample sets), and by the extractability and accountability obtained by the analysis of weathered radioactive substrates using Analytical Method AG-671A. Validation results are reported in ABR-97117, Final Report Amendment Number 1⁴.

V. CERTIFICATION

The reports and experimental results included in this study, Laboratory Project I.D. Analytical Method AG-671A, are certified to be authentic accounts of the experiments.

Carol G Hayworth
Carol G Hayworth, B.S.
Study Director
Scientist I, Residue Chemistry
Human Safety Department
(336-632-6000)

3/25/78
Date

VI. TABLES AND FIGURES

TABLE I LIQUID CHROMATOGRAPHIC OPERATING CONDITIONS FOR ANALYSIS OF CGA-210007

Instrument Conditions ^{1 2}	
Column 1/Typical Retention Time	Phenomenex Spherisorb [®] SAX/SB, 150 x 3.2 mm I.D., 5 µm, HPLC Column Cat. #00F-0110-R0 Retention Time ~7 min.
Mobile Phase 1/Flow Rate	80:20 0.02M H ₃ PO ₄ /ACN Flowrate 0.30 mL/min
Column 2	YMC-Pack C8 250 x 4.6 mm I.D., S-5µm, 120A, HPLC Column, Cat #OC12S05-2546WT A-203-5
Mobile Phase 2/Flow Rate	75:25 0.02M H ₃ PO ₄ /ACN Flowrate 1.0 mL/min
Typical Retention Time Through Column 1 & Column 2. Column Oven/Temperature:	Retention Time ~22 min Eppendorf column heater CH-30 Cat. #4301-000-028 and gauge TC-50 Cat. #4300-000-008 or equivalent Set at 30 °C
Detection/UV Wavelength	Perkin Elmer LC-95 UV/Visible Spectrophotometer Detector, or equivalent Wavelength, 235 nm
Pumps(2):	Perkin Elmer Isocratic LC Pump Model 250, or equivalent
Column Switching Valve:	VICI EQ60 HPLC Switching Valve, or equivalent
Auto Injector:	Perkin Elmer Model ISS 100, or equivalent
Volume Injected	50 µL
Data Acquisition.	MicroVAXII (Q) Operating System, VMS Version 5.3-1 Application Software VG Multichrom Version 2.0 Worksheet Version. 1.6.0, or equivalent

¹ HPLC columns should be converted from the storage solvent to the appropriate mobile phase according to the manufacturer's instructions. HPLC columns should be equilibrated until the retention time of the standard (CGA-210007) is stabilized.

² Retention times listed are approximate and will vary depending on instrument setup.

TABLE II. SOURCE OF ACCURACY AND PRECISION DATA

Project #	Study #	Commodity	Field Trial # or Analytical Set Name	Project #	Study #	Commodity	Field Trial # or Analytical Set Name
436001	491-97	¹⁴ C Tobacco	Reference 5	436003	92-96	Celery	OW FR 539-96
436001	491-97	¹⁴ C Lettuce	Reference 6	436003	92-98	Celery	OW FR 540-96
436001	491-97	¹⁴ C Tomato	Reference 7	436003	92-96	Celery	OW FR 541-96
436001	91-96	Tobacco	OS FR 607-96	436003	92-96	Celery	OW FR 542-96
436001	91-96	Reagent Blank	91-96 BBR	436003	92-96	Celery	TL FR 401-96
436001	91-96	Tobacco	NE FR 211-96	436003	92-96	Celery	NE FR 727-96
436001	91-96	Tobacco	NE FR 212-96	436004	93-98	Tomato, Whole	02 FR 030-96
436003	92-96	Head Lettuce	02 FR-028-96	436004	93-98	Tomato, Whole	OW FR 412-96
436003	92-96	Reagent Blank	92-96-BBG	436004	93-96	Tomato, Whole	TL FR 404-96
436003	02-96	Head Lettuce	05 FR-008-96	436004	93-96	Reagent Blank	Contract Method Trial Data ¹ FL FR 404-96
436003	92-96	Head Lettuce	04 FR-533-96	436004	93-96	Tomato, Wei Pomato	Contract Method Trial Data ¹ 02 FR-030-96
436003	92-96	Head Lettuce	04 FR-534-96	436004	93-96	Reagent Blank	93-96-BA2
436003	92-96	Head Lettuce	04 FR-535-96	436004	93-96	Tomato, Wei Pomato	OW FR 412-96
436003	92-96	Head Lettuce	FL FR-02-96	436004	93-96	Reagent Blank	93-96-031
436003	92-96	Leaf Lettuce	02-FR-027-96	436004	93-96	Tomato Juice	02 FR-030-96
436003	92-96	Reagent Blank	92-96-BA0	436004	93-96	Tomato Juice	OW FR-412-96
436003	92-96	Leaf Lettuce	05-FR-010-96	436004	93-96	Tomato Puree	02 FR-030-96
436003	92-96	Leaf Lettuce	07-FR-006-96	436004	93-96	Tomato Puree	OW FR-412-96
436003	92-96	Leaf Lettuce	04 FR-536-96	436004	93-96	Tomato Paste	02-FR-030-96
436003	92-96	Leaf Lettuce	04 FR-537-96	436004	93-96	Tomato Paste	OW FR 412-96
436003	92-96	Leaf Lettuce	04 FR-538-96	436004	93-96	Hot Pepper	02-FR-071-96
436003	92-96	Spinach	02 FR-028-96	436004	93-96	Hot Pepper	Contract Method Trial Data ¹ 05 FR 323-96
436003	92-96	Reagent Blank	92-96-BA2				Contract Method Trial Data ¹
436003	92-96	Spinach	05 FR 318-96				
436003	92-96	Spinach	04 FR 543-96				
436003	92-96	Spinach	04 FR 317-96				
436003	92-96	Spinach	NE-TR 304-96				
436003	92-96	Spinach	NE FR 501-96				

¹ Contract method trial data (Reference 10) was generated by Rizeema, Inc.

TABLE III. ¹⁴C-TOBACCO MULTICHROM WORKSHEET

Worksheet	491-97-044	Unit of Duplicate	0.010 ppm	Extracted Analyzed Channel	4-Sep-1987	Product	CON-245/04		
Area Scan Name	4919/10	Corr Factor A	1.000	Channel	22-Sep-1987	Sample	CON-210007		
Analysis	TOMATO/TOBACCO/LEAF	Corr Factor B	1.000				Project		
Subject	TOMATO/TOBACCO/LEAF						491-97		
1	TC	Field Test	Sap Crude	Moist	Smp Wt	Final	Final	Inf	Smp Wt
2	S			%	g	ml	ml	ul	mg
3	1	1	0.00271	0.00	4.00	172.00	172.00	50.00	50.00
4	1	1	0.0000	0.00	4.00	172.00	172.00	50.00	50.00
5	1	1	0.0000	0.00	4.00	172.00	172.00	50.00	50.00
6	1	1	0.0000	0.00	4.00	172.00	172.00	50.00	50.00
7	1	1	0.0000	0.00	4.00	172.00	172.00	50.00	50.00
8	1	1	0.0000	0.00	4.00	172.00	172.00	50.00	50.00
9	1	1	0.0000	0.00	4.00	172.00	172.00	50.00	50.00
10	1	1	0.0000	0.00	4.00	172.00	172.00	50.00	50.00
11	1	1	0.0000	0.00	4.00	172.00	172.00	50.00	50.00
12	1	1	0.0000	0.00	4.00	172.00	172.00	50.00	50.00
13	1	1	0.0000	0.00	4.00	172.00	172.00	50.00	50.00

Std Wt	Std Wt	Peak Weight	Analyte	Residue	Value	Added	Avg	Correct	Comments
mg	mg	μg	ng	ppm	control	ppm	Recovery	Residue	
1	0.20000	1.00000	1514	0.001			98.21	<0.020	STANDARD
2	0.04000	2.00000	295	-0.042					CONTROL TUL-C1
3	0.50000	25.00000	3753	0.021	0.021	1.500	105.76		STANDARD
4	0.00700	0.35000	56	0.021	1.360	1.500	90.67		STANDARD
5	0.10000	5.00000	167	1.088	1.088				CONTROL #0 02PPM TUL-R1
6	0.00000	0.00000	177	1.088	1.088				CONTROL #1 50PPM TUL-R2
7	0.00000	0.00000	182	1.088	1.088				STANDARD
8	0.00000	0.00000	1872	1.088	1.088				STANDARD
9	0.00000	0.00000	1799	1.088	1.088				STANDARD
10	0.00000	0.00000	1803	1.088	1.088				STANDARD
11	0.00000	0.00000	94	0.458	0.458				STANDARD
12	0.00000	0.00000	1657	10.827	10.827				STANDARD

Notes: Study 401-97, Radiovalidation of AG-671
 Substrate tobacco plant leaves, uncurd
 Source Single Metabolism Study 94M04, Harvest Interval 5

Blank Control F Freezer R Recovery S Standard V Solvent X Sample
 Control used for this group Negative controls are treated as 0 Underlined, Strikethrough Manual counting rejected

TABLE IV ¹⁴C-LETTUCE MULTICHROM WORKSHEET

Worksheet No	491 97-11-2	Limit of Detect	0.020 ppm	Extraction Channel	24 Sep-1991	Product	CO ₂ -245704	
Anal. Sub Name	491 97-11-2	Corr Factor A	1.0000	Analyzed Channel	24 Sep-1991	Analyte	CO ₂ -210007	
Subject	TORNADO TOMATO/LETTUCE	Corr Factor B	1.1100		24	Project	491-97	
1	1 g Fluid Test	Moist %	wt ml	Integra ml	Aliq 2 ml	Final vol ml	Inj vol ul	Sep Wt mg
2	S	0.00	4.00	174.00	174.00	4.00	4.0	50.000
3	C	0.00	4.00	174.00	174.00	4.00	4.0	50.000
4	R	0.00	4.00	174.00	174.00	4.00	4.0	50.000
5	R	0.00	4.00	174.00	174.00	4.00	4.0	50.000
6	R	0.00	4.00	174.00	174.00	4.00	10.0	20.000
7	S	0.00	4.00	174.00	172.00	3.954	10.0	19.770
8	S	0.00	4.00	174.00	172.00	3.954	10.0	19.770
9	X	0.00	4.00	174.00	172.00	3.954	10.0	19.770
10	X	0.00	4.00	174.00	172.00	3.954	10.0	19.770
11	X	0.00	4.00	174.00	172.00	3.954	10.0	19.770
12	S	0.00	4.00	174.00	172.00	3.954	10.0	19.770
13	S	0.00	4.00	174.00	172.00	3.954	10.0	19.770

Substrate	Std	Conc	Peak Height	Analyte Found	Value	Added	Recovery	Avg Recov	Corr	Dry Wt	Comment
1	g	mg/ml	μV	ppm	correct for control	ppm	%	%	ppm	ppm	
1	S	0.20000	1513	0.703	0.000	0.020	95.23	95.23	-0.020	50.000	STANDARD
2	S	0.01000	86	0.311	0.000	0.000				50.000	CONTROL 17-C1
3	C	0.50000	3909	25.131	0.019	0.020	91.68	91.68		50.000	CONTROL 02 RUN 17-B1
4	S	0.01000	131	0.801	0.968	1.000	96.77	96.77		50.000	CONTROL 01 OFFN 17-B2
5	R	0.01000	2575	16.542	0.329					50.000	STANDARD
6	R	0.01000	185	1.709	0.301					50.000	STANDARD
7	S	0.04000	796	0.871	0.270					50.000	INCURRED SAMPLE 17-T1
8	S	0.07000	715	0.885	0.316					50.000	INCURRED SAMPLE 17-T12
9	X	0.10000	783	0.562	0.284					50.000	INCURRED SAMPLE 17-T13
10	X	0.10000	1559	5.998						50.000	STANDARD
11	X	0.20000								50.000	STANDARD
12	S	0.10000								50.000	STANDARD
13	S	0.20000								50.000	STANDARD

Notes: Study 491-97, Radiovalidation of AG-671
 Solvent: V
 Source: Basile Metabolism Study 96N02.1, Intv 2
 F: Freezer
 R: Recovery
 S: Standard
 X: Sample
 Negative controls are treated as 0.0 Underline Strikethrough: Admin pyridin 1a3etm

TABLE V. ¹⁴C-TOMATO MULTICHROM WORKSHEET

Run	Field Test	Sample Code	Moist %	Net wt g	Ext g	Net wt g	Ext g	Exl vol ml	Aliq 1 vol ml	Interm vol ml	Aliq 2 vol ml	Final wt g	Final vol ml	Inj vol ul	Spk wt mg
1 S															
2 S															
3 C	1	02HRO3D96	0.00	4.00	1.72	0.00	1.72	0.00	1.72	0.00	4.00	4.00	4.0	50.00	50.000
4 S															
5 R															
6 R															
7 R															
8 S															
9 X	1	RHW01697	0.00	4.00	1.72	0.00	1.72	0.00	1.72	0.00	3.953	4.0	50.00	49.412	
10 X	1		0.00	4.00	1.72	0.00	1.72	0.00	1.72	0.00	3.953	4.0	50.00	49.419	
11 X	1		0.00	4.00	1.72	0.00	1.72	0.00	1.72	0.00	3.953	4.0	50.00	49.419	
12 S															
13 S															

DI	Std vol	Std wt	Analyte Found	Residue Found	Value corr'd for control	Added Amount ppm	Recov %	Net Recov %	Corr'd Residue ppm	Comments
1 S	0.20000	10.0000	10.197	0.514	0.003	0.020	96.7	92.24	0.020	STANDARD CONTROL TO C-1
2 S	0.01000	0.5000	1.77	0.147	0.003	0.020	96.7	92.24	0.020	STANDARD CONTROL TO C-1
3 C	0.04000	2.0000	5.5	0.441	0.003	0.500	88.21	92.24	0.020	STANDARD CONTROL TO C-1
4 S			31.5	2.041	0.023	0.500	88.21	92.24	0.020	STANDARD CONTROL TO C-1
5 S			59.96	18.995	0.444	0.500	88.21	92.24	0.020	STANDARD CONTROL TO C-1
6 R			78.77	24.961	0.441	0.500	88.21	92.24	0.020	STANDARD CONTROL TO C-1
7 C	0.50000	25.0000	140.0	8.289	0.181					INCURRED SAMPLE TO-11
8 S	0.00100	0.3500	24.30	1.950	0.084					INCURRED SAMPLE TO-12
9 X			27.80	1.950	0.222					INCURRED SAMPLE TO-13
10 X			29.65	1.971	0.222					STANDARD
11 X			29.65	1.971	0.222					STANDARD
12 S	0.10000	5.0000	15.44	4.671						STANDARD
13 S	0.20000	10.0000	31.47	9.957						STANDARD

Substrate: UNWRAPPED TOMATOES
TREATMENT: COA-245/04

Motor: Study 491-97, Revalidation of AG-671
Substrate: Tomato, after washing
Source: Bvto Macabaita Study 51MA05, Interval 4

D Blank
C Control
F Freezer
R Recovery controls are treated as 0
S Standard
V Solvent
X Sample
Negative controls are treated as 0
Underline, Strike through: Manual override, Rejected

TABLE VI. RECOVERY DATA FROM CONTROL SUBSTRATES FORTIFIED WITH CGA-245704 AND ANALYZED BY ANALYTICAL METHOD AG-671A

Commodity	Substrate	Sample Code	Control (µg/ml)	CGA 245704 Added (µg/ml)	CGA-245704 Found (µg/ml)	Recovery CGA-245704 Corrected (%) ^{1,2}	Mean Recovery (%)	Standard Deviation	Mean CV (%)		
Tobacco	Resonant Blank Lower Leaves, Uncured Green	TUL RI, TUL R2	-0.02	1.5	<0.001	-0.02	91	84	9.1 (n=17)	11	
		12-A	0.02	0.02	0.021	1.350	106				
		14-A	-0.02	0.50	0.452	-	90				
		14-A	-0.02	0.10	0.074	0.015	74				
		18-A	-0.02	1.5	1.330	2.892	89				
		18-A	-0.02 (0.001)	0.02	0.10	0.063	75				
		18-A	-0.02	0.50	0.402	0.80	80				
		18-A	-0.02	0.50	0.382	0.76	76				
		14-A	-0.02 (0.003)	0.02	0.029	0.178	84				
		14-A	-0.02	1.0	0.859	0.88	88				
		14-A	-0.02	0.02	0.10	0.014	0.076				70
		14-A	-0.02	0.50	0.453	0.91	91				
		13-A	-0.02 (0.002)	5.0	3.917	78					
		15-A	-0.02	4.0	3.036	76					
Head Lettuce	Resonant Blank With Wrapper Leaves	LT RI, LT R2	-0.02	1.0	<0.001	-0.02	88	88	11 (n=7)	12	
		15-A	0.02	0.02	0.019	0.669	67				
		15-A	-0.02	0.02	0.023	0.083	117				
		15-A	-0.02	0.02	0.022	0.114	104				
		15-A	-0.02	0.02	0.018	0.018	91				
		15-A	-0.02	0.10	0.089	0.093	88				
15-A	-0.02	0.10	0.092	0.093	92						
15-A	-0.02	0.50	0.459	0.458	83						
15-A	-0.02	0.50	0.454	0.454	91						

¹ Minor discrepancies in residue recoveries (%) between raw data and that reported in this table are a function of round-off differences
² Results are corrected for control residues, if detected

TABLE VI. RECOVERY DATA FROM CONTROL SUBSTRATES FORTIFIED WITH CGA-245704 AND ANALYZED BY ANALYTICAL METHOD AG-671A (Continued)

Spotting	Substrate	Sample Code	Control (nm)	CGA-245704 Added (nm)	CGA-245704 Found (nm)	Recovery CGA-245704 Corrected (%) ¹	Mean Recovery (%)	Standard Deviation	Mean CV (%)
Leaf Lettuce	Reagent Blank Leaves	1-2 A	-0.02 (0.005)	0.20	0.169	-0.001	84		
		1-3 A	-0.02 (0.009)	0.10	0.089	0.000	81		
		1-4 A	-0.02	0.50	0.423	0.453	115		
		1-5 A	-0.02	0.02	0.020	0.004	93		
		1-6 A	-0.02	0.02	0.017	0.030	84		
		1-7 A	-0.02	0.02	0.017	0.017	80		
		1-8 A	-0.02 (0.006)	0.02	0.023	0.023	78		
		1-9 A	-0.02 (0.001)	0.10	0.091	0.091	90		
		1-10 A	-0.02	0.02	0.049	0.049	84		
		1-11 A	-0.02	0.02	0.018	0.018	81		9.0 (n=14)
Spinach	Reagent Blank Leaves	1-2 A	-0.02	0.02	0.020	-0.001	55		
		1-3 A	-0.02	0.10	0.095	0.501	85		
		1-4 A	-0.02 (0.003)	0.02	0.018	0.500	78		
		1-5 A	-0.02 (0.003)	0.10	0.095	1.365	82		
		1-6 A	-0.02	0.02	0.015	0.084	73		
		1-7 A	-0.02	0.02	0.018	0.018	80		
		1-8 A	-0.02 (0.005)	0.02	0.020	0.020	78		
		1-9 A	-0.02	0.02	0.017	0.017	85		
		1-10 A	-0.02 (0.008)	0.02	0.024	0.024	84		
		1-11 A	-0.02	0.02	0.024	0.024	84		8.5 (n=13)

¹ Minor discrepancies in residue recoveries (%) between raw data and that reported in this table are a function of round-off differences
² Results are corrected for control residues if detected

TABLE VI. RECOVERY DATA FROM CONTROL SUBSTRATES FORTIFIED WITH CGA-245704 AND ANALYZED BY ANALYTICAL METHOD AG-671A (Continued)

Commodity	Substrate	Sample Code	Control Level ¹	CGA-245704 Added (µg/ml)	CGA-245704 Found (µg/ml)	Recovery CGA-245704 Corrected (%) ²	Mean Recovery (%)	Standard Deviation	Mean CV (%)	
Celery	Leaf Steak	1-2-B	-0.02	0.02	0.030	98	95			
		1-3-A	-0.02 (0.002)	0.10	0.092	90	92			
		1-4-A	-0.02 (0.01)	0.02	0.017	80	92			
		1-5-A	-0.02	0.10	0.094	84	97			
		1-6-A	-0.02 (0.002)	0.02	0.013	85	87			
		1-7-A	-0.02	0.02	0.020	72	72			
		1-8-A	-0.02	0.02	0.014	72	82			
		1-9-A	-0.02	0.10	0.098	80	105			
		1-10-A	-0.02	0.02	0.018	80	80		9.2 (n=17)	10
		1-11-A	-0.02 (0.003)	0.02	0.023	0.444	86	88		
Tomato	Fruit	TO-R1, TO-R2	-0.02 (0.007)	0.10	0.102	94	73			
		1-8-AB-1	-0.02 (0.007)	0.50	0.501	99	99			
		1-9-AB-1	-0.02 (0.004)	0.50	0.501	99	99			
		1-4-A	-0.02 (0.009)	0.02	0.021	92	90			
		1-5-A	-0.02 (0.001)	0.20	0.209	84	71			
		1-6-A	-0.02	1.0	0.974	81	87			
		1-7-A	-0.02	1.0	0.974	83	97			
		1-8-A	-0.02	1.0	0.974	87	87			
		1-9-A	-0.02	1.0	0.974	87	87			
		1-10-A	-0.02	1.0	0.974	87	87			
Fruit Machine, Unwashed	Fruit Machine, Unwashed	1-8-AB-4	-0.02	0.02	0.001	-0.02	-0.02			
		1-9-AB-4	-0.02	0.10	0.025	0.078	83	87		
		1-10-AB-4	-0.02 (0.007)	0.10	0.099	1.078	99	98		
		1-11-AB-4	-0.02 (0.004)	0.10	0.110	0.760	100	75		
		1-12-AB-4	-0.02	0.10	0.062	0.310	98	102		
		1-13-AB-4	-0.02 (0.011)	0.02	0.038	0.114	111	103		
		1-14-AB-4	-0.02	0.02	0.031	1.009	94	93		
		1-15-AB-4	-0.02	0.02	0.018	0.319	92	104		
		1-16-AB-4	-0.02	0.10	0.14	1.004	91	107		
		1-17-AB-4	-0.02	0.02	0.040	0.849	90	100		
Fruit Machine, Washed	Fruit Machine, Washed	1-8-AB-7	-0.02 (0.002)	0.02	0.019	85	104			
		1-9-AB-7	-0.02	0.02	0.022	101	101			
		1-10-AB-7	-0.02	0.02	0.022	101	101			
		1-11-AB-7	-0.02	0.02	0.022	101	101			
		1-12-AB-7	-0.02	0.02	0.022	101	101			
		1-13-AB-7	-0.02	0.02	0.022	101	101			
		1-14-AB-7	-0.02	0.02	0.022	101	101			
		1-15-AB-7	-0.02	0.02	0.022	101	101			
		1-16-AB-7	-0.02	0.02	0.022	101	101			
		1-17-AB-7	-0.02	0.02	0.022	101	101			
Reagent/Blank	Reagent/Blank	1-8-AB-1	-0.02	0.02	0.001	-0.02	-0.02			
		1-9-AB-1	-0.02	0.02	0.001	-0.02	-0.02			
		1-10-AB-1	-0.02	0.02	0.001	-0.02	-0.02			
		1-11-AB-1	-0.02	0.02	0.001	-0.02	-0.02			
		1-12-AB-1	-0.02	0.02	0.001	-0.02	-0.02			
		1-13-AB-1	-0.02	0.02	0.001	-0.02	-0.02			
		1-14-AB-1	-0.02	0.02	0.001	-0.02	-0.02			
		1-15-AB-1	-0.02	0.02	0.001	-0.02	-0.02			
		1-16-AB-1	-0.02	0.02	0.001	-0.02	-0.02			
		1-17-AB-1	-0.02	0.02	0.001	-0.02	-0.02			
Wet Pomace	Wet Pomace	1-8-AB-4	-0.02	0.02	0.001	-0.02	-0.02			
		1-9-AB-4	-0.02	0.10	0.025	0.078	83	87		
		1-10-AB-4	-0.02 (0.007)	0.10	0.099	1.078	99	98		
		1-11-AB-4	-0.02 (0.004)	0.10	0.110	0.760	100	75		
		1-12-AB-4	-0.02	0.10	0.062	0.310	98	102		
		1-13-AB-4	-0.02	0.02	0.038	0.114	111	103		
		1-14-AB-4	-0.02	0.02	0.031	1.009	94	93		
		1-15-AB-4	-0.02	0.02	0.018	0.319	92	104		
		1-16-AB-4	-0.02	0.10	0.14	1.004	91	107		
		1-17-AB-4	-0.02	0.02	0.040	0.849	90	100		
Fresh Juice	Fresh Juice	1-8-AB-7	-0.02 (0.002)	0.02	0.019	85	104			
		1-9-AB-7	-0.02	0.02	0.022	101	101			
		1-10-AB-7	-0.02	0.02	0.022	101	101			
		1-11-AB-7	-0.02	0.02	0.022	101	101			
		1-12-AB-7	-0.02	0.02	0.022	101	101			
		1-13-AB-7	-0.02	0.02	0.022	101	101			
		1-14-AB-7	-0.02	0.02	0.022	101	101			
		1-15-AB-7	-0.02	0.02	0.022	101	101			
		1-16-AB-7	-0.02	0.02	0.022	101	101			
		1-17-AB-7	-0.02	0.02	0.022	101	101			
Purée	Purée	1-8-AB-4	-0.02	0.02	0.001	-0.02	-0.02			
		1-9-AB-4	-0.02	0.10	0.025	0.078	83	87		
		1-10-AB-4	-0.02 (0.007)	0.10	0.099	1.078	99	98		
		1-11-AB-4	-0.02 (0.004)	0.10	0.110	0.760	100	75		
		1-12-AB-4	-0.02	0.10	0.062	0.310	98	102		
		1-13-AB-4	-0.02	0.02	0.038	0.114	111	103		
		1-14-AB-4	-0.02	0.02	0.031	1.009	94	93		
		1-15-AB-4	-0.02	0.02	0.018	0.319	92	104		
		1-16-AB-4	-0.02	0.10	0.14	1.004	91	107		
		1-17-AB-4	-0.02	0.02	0.040	0.849	90	100		
Pulp	Pulp	1-8-AB-7	-0.02 (0.002)	0.02	0.019	85	104			
		1-9-AB-7	-0.02	0.02	0.022	101	101			
		1-10-AB-7	-0.02	0.02	0.022	101	101			
		1-11-AB-7	-0.02	0.02	0.022	101	101			
		1-12-AB-7	-0.02	0.02	0.022	101	101			
		1-13-AB-7	-0.02	0.02	0.022	101	101			
		1-14-AB-7	-0.02	0.02	0.022	101	101			
		1-15-AB-7	-0.02	0.02	0.022	101	101			
		1-16-AB-7	-0.02	0.02	0.022	101	101			
		1-17-AB-7	-0.02	0.02	0.022	101	101			

¹ Minor discrepancies in residue recoveries (%) between raw data and that reported in this table are a function of round-off differences.

² Results are corrected for control residues, if detected.

TABLE VI RECOVERY DATA FROM CONTROL SUBSTRATES FORTIFIED WITH CGA-245704 AND ANALYZED BY ANALYTICAL METHOD AG-671A (Continued)

Commodity	Substrate	Sample Code	Control (ppm)	CGA-245704 Added (ppm)	CGA-245704 Found (ppm)	Recovery CGA 245704 Corrected (%) ^{1,2}	Mean Recovery (%)	Standard Deviation	Mean CV (%)
Red Pepper	Pepper	14A	-0.02 (0.003)	0.02	0.019	83	87	4.7 (n=8)	5.6
		"	-0.02 (0.003)	0.02	0.019	84	84		
		"	0.20	0.20	0.192	95	79		
		"	1.0	1.0	0.813	81	82		
Hot Pepper	Pepper	14A	-0.02 (0.003)	0.02	0.025	88	114	12 (n=9)	13
		"	-0.02 (0.003)	0.02	0.032	116	83		
		"	0.20	0.20	0.189	97	95		
		"	1.0	1.0	0.811	82	85		
							84		
							86		
							96		

n 142
Overall Range 79-141

¹ Minor discrepancies in residue recoveries (%) between raw data and that reported in this table are a function of round-off differences
² Results are corrected for control residues, if detected

TABLE VII. STATISTICAL RESULTS FROM THE ANALYSIS OF FORTIFIED CONTROLS BASED ON FORTIFICATION LEVEL

Fortification Level (ppm)	Tobacco Green		Tobacco Cured		Hard Lettuce		Loaf Lettuce		Spinach		Cajery		Tomato		Tomato Processing		Bell Pepper		Hot Pepper		Overall			
	mean	std dev	mean	std dev	mean	std dev	mean	std dev	mean	std dev	mean	std dev	mean	std dev	mean	std dev	mean	std dev	mean	std dev	Mean	Standard Dev.	Mean CV	n
0.02	14	1.4	11	1.1	16	1.6	12	1.2	8.8	0.88	9.1	0.91	2.0	0.20	2.8	0.28	2.1	0.21	1.8	0.18	92	14	15	52
	5	0.5	2	0.2	6	0.6	7	0.7	8	0.8	8	0.8	4	0.4	5	0.5	3	0.3	3	0.3				
	81	8.1	10	1.0	24	2.4	5.5	0.55	5.2	0.52	9.1	0.91	9.4	0.94	102	10.2					92	7.6	6.2	28
0.10	12	1.2	3	0.3	4	0.4	4	0.4	5.7	0.57	5.3	0.53	1	0.1	6	0.6								
	87	8.7																						
	81	8.1																			86	8.1	9.4	11
0.50	84	8.4	85	8.5	92	9.2	1.0	1.0	0	0	0	0	0	0	0	0	0	0	0	0	93	9.2	9.9	21
	74	7.4																						
	8.8	0.88	1.1	0.11	3	0.3	1	0.1	2	0.2	2	0.2	3	0.3	5	0.5								
1.0	85	8.5			97	9.7															89	8.9	9.9	19
	1.4	0.14	1.6	0.16	2	0.2	1	0.1	1	0.1	1	0.1	1	0.1	2	0.2								
1.5	89	8.9			84	8.4															93	9.3	8.6	6
	1.4	0.14	1.6	0.16	2	0.2	1	0.1	1	0.1	1	0.1	1	0.1	2	0.2								
2.0	89	8.9			1	0.1															89	8.9	8.6	1
	1.4	0.14	1.6	0.16	2	0.2	1	0.1	1	0.1	1	0.1	1	0.1	2	0.2								
3.0	90	9.0			90	9.0															90	9.0	8.6	1
	1.4	0.14	1.6	0.16	2	0.2	1	0.1	1	0.1	1	0.1	1	0.1	2	0.2								

1. Overall statistics are a compilation of individual data points from Table VI

TABLE VII STATISTICAL RESULTS FROM THE ANALYSIS OF FORTIFIED CONTROLS BASED ON FORTIFICATION LEVEL (Continued)

Fertilization Level (ppm)	Tobacco Green		Tobacco Cured		Head Lettuce		Leaf Lettuce		Spinach		Celery		Tomato Processing		Bell Pepper		Hot Pepper		Overall	
	mean	standard dev	mean	standard dev	mean	standard dev	mean	standard dev	mean	standard dev	mean	standard dev	mean	standard dev	mean	standard dev	mean	standard dev	Mean	CV
4.0	76	1	76	1	76	1	76	1	76	1	76	1	76	1	76	1	76	1	76	1
5.0	76	1	76	1	76	1	76	1	76	1	76	1	76	1	76	1	76	1	76	1
20	86	1	86	1	86	1	86	1	86	1	86	1	86	1	86	1	86	1	86	1
Overall	84	84	88	84	90	90	90	90	89	89	90	90	87	100	84	84	85	85	91	11
mean	8.1	11	11	10	9.3	10	9.3	9.2	8.1	8.1	8.1	12	12	4.7	4.7	12.4	12.4	12.4	12	12
standard dev	11	12	12	11	10	10	9.8	10	9.3	9.3	10	10	9.3	12	5.6	5.6	13	13	13	13
mean CV	17	17	7	7	14	14	14	14	15	15	17	17	15	26	9	9	9	9	9	9

1 Overall statistics are a compilation of individual data points from Table VI

TABLE VIII. EXTRACTABILITY AND ACCOUNTABILITY DATA FROM ANALYSES OF RADIOLABELLED TOBACCO, LETTUCE, AND TOMATO ANALYZED BY ANALYTICAL METHOD AG-671A

Community	Substrate	Sample Code	Extraction Date	Analyst Date	Total Radioactive Residue (TRR) (ppm)	¹⁴ C Liquid Scintillation Analyses ¹				HPLC Analyses ²		Accountability (%)
						Extractability ³ Initial Extract (ppm)	Extractability ³ Final Fraction (%)	Extractability ³ Final Fraction (ppm)	Extractability ³ Final Fraction (%)	Final Fraction ⁴ (ppm Formol)	Final Fraction ⁴ (ppm Corrected)	
Tobacco	Lower Leaves, uncut	TUL-11	9/4/67	9/22/67	1.388	1.500	100	1.094	78	1.088	1.108	80
		TUL-12	9/4/67	9/22/67	1.388	1.744	111	1.095	78	1.045	1.064	77
		TUL-13	9/4/67	9/22/67	1.388	1.478	106	1.000	78	1.048	1.067	77
						mean = 108 std.dev = 2.5 Mean CV = 2.3 n = 3				mean = 108 std.dev = 1.7 Mean CV = 2.2 n = 3		
Head Lettuce	With Wrapper Leaves	L1-11	9/24/67	9/24/67	0.980	0.850	90	0.505	52	0.729	0.346	36
		L1-12	9/24/67	9/24/67	0.880	0.649	88	0.497	52	0.301	0.316	33
		L1-13	9/24/67	9/24/67	0.960	0.784	82	0.439	46	0.270	0.281	30
						mean = 87 std.dev = 4.2 Mean CV = 4.9 n = 3				mean = 87 std.dev = 3.0 Mean CV = 6.1 n = 3		
Tomato	Fruit	TC-11	9/28/67	9/10/67	0.312	0.287	92	0.208	67	0.181	0.186	63
		TC-12	9/28/67	9/10/67	0.312	0.321	103	0.233	74	0.208	0.228	72
		TC-13	9/28/67	9/10/67	0.312	0.317	102	0.241	77	0.222	0.241	77
						mean = 99 std.dev = 6.1 Mean CV = 6.1 n = 3				mean = 99 std.dev = 7.1 Mean CV = 7.1 n = 3		

1 Minor discrepancies in residue recoveries (%) between raw data and that reported in this table are a function of round-off differences.
 2 Total radioactive residues (TRR) as determined by combustion (LSC) analyses in the metabolism study (Reference 5, 6, 7).
 3 ¹⁴C CGA-245704 residues determined by liquid scintillation counting of aliquots of the sample extract (AG-671A Section ILC.2.2.3).
 4 Extractability of initial extract (%) = (Total ¹⁴C CGA-245704 residues in sample initial extract (ppm) ÷ sample TRR (ppm)) x 100
 5 ¹⁴C CGA-245704 residues determined by liquid scintillation counting of aliquots of the sample final fraction (AG-671A Section ILC.3.4.2).
 6 Extractability of final fraction (%) = (¹⁴C CGA-245704 residues in sample final fraction (ppm) ÷ sample TRR (ppm)) x 100
 7 Residue (ppm) determined by High Performance Liquid Chromatography (HPLC).
 8 Residue (ppm) determined by HPLC corrected by mean recovery (%) from fortified controls in respective sample set (Tables III-V).
 9 HPLC Accountability (%) = (Total CGA-245704 residues in sample final fraction (ppm) ÷ sample TRR (ppm)) x 100

TABLE 1X. COMPARISON OF VALIDATION EXTRACTABILITY AND ACCOUNTABILITY AND ACCOUNTABILITY WITH QUANTITATION OF METABOLITE FRACTIONS USING 2-DIMENSIONAL THIN LAYER CHROMATOGRAPHY

Commodity	Substrate	Sample Code	Total Radioactive Residue (TRR) ¹ (ppm)	¹⁴ C Liquid Scintillation Analysis Extractability ² (%)	¹⁴ C Liquid Scintillation Analysis 2-D TLC Quantitation ³ (%)	HPLC Analysis Accountability ⁴ (%)	¹⁴ C Liquid Scintillation Analysis 2-D TLC Accountability ⁵ (%)
Tobacco	Lower Leaves, Un cured	TUL-41	1,368	108	-	80	-
		TUL-42	1,368	111	-	77	-
		TUL-43	1,368	106	-	77	-
		TUL Average 2-D TLC	1,368	108	93	78	67
Head Lettuce	With Wrapper Leaves	LT-4	0.960	80	-	36	-
		LT-2	0.960	85	-	33	-
		LT-3	0.960	82	-	30	-
		LT Average 2-D TLC	0.960	87	98	33	37 ⁶
Tomato	Fruit	TO-11	0.312	92	-	63	-
		TO-12	0.312	103	-	72	-
		TO-13	0.312	102	-	77	-
		TO Average 2-D TLC	0.312	99	97	71	64

¹ Total radioactive residues (TRR) as determined by combustion (SC) analysis in the metabolism study (Reference 5, 6, 7)
² Extractability of initial extract (%) = (Total ¹⁴C-GGA-245704 residues in sample initial extract (ppm) ÷ sample TRR (ppm)) × 100
³ Extractability data obtained in metabolism study (References 5, 6, 7)
⁴ HPLC Accountability (%) = (Total ¹⁴C-GGA 245704 residues in sample final fraction (ppm) ÷ sample TRR (ppm)) × 100
⁵ CGA-210007 quantitation results from metabolism study (References 5 & 7) after treatment with 0.1N NaOH (1 hr @ room temperature) radioactivity after treatment with 0.1N NaOH and cleavage with cellulase (Reference 6) Extractability and accountability residue results reported herein were determined on the unwashed lettuce substrate (combined surface & penetrator radioactivity). The metabolism 2-D TLC surface and penetrator radioactivity results were combined for comparison with the residue accountability results. The TRR of the unwashed lettuce was 0.960 ppm.
⁶ Compare metabolism hydrolysis conditions were not available for the head lettuce. Metabolism 2-D-TLC quantitation for head lettuce was conducted on the penetrator

TABLE X. COMPARATIVE RESULTS OF LC/MS DETECTION VERSUS HPLC/UV DETECTION IN TOBACCO, LETTUCE, AND TOMATO SUBSTRATES

Commodity	Substrate	Sample Type	LC-MS Detection		HPLC - UV Detection	
			CGA-210007, ppm	Recovery %	CGA-210007, ppm	Recovery %
TOBACCO	FLUE-CURED	CTRL	<0.05		<0.10	
		REC + 1.0 PPM		76		83
		REC + 1.0 PPM		105		85
		FREEZER	0.95		0.96	
		FREEZER	0.97		0.94	
			CGA-245704, ppm	Recovery %	CGA-245704, ppm	Recovery %
LETTUCE	HEAD	CTRL	<0.02		<0.02	
		REC + 0.50 PPM		92		95
		REC + 0.50 PPM		94		87
		1X, 0-DAY PHI	0.47		0.44	
	1X, 0-DAY PHI	0.50		0.48		
	HEAD MATURE	CTRL	<0.02		<0.02	
		REC + 0.50 PPM		91		94
		REC + 0.50 PPM		92		91
1X 7-DAY PHI		0.08		0.09		
1X 7-DAY PHI	0.09		0.09			
TOMATOES	FRUIT (UNWASHED)	CTRL	<0.02		<0.02	
		REC + 0.10 PPM		94		98
		REC + 0.50 PPM		101		103
		1X, 14-DAY PHI	0.28		0.27	
		3X 14-DAY PHI	0.64		0.64	
	5X, 14-DAY PHI	0.75		0.71		
	WET POMACE	CTRL	<0.02		<0.02	
		REC + 0.50 PPM		79		74
		REC + 0.50 PPM		91		87
		1X 14-DAY PHI	0.20		0.19	
		3X 14-DAY PHI	0.50		0.48	
	5X 14-DAY PHI	0.75		0.79		
	PASTE	CTRL	0.05		0.03	
		REC + 0.02 PPM		81		93
REC + 0.10 PPM			74		94	
REC + 0.50 PPM			72		85	
REC + 1.0 PPM			84		91	
1X, 14-DAY PHI		1.5		1.3		
3X, 14-DAY PHI		3.9		3.4		
5X, 14-DAY PHI	5.8		5.2			

* Field-incurred residues of CGA-245704 from a representative sample of flue-cured tobacco were not available. Instead, this data shows the results of analysis of flue-cured tobacco that has been lab-fortified with 1.00 ppm of CGA-210007 as part of an in-progress storage stability study.

TABLE XI. OPERATING CONDITIONS FOR THE CONFIRMATION OF
CGA-245704 (AS CGA-210007) BY HIGH PERFORMANCE
LIQUID CHROMATOGRAPHY WITH MASS SPECTROMETRIC
DETECTION

System Description

ESI/APCI quadrupole LC/MS, Platform-LC™ model, equipped with HP-1100 high-performance liquid chromatography inlet system (autosampler, quaternary HPLC pump, column oven, and in-line UV detector).

- Micromass UK Ltd., Wythenshawe, United Kingdom

Liquid Chromatographic Conditions

Analytical HPLC Column: ZORBAX LC-MS SB-C18,
50 mm x 2.1 mm I.D., 5.0 µm,
Part No. 86095.902

Mobile Phase Composition: 80/20 0.5% formic acid in
water/acetonitrile

Mobile Phase Flow Rate: 0.20 mL/min

Typical Retention Time: ~5.0 min

Column Oven Temperature: 25 °C

Volume Injected: 50 µL

Column Effluent/Detector Split ~0.05 mL/min to waste
Ratio: ~0.18 mL/min to detector

MS Detector Setpoints


Ionization Mode	electrospray	Source Heater	90°C
Polarity	negative ion	Low Mass Res.	15.0 amu
Monitor Ion	179 m/z [M-H] ⁻	High Mass Res.	±5.0 amu
Capillary/Needle/Probe	3.6 kV	Multiplier	650 V
Sample Cone	30 V	Analyzer Vacuum	1.5e-5 torr
Skimmer Lens Offset	5 V	N ₂ Drying Gas Flow	410 L/hr.
Ion Energy	1.1 V		

Software Version

Acquisition, Control, and
Quantitation: Masslynx v2.3, Build 007


Operating System: Windows NT, v4.0, Build 1381

FIGURE 1. STRUCTURES AND CHEMICAL NAMES OF CGA-245704 AND CGA-210007



CGA-245704

Benzo[1,2,3]thiadiazole-7-carbothioic acid S-methyl ester
CAS Registry Number 135158-54-2
Molecular Weight = 210.28 g/mol



CGA-210007

Benzo[1,2,3]thiadiazole-7-carboxylic acid
CAS Registry Number 35272-27-6
Molecular Weight = 180.18 g/mol

FIGURE 2. FLOW DIAGRAM FOR ANALYTICAL METHOD AG-671A




FIGURE 3. SCHEMATIC DIAGRAM OF THE HPLC COLUMN SWITCHING SYSTEM AND REPRESENTATIVE CHROMATOGRAM OF STANDARD PROFILE FROM COLUMN 1






FIGURE 3. SCHEMATIC DIAGRAM OF THE HPLC COLUMN SWITCHING SYSTEM AND REPRESENTATIVE CHROMATOGRAM OF STANDARD PROFILE FROM COLUMN 1 (Continued)




Representative chromatogram of CGA-210007 (25 ng/ μ l, 24.904 μ V) standard elution from column 1 (Spherisorb[®] SAX/S3). The valve switch (cuttime) was set at 6.0-6.8 minutes

FIGURE 4. REPRESENTATIVE HPLC CHROMATOGRAMS OF CGA-210007
STANDARDS (¹⁴C-Tomato Analytical Set)




1. Standard, CGA-210007, 0.35 ng, 100 μ V




2. Standard, CGA-210007, 0.50 ng, 177 μ V

FIGURE 4. REPRESENTATIVE HPLC CHROMATOGRAMS OF CGA-210007
STANDARDS (¹⁴C-Tomato Analytical Set) (Continued)




3. Standard, CGA-210007, 2.0 ng, 652 μ V




4. Standard, CGA-210007, 5 ng, 1544 μ V

FIGURE 4. REPRESENTATIVE HPLC CHROMATOGRAMS OF CGA-210007 STANDARDS (¹⁴C-Tomato Analytical Set) (Continued)




5. Standard, CGA-210007, 10 ng, 3223 μ V




6. Standard, CGA-210007, 10 ng, 3147 μ V

FIGURE 4. REPRESENTATIVE HPLC CHROMATOGRAMS OF CGA-210007
STANDARDS (¹⁴C-Tomato Analytical Set) (Continued)



7 Standard, CGA-210007, 25 ng, 7877 μ V


FIGURE 5. REPRESENTATIVE CALIBRATION CURVE USING THE CGA-210007 STANDARDS FROM FIGURE 4




CGA-210007 Calibration Data

ng injected	peak height, µV
0.35	100
0.50	177
2.0	652
5.0	1544
10	3223
10	3147
25	7877

FIGURE 6. REAGENT BLANK CHROMATOGRAMS




1. Reagent Blank, Tobacco Analytical Set 91-96-BBR,
Protocol 91-96




2. Reagent Blank, Head Lettuce Analytical Set 92-96-BBG,
Protocol 92-96

FIGURE 6 REAGENT BLANK CHROMATOGRAMS (Continued)




3. Reagent Blank, Leaf Lettuce Analytical Set 92-96-BA0,
Protocol 92-96



4. Reagent Blank, Spinach Analytical Set 92-96-BAG2,
Protocol 92-96


FIGURE 6. REAGENT BLANK CHROMATOGRAMS (Continued)

Sample Name reagent blank
Data File C:\TC4\NOVA\NOV7019A.RAW Date 11/6/97 01:30 AM
Sequence File: C:\TC4\NOVA\NOV1104.SEG




CGA-210007 - TOMATOES - METHOD VALIDATION

5. Reagent Blank, Whole Tomato, Ricerca, Inc. Method Trial Data, Protocol 93-96




5. Reagent Blank, Tomato Pomace Analytical Set 93-96-BA2, Protocol 93-96

FIGURE 6. REAGENT BLANK CHROMATOGRAMS (Continued)




7. Reagent Blank, Tomato Pomace Analytical Set 93-96-BH,
Protocol 93-96

FIGURE 7. HPLC CHROMATOGRAMS FROM THE ANALYSIS OF ¹⁴C-TOBACCO (Control, Procedural Recovery, and Radiolabeled Samples)




1. TUL-C1, control, 50 000 mg injected, <0.35 ng CGA-210007 quantitated, <0.02 ppm CGA-245704 determined




2. TUL-R1, control + 0.02 ppm CGA-245704, 50,000 mg injected, 0.904 ng CGA-210007 quantitated, 0.021 ppm CGA-245704 determined, 106% CGA-245704 recovered

FIGURE 7. HPLC CHROMATOGRAMS FROM THE ANALYSIS OF
¹⁴C-TOBACCO (Control, Procedural Recovery, and
Radiolabeled Samples) (Continued)




3. TUL-R2, control + 1.5 ppm CGA-245704, 10.000 mg injected,
11.624 ng CGA-210007 quantitated, 1.360 ppm CGA-245704
determined, 91% CGA-245704 recovered




4. TUL-I1, radiolabeled incurred substrate, 13.178 mg injected,
12.253 ng CGA-210007 quantitated, 1.088 ppm CGA-245704
determined, 1.1 ppm CGA-245704 (corrected for average
recovery) determined

FIGURE 7. HPLC CHROMATOGRAMS FROM THE ANALYSIS OF ¹⁴C-TOBACCO (Control, Procedural Recovery, and Radiolabeled Samples) (Continued)




5. TUL-I2, radiolabeled incurred substrate, 13.178 mg injected, 11 770 ng CGA-210007 quantitated, 1 045 ppm CGA-245704 determined, 1.1 ppm CGA-245704 (corrected for average recovery) determined




6. TUL-I3, radiolabeled incurred substrate, 13.178 mg injected, 11.801 ng CGA-210007 quantitated, 1.048 ppm CGA-245704 determined, 1.1 ppm CGA-245704 (corrected for average recovery) determined

FIGURE 8. HPLC CHROMATOGRAMS FROM THE ANALYSIS OF ¹⁴C-LETTUCE (Control, Procedural Recovery, and Radiolabeled Samples)




1. LT-C1, control, 50.000 mg injected, <0.35 ng CGA-210007 quantitated, <0.02 ppm CGA-245704 determined




2. LT-R1, control + 0.02 ppm CGA-245704, 50.000 mg injected, 0.801 ng CGA-210007 quantitated, 0.019 ppm CGA-245704 determined, 94% CGA-245704 recovered

FIGURE 8. HPLC CHROMATOGRAMS FROM THE ANALYSIS OF ¹⁴C-LETTUCE (Control, Procedural Recovery, and Radiolabeled Samples) (Continued)




3. LT-R2, control + 1.0 ppm CGA-245704, 20.000 mg injected, 16.542 ng CGA-210007 quantitated, 0.968 ppm CGA-245704 determined, 97% CGA-245704 recovered




4. LT-I1, radiolabeled incurred substrate, 19.770 mg injected, 5.565 ng CGA-210007 quantitated, 0.329 ppm CGA-245704 determined, 0.346 ppm CGA-245704 (corrected for average recovery) determined

FIGURE 8. HPLC CHROMATOGRAMS FROM THE ANALYSIS OF ¹⁴C-LETTUCE (Control, Procedural Recovery, and Radiolabeled Samples) (Continued)




5. LT-12, radiolabeled incurred substrate, 19.770 mg injected, 5 085 ng CGA-210007 quantitated, 0.301 ppm CGA-245704 determined, 0.316 ppm CGA-245704 (corrected for average recovery) determined




6. LT-13, radiolabeled incurred substrate, 19.770 mg injected, 4 562 ng CGA-210007 quantitated, 0.270 ppm CGA-245704 determined, 0.284 ppm CGA-245704 (corrected for average recovery) determined

FIGURE 9. HPLC CHROMATOGRAMS FROM THE ANALYSIS OF ¹⁴C-TOMATO
(Control, Procedural Recovery, and Radiolabeled Samples)




1. TO-C1, control, 50 000 mg injected, <0.35 ng (0.147 ng) CGA-210007 quantitated, <0.02 ppm (0.003) CGA-245704 determined




2. TO-R1, control + 0.02 ppm CGA-245704, 50.000 mg injected, 0.970 ng CGA-210007 quantitated, 0.023 ppm CGA-245704 determined, 0.019 ppm CGA-245704 (corrected for control) determined, 96% CGA-245704 recovered

FIGURE 9. HPLC CHROMATOGRAMS FROM THE ANALYSIS OF ¹⁴C-TOMATO (Control, Procedural Recovery, and Radiolabeled Samples) (Continued)




3. TO-R2, control - 0.50 ppm CGA-245704, 50.000 mg injected, 18.995 ng CGA-210007 quantitated, 0.444 ppm CGA-245704 determined, 0.441 ppm CGA-245704 (corrected for control) determined, 88% CGA-245704 recovered




4. TO-I1, radiolabeled incurred substrate, 49.419 mg injected, 7.650 mg CGA-210007 quantitated, 0.181 ppm CGA-245704 determined, 0.196 ppm CGA-245704 (corrected for average recovery) determined

FIGURE 9. HPLC CHROMATOGRAMS FROM THE ANALYSIS OF ¹⁴C-TOMATO
(Control, Procedural Recovery, and Radiolabeled
Samples) (Continued)




- 5 TO-12, radiolabeled incurred substrate, 49 419 mg injected,
8 791 ng CGA-210007 quantitated, 0.208 ppm CGA-245704
determined, 0 226 ppm CGA-245704 (corrected for average
recovery) determined




- 6 TO-13, radiolabeled incurred substrate, 49 419 mg injected,
9 377 ng CGA-210007 quantitated, 0 222 ppm CGA-245704
determined, 0.241 ppm CGA-245704 (corrected for average
recovery) determined

FIGURE 10. REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF FLUE-CURED TOBACCO




1. Calibration Standard, 2.00 ng injected.




2. 3-TF-21-C, control, 10.0 mg injected, <0.50 ng CGA-210007 quantitated, <0.05 ppm CGA-210007 determined.

* Samples of flue-cured tobacco analyzed during development of the LC/MS confirmatory procedure were obtained from an in-progress storage stability study and were fortified with CGA-210007 rather than CGA-2-5704

FIGURE 10 REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF FLUE-CURED TOBACCO
(Continued)




3. 3-TF-21-R-A, control - 1.00 ppm CGA-210007, 10.0 mg injected, 7.593 ng CGA-210007 quantitated, 0.759 ppm CGA-210007 determined, 76% CGA-210007 recovered.




4. 3-TF-21-F-A, 1.00 PPM 2-Month freezer-stored sample, 10.0 mg injected, 8.599 ng CGA-210007 quantitated, 0.860 ppm CGA-245704, 0.95 ppm CGA-210007 determined (corrected for 91% average CGA-210007 recovery)

* Samples of flue-cured tobacco analyzed during development of the GC/MS confirmatory procedure were obtained from an in-progress storage stability study and were fortified with CGA-210007 rather than CGA-245704

FIGURE 11. REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF HEAD LETTUCE (0-Day PHI)




1. Calibration Standard, 10.00 ng injected.




2 1-2-A, control, 50.0 mg injected, <0.35 ng CGA-210007 quantitated, <0.02 ppm CGA-245704 equivalent determined.

FIGURE 11. REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF HEAD LETTUCE (0-Day PHI) (Continued)




3. 1-2-A, control + 0.50 ppm CGA-245704, 20.0 mg injected, 7.841 ng CGA-210007 quantitated, 0.459 ppm CGA-245704 equivalent determined, 92% CGA-245704 recovered.




4. 2-2-A, 1X treated head lettuce, 20.0 mg injected, 7.491 ng CGA-210007 quantitated, 0.438 ppm CGA-245704, 0.47 ppm CGA-245704 equivalent determined (corrected for 93% average recovery)

FIGURE 12. REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF UNWASHED TOMATOES




1. Calibration Standard, 5.00 ng injected.




2. 1-8-AB-1, control, 50.0 mg injected, <0.50 ng CGA-210007 quantitated, <0.02 ppm CGA-2457C4 equivalent determined

FIGURE 12. REPRESENTATIVE ION CHROMATOGRAMS FROM THE
CONFIRMATORY ANALYSIS OF UNWASHED TOMATOES
(Continued)




3. 1-8-AB-1, control + 0.10 ppm CGA-245704, 50.0 mg injected,
4.034 ng CGA-210007 quantitated, 0.094 ppm CGA-245704
equivalent determined, 94% CGA-245704 recovered.




4. 2-8-AB-1, 1X treated unwashed tomato, 40.0 mg injected,
9.178 ng CGA-210007 quantitated, 0.268 ppm CGA-245704, 0.28
ppm CGA-245704 equivalent determined (corrected for 98%
average recovery)

FIGURE 13. REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF TOMATO WET POMACE




1. Calibration Standard, 10.0 ng injected.




2. 1-6-AB-4, control, 50.0 mg injected, <0.50 ng (0.243 ng) CGA-210007 quantitated, <0.02 ppm (0.006 ppm) CGA-245704 equivalent determined.

FIGURE 13. REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF TOMATO WET POMACE
(Continued)




3. 1-6-AB-4, control + 0.50 ppm CGA-245704, 25.0 mg injected, 8.561 ng CGA-210007 quantitated, 0.395 ppm CGA-245704 equivalent determined (corrected for control residue), 79% CGA-245704 recovered




4. 2-6-AB-4, 1X treated tomato wet pomace, 50.0 mg injected, 7.377 ng CGA-210007 quantitated, 0.173 ppm CGA-245704, 0.20 ppm CGA-245704 equivalent determined (corrected for 85% average recovery)

FIGURE 14. REPRESENTATIVE ION CHROMATOGRAMS FROM THE CONFIRMATORY ANALYSIS OF TOMATO PASTE




1. Calibration Standard, 50 ng injected.




2. 1-8-AB-7, control, 50.0 mg injected, 1.649 ng CGA-210007 quantitated, 0.039 ppm CGA-245704, 0.05 ppm CGA-245704 equivalent determined (corrected for 78% average recovery)

FIGURE 14. REPRESENTATIVE ION CHROMATOGRAMS FROM THE
CONFIRMATORY ANALYSIS OF TOMATO PASTE
(Continued)



- 3 1-8-AB-7, control + 0.02 ppm CGA-245704, 50.0 mg injected, 2.342 ng CGA-210007 quantitated, 0.055 ppm CGA-245704 equivalent determined, 0.016 ppm CGA-245704 (corrected for control), 81% CGA-245704 recovered.



4. 2-8-AB-7, 1X treated tomato paste, 10.0 mg injected, 9.640 ng CGA-210007 quantitated, 1.128 ppm CGA-245704, 1.5 ppm CGA-245704 equivalent determined (corrected for 78% average recovery)

VII. REFERENCES

1. Formica, G., REM 172.02, "Determination of Total Residues of CGA-245704 as CGA-210007 by High Performance Liquid Chromatography (HPLC)," Cereals (Grain, Straw, and Green Aerial Part), December 6, 1993.
2. Formica, G., REM 172.11, "Determination of Total Residues of CGA-245704 as CGA-210007 by High Performance Liquid Chromatography," Plant Material and Milk, October 20, 1995.
3. Moore, M., "Preparation of Crop Samples for Residue Analysis," Product Safety Group SOP Number 7.21, Revision 2, Novartis Crop Protection, Inc., March 5, 1997.
4. Hayworth, C. and Grunenwald, M., ABR-97117, Final Report Amendment Number 1, "Validation of Draft Analytical Method AG-671A for the Determination of Total Residues of CGA-245704 as CGA-210007 in Tobacco, Leafy Vegetables, and Fruiting Vegetables by Column Switching High Performance Liquid Chromatography."
5. Nicollier, G., CIBA GEIGY Ltd., Basel, Switzerland, Plant Metabolism, "Behavior and Metabolism of CGA 245704 in Greenhouse Grown Tobacco after Foliar Spray Application of [U-¹⁴C-Phenyl]-CGA 245704," Study 94WA04.
6. Kiffe, M., CIBA GEIGY Ltd., Basel, Switzerland, Plant Metabolism, "Metabolism of CGA 245704 in Greenhouse Grown Lettuce after Treatment with [Phenyl-U-¹⁴C] Labelled Material, Study 96MK02.1," In-progress.
7. Nicollier, G., CIBA GEIGY Ltd., Basel, Switzerland, Plant Metabolism, "Metabolism of CGA 245704 in Greenhouse Grown Tomatoes after Treatment [Phenyl-U-¹⁴C] Labeled Material," Study 94WA05.
8. Hayworth, C., Protocol 91-96, "CGA-245704 - Magnitude of the Residues in or on Tobacco," Novartis Crop Protection, Inc., April 17, 1996

VII. REFERENCES (Continued)

9. Hayworth, C , Protocol 92-96, "CGA-245704 - Magnitude of the Residues in or on Crop Group 4: Leafy Vegetables," Novartis Crop Protection, Inc , February 16, 1996.
10. Hayworth, C., Protocol 93-96, "CGA-245704 - Magnitude of the Residues in or on Crop Group 8: Fruiting Vegetables," Novartis Crop Protection, Inc , March 8, 1996.

APPENDIX I

SEPARATE DOCUMENTS ACCOMPANYING THIS REPORT

Formica, G., REM 172.02, "Determination of Total Residues of CGA-245704 as CGA-210007 by High Performance Liquid Chromatography (HPLC)," Cereals (Grain, Straw, and Green Aerial Part), December 6, 1993

Formica, G., REM 172.11, "Determination of Total Residues of CGA-245704 as CGA-210007 by High Performance Liquid Chromatography," Plant Material and Milk, October 20, 1995

SUBMITTER/SPONSOR: Novartis Crop Protection, Inc., P. O. Box 18300, Greensboro, NC 27419-8300


CGA 245704	RESIDUE METHOD VALIDATED	REM 172.02
DETERMINATION OF TOTAL RESIDUES OF CGA 245704 AS CGA 210007 BY HIGH PER- FORMANCE LIQUID CHROMATOGRAPHY (HPLC)	CEREALS (GRAIN, STRAW, AND GREEN AERIAL PART	DEC. 06, 1993 PP 2.53/Gi

TABLE OF CONTENTS


1. INTRODUCTION	3. RESULTS AND DISCUSSION
1.1 Scope of the Method	3.1 Accuracy and Precision
1.2 Principle of the Method	3.2 Limit of Determination
2. MATERIALS AND METHODS	3.3 Ruggedness
2.1 Equipment	3.4 Limitations
2.2 Reagents and Standards	3.5 Method Validation
2.3 Analytical Procedure	3.6 Correction for the Increase of the Effective Extraction Volume by Soluble Sample Constituents.
2.4 Instrumentation	4. CERTIFICATION
2.5 Interferences	5. REFERENCES
2.6 Confirmatory Techniques	6. FIGURES AND TABLES
2.7 Time Required for Analysis	6.1 Figures
2.8 Modifications and Potential Problems	6.2 Representative Chromatograms
	6.3 Recoveries

CHEMICAL STRUCTURES

CGA 245704



CGA 210007



Proprietary Information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

1. INTRODUCTION

1.1 Scope of the Method

The method allows the quantitative determination of the combined residues of CGA 245704 and its metabolites which contain the CGA 210007 moiety in cereal (grain, straw and green aerial part or abbreviated simply "plant"). (See section 6.1, Figure 1 for structures and chemical names).

The lower practical level of quantitation by this method is 0.019 mg/kg for grain and 0.047 mg/kg for straw and plant calculated as CGA 245704 equivalents.

1.2 Principle of the Method

CGA 245704 and possible metabolites containing the CGA 210007 moiety are hydrolyzed to CGA 210007 by heating a homogenized subsample in the presence of sodium hydroxide. Methanol is added and the mixture shaken. An aliquot of the extract is cleaned up by acidic/alkaline/acidic partition. The reextract is cleaned up by a silica column. Additional clean up and final determination are performed by a two column HPLC switching system with UV-detection.

2. MATERIALS AND METHODS

Standard laboratory equipment is not listed. All equipment and chemicals mentioned herein can be substituted by suitable products of any origin. Prove suitability of reagents by analyzing reagent blanks.

2.1 Equipment

- 2.1.1 Rotating evaporator, Büchi, Rotavapor RE, Büchi AG, Flawil, CH.
- 2.1.2 Vacuum system, Terno Duo 500 S, Terno AG, Köblis, CH or commercial water aspirator.
- 2.1.3 Circulation cooler, mgw Lauda WK 450, Meßgeräte-Werke-Lauda, Dr. R. Wobser KG, Lauda-Königshofen, FRG.
- 2.1.4 Lab-shaker, A. Kühner AG, 4000 Basel, CH.
- 2.1.5 Folded filter paper, 15 cm diameter, Macherey-Nagel, 5160 Düren, FRG, Cat. No. MN 713 1/4.
- 2.1.6 Glass wool.
- 2.1.7 Ultrasonic bath, Bransonic 220, Branson Cleaning Equipment Co., Parrot Drive, Shelton, CT 06484-0768, USA.
- 2.1.8 Vortex mixer, Model Vortex-Genie, Scientific Industries Inc., Springfield, Massachusetts, USA.
- 2.1.9 Empty borosilicate glass cartridge, 8 ml volume, J. T. Baker Inc., Phillipsburg, N.J. 08865, USA, represented in Switzerland by P. H. Stehelin and Cie AG, 4003 Basel, Cat. No. 7328-06.
- 2.1.10 PTFE frits for 8 ml borosilicate glass cartridge, J. T. Baker Inc., Cat. No. 7329-06.
- 2.1.11 Adapter for 8 ml glass cartridge, J. T. Baker Inc., Cat. No. 4528.
- 2.1.12 Glass reservoirs built in-house of about 25 ml content to be attached to the cartridges or equivalent (e.g. 25 ml Bond-Elut reservoirs, Varian, Harbor City CA, USA, Cat. No. A-1213-1011).
- 2.1.13 Vacuum manifold to accommodate solid-phase extraction cartridges, built in-house or commercial equivalent, e.g. VISIPREP, Supelco Inc., Bellefonte, Pennsylvania, USA, represented in Switzerland by Supelco SA, Gland, CH, Cat. No. 5-7030.
- 2.1.14 High Performance Liquid Chromatograph: refer to section 2.4.

Proprietary information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG.

2.2 Reagents and Standard

Main suppliers' addresses: E. Merck AG, 6100 Darmstadt, FRG.
Fluka Chemie AG, 9470 Buchs, CH.
J. T. Baker Inc., Phillipsburg, N.J. 08864, USA.

- 2.2.1 Water, demineralized and HPLC-grade, prepared in house, Ciba or HPLC-grade, J. T. Baker Inc., Cat. No. 4218.
- 2.2.2 Methanol for extraction, distilled in house, Ciba or analytical grade, E. Merck, Cat. No. 6009.
- 2.2.3 n-Hexane for residue analysis, Merck, Cat. No. 4371.
- 2.2.4 tert.-Butyl methyl ether (TBME), puriss., analytical grade, Fluka, Cat. No. 20249.
- 2.2.5 Methanol, LiChrosolv, HPLC-grade, Merck, Cat. No. 6018.
- 2.2.6 Acetonitrile, LiChrosolv, chromatography grade, Merck, Cat. No. 14291.
- 2.2.7 Hydrochloric acid 1 M, Ciba, Cat. No. 24 or J. T. Baker Inc., Cat. No. 7088 or about 1 M prepared by diluting of an analytical grade concentrated hydrochloric acid with HPLC-water.
- 2.2.8 ortho-Phosphoric acid, 85%, analytical grade, E. Merck, Cat. No. 573. Prepare a solution containing about 0.02 mole of phosphoric acid per liter HPLC water by diluting 2.3 g concentrated acid to 1 liter with water for HPLC.
- 2.2.9 Sodium hydroxide 1 M, Ciba, Internal Cat. No. 17 or J. T. Baker Inc., Cat. No. 7097.
- 2.2.10 Silica Gel for flash chromatography, J. T. Baker Inc., Cat. No. 7024-01
- 2.2.11 Calcium chloride dihydrate, analytical grade, Merck, Cat. No. 2382.
- 2.2.12 Sodium chloride, analytical grade, Merck, Cat. No. 6404. Prepare a saturated solution of sodium chloride in water for HPLC (brine).
- 2.2.13 Ammonium chloride, analytical grade, Merck, Cat. No. 1145. Prepare a 0.05 M solution of ammonium chloride by dissolving 2.7 g salt in one liter of 0.02 M phosphoric acid (see above step 2.2.8).
- 2.2.14 CGA 245704 reference substance for recovery experiments. Prepare a stock solution containing 200 µg CGA 245704/mL acetonitrile LiChrosolv.
- 2.2.15 CGA 210007 reference substance for standardization. Prepare a stock solution containing 200 µg CGA 210007/mL acetonitrile LiChrosolv.

2.3 Analytical Procedure

2.3.1 Preparation of samples and subsamples

Prepare representative laboratory samples from the field samples (e.g. mix grain with dry ice and mill in a cross beater mill. Mix straw and plant with dry ice and chop using a suitable cutter). Analyze the samples immediately after preparation or store at about -20 °C until analysis.

For analysis weigh subsamples of e.g. 5 g of grain or e.g. 2 g of straw or plant into a 250 mL wide mouth jar. Avoid thawing of the frozen laboratory samples but preferably scratch off or drill out the subsample needed for analysis.

2.3.2 Fortification

To regularly check the performance of the method, analyze also at least two fortified control samples with each series of analyses. To prepare these samples, add known amounts of CGA 245704 to control samples prior to extraction.

Proprietary information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

Select fortification levels to be one and ten times the lower practical level of determination or in the range of the expected residue levels. Make sure that control samples neither are contaminated nor show interfering signals.

To fortify samples with 0.02 and 0.2 mg/kg for grain or 0.05 and 0.5 mg/kg for straw and plant, prepare solutions of CGA 245704 containing 0.1 and 1.0 µg per mL acetonitrile LiChrosolv by appropriate dilution of the stock solution (cf section 2.2.14). To prepare fortified samples, add 1 mL of one of the solutions to an untreated subsample. Proceed as described in section 2.3.3.

2.3.3 Hydrolysis/Extraction

Add as much water to the weighed subsample (cf. section 2.3.1) as to achieve a total water volume of 18 mL taking into account the crop solubles content being 1 mL for grain, 0.4 mL for straw and 1.6 mL for plant (see section 3.6 for moisture content). Add 2 mL of sodium hydroxide 1 M, seal the jar tightly, swirl the mixture and allow the jar to stand for about 30 minutes in an oven preheated to 50 - 55 °C. Swirl the jar occasionally. Allow the mixture to cool, add 80 mL (79 mL to fortified samples) of distilled methanol and 2 g of calcium chloride (the addition of the salt is necessary for grain and optionally for plant and straw. The salt facilitates the following filtration and reduces the formation of emulsion during the acidic partition, section 2.3.4). Shake the tightly sealed jar for about 30 minutes. Total volume of the extract is 100 mL. Filter immediately about 50 mL extract through a folded filter paper of 15 cm diameter into a 100 mL Erlenmeyer flask and continue according to section 2.3.4.

Remark: To avoid formation of emulsion during the following partition step, especially for wheat grain samples, filtration and subsequent partition cleanup (section 2.3.4) have to be performed without interruption immediately after shaking with methanol. If this is not possible, the extraction mixtures are to be immediately stored at about -18 °C until filtration and partition can be done. After aliquoting for analysis, the remaining filtrates are to be immediately stored at about -18 °C, as reserve for potential further analysis.

2.3.4 Cleanup by Acidic Partition

Transfer 20 mL of filtrate, corresponding to 1 g grain or 0.4 g straw and plant material to a 250 mL separatory funnel. Add 50 mL water, 40 mL saturated sodium chloride solution and 8 mL 1 M hydrochloric acid. Shake the mixture to homogenize. Add 50 mL of a solvent mixture of hexane + tert.-butyl methyl ether 7 vol. + 3 vol. solvent mixture A). Shake and allow phases to separate. Drain the aqueous phase into a 250 mL Erlenmeyer flask. Filter the organic phase including emulsion, if present, into a clean 250 mL separatory funnel through a funnel plugged up with a loose plug of glass wool. Transfer the aqueous phase back to the separatory funnel. Add 50 mL of the solvent mixture A and shake. Allow phases to separate. Drain and discard the aqueous phase. Filter the organic phase into the separatory funnel containing the first reextract through the funnel plugged with glass wool. Observe if water is present below the organic phase. In this case, drain and discard the water.

Proprietary information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

2.3.5 Cleanup by Alkaline/Acidic Partition

Add 20.0 mL of 1 M sodium hydroxide to the organic phase of the acidic partition clean up (section 2.3.4), stopper the separatory funnel and shake well. Allow phases to separate. Drain 15.0 mL of the alkaline phase, corresponding to 0.75 g of grain or to 0.3 g plant/straw, into a graduated 25 mL test tube and discard the rest. Add 2 mL of 85% phosphoric acid and 10 mL of the solvent mixture A (see section 2.3.4). Shake and allow phases to separate. Transfer the organic phase into a 100 mL round bottom flask by means of a Pasteur pipette. Repeat the partition two times more with portions of 10 mL of the same solvent mixture, each. Transfer the organic phase each time into the round bottom flask. Discard the aqueous phase. Evaporate the combined organic phases to dryness under reduced pressure using a rotating evaporator (water bath temp.: about 50 °C).

2.3.6 Cleanup by Silica Cartridge

Remark: To pack the cartridges used for this cleanup proceed as follows:

Insert a PTFE frit on the bottom of a borosilicate glass tube of 8 mL content, fill 1 g of silica gel into the tube and place a second frit on the top of the filling. Vibrate for ~1 minute for the material to settle, using a Vortex mixer and press the upper frit on the top of the adsorbent. Add 2 g of anhydrous sodium sulfate. The cartridge is ready for use.

Attach the cartridge packed as described above, to the vacuum manifold. Tightly connect a 25 mL glass reservoir to the top of the cartridge. Precondition the cartridge by passing through 5 mL of a solvent mixture of TBME + methanol LiChrosolv 85 vol. + 15 vol. by suction.

Dissolve the residue of the alkaline/acidic partition (section 2.3.5) in 5 mL of a solvent mixture of TBME + methanol LiChrosolv 85 vol. + 15 vol. and transfer the solution to the reservoir. Pass the mixture through the cartridge at a rate of 5 - 7 mL per minute by suction. Discard the eluate. Wash the flask with 5 mL of the same solvent mixture, transfer the wash to the reservoir and pass it through the cartridge. Discard the eluate. Elute CGA 210007 with 10 mL of a solvent mixture of TBME + methanol LiChrosolv 6 vol. + 4 vol. into a 25 mL round bottom flask by suction.

2.3.7 Preparation of the Final Solution

Evaporate the eluate of the silica cartridge cleanup (section 2.3.6) to dryness under reduced pressure using a rotating evaporator (water bath temp.: about 50 °C). Dissolve the residue in 3 mL of 0.02 M phosphoric acid in HPLC water using an ultra sonic bath (final solution).

2.4 Instrumentation

2.4.1 High Performance Liquid Chromatographic Systems (HPLC)

For the determination of CGA 210007 use a HPLC two-column switching system with UV-detector, pumps, autosampler-injector, columns and conditions as follows:

Proprietary information of CIBA-GEIGY AG.
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG.

Detector: Spectroflow 783 (ABI Analytical Kratos Division, Ramsey, NJ07446, USA) or equivalent.

Pumps: Two Shimadzu Solvent Delivery Modules LC-9A, (Shimadzu Corporation, Analytical Instruments Division, Kyoto, Japan) or equivalent.

Valves: Multiport Streamswitch "MUST" equipped with two low dead volume valves for column switching (Spark Holland, 7800 AJ Emmen, Holland) in connection with an automatic injection unit PROMIS II (Spark Holland) or equivalent.

Recorder: Strip chart SE 120, dual channel (ABB, Goerz Metrawatt, 1101 Vienna, A) or equivalent.

Column oven: Modell STH 585. Temperature range: 5 - 80 °C (Gynkotec GmbH) or equivalent.

REMARK: Column oven is absolutely necessary in non thermostated laboratories to ensure constant elution times and reproducible cut of the a.i. from column to column (see below "Determination of Switching Interval").

Optional (for system automation):

Control and data collection unit: HP 3350A Laboratory Automation System (Hewlett-Packard, Palo Alto, CA 94304, USA).

Column 1: Stainless steel tube, 10 cm length, 2 mm i.d., packed with Nucleosil SB (strongly basic anion exchanger), particle size 5 µm (Dr. H. Knauer KG, 6370 Oberursel, FRG, Cat. No.: B45 - Y91) or equivalent.

Remark: A new column has to be conditioned with the mobile phase 1 (see below) for about 8 hours or overnight at a flow rate of 0.2 mL/min to obtain reproducible chromatograms.

Column 2: Stainless steel tube, 10 cm length, 2 mm i.d., packed with Nucleosil 100 C₁₈, particle size 5 µm (Dr. H. Knauer KG, Cat. No.: B45 - Y76) or equivalent.

For the HPLC system, connect pumps, injector, columns, valves and detector according to Figure 3 (section 6.1).

Mobile phase 1: 0.05 M ammonium chloride in 0.02 M phosphoric acid + acetonitrile Lichrosolv 85 vol. + 15 vol.

Mobile phase 2: 0.02 M phosphoric acid in Water for HPLC + acetonitrile Lichrosolv 76 vol. + 24 vol.

Flow rate 1: 0.2 mL/min

Retention time: about 13 min (column 1)
 Flow rate 2: 0.2 ml/min
 Retention time: about 22 min (columns 1 and 2)
 Injection vol.: 50 µL
 Recorder: 10 mV full scale
 Column oven: 23 °C
 Chart speed: 0.5 - 1.0 cm/min
 Detector: Wave length: 235 nm
 sensitivity: 0.005 aufs

Determination of Switching Interval

Determine actual switching times each time a series of samples is to be quantitated. Connect the outlet of column 1 directly to the UV-detector. Inject 50 µL of a 0.08 µg CGA 210007/mL standard solution (see below section 2.4.2) and measure the begin and the end times of the peak of CGA 210007. Thereafter, connect column 2 to the detector and switch that portion of eluate from column 1 which contains CGA 210007, e.g. portion between 12 min 00 s and 13 min 24 s to column 2 (see section 6.2.7, chromatogram A for an example of measurement of switching times).

2.4.2 Preparation of Standard Solution

Standardize the chromatographic system each time a series of samples is to be quantitated. The range of the concentrations is depending on the range of residues to be determined, in particular, the lowest standard concentration is depending on the lower practical level.

Calculate the lowest standard concentration (C) as follows:

$$C = \frac{L \cdot F \cdot A}{V_f} \quad [\mu\text{g/mL}]$$

L = lower practical level [mg/kg] or [µg/g]
 F = conversion factor from CGA 245704 to CGA 210007
 A = aliquot of crop cleaned up through silica [g]
 V_f = volume of the final solution [mL]

With the values proposed in this method, the lowest standard concentration is 0.004 µg CGA 210007/mL as presented below:

$$C = \frac{0.019 \times 0.857 \times 0.75}{3} = 0.004 \mu\text{g/mL} \quad \text{for grain and}$$

$$C = \frac{0.047 \times 0.857 \times 0.30}{3} = 0.004 \mu\text{g/mL} \quad \text{for plant/straw}$$

Prepare at least four standard solutions of different concentrations by appropriately diluting the stock solution of CGA 210007 (section 2.2.15) with 0.02 M phosphoric acid in water for HPLC.

Select the concentrations as required; typical values are 0.08, 0.04, 0.01 and 0.004 µg/mL.

The conversion factor from CGA 210007 to CGA 245704 is: 1.17

2.4.3 Quantitation of Residues

As the detector response may change on injection of coextractives, it is recommended to inject standards and final solutions alternately. Inject 50 µL of standard solutions (corresponding to 4.0, 2.0, 0.5 and 0.2 ng CGA 210007, respectively) and final solutions. Measure the response of the analyte at the characteristic retention time and calculate response function and residues as detailed in REM 119.04 (Ciba, Plant Protection Division, Residue Analysis, 1991).

2.5 Interferences

No interferences which disturbed the determination of the CGA 210007 were encountered so far.

2.6 Confirmatory Techniques

No specific conditions have been worked out so far.

2.7 Time Required for Analysis

A total of about 15 hours is required to work up a set of 12 samples to the point of HPLC injection. Automated HPLC analysis can be performed overnight.

2.8 Modifications and Potential Problems

- a) The filtration of the extraction mixtures (section 2.3.3) may be slow if calcium chloride is not added. Furthermore, the filtrate can be turbid.
- b) Extracts and filtrates (see section 2.3.3) must not remain at room temperature before partition, otherwise an emulsion appears during the partition cleanup.

3. RESULTS AND DISCUSSION

3.1 Accuracy and Precision

The overall accuracy and precision of the method is: (status as of Nov. 1993)

$\bar{X} = 91\%$, $s_{\text{abs}} = 6.6\%$, $n = 50$.

Find typical recovery percentages in section 6.3.

3.2 Limit of Determination

The ultimate limit of determination of this method was not determined. The lower practical level of CGA 245704 quantitated by this method is: 0.019 mg/kg for grain and 0.047 mg/kg for plant and straw.

Proprietary information of CIBA-GEIGY AG.
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

3.3 Ruggedness

The method was tested by three technicians using fortified samples. The analyses were performed on different days with different instrumentations.

Results obtained were comparable with regard to recovery percentages and chromatographic profile.

3.4 Limitations

None known, so far.

3.5 Method Validation

The performance of the method has been checked with recovery experiments, section 6.3.

3.6 Correction for the Increase of the Effective Extraction Volume by Soluble Sample Constituents

The soluble part of the sample (i.e. water and low molecular plant constituents) has to be considered for correct results. Exact values for this contribution are not easily available. In view of the relatively high imprecision of residue analytical results, rough estimates of the volume of the soluble portion of the samples are acceptable, as long as the total volume is known with reasonable accuracy.

In view of an overall precision of $\pm 2\%$ for the total volume of the extract, it is sufficient to know the volume contribution of the solubles by $\pm 20\%$, provided that the ratio sample size/solvent added is 0.1, e.g. sample size 10 g, solvent added 100 mL (specific weight of the soluble = 1 is supposed, the volume contraction on mixing is neglected).

For the calculation procedure of this method, the unknown exact volume of the solubles are replaced by rough estimates:

The various crop materials are classified into three groups. The volumes of soluble constituents are estimated using group factors (G) of 0.8, 0.5 and 0.2 mL/g, respectively. The following classification bases on average values of solubles.

Group I: Volume contribution 0.8 mL per gram matrix.
With the limits $\pm 20\%$ as given above, the group value of 0.8 mL/g covers "true values" between 0.6 and 1.0 mL/g with sufficient precision.
e.g.: - Fruit
- Vegetables
- Green parts
- Fodder
- Meat

Group II: Volume contribution 0.5 mL per gram matrix.
"True values" between 0.3 and 0.7 mL/g are covered with sufficient precision.
e.g.: - Nuts
- Oil seeds

Proprietary information of CIBA-GEIGY AG.
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

Group III: Volume contribution 0.2 mL per g matrix.
"True values" between 0 and 0.4 mL/g are covered with sufficient precision.
e.g.: - Cereals
- Straw
- Peas and beans
- Other dried crops
- Soil

4. CERTIFICATION

The experimental results included in this method REM 172.02 are certified to be authentic accounts of the experiments. All raw data of results mentioned have been audited by the internal Quality Assurance Unit (QUA) and are stored in archives of Residue Analysis PP 2.53.

Statement of compliance with Good Laboratory Practice (GLP).

The recovery studies were performed in compliance with "Procedure and Principles of GLP in Switzerland", issued by the Department of the Interior, Bern, 1986, equivalent to "OECD Principles of Good Laboratory Practice", Paris, 1981. The facility was accredited by Swiss authorities to be in compliance with OECD GLP-principles.

5. REFERENCES

- REM 119.04: General Method for Calculation of Residue Analytical Results, Ciba, Plant Protection Division, Residue Analysis, 4002 Basel, CH.
- REM 172.03: "CGA 245704, DETERMINATION OF TOTAL RESIDUES AS CGA 210007 BY HPLC IN SOIL", Dr. G. Formica and C. Giannone, Sep. 10, 1993, CIBA, Plant Protection Division, Residue Analysis, 4002 Basel, CH.

Collaborator: H Borell


Formica

Dr. G. Formica
Ciba
Plant Protection Division
Residue Analysis (PP 2.53)
R-1098.P.09
CH-4002 Basel
Tel. (061) 697 27 85

C. Giannone


C. Giannone

Proprietary information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

6. FIGURES AND TABLES6.1 FiguresFigure 1: Structures and Chemical NamesCGA 245704 $C_8H_8N_2OS_2$


Molecular mass: 210.28

7-S-methyl-1,2,3-benzothiadiazolethioate

CGA 210007 $C_7H_6N_2O_2S$


Molecular mass: 180.18

1,2,3-Benzothiadiazol-7-carboxylic acid

Figure 2: Procedure Flow Diagram

* The moisture content of crop is to be taken into consideration (see section 2.3.3 and section 3.6)


Proprietary information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

Figure 3: Chromatographic System for Column Switching

Proprietary information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

6.2 Representative Chromatograms


6.2.1 Representative chromatograms for wheat grain. File: E. Wilhelm; book 10; Nov. 17, 1993



- A: Standard solution of CGA 210007. Volume injected corresponds to 2.0 ng a.i.
B: Standard solution of CGA 210007. Volume injected corresponds to 0.2 ng a.i.
C: Control wheat grain. Aliquot injected corresponds to 12.5 mg grain.

Proprietary information of CIBA-GEIGY AG.
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG.

6.2.1 (cont.)




D: Control wheat grain fortified with 0.02 mg CGA 245704 per kg.
 Aliquot injected corresponds to 12.5 mg grain or to 0.214 ng CGA 210007.
 Value found: 0.191 ng CGA 210007, corresponding to 90%.

E: Control wheat grain fortified with 0.20 mg CGA 245704 per kg.
 Aliquot injected corresponds to 12.5 mg grain or to 2.14 ng CGA 210007.
 Value found: 1.86 ng CGA 210007, corresponding to 87%.

F: Sample from treated field. Aliquot injected corresponds to 12.5 mg grain.

Proprietary information of CIBA-GEIGY AG
 Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG.

6.2.2 (cont.)




D: Control wheat straw fortified with 0.05 mg CGA 245704 per kg.
Aliquot injected corresponds to 5 mg straw or to 0.214 ng CGA 210007.
Value found: 0.192 ng CGA 210007, corresponding to 90%.

E: Control wheat straw fortified with 0.50 mg CGA 245704 per kg.
Aliquot injected corresponds to 5 mg straw or to 2.14 ng CGA 210007.
Value found: 1.73 ng CGA 210007, corresponding to 81%.

F: Sample from treated field. Aliquot injected corresponds to 5 mg straw.

Proprietary information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG


6.2.3 Representative chromatograms for wheat plant.
File: E. Wilhelm; book 10; Nov. 26, 1993



- A: Standard solution of CGA 210007. Volume injected corresponds to 2.0 ng a.i.
B: Standard solution of CGA 210007. Volume injected corresponds to 0.2 ng a.i.
C: Control wheat plant. Aliquot injected corresponds to 5 ng plant.

Proprietary information of CIBA-GEIGY AS
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AS

6.2.3 (cont.)



D: Control wheat plant fortified with 0.05 mg CGA 245704 per kg.
 Aliquot injected corresponds to 5 mg plant or to 0.214 ng CGA 210007.
 Value found: 0.23 ng CGA 210007, corresponding to 107%.


E: Control wheat plant fortified with 0.50 mg CGA 245704 per kg.
 Aliquot injected corresponds to 5 mg plant or to 2.14 ng CGA 210007.
 Value found: 1.92 ng CGA 210007, corresponding to 90%.

F: Sample from treated field. Aliquot injected corresponds to 5 mg plant.
 Value found: 1.27 ng CGA 210007, corresponding to 1.49 mg CGA 245704/kg.

Proprietary information of CIBA-GEIGY AG

Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG


6.2.4 Representative chromatograms for barley grain.
File: H. Borell; book 12; Nov. 18, 1993



A: Standard solution of CGA 210007. Volume injected corresponds to 2.0 ng a.i.
B: Standard solution of CGA 210007. Volume injected corresponds to 0.2 ng a.i.
C: Control barley grain. Aliquot injected corresponds to 12.5 mg grain.

Proprietary Information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

6.2.4 (cont.)




D: Control barley grain fortified with 0.02 mg CGA 245704 per kg.
Aliquot injected corresponds to 12.5 mg grain or to 0.214 ng CGA 210007.
Value found: 0.206 ng CGA 210007, corresponding to 97%.

E: Control barley grain fortified with 0.20 mg CGA 245704 per kg.
Aliquot injected corresponds to 12.5 mg grain or to 2.14 ng CGA 210007.
Value found: 1.95 ng CGA 210007, corresponding to 91%.

F: Sample from treated field. Aliquot injected corresponds to 12.5 mg grain.

Proprietary information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG


6.2.5 Representative chromatograms for barley straw.
File: H. Borell; book 12; Nov. 16, 1993



A: Standard solution of CGA 210007. Volume injected corresponds to 2.0 ng a.f.
B: Standard solution of CGA 210007. Volume injected corresponds to 0.2 ng a.f.
C: Control barley straw. Aliquot injected corresponds to 5 mg straw.

Proprietary information of CIBA-GEIGY AG.
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG.

6.2.5 (cont.)




D: Control barley straw fortified with 0.05 mg CGA 245704 per kg.
Aliquot injected corresponds to 5 mg straw or to 0.214 ng CGA 210007.
Value found: 0.203 ng CGA 210007, corresponding to 95%.

E: Control barley straw fortified with 0.50 mg CGA 245704 per kg.
Aliquot injected corresponds to 5 mg straw or to 2.14 ng CGA 210007.
Value found: 1.93 ng CGA 210007, corresponding to 90%.

F: Sample from treated field. Aliquot injected corresponds to 5 mg straw.

Proprietary information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG.

6.2.6 Representative chromatograms for barley plant.
File: H. Borell; book 12; Nov. 24, 1993




A: Standard solution of CGA 210007. Volume injected corresponds to 2.0 ng a.i.

B: Standard solution of CGA 210007. Volume injected corresponds to 0.2 ng a.i.

C: Control barley plant. Aliquot injected corresponds to 5 ng plant.

Proprietary information of CIBA-GEIGY AG.
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

6.2.6 (cont.)




D: Control barley plant fortified with 0.05 mg CGA 245704 per kg.
 Aliquot injected corresponds to 5 mg plant or to 0.214 ng CGA 210007.
 Value found: 0.215 ng CGA 210007, corresponding to 100%.

E: Control barley plant fortified with 0.50 mg CGA 245704 per kg.
 Aliquot injected corresponds to 5 mg plant or to 2.14 ng CGA 210007.
 Value found: 1.99 ng CGA 210007, corresponding to 93%.

F: Sample from treated field. Aliquot injected corresponds to 5 mg plant.
 Value found: 1.46 ng CGA 210007, corresponding to 0.34 mg CGA 245704/kg.

Proprietary information of CIBA-GEIGY AG.
 Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG.

6.2.7 Typical chromatogram for the determination of the switching interval.



A: Example of standard solution chromatographed only through column 1 to determine the switching interval (cut). In this case the "cut" is between 12 min 00 s and 13 min 24 s. Amount of CGA 210007 injected: 4 ng.


Proprietary information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

6.2.8 Typical calibration curves of CGA 210007.

a) File: H. Borelli; Book 12; Oct. 27, 1993.
Wheat grain; Sequence: HBU_4393.

b) File: H. Borelli; Book 12; Nov. 18, 1993.
Barley grain; Sequence: HBA_4693.

a)




CGA 210007 ng	Peak height	
	x mm	\bar{x}^* mm
4.0	121/114/109	113
2.0	60/56/63	56.8
0.5	15/14/13.5	14.3
0.2	6.0/5.8/5.5	5.75

* Mean values calculated according to REM 119.04 (see References, section 5).

Absolute standard deviation: 4.6%.

b)



CGA 210007 ng	Peak height	
	x mm	\bar{x}^* mm
4.0	120/122	120
2.0	60.4/60	60.0
0.5	14.7/14.9	15.1
0.2	6.30/6.16	6.18

* Mean values calculated according to REM 119.04 (see References, section 5).

Absolute standard deviation: 1.7%.

Proprietary information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG.

6.3 RECOVERIES

Status: As of Nov. 1993

TYPICAL PERCENT RECOVERIES USING METHOD REM 172.02

Substrate	Fortification level [mg/kg]	Percent recovered
Wheat grain	0.02	98/98/97/90/93
	0.04	98/86/93/88/92
	0.20	94/89/88/87/84
Wheat plant	0.05	75/94/107
	0.10	80/97/95
	0.50	93/98/90
Wheat straw	0.05	98/78/88/90/92
	0.10	93/86/87/80/76
	0.50	93/80/89/81/82
Barley grain	0.02	97/97
	0.04	92/93
	0.20	91/89
Barley plant	0.05	100
	0.50	93
Barley straw	0.05	95
	0.10	91
	0.50	90

Average recovery: $\bar{X} = 91\%$, $s_{obs} = 6.6\%$, $n = 50$

Acceptable recovery percentages: 71 - 111%

Further recovery values may be added and the range of acceptable recoveries may be updated accordingly.

Proprietary information of CIBA-GEIGY AG
 Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

The above recoveries are documented in the following files:

H. Borell; Book 12; Oct. 27, 1993. REPORT ON RESIDUE STUDY: 2067/92
H. Borell; Book 12; Oct. 27, 1993. REPORT ON RESIDUE STUDY: 2068/92
H. Borell; Book 12; Nov. 02, 1993. REPORT ON RESIDUE STUDY: 2067/92
H. Borell; Book 12; Nov. 08, 1993. REPORT ON RESIDUE STUDY: 2068/92
H. Borell; Book 12; Nov. 24, 1993. REPORT ON RESIDUE STUDY: 2027/93
H. Borell; Book 12; Nov. 18, 1993. REPORT ON RESIDUE STUDY: 2027/93
H. Borell; Book 12; Nov. 16, 1993. REPORT ON RESIDUE STUDY: 2027/93
H. Borell; Book 12; Nov. 16, 1993. REPORT ON RESIDUE STUDY: 2026/93
H. Borell; Book 12; Nov. 22, 1993. REPORT ON RESIDUE STUDY: 2026/93
H. Borell; Book 12; Nov. 22, 1993. REPORT ON RESIDUE STUDY: 2028/93

C. Giannone; Book 3R; Oct. 20, 1993. REPORT ON RESIDUE STUDY: 2067/92
C. Giannone; Book 3R; Oct. 20, 1993. REPORT ON RESIDUE STUDY: 2068/92

E. Wilhelm; Book 10; Nov. 26, 1993. REPORT ON RESIDUE STUDY: 2023/93
E. Wilhelm; Book 10; Nov. 24, 1993. REPORT ON RESIDUE STUDY: 2023/93
E. Wilhelm; Book 10; Nov. 24, 1993. REPORT ON RESIDUE STUDY: 2024/93
E. Wilhelm; Book 10; Nov. 17, 1993. REPORT ON RESIDUE STUDY: 2024/93
E. Wilhelm; Book 10; Nov. 17, 1993. REPORT ON RESIDUE STUDY: 2023/93

For all experiments mentioned, a Quality Assurance Statement was released by the internal QUA.

FILE /REP/CG/REN172.02

Proprietary information of CIBA-GEIGY AG
Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG

Dr. G. Formica
Residue Analysis
R-1098 P.09
Tel. 72785/72795

15 June 1994

AMENDMENT 1 *

TO CGA 245704, REM 172.02

*also filed
Remo 1993/HW*

Page 4 of the above Method has been revised.

Reason

Special attention must be paid to the word "preated". In the revised page, the word "preheated" has been emphasized.

15. Jun 1994 Formica
date Dr. G. Formica

* Revised page enclosed

File Am to REM172 02

Select fortification levels to be one and ten times the lower practical level of determination or in the range of the expected residue levels. Make sure that control samples neither are contaminated nor show interfering signals.

To fortify samples with 0.02 and 0.2 mg/kg for grain or 0.05 and 0.5 mg/kg for straw and plant, prepare solutions of CGA 245704 containing 0.1 and 1.0 µg per mL acetonitrile LiChrosolv by appropriate dilution of the stock solution (cf section 2.2.14). To prepare fortified samples, add 1 mL of one of the solutions to an untreated subsample. Proceed as described in section 2.3.3.

2.3.3 Hydrolysis/Extraction

Add as much water to the weighed subsample (cf. section 2.3.1) as to achieve a total water volume of 18 mL taking into account the crop solubles content being 1 mL for grain, 0.4 mL for straw and 1.6 mL for plant (see section 3.6 for moisture content). Add 2 mL of sodium hydroxide 1 M, seal the jar tightly, swirl the mixture and allow the jar to stand for about 30 minutes in an oven preheated to 50-55 °C. Swirl the jar occasionally. Allow the mixture to cool, add 80 mL (79 mL to fortified samples) of distilled methanol and 2 g of calcium chloride (the addition of the salt is necessary for grain and optionally for plant and straw. The salt facilitates the following filtration and reduces the formation of emulsion during the acidic partition, section 2.3.4). Shake the tightly sealed jar for about 30 minutes. Total volume of the extract is 100 mL. Filter immediately about 50 mL extract through a folded filter paper of 15 cm diameter into a 100 mL Erlenmeyer flask and continue according to section 2.3.4.

Remark: To avoid formation of emulsion during the following partition step, especially for wheat grain samples, filtration and subsequent partition cleanup (section 2.3.4) have to be performed without interruption immediately after shaking with methanol. If this is not possible, the extraction mixtures are to be immediately stored at about -18 °C until filtration and partition can be done. After aliquoting for analysis, the remaining filtrates are to be immediately stored at about -18 °C, as reserve for potential further analysis.

2.3.4 Cleanup by Acidic Partition

Transfer 20 mL of filtrate, corresponding to 1 g grain or 0.4 g straw and plant material to a 250 mL separatory funnel. Add 50 mL water, 40 mL saturated sodium chloride solution and 8 mL 1 M hydrochloric acid. Shake the mixture to homogenize. Add 50 mL of a solvent mixture of hexane + tert.-butyl methyl ether 7 vol + 3 vol. solvent mixture A). Shake and allow phases to separate. Drain the aqueous phase into a 250 mL Erlenmeyer flask. Filter the organic phase including emulsion, if present, into a clean 250 mL separatory funnel through a funnel plugged up with a loose plug of glass wool. Transfer the aqueous phase back to the separatory funnel. Add 50 mL of the solvent mixture A and shake. Allow phases to separate. Drain and discard the aqueous phase. Filter the organic phase into the separatory funnel containing the first reextract through the funnel plugged with glass wool. Observe if water is present below the organic phase. In this case, drain and discard the water.

* Revised page according to Amendment 1

Proprietary information of CIBA-GEIGY AG

Not to be disclosed to third parties without previous consent of CIBA-GEIGY AG


CGA 245704	RESIDUE METHOD VALIDATED	REM 172.11
DETERMINATION OF TOTAL RESIDUES OF CGA 245704 AS CGA 210007 BY HIGH PERFORMANCE LIQUID CHROMATOGRAPHY	PLANT MATERIAL AND MILK	Oct. 20, 1995 PP 2 533/Fo

TABLE OF CONTENTS


1	INTRODUCTION	3	RESULTS AND DISCUSSION
1 1	Scope of the Method	3 1	Accuracy and Precision
1 2	Principle of the Method	3 2	Limit of Quantitation
2	MATERIALS AND METHODS	3 3	Ruggedness
2 1	Equipment	3 4	Limitations
2 2	Reagents and Standards	3 5	Method Validation
2 3	Analytical Procedure	4	CERTIFICATION
2 4	Instrumentation	5	REFERENCES
2 5	Interferences	6	FIGURES AND TABLES
2 6	Confirmatory Techniques	6 1	Figures
2 7	Time Required for Analysis	6 2	Representative Chromatograms
2 8	Modifications and Potential Problems	7	APPENDICES
2 9	Calculation procedure	7.1	Report on Method Validation (Study 130/95)
2 10	Estimation of the Effective Extraction Volume		

CHEMICAL STRUCTURES

CGA 245704



CGA 210007



1. INTRODUCTION

1.1 Scope of the Method

The method allows the quantitative determination of the combined residues of CGA 245704 and its metabolites which contain the CGA 210007 moiety in plant material and milk (see section 6.1, Figure 1 for structures and chemical names)

The limits of quantitation of this method are 0.02 mg/kg for plant material and 0.005 mg/kg for milk, calculated as CGA 245704 equivalents

1.2 Principle of the Method

CGA 245704 and possible metabolites containing the CGA 210007 moiety are hydrolyzed to CGA 210007 by heating a homogenized subsample in the presence of sodium hydroxide. Methanol is added and the mixture shaken. An aliquot of the extract is cleaned up by acidic/alkaline/acidic partition. A further cleanup by a silica column (see section 2.3.6) may be necessary for some crop materials. Final determination is performed by a two column HPLC switching system with UV-detection. The principle of this method and the analytical procedure described are based on REM 172.02 (see reference 5.1)

2. MATERIALS AND METHODS

Standard laboratory equipment is not listed. All equipment and chemicals mentioned herein can be substituted by suitable products of any origin. Prove suitability of reagents by analyzing reagent blanks

2.1 Equipment

- 2.1.1 Rotating evaporator, Büchi, Rotavapor RE, Büchi AG, Flawil, CH
- 2.1.2 Vacuum system, Temo Duo 500 S, Temo AG, KÖblis, CH or commercial water aspirator
- 2.1.3 Circulation cooler, mgw Lauda WK 450, Messgeräte-Werke-Lauda, Dr. R. Wobser KG, Lauda-Königshofen, FRG, Cat. No. 7329-06.
- 2.1.4 Lab-shaker, A. Köhner AG, 4000 Basel, CH
- 2.1.5 Folded filter paper, 15 cm diameter, Macherey-Nagel, 5160 Düren, FRG, Cat. No. MN 713 1/4
- 2.1.6 Ultrasonic bath, Branson 220, Branson Cleaning Equipment Co., Parrot Drive, Shelton, CT 06484-0768, USA.
- 2.1.7 High Performance Liquid Chromatograph, refer to section 2.4

The following equipment is needed only when the silica column cleanup (see section 2.3.6) is necessary

- 2.1.8 Vacuum manifold to accommodate solid-phase extraction cartridges, built in-house or commercial equivalent, e.g. VISIPREP, Supelco, Bellefonte, Pennsylvania, USA, represented in Switzerland by Supelco SA, 9471 Buchs, CH, Cat. No. 5-7030
- 2.1.9 Vortex mixer, Model Vortex-Genie, Scientific Industries Inc., Springfield, Massachusetts, USA.
- 2.1.10 Empty borosilicate glass tube, 8 mL volume, J. T. Baker Inc., Phillipsburg, N.J. 08865, USA, represented in Switzerland by P. H. Stehelin and Cie AG, 4003 Basel, Cat. No. 7328-06.
- 2.1.11 PTFE frits for 8 mL borosilicate glass cartridge, J. T. Baker Inc., Cat. No. 7329-06
- 2.1.12 Adapter for 8 mL glass tube, J. T. Baker Inc., Cat. No. 4528.
- 2.1.13 Glass reservoirs built in-house of about 25 mL content to be attached to the cartridges or equivalent (e.g. 25 mL Bond-Ekret reservoirs, Varian, Harbor City CA, USA, Cat. No. A-1213-1011)

2.2 Reagents and Standard

Main suppliers' addresses. E. Merck AG, 6100 Darmstadt, FRG.
Fluka Chemie AG, 9470 Buchs, CH.
J. T. Baker Inc., Phillipsburg, N.J. 08864, USA

- 2.2.1 Water, HPLC-grade, prepared in house, Ciba or J. T. Baker Inc., Cat. No. 4218
- 2.2.2 Methanol for extraction, distilled in house, Ciba or analytical grade, E. Merck, Cat. No. 6009
- 2.2.3 n-Hexane for residue analysis, Merck, Cat. No. 4371
- 2.2.4 tert-Butyl methyl ether (TBME), puriss., analytical grade, Fluka, Cat. No. 20249.

- 2.2.5 Methanol, LiChrosolv, HPLC-grade, Merck, Cat. No. 6018
2.2.6 Acetonitrile, LiChrosolv, chromatography grade, Merck, Cat. No. 14291.
2.2.7 Hydrochloric acid 1 M, Ciba, Cat. No. 24 or J. T. Baker Inc., Cat. No. 7088 or about 1 M prepared by diluting analytical grade concentrated hydrochloric acid with HPLC-water
2.2.8 ortho-Phosphoric acid, 85%, analytical grade, E. Merck, Cat. No. 573. Prepare a solution containing about 0.02 mole of phosphoric acid per liter HPLC water by diluting 2.3 g concentrated acid to 1 liter with water for HPLC
2.2.9 Sodium hydroxide 1 M, Ciba, Internal Cat. No. 17 or J. T. Baker Inc., Cat. No. 7087.
2.2.10 Silica Gel for flash chromatography, J. T. Baker Inc., Cat. No. 7024-01
2.2.11 Calcium chloride dihydrate, analytical grade, Merck, Cat. No. 2382.
2.2.12 Sodium chloride, analytical grade, Merck, Cat. No. 6404. Prepare a saturated solution of sodium chloride in water for HPLC (brine).
2.2.13 Ammonium sulfate, analytical grade, Merck, Cat. No. 1217. Prepare a 0.05 M solution of ammonium sulfate by dissolving 6.61 g salt in one liter of 0.02 M phosphoric acid (see section 2.2.8)
2.2.14 CGA 245704 reference substance for recovery experiments. Prepare a stock solution containing 200 µg CGA 245704/mL acetonitrile LiChrosolv.
2.2.15 CGA 210007 reference substance for standardization. Prepare a stock solution containing 200 µg CGA 210007/mL acetonitrile LiChrosolv.

REMARK

CGA 245704 solutions in acetonitrile have been demonstrated to rapidly degrade if exposed to day light at room temperature. Therefore prepare solutions of CGA 245704 in brown bottles. However, if a stock solution of 200 µg/mL acetonitrile is stored in the dark at -20 °C, the a.i. is stable at least for 6 months.

2.3 Analytical Procedure

2.3.1 Preparation of Specimens and Subspecimens

Prepare representative laboratory specimens from the field specimens (e.g. mix plant material with dry ice and chop using a suitable cutter).

Analyze the specimens immediately after preparation or store them at about -20 °C until analysis.

For analysis let milk specimens thaw and shake vigorously. Weigh subspecimens of a g 10 g of plant material or 20 g milk into a 250 mL wide mouth jar. Avoid thawing of the frozen laboratory plant material specimens, but preferably scratch off or drill out the subspecimens needed for analysis.

2.3.2 Fortification

To regularly check the performance of the method, analyze also at least two fortified control specimens with each series of analyses. To prepare these specimens, add known amounts of CGA 245704 to control specimens prior to extraction. Select fortification levels to be one and ten times the limit of quantitation (LOQ) or in the range of the expected residue levels. Make sure that control specimens neither are contaminated nor show interfering signals.

To fortify specimens with 0.02 and 0.2 mg/kg (plant material) or 0.005 and 0.05 mg/kg (milk), prepare standard solutions of CGA 245704 containing 0.2 and 2.0 µg/mL acetonitrile LiChrosolv by appropriate dilution of the stock solution (see section 2.2.14). To prepare 0.02 and 0.2 mg/kg fortified specimens (for plant material) add 1 mL of one of the standard solutions to an uncontaminated plant material subspecimen. To prepare 0.005 and 0.05 mg/kg fortified specimens (for milk), add 0.5 mL of one of the standard solutions to an uncontaminated milk subspecimen. Proceed as described in section 2.3.3.

REMARK

The standard solutions used for the fortification were demonstrated to be chemically stable for at least 24 days, if they are kept in brown volumetric flasks and stored in a refrigerator at 4 °C.

2.3.3 Hydrolysis/Extraction**2.3.3.1 Plant Material**

Add as much water to the weighed plant material subspecimens (cf. section 2.3.1) as to achieve a total water volume of 18 mL, taking into account the soluble part of the subspecimens (see section 2.10). Add 2 mL of sodium hydroxide 1 M, seal the jar tightly, swirl the mixture and allow the jar to stand for about 30 minutes in an oven preheated to 55-60 °C. Swirl the jar occasionally. Allow the mixture to cool, add 80 mL (79 mL to fortified subspecimens) of distilled methanol and 2 g of calcium chloride (the addition of the salt may improve, for some materials, the following filtration and reduces the formation of an emulsion during the acidic partition, see section 2.3.4). Shake the tightly sealed jar for about 30 minutes. The total volume of the extract is 100 mL. Filter immediately about 50 mL extract through a folded filter paper of 15 cm diameter into a 100 mL Erlenmeyer flask and continue according to section 2.3.4.

2.3.3.2 Milk

Add 4 mL of sodium hydroxide 1 M, seal the jar tightly, swirl the mixture and allow the jar to stand for about 30 minutes in an oven preheated to 55-60 °C. Swirl the jar occasionally. Allow the mixture to cool, add 101 mL of distilled methanol (100.5 mL to fortified milk subspecimens) and 2 g of calcium chloride. Shake the tightly sealed jar for about 30 minutes. The total volume of the extract is 125 mL. Filter immediately about 50 mL extract through a folded filter paper of 15 cm diameter into a 100 mL Erlenmeyer flask and continue according to section 2.3.4.

2.3.4 Cleanup by Acidic Partition

Transfer 10 mL (plant material) and 25 mL (milk) of filtrate, corresponding to 1 g plant material and 4 g milk, respectively, to a 500 mL separatory funnel. Add 50 mL water, 40 mL saturated sodium chloride solution and 8 mL 1 M hydrochloric acid. Shake the mixture to homogenize. Add 50 mL of a solvent mixture of hexane + tert-butyl-methyl ether 7 vol. + 3 vol (solvent mixture A). Shake and allow phases to separate. Drain the aqueous phase and the emulsion between both phases, if there is any, into a 250 mL Erlenmeyer flask. Drain the organic phase into a clean 250 mL separatory funnel. Transfer the content of the Erlenmeyer flask back to the separatory funnel. Add 50 mL of the solvent mixture A and shake. Allow phases to separate. Drain and discard the aqueous phase. Drain the organic phase into the separatory funnel containing the first reextract. Observe if water is present below the organic phase. If this is the case, drain and discard the water.

2.3.5 Cleanup by Alkaline/Acidic Partition

Add 20.0 mL of 1 M sodium hydroxide to the organic phase of the acidic partition clean up (section 2.3.4), stopper the separatory funnel and shake well. Allow phases to separate. Drain 15.0 mL of the alkaline phase, corresponding to 0.75 g of plant material or 3 g of milk, into a graduated 25 mL test tube and discard the rest. Add 2 mL of 85% phosphoric acid and 10 mL of the solvent mixture A (see section 2.3.4). Shake and allow phases to separate. Transfer the organic phase into a 100 mL round bottom flask by means of a Pasteur pipette, taking care not to take along water drops! Repeat the partition two times more with portions of 10 mL of the same solvent mixture, each. Transfer the organic phase each time into the round bottom flask and discard the aqueous phase. Observe again the combined organic phase to make sure that no water drops have been taken along. Evaporate the organic phase to dryness under reduced pressure using a rotating evaporator (water bath temp. about 50 °C).

REMARK

The following described cleanup by silica is not necessary for the crops tested in this method and for milk. For these materials, it can be omitted. However, for other, not yet analyzed plant materials, this additional silica cleanup step may be necessary.

2.3.6 Cleanup by Silica Cartridge

To pack the cartridges used for this cleanup, proceed as follows: Put down a PTFE frit onto the bottom of a borosilicate glass tube of 8 mL content, fill 1 g of silica gel into the tube and place a second frit onto the top of the filling. Using a Vortex mixer, vibrate for approx. 1 minute for the material to settle. Press the upper frit onto the top of the adsorbent. Add 2 g of anhydrous sodium sulfate and place a third frit on the top. Vibrate as before for about 30 s and press the upper frit on the top. The cartridge is ready for use.

Attach the cartridge packed as described above, to the vacuum manifold. Tightly connect a 25 mL glass reservoir to the top of the cartridge. Precondition the cartridge by sucking through 5 mL of a solvent mixture of TBME + methanol LiChrosolv 85 vol. + 15 vol.

Dissolve the residue of the alkaline/acidic partition (see section 2.3.5) in 5 mL of a solvent mixture of TBME + methanol LiChrosolv 85 vol. + 15 vol. and transfer the solution to the reservoir. Pass the mixture through the cartridge at a rate of 5 - 7 mL per minute by suction. Discard the eluate. Wash the flask with 5 mL of the same solvent mixture, transfer the wash to the reservoir and pass it through the cartridge, too. Discard the eluate. Elute CGA 210007 with 10 mL of a solvent mixture of TBME + methanol LiChrosolv 6 vol. + 4 vol. into a 25 mL round bottom flask by suction. Evaporate the eluate to dryness under reduced pressure using a rotating evaporator (water bath temp.: about 50 °C).

2.3.7 Preparation of the Final Solution

Dissolve the residue obtained according to the procedure described in section 2.3.5 or 2.3.6 in 3 mL of 0.02 M phosphoric acid in HPLC water using an ultra sonic bath (final solution)

2.4 Instrumentation**2.4.1 High Performance Liquid Chromatographic System (HPLC)**

For the determination of CGA 210007 use a HPLC two-column switching system with UV-detector, pumps, autosampler-injector, columns and conditions as follows

Detector:	Spectroflow 783 (ABI Analytical, Kratos Division, Ramsey, J07446, USA) or equivalent.
Pumps:	Pumps must allow constant, pulsefree flowrate of 0.2 mL/min. Two Shimadzu Solvent Delivery Modules LC-9A (Shimadzu Corporation, Analytical Instruments Division, Kyoto, Japan) or equivalent.
Valves:	Multiport Streamswitch "MUST" equipped with two low dead volume valves for column switching (Spark Holland, 7800 AJ Emmen, Holland) in connection with an automatic injection unit PROMIS II (Spark Holland) or equivalent.
Recorder:	Strip chart SE 120, dual channel (ABB, Goerz Metrawatt, 1101 Vienna, A) or equivalent.
Column oven:	Temperature range 5 - 80 °C. Modell STH 585 (Gynkotek GmbH) or equivalent.

REMARK

A column oven is absolutely necessary in non thermostated laboratories to ensure constant elution times and a reproducible cut of the a.s. from column to column (see below "Determination of Switching Interval")

Optional (for system automation)
 Control and collection unit HP 3350A Laboratory Automation System (Hewlett-Packard, Palo Alto, data CA 94304, USA).
 Column 1: Stainless steel tube, 10 cm length, 2 mm i.d., packed with Nucleosil SB, strongly basic anion exchanger, particle size 5 μm , Cat. No : B45 - Y81, ("Wissenschaftliche Gerätebau" Dr. H. Knauer GmbH, 14163 Berlin, FRG) or equivalent.

REMARK
 A new column 1 has to be conditioned with the mobile phase 1 (see below) overnight at a flow rate of 0.2 mL/min. The retention time of CGA 210007 has to range between 10 and 15 min. If the retention time is above this range, the elution strength of mobile phase 1 has to be adjusted accordingly. For example, if the measured retention time is 20 min, the amount of acetonitrile in the mobile phase 1 has to be gradually increased to a maximum of 20%. (Remember: peak broadening occurs with an increasing amount of acetonitrile in the mobile phase 1, which results in a loss of overall sensitivity. Therefore keep the amount of acetonitrile in the mobile phase 1 as low as possible. The retention time on the column 1 has to be checked at each start of a new sequence and, if necessary, the concentration of acetonitrile has to be adjusted accordingly.)

Column 2: Stainless steel tube, 10 cm length, 2 mm i.d., packed with Nucleosil 100 C₁₈, particle size 5 μm , Cat. No : B45 - Y76 ("Wissenschaftliche Gerätebau" Dr. H. Knauer GmbH, 14163 Berlin, FRG) or equivalent

For the HPLC system, connect pumps, injector, columns, valves and detector according to Figure 3 (see section 6.1).

Mobile phase 1: 0.05 M ammonium sulfate in 0.02 M phosphoric acid + acetonitrile LiChrosolv 85 vol. + 15 vol. The composition of this mobile phase has to be adjusted as explained in the remark above

Mobile phase 2: 0.02 M phosphoric acid in water for HPLC + acetonitrile LiChrosolv 76 vol + 24 vol

Flow rate 1: 0.2 ml/min.
 Retention time: 10-15 min (column 1)
 Flow rate 2: 0.2 ml/min.
 Retention time: 15-19 min. (columns 1 and 2)
 Injection vol.: 50 μL
 Recorder: 10 mV full scale
 Column oven: 23 °C
 Chart speed: 0.5 - 1.0 cm/min
 Detector: Wave length: 235 nm; sensitivity: 0.004 aufs

Determination of Switching Interval

Determine actual switching times each time a series of specimens is to be quantitated. Connect the outlet of column 1 directly to the UV-detector. Inject 50 μL of a 0.08 μg CGA 210007/mL standard solution (see section 2.4.2) and measure the begin and the end times of the CGA 210007 peak. Thereafter, connect column 2 to the detector and switch that portion eluate 1 which contains CGA 210007, e.g. the portion between 11 min 30 s and 13 min 30 s from column 1 to column 2 (in this case the cut time is 2 min.). See section 6.2.2 for an example of the measurement of switching times.

REMARK

If an interference peak, which disturbs the quantitation of CGA 210007, appears in the control sample, reduce the "cut time" in portions of about 15 s each, until the interference disappears or is reduced to an insignificant height (below the limit of quantitation)

2.4.2 Calibration of the Chromatographic System by External Standards

Standardize the chromatographic system each time a series of specimens is to be quantitated. The range of the concentrations is depending on the range of residues to be determined, in particular, the lowest standard concentration is depending on the limit of quantitation.

Calculate the lowest standard concentration (C) as follows:

$$C = \frac{LOQ \cdot F \cdot A}{VF} \text{ [}\mu\text{g/mL]}$$

LOQ Limit of quantitation [mg/kg] (0.02 for plant material or 0.005 for milk)
F Conversion factor from CGA 245704 to CGA 210007 (0.857)
VF Volume of the final solution [mL]
A Aliquot of crop cleaned up in g (0.76 g for plant or 3.0 g for milk)

With the values proposed in this method, the lowest standard concentration is 0.004 μg CGA 210007/mL, as calculated below:

$$C = \frac{0.02 (0.005) \times 0.857 \times 0.76 (3.0)}{3} = 0.004 \text{ [}\mu\text{g/mL]}$$

Prepare at least four standard solutions of different concentrations by appropriately diluting the stock solution of CGA 210007 (see section 2.2.15) with 0.02 M phosphoric acid in water for HPLC.

Select the concentrations as required, typical values are 0.08, 0.04, 0.01 and 0.004 $\mu\text{g/mL}$. As the detector response may change on injection of coextractives, it is recommended to inject standards and final solutions alternately. Inject 50 μL of standard solutions (corresponding to 4.0, 2.0, 0.5 and 0.2 ng CGA 210007, respectively) and final solutions.

The conversion factor from CGA 210007 to CGA 245704 is 1.17.

2.5 Interferences

No interferences which hindered the determination of the CGA 210007 have been encountered so far.

2.6 Confirmatory Techniques

No specific conditions have been worked out so far.

2.7 Time Required for Analysis

A total of about 2 working days is required to work up a set of 12 specimens to the point of HPLC injection. Automated HPLC analysis can be performed overnight.

2.8 Modifications and Potential Problems

See last remark on page 6.

2.9 Calculation Procedure

Measure the response of the analyte at the characteristic retention time and calculate response function and residues as detailed in REM 119 06 (reference 5.2).

2.10 Estimation of the Effective Extraction Volume

The soluble part of the specimen (i.e. water and low molecular plant constituents) has to be considered for correct results. Exact values for this contribution are not easily available. In view of the relatively high imprecision of residue analytical results, rough estimates for the volume of the solubles are acceptable as long as the total volume is known with reasonable accuracy

In view of an overall accuracy of $\pm 2\%$ for the total extraction volume it is sufficient to know the volume contribution of the solubles by $\pm 20\%$, provided that the ratio specimen size/solvent added is 0.1, e.g. specimen size 10 g and solvent added 100 mL (supposing the specific weight of the solubles to be 1 g/mL and the volume contraction on mixing to be neglected)

For the calculation procedure of this method, the unknown exact volumes of the solubles are replaced by rough estimates:

The various crop materials are classified into three groups. The volumes of soluble constituents are estimated using group-factors (G) of 0.8, 0.5 and 0.2 mL/g. The following classification is based on average values of solubles (5.3):

Group I: Volume contribution of 0.8 mL per gram matrix
With the limits $\pm 20\%$ as given above, the group value of 0.8 mL/g covers "true values" between 0.6 and 1.0 mL/g with sufficient accuracy
E.g.: Fruits, vegetables, green parts, fodder, meat

Group II: Volume contribution of 0.5 mL per gram matrix.
"True values" between 0.3 and 0.7 mL/g are covered with sufficient accuracy
E.g.: Nuts, oil seed.

Group III: Volume contribution of 0.2 mL per gram matrix.
"True values" between 0 and 0.4 mL/g are covered with sufficient accuracy
E.g.: Cereals, straw, peas and beans, other dried crops, often also soil but check water content first.

3. RESULTS AND DISCUSSION

3.1 Accuracy and Precision

The overall accuracy and precision of the method is (status as of Sep.28, 1995)

$$\bar{x} = 93.5\% \quad \text{Sabs} = 7.2\% \quad n = 76$$

Acceptable recovery percentages 73 - 113 %

3.2 Limit of Quantitation

The limit of quantitation LOQ of this method is 0.02 mg/kg for plant material and 0.005 mg/kg for milk.

3.3 Ruggedness

The method was successfully tested by two technicians using fortified specimens. The analyses were performed on different days with different instrumentation. The results obtained were comparable with regard to recovery percentages and chromatographic profile.

3.4 Limitations

None known so far.

3.5 Method Validation

The method has been validated by performing recovery experiments. The report on this validation study (130/95) including QA-statement is attached to this method as section 7.1

4. CERTIFICATION

The experimental results included in this method REM 172.11 are certified to be authentic accounts of the experiments.

Origin of data: Special Study 130/95.

5. REFERENCES

- 5.1 REM 172.02, "CGA 245704, DETERMINATION OF TOTAL RESIDUES AS CGA 210007 BY HPLC IN CEREALS", Dr. G. Formica and C. Giannone, Dec. 08, 1993, CIBA, Crop Protection Division, Residue Analysts, 4002 Basel, CH.
- 5.2 REM 119.06: General Method for Calculation of Residue Analytical Results, Ciba, Crop Protection Division, Residue Analysis, 4002 Basel, CH
- 5.3 Documenta Geigy, Wissenschaftliche Tabelle, Ed. J. R. Geigy AG, CH-Basel, 7th ed. 1968, p. 494-511

Collaborator: H. Borell


Formica

Dr. G. Formica
Ciba
Crop Protection Division
Residue Analysis (PP 2 533)
R-1098.P 09
CH-4002 Basel
Tel. (061) 697 27 85


H. Borell

H. Borell

File: REM17211

6. FIGURES AND TABLES**6.1 Figures****Figure 1. Structures and Chemical Names**CGA 245704

Benzo[1,2,3]thiadiazole-7-carbothioic acid S-methyl ester

 $C_8H_6N_2OS_2$ Molecular mass: 210.28CGA 210007

Benzo[1,2,3]thiadiazole-7-carboxylic acid

 $C_8H_6N_2OS$ Molecular mass: 180.18

Figure 2: Procedure Flow Diagram






Figure 3. Setup of the Column Switching System




6.2 Representative Chromatograms

All chromatograms shown were recorded during the validation study (see section 7.1) Peak responses are expressed in Peak Height Units (PHU). However, all peak heights in the following Tables are given in mm as measured on the chart recorder


6.2.1 CGA 210007 in banana peel




Standard solution 0.08 µg/mL CGA 210007
(4 ng injected), peak 120 mm



Standard solution 0.004 µg/mL CGA 210007
(0.2 ng injected), peak 6 mm




Untreated control specimen (banana peel),
formal specimen size injected: 12.5 mg,
peak 0 mm




Control specimen fortified with 0.02 mg/kg
CGA 245704, formal specimen size injected
12.5 mg, corresponding to 0.214 ng CGA
210007 Peak found: 6.5 mm, corresponding
to 0.22 ng CGA 210007 or to 102%


6.2.1 CGA 210007 in banana peel - Continued



Control specimen fortified with 0.02 mg/kg CGA 245704, formal specimen size injected 12.5 mg, corresponding to 0.214 ng CGA 210007. Peak found: 7.3 min, corresponding to 0.24 ng CGA 210007 or to 114%



Control specimen fortified with 0.02 mg/kg CGA 245704, formal specimen size injected 12.5 mg, corresponding to 0.214 ng CGA 210007. Peak found: 7.0 min, corresponding to 0.23 ng CGA 210007 or to 109%




Control specimen fortified with 0.2 mg/kg CGA 245704, formal specimen size injected 12.5 mg, corresponding to 2.14 ng CGA 210007. Peak found: 12.9 min, corresponding to 2.10 ng CGA 210007 or to 98%.

Further details are exactly described in this method and are given again below here. They are valid for this sequence as well as for all mentioned sequences in this method (except for milk, where further details are given in the pages containing the corresponding chromatograms):

Weight of plant: 10 g
Total volume for extraction: 100 mL
Aliquot: 0.75 g plant. Final volume: 3 mL
Volume injected: 50 μ L

Details on the calibration curve are given on the section 7

6.2.2 Chromatogram of the Analyte on the First Column of the HPLC System, used for the determination of the switching interval




Standard solution 0.08 $\mu\text{g}/\text{mL}$ CGA 210007,
50 μL injected, corresponding to 4 ng


The two vertical bars on the peak of the analyte, mark a possible choice for the cutting interval (CUT). In this case the interval is between 11 min. 30 s and 13 min 30 s. The peak fraction within this interval is to be switched to the second column, as described in the section 2.4.1

The given chromatogram was recorded just before the start of the sequence PGZ_3595 for the validation of the method for tomatoes (see also the chromatograms recorded of this sequence on the following page).


6.2.3 Analysis for CGA 210007 in tomatoes




Standard solution 0.08 µg/mL CGA 210007
(4 ng injected), peak: 12.3 min



Standard solution 0.004 µg/mL CGA 210007
(0.2 ng injected), peak: 6 min




Untreated control specimen (tomatoes), formal
specimen size injected. 12.5 mg,
peak: 0 min




Control specimen fortified with 0.02 mg/kg
CGA 245704, formal specimen size injected.
12.5 mg, corresponding to 0.214 ng CGA
210007. Peak found 4.8 min, corresponding
to 0.20 ng CGA 210007 or to 92 %


6.2.3 CGA 210007 in tomatoes - Continued



Control specimen fortified with 0.02 mg/kg
CGA 245704, formal specimen size injected
12.5 mg, corresponding to 0.214 ng CGA
210007. Peak found: 5.5 min, corresponding
to 0.22 ng CGA 210007 or to 102 %




Control specimen fortified with 0.02 mg/kg
CGA 245704, formal specimen size injected
12.5 mg, corresponding to 0.214 ng CGA
210007. Peak found: 5.5 min, corresponding
to 0.22 ng CGA 210007 or to 102 %




Control specimen fortified with 0.2 mg/kg
CGA 245704, formal specimen size injected
12.5 mg, corresponding to 2.14 ng CGA
210007. Peak found: 5.8 min, corresponding to
1.9 ng CGA 210007 or to 89% recovery


6.2.4 CGA 210007 in cucumbers




Standard solution 0.08 µg/mL CGA 210007
(4 ng injected), peak 10.6 min



Standard solution 0.004 µg/mL CGA 210007
(0.2 ng injected), peak 4.7 min




Untreated control specimen (cucumbers),
formal specimen size injected 12.5 mg,
peak 0 min




Control specimen fortified with 0.02 mg/kg
CGA 245704, formal specimen size injected
12.5 mg, corresponding to 0.214 ng CGA
210007. Peak found 4.9 min, corresponding
to 0.20 ng CGA 210007 or to 95%


6.2.4 CGA 210007 in cucumbers - Continued



Control specimen fortified with 0.02 mg/kg CGA 245704, formal specimen size injected: 12.5 mg, corresponding to 0.214 ng CGA 210007. Peak found: 5.3 min, corresponding to 0.22 ng CGA 210007 or to 102%.




Control specimen fortified with 0.02 mg/kg CGA 245704, formal specimen size injected: 12.5 mg, corresponding to 0.214 ng CGA 210007. Peak found: 4.2 min, corresponding to 0.18 ng CGA 210007 or to 83%.




Control specimen fortified with 0.2 mg/kg CGA 245704, formal specimen size injected: 12.5 mg, corresponding to 2.14 ng CGA 210007. Peak found: 5.2 min, corresponding to 2.0 ng CGA 210007 or to 91% recovery.


6.2.5 CGA 210007 in milk




Standard solution 0.08 µg/mL CGA 210007
(4 ng injected), peak: 10.4 min



Standard solution 0.004 µg/mL CGA 210007
(0.2 ng injected), peak: 5 min




Untreated control specimen (milk), formal
specimen size injected 50 mg,
peak: 0 min




Control specimen fortified with 0.005 mg/kg
CGA 245704, formal specimen size injected
50 mg, corresponding to 0.214 ng CGA
210007. Peak found: 5 min, corresponding to
0.20 ng CGA 210007 or to 95%


6.2.5 CGA 210007 in Milk - Continued



Control specimen fortified with 0.005 mg/kg
CGA 245704, formal specimen size injected:
50 mg, corresponding to 0.214 ng CGA
210007. Peak found: 5 min, corresponding to
0.20 ng CGA 210007 or to 95 %



Control specimen fortified with 0.005 mg/kg
CGA 245704, formal specimen size injected:
50 mg, corresponding to 0.214 ng CGA
210007. Peak found: 5 min, corresponding to
0.20 ng CGA 210007 or to 95 %



Control specimen fortified with 0.05 mg/kg
CGA 245704, formal specimen size injected:
50 mg, corresponding to 2.14 ng CGA
210007. Peak found: 5.4 min, corresponding to
2.10 ng CGA 210007 or to 98 %

Further details for milk sequences:
Weight of milk: 20 g
Total volume for extraction: 125 mL
Aliquot: 3 g Final volume: 3 mL
Volume injected: 50 µL

7 APPENDICES**7.1 Report on Validation Study**

Report on Special Study 130/95	VALIDATION of Method REM 172.11 Validation by Analysis of Fortified Specimens and Determination of Recoveries
--------------------------------------	---------------------------------------------------------------------------------------------------------------------

STUDY DIRECTOR Dr G Formica, R-1098.P.09, Ciba, Crop Protection
Product Safety/Dietary Exposure/Residue Analysis
Date *Oct 20, 95* Signature *Formica*

TEST SUBSTANCES CGA 245704 and CGA 210007
TEST SYSTEMS bananas; tomatoes; cucumbers and milk

1 GENERAL

testing facility Ciba-Geigy AG, Crop Protection/Residue Analysis CH-4002 Basel
sponsor Ciba-Geigy AG, Crop Protection, Dietary Exposure Assessment, CH-4002
Basel, Dr S. Dennis
quality assurances Ciba-Geigy AG, Crop Protection/Product Safety, CH-4002 Basel
test substance Ciba code CGA 245704
Chemical name Benzo[1,2,3]thiadiazole-7-carboxylic acid S-methyl ester
Ciba code: CGA 210007
Chemical name: 1,2,3-Benzothiadiazol-7-carboxylic acid
study dates Protocol issued 12 Jul 1995
Experimental Date 25 July, 95 - 9 Sep. 85
archives Ciba, Crop Protection/Residue Analysis.
All raw data (including draft method) and original report.
guidelines - Procedures and Principles of GLP (OECD, Paris 1981, CH-Bern 1986)
- European Directive EEC 91/414, European Uniform Principles

2 OBJECTIVE

- The objective of this study is to establish the performance of method REM 172.11 by determining recovery values obtained by analyses of fortified specimens in two Laboratories according to the procedure described in the method.
- The report on validation including the results is issued as part of the final method REM 172.11.

3 RESULTS

See Table on section 7 of this validation
Average recovery (Banana, Tomatoes, Cucumbers and Milk) 93.5%, Sabs=7.2%, n = 76
Acceptable recovery percentage: 73.5 - 113.5
Limit of Quantitation (LOQ) 0.02 mg/Kg for Plant and 0.005 mg/Kg for Milk

4 EXPERIMENTAL DESIGN FOR RECOVERIES

Determination of residues of fortified specimens under repeatability and reproducibility conditions, including checking of performance criteria according to validation requirements of EEC 91/414. Details are as follows.

- **Laboratory 1** (P. Gänswein, R-1098 P 05), **substrate group 1** (banana pulp)
8 specimens at limit of quantitation (LOQ) of 0.02 mg/kg
8 specimens at 10 times LOQ
- **Laboratory 2** (H. Borell, R-1098 P 05), **substrate group 1** (banana pulp)
3 specimens at limit of quantitation (LOQ) of 0.02 mg/kg
- **Laboratory 1** (P. Gänswein), **substrate group 2** (tomatoes)
3 specimens at limit of quantitation (LOQ) of 0.02 mg/kg
3 specimens at 10 times LOQ
- **Laboratory 1** (P. Gänswein), **substrate group 3** (cucumbers)
3 specimens at limit of quantitation (LOQ) of 0.02 mg/kg
3 specimens at 10 times LOQ
- **Laboratory 1** (P. Gänswein), **substrate group 4** (milk)
3 specimens at limit of quantitation (LOQ) of 0.005 mg/kg
3 specimens at 10 times LOQ


5 ANALYTICAL METHOD

Method REM 172.11, draft issued before start of validation, no modifications to this draft were made during the study.

6 Preparation of the Calibration Curve

The chromatograms given in section 6.2.1 and originate from the same sequence: HBT_3195. All the responses of the standard solutions injected for the preparation of the calibration curve are summarized in the following table

Concentration of standard solution [$\mu\text{g/mL}$]	Amount of standard injected [ng]	Peak response calculated [mm]	Peak responses measured [mm]
0.08	4	120	120, 119; 119
0.04	2	59.9	60; 60; 61
0.01	0.5	15.0	15, 15, 15
0.004	0.2	6.0	6, 6; 6



Standard deviation

Quality criterion for standard curve, standard deviation s (square root of the variance of all relative deviations of the measured response from the corresponding value of the response function), calculated according to REM 119.08 (see reference given in section 5.2)

7 RESULTS (TABLE OF ALL RECOVERIES)

The following table and section 6 2 1 and 6 2.5 contain all information, which is necessary to calculate the individual recoveries.

Substrate	Fortificat. level (mg/kg)	Laboratory 1 Person 1: P. Gärsvén (not involved in the Method development)				Laboratory 2 Person 2: H. Bonell			
		recovery %	peak height mm	sequence * calibration curve $y=A + Bx$	average and standard deviation	recovery %	peak height mm	sequence, calibration curve $y=A + Bx$	average and standard deviation
banana peel	0.02	95	6.0	PGJ_3465 A=0.5625 B=28.62	$\bar{x} = 88.5\%$ $s = 7.1$ n=6	102	6.5	HEJ_3195 A=0.008553 B=29.98	$\bar{x} = 103.7\%$ $s = 6.5\%$ n=6
		78	5.0			114	7.3		
		82	6.2			109	7.0		
0.2	94	54	$s = 8.1$	98	63	$s = 1.2\%$	98	63	
	89	51		98	63				
	82	53		101	68				
banana pulp	0.02	98	5.0	PGP_3395 A=0.6960 B=28.15	$\bar{x} = 92.3\%$ $s = 7.2\%$ n=6	89	6.2	HEJ_3195 A=0.008553 B=29.98	$\bar{x} = 91.3\%$ $s = 3.3\%$ n=6
		95	4.8			86	5.0		
		98	5.0			94	5.5		
0.2	91	54	$s = 2.7\%$	85	59	$s = 2.1\%$	85	59	
	101	60		92	68				
	87	54		92	57				
0.02	86	5.5	PGW_3465 A=0.8590 B=25.33	$s = 6.2\%$	86	5.5			
	76	5.0			86	5.5			
	85	6.5			86	5.5			
0.2	89	5.7	$s = 6.2\%$	102	55				
	91	5.8		100	55				
	100	52		94	52				
tomatoes	0.02	82	4.8	PGZ_3595 A=1.309 B=31.12	$\bar{x} = 93.7\%$ $s = 6.6\%$ n=6	112	6.0	HEJ_3655 A=0.4801 B=26.92	$\bar{x} = 93.0\%$ $s = 13.0\%$ n=6
		102	5.5			95	5.0		
		102	5.5			95	5.0		
0.2	89	58	$s = 3.5\%$	72	41	$s = 3.0\%$	88	58	
	88	58		88	58				
	89	58		88	58				
cucumbers	0.02	98	5.0	PGE_3695 A=0.8046 B=27.34	$\bar{x} = 95.0\%$ $s = 5.1\%$ n=6	95	4.9	HEJ_3695 A=0.5973 B=26.93	$\bar{x} = 93.2\%$ $s = 6.2\%$ n=6
		98	5.0			102	5.3		
		104	5.5			83	4.2		
0.2	90	52	$s = 3.8\%$	91	52	$s = 1.7\%$	85	54	
	91	53		85	54				
	93	54		85	53				
milk	0.005	91	5.0	PGG_3695 A=0.2135 B=29.92	$\bar{x} = 90.2\%$ $s = 6.6\%$ n=6	103	5.5	HEJ_3495 A=0.5383 B=27.39	$\bar{x} = 95.3\%$ $s = 3.0\%$ n=12
		82	4.5			85	5.0		
		82	4.5			85	5.0		
0.06	96	55	$s = 4.6\%$	97	56	$s = 2.4$	97	56	
	95	54		93	54				
	95	54		93	54				
0.005				94	5.0	HEJ_3695 A=0.1883 B=25.70	$s = 1.2\%$	95	52
				95	52				
				95	52				
0.05				98	5.0	HEJ_3595 A=0.2354 B=25.83	$s = 2.4\%$	98	5.0
				98	5.0				
				98	5.0				
0.05				98	54			98	54
				92	51				
				92	51				

* Name of the file containing the chromatograms in the related substrate
No residues were found in all control specimens

8 CRITERIA FOR METHOD VALIDATION

To check the criteria described in the Annex VI of European Directive 91/414/EEC for validation of this method only the necessary recovery values have been considered in the following tables:

8.1 FIRST SUBSTRATE GROUP (banana pulp)**Recovery Table (Results)**

Lab 1	con-	recoveries (%) found (individual data)							
Substrate 1	trol	1	2	3	4	5	6	7	8
0.02 mg/kg	<0.02	98	95	98	86	78	95	89	91
0.20 mg/kg	<0.02	91	101	91	87	82	102	100	94

Recovery-check

Lab 1	N	average	EEC-criterion	standard deviation	EEC-criterion	Decision/Comments
Substrate 1						
0.02 mg/kg	8	91 %	70-110 %	7.4 %	≤20 %	EEC criterion fulfilled
0.20 mg/kg	8	94 %	70-110 %	7.1 %	≤20 %	EEC criterion fulfilled

Repeatability Table at LOQ

Lab 1	N	Recoveries found at LOQ			Decision/Comments criterion: difference ≤ 50 %
Substrate 1		minimum	maximum	difference	
0.02 mg/kg (LOQ)	8	76 %	98 %	22 %	EEC criterion fulfilled

Generation of results under repeatability conditions (same lab, operator, within short interval of time). Exception 2 different equipments were used, i.e. the results were obtained under worst conditions and nevertheless the EEC criterion was fulfilled

Recovery Table (Results)

Substrate 1	con-	recoveries (%) found at LOD (individual data)		
at LOQ	trol	1	2	3
Lab 2	<0.02	89	86	94

Generation of results under reproducibility conditions (different lab, operator, equipment with exception of 3 results)

Reproducibility Table at LOQ

Lab 1 + 2	N	Combined recoveries at LOQ			Comment criterion: difference ≤ 100 %
Substrate 1		minimum	maximum	difference	
0.02 mg/kg (LOQ)	11	76 %	98 %	22 %	EEC criterion fulfilled

8.2 SECOND SUBSTRATE GROUP (tomatoes)**Recovery Table (Results)**

Lab 1	con-	recoveries (%) found (individual data)		
Substrate 2	trol	1	2	3
0.02 mg/kg	<0.02	92	102	102
0.20 mg/kg	<0.02	89	88	89

Recovery-check

Lab 1	N	average	EEC-criterion	standard deviation	EEC-criterion	Decision/Comments
Substrate 2						
0.02 mg/kg	3	99 %	70-110 %	6 %	≤20 %	EEC criterion fulfilled
0.20 mg/kg	3	89 %	70-110 %	1 %	≤20 %	EEC criterion fulfilled

8.3 THIRD SUBSTRATE GROUP (cucumbers)

Recovery Table (Results)

Lab 1 Substrate 3	con- trol	recoveries [%] found (individual data)		
		1	2	3
0.02 mg/kg	<0.02	96	88	104
0.20 mg/kg	<0.02	90	91	93

Recovery-check

Lab 1 Substrate 3	N	average	EEC- criterion	standard deviation	EEC- criterion	Decision/Comments
0.02 mg/kg	3	99 %	70-110 %	5 %	≤ 20 %	EEC criterion fulfilled
0.20 mg/kg	3	91 %	70-110 %	2 %	≤ 20 %	EEC criterion fulfilled

8.4 FOURTH SUBSTRATE GROUP (milk)

Recovery Table (Results)

Lab 1 Substrate 4	con- trol	recoveries [%] found (individual data)		
		1	2	3
0.005 mg/kg	< 0.05	91	82	82
0.05 mg/kg	<0.05	96	95	95

Recovery-check

Lab 1 Substrate 4	N	average	EEC- criterion	standard deviation	EEC- criterion	Decision/Comments
0.005 mg/kg	3	85 %	70-110 %	5 %	≤ 20 %	EEC criterion fulfilled
0.05 mg/kg	3	95 %	70-110 %	1 %	≤ 20 %	EEC criterion fulfilled

8.5 OVERALL RECOVERY TABLE (all fortification levels, four substrate groups, two labs)

Labs 1 + 2	N	data of recoveries from above			average criterion 70-110 %	standard deviation
		minimum	maximum	difference		
Substrate 1	19	76 %	102 %	26 %	92 %	7 %
Substrate 2	6	88 %	102 %	14 %	94 %	7 %
Substrate 3	6	90 %	104 %	14 %	95 %	5 %
Substrate 4	6	82 %	96 %	14 %	90 %	7 %
Substrate 1 + 2 + 3 + 4	37	76 %	104 %	28 %	92 %	6 %

9 Collaborators

H Borell, P Gänswain

10 Conversion Factor

The conversion factor from CGA 210007 to CGA 245704 is: 1.17

11 CONCLUSION

According to the criteria described in the Annex VI of European Directive 91/414/EEC this method has been successfully validated.

Statement of Compliance

This study was performed in a facility which was certified to be in compliance with GLP. The study was conducted in full compliance with procedures and principles of GLP in Switzerland, and OECD principles of GLP.

Quality Assurance Statement

See next page.

Quality Assurance Statement
Ciba-Geigy Ltd , GLP Quality Assurance Product Safety, 4002 Basel

Project 130/95
Test Substance CGA 245704
Study Title Recovery Validation of Method REM 172.11 "CGA 245704,
Determination of Total Residues of CGA 245704 as CGA 210007
in Plant Material and Milk"
Study Director Dr Giuseppe Formica
QA-Inspector Pierre Fontana

I hereby certify that the following Quality Assurance activities were performed:

Activity	Performed	Reported
Facility Inspection	May 11, 1995	May 29, 1995
Protocol Audit	July 18, 1995	July 18, 1995
Final Report Audit	October 19, 1995	October 20, 1995

October 20, 1995
Date
Form QES/TA/11


.....
Pierre Fontana
Inspector Quality Assurance

APPENDIX II

SEPARATE DOCUMENTS ACCOMPANYING THIS SUBMISSION

Hayworth, C. and Grunenwald, M , ABR-97117, Final Report Amendment Number 1, "Validation of Draft Analytical Method AG-671A for the Determination of Total Residues of CGA-245704 as CGA-210007 in Tobacco, Leafy Vegetables, and Fruiting Vegetables by Column Switching High Performance Liquid chromatography."

Kiffe, M , CIBA GEIGY Ltd , Basel, Switzerland, Plant Metabolism, "Metabolism of CGA 245704 in Greenhouse Grown Lettuce after Treatment with [Phenyl-U-¹⁴C] Labelled Material, Study 96MK02 1," In-progress.

Nicollier, G., CIBA GEIGY Ltd., Basel, Switzerland, Plant Metabolism, "Metabolism of CGA 245704 in Greenhouse Grown Tomatoes after Treatment [Phenyl-U-¹⁴C] Labeled Material," Study 94WA05.

SUBMITTER/SPONSOR: Novartis Crop Protection, Inc., P. O. Box 18308, Greensboro, NC 27419-8380

APPENDIX III

PREVIOUSLY SUBMITTED DOCUMENTS WITH EPA MRID NUMBERS

Nicollier, G , CIBA GEIGY Ltd , Basel, Switzerland, Plant Metabolism, "Behavior and Metabolism of CGA 245704 in Greenhouse Grown Tobacco after Foliar Spray Application of [U-¹⁴C-Phenyl]-CGA-245704," Study 94WA04 MRID No. 44014255

SUBMITTER/SPONSOR: Novartis Crop Protection, Inc., P O Box 18308, Greensboro, NC 27419-8300

EPA ADDENDUM

PP#8F4974

Acibenzolar-s-methyl on Bananas, Spinach, Tomatoes, and Tomato Paste

Conditions and Substitutions used by ACB

1 ACB used a Hewlett-Packard (HP) 1100 Series High Performance Liquid Chromatograph (HPLC) equipped with a HP 1100 Series variable wavelength ultraviolet/visible (UV/VIS) detector (HPLC-UV). The detector was set at 235 nm. The detector was interfaced to a TurboChrom (ver 6.1.1) data system. Two analytical columns were used. The columns were connected to a column switching valve on the HPLC, therefore, only one instrument (pump) was used. The column temperatures were programmed according to the method (30 °C). Column 1 was a Phenomenex® Spherclone SAX, 150 x 4.6 mm i.d., 5 µm, Cat. #OOF-4149-EO. The retention time for this column was ~8 minutes. Column 2 was a Phenomenex® Luna C₈, 250 x 4.6 mm i.d., 5 µm, Cat. #OOG-4040-EO. The retention time for this column (Col.1 + Col. 2) was ~42 minutes. The total analysis time was 65 minutes. The flow rate through both columns was 0.6 mL/min using mobile phase 1 and 2 as prepared in the method. However, when the flow was switched to Col. 2, 2% acetonitrile was added from another reservoir into the mixing chamber (2% acetonitrile + 98% mobile phase 2) to bring the analyte into a retention time window free from endogenous peaks. A 100 µL injection volume was used instead of 50 µL as suggested in the method. The method called for Column 1 to be a Phenomenex® Spherisorb SAX/SB, 150 x 3.2 mm i.d., 5 µm with a flow rate of 0.3 mL/min. The method called for Column 2 to be a YMC-Pack C₈, 250 x 4.6 mm i.d., 5 µm with a flow rate of 1.0 mL/min. By using two columns of the same i.d. (inside diameter) on one instrument, the same flow rate could be used with good chromatographic results. With a flow rate of 0.6 mL/min, back pressure was controlled and the amount of mobile phase used was minimized (hence hazardous waste reduced).

2. The TurboVap®II (Zymark Corporation, Hopkinton, MA) operated at 60 °C under a nitrogen stream was used in place of the rotary vacuum evaporators to concentrate samples during the methanol evaporation. ACB highly recommends this substitution be incorporated into the procedure. This substitution eliminated the need to monitor the samples during this step against any bumping and/or foaming that was forewarned in the method had the rotary evaporators been used. In addition, TurboVap® glassware has graduated volumes. This was important because the final volume was recommended to be 8-10 mL. An 8 mL aliquot of acid is added to the 8-10 mL remaining volume and the mixture is loaded onto a Chem Elute™ column. This was a critical step because the Chem Elute™ columns have a finite capacity of 20 mL.

EPA ADDENDUM

PP#8F4974

Acibenzolar-s-methyl on Bananas, Spinach, Tomatoes, and Tomato Paste

Page 2

3a. ACB found no information in the method regarding the time that the separatory funnels were to be shaken during the partitioning steps. ACB shook each sample for one minute following venting procedures.

3b ACB found no information concerning stopping points in the procedure. ACB found that storing the extracts stoppered in the refrigerator overnight following the elution from the Chem Elute™ columns (100 mL 9:1 Dichloromethane/ethyl acetate) did not diminish recoveries. ACB found that storing the extracts stoppered in the refrigerator overnight following the final partitioning with ethyl acetate prior to evaporating to dryness did not diminish recoveries.

3c. ACB found no information concerning the stability of the final extracts contained in auto-sampler vials awaiting analysis. ACB found that the extracts were stable at room temperature for a least 48 hours in capped auto-sampler vials.

4. The method suggested the use of calibration curves to determine sample concentrations (peak height vs concentration). ACB ran calibration curves to demonstrate linearity over the dynamic range of the fortified samples. ACB determined sample concentrations from a ratio of sample responses to the average of standard responses that bracketed the samples during analysis. The bracketing standards were made to be the same concentration as the expected concentration of the fortified samples if 100% recovery were obtained. Prior to analyzing a sample set, the highest concentration of standard was injected onto column 1 to determine the retention time of the standard. After the retention time was established, the timetable of the method was revised to optimize the column switching event.