

US EPA ARCHIVE DOCUMENT



Report No. AV96R001

Page 155

Appendix IV (continued)

**ANALYTICAL PROCEDURE FOR THE DETERMINATION OF
PROPAMOCARB IN POTATOES**

ANALYTICAL METHOD: XAM-34

Issued: October 15, 1995

Author: Nicole Bertrand

Developed by: Xenos Laboratories Inc.
2319 St-Laurent Blvd., Unit 100
Ottawa, Ontario
K1G 4K6, Canada

Approved by: R. Greenhalgh
Xenos Laboratories Inc.
2319 St-Laurent Blvd., Unit 100
Ottawa, ON
K1G 4K6, Canada

Appendix IV (continued)

Xeno Method XAM-34	1
TABLE OF CONTENTS	
TITLE PAGE	1
TABLE OF CONTENTS	2
1.0 INTRODUCTION.....	4
2.0 REFERENCE COMPOUNDS	4
2.1 Hoe 074189, Propamocarb	4
3.0 PRINCIPLE OF THE METHOD	5
4.0 EQUIPMENT.....	5
5.0 CHEMICALS/REAGENTS	6
6.0 PREPARATION OF STANDARD SOLUTIONS	6
6.1 Propamocarb Fortification Solutions	6
6.2 Propamocarb Calibration Solutions	7
7.0 ANALYTICAL PROCEDURE	7
7.1 Sample Preparation	7
7.2 Extraction	7
7.3 Liquid/Liquid Partition	8
7.4 Chem-But Clean-Up	8
8.0 QUANTITATION GC/MS.....	9
8.1 Instrumentation	9
8.2 Gas Chromatographic Conditions	9
8.3 Mass Spectrometer Conditions	9
8.3.1 Mass spectrometer Acquisition Segments	9
8.4 GC/MS Calibration	10
8.5 Sample Analysis	10
9.0 QUANTITATION GC/TSD	10
9.1 Instrumentation	10
9.2 Gas Chromatographic Conditions	11
9.3 TSD Calibration	11
9.4 Sample Analysis	12
10.0 CALCULATION OF THE RESIDUE	12
10.1 Calculation	12
10.2 Example Calculation of Procedural Recovery-GC/MS	13
10.3 Example Calculation of Procedural Recovery-GC/TSD	13

Appendix IV (continued)

Xenos Method XAM-34 3

Table of Contents continued

11.0 QUALITY CONTROL PROCEDURES.....	14
11.1 Laboratory Purifications and Control Samples	14
11.2 GC Analysis	14
11.3 Sample Storage	14
12.0 RESULTS.....	14
Table I: Recoveries of Propamocarb From Fortified Potato - GC/MS	15
Table II: Recoveries of Propamocarb From Fortified Potato - GC/TSD.....	15
13.0 REFERENCES.....	16

APPENDICES:

- Appendix I: Flow Diagram of Method XAM-34
Appendix II: GC/MS; Representative Ion Chromatograms, Typical Calibration Curve and Chromatographic Data Sheet
Appendix III: GC/TSD; Representative Chromatograms, Typical Calibration Curve and Chromatographic Data Sheet

Appendix IV (continued)**Xenex Method XAM-34**

4

1.0 INTRODUCTION


Propamocarb is a systemic fungicide when applied to soil but it can also be used as a foliar spray. Propamocarb is specific against *Phycomyces* and is recommended as a preventative treatment.

Analytical procedures for fruit, vegetable and sandy loam soils have been reported. These procedures generally used gas chromatography (GC) with a nitrogen/phosphorus specific detector. A procedure using mass spectrometry (MS) in the electron impact (EI) mode for quantitation has also been reported (Chambers and Charters, 1994). The analytical procedure described in this report utilizes both a thermionic selective detector (TSD) and an ion trap MS for quantitation. The latter operates in the Mass Instability Chemical Ionization (MICI) mode.

This analytical method is suitable for the determination of extractable residues of the fungicide propamocarb in potato. The limit of quantitation (LOQ) for this procedure is 0.05 ppm.

2.0 REFERENCE COMPOUND**2.1 Hot 074189: Propamocarb**

CAS Name: propyl [3-(dimethylamino)propyl]carbamate
IUPAC Name: propyl (3-dimethylaminopropyl)carbamate
Structural formula:



Appendix IV (continued)

Xenos Method XAM-34

5

3.0 PRINCIPLE OF THE METHOD

Extractable residues of propanocarb are removed from potatoes by homogenizing with acidified methanol. The extract is filtered and the methanol is removed using a rotary evaporator. The acidic extract is partitioned with dichloromethane and diisopropyl ether. The aqueous extract is basified by addition of 10 N NaOH and the solution is applied to a Chem-Elut column. Propanocarb is eluted from the Chem-Elut, with diisopropyl ether. The eluate is evaporated using a rotary evaporator. The residue is made up to final volume with diisopropyl ether and the extract is analyzed by GC/TSD or GC/Ion Trap MS operating in the SECI mode. The latter uses the 189 m/z for quantitation of the propanocarb residues.

A flow diagram of the analytical method is presented in Appendix I.

4.0 EQUIPMENT

Unless otherwise indicated, the equipment listed below may be substituted with functionally equivalent equipment as may be available.

- Balance: Mettler model BB-2440, for samples
- Balance: Mettler model AE 240, for standards
- Commercial Food Processor: Rubo-Coupe
- Microliter syringes: 100, 250, 500, 1000 microliter, Hamilton
- Volumetric pipettes, glass: 1.0-mL, 5.0-mL and 10.0-mL
- Polytron Homogenizer, equipped with PT-35 generator
- Beakers, glass: 250-mL
- Büchner Funnels, porcelain
- Adapter, glass: to filter directly into round bottom flask
- Filter papers: glass fiber 934-AH
- Round bottom flask, glass: 500-mL, 24/40 joint
- Separatory funnel, glass: 125-mL
- Round bottom flask, glass: 250-mL
- Rotary Evaporator equipped with water bath
- Disposable Pasteur pipettes, glass: 22.5 cm x 0.7 mm O.D.
- Graduated cylinder, glass: 10-mL, 25-mL, 100-mL
- Volumetric flasks, glass: 5-mL, 10-mL, 25-mL, 50-mL
- pH paper

Appendix IV (continued)

Xenon Method XAM-34

6

5.0 CHEMICALS/REAGENTS

Alternate suppliers/brands of reagents having comparable specifications may be used

- Methanol: (pesticide quality)
- Diisopropyl ether: ACS
- Water: distilled
- Sodium hydroxide: (Reagent Grade)
- Hydrochloric acid: concentrated
- Chem-Eut Cartridge, CE 1020: Analytichem/Varian
- Celite 545

6.0 PREPARATION OF STANDARD SOLUTIONS

The analytical standard is supplied by the sponsor: AgrEvo USA Company or AgrEvo Canada Inc. It should be stored in a secured freezer when not in use, at about -13°C.

6.1 Propamocarb Fortification Solutions**Solution A:**

Weigh 50 mg (± 0.1 mg) of Hoe 074189 and quantitatively transfer into a 50-mL volumetric flask. The volume is made up with acidified methanol (1 mL of 1N HCl added to 100 mL of methanol) to give a concentration of 1.0 mg/mL of propamocarb. This stock solution will be used to prepare fortification solutions. The solution is stored in a freezer and must be prepared anew every six months.

Solution B:

Transfer 500 μ L of primary stock solution A to a 50-mL volumetric flask. The volume is made up with methanol to give a concentration of 10.0 μ g/mL of Propamocarb. This solution is stored in a freezer and is prepared anew every three months.

Solution C:

Transfer 5.0 mL of solution B to a 50-mL volumetric flask. The volume is made up with methanol to give a concentration of 1.0 μ g/mL of Propamocarb. This solution is stored in a freezer and is prepared anew every three months. Additional fortification solutions may be prepared if necessary.

6.2 Propamocarb Calibration Solutions**Solution D:**

Weigh 50 mg (± 1 mg) of Hoe 074189 and quantitatively transfer into a 50-mL volumetric flask. The volume is made up with diisopropyl ether to give a concentration of 1.0 mg/mL of propamocarb. This stock solution will be used to prepare GC calibration standards. The solution is stored in a freezer and must be prepared anew every six months.

Appendix IV (continued)

Xenos Method XAM-34

7

6.3 Propamocarb Calibration Solutions, continued

Solution E

Transfer 500 μ L of solution D to a 50-mL volumetric flask. The volume is made up with diisopropyl ether to give a concentration of 10.0 μ g/mL of Propamocarb. This solution is stored in a freezer and is prepared anew every month.

The calibration solutions are prepared by dilution of solution E. Typical dilutions are shown below:

Final concentration (μ g/mL)	Final volume (mL)	Volume of solution E used (μ L)
0.10	50.0	500
0.15	50.0	750
0.25	50.0	1250
0.50	50.0	2500
0.75	50.0	3750
1.0	50.0	5000
1.5	50.0	7500

The calibration standard solutions are prepared anew every month. They are kept in a refrigerator (about 4°C) when not in use.

7.0 ANALYTICAL PROCEDURE

7.1 Sample Preparation:

The potato sample is ground using a food processor (Robo-Coupe) and a representative analytical sample obtained.

7.2 Extraction:

A representative analytical sample (20 ± 0.1 g) is weighed into a 250-mL glass beaker. HCl (1 N, 1.0 mL) is added, followed by methanol (100 mL). The mixture is homogenized, using the Polytron at speed 4, for one minute. The extract is filtered through a glass fiber filter covered with a thin layer of Celite 545 (approximately 1.5 g) into a 500-mL round bottom flask. The beaker is rinsed with 50 mL of methanol which is poured through the filter pad and collected into the same round bottom flask. The methanol is rotary evaporated with a water bath temperature set at approximately 45 °C, leaving only the acidic aqueous extract (approximately 5 mL).

Appendix IV (continued)

Xenos Method XAM-34

8

7.3 Liquid/Liquid Partition:

Additional HCl (5mL, 1N) is added to the aqueous extract. The pH is verified at this stage, using pH paper to assure that it is acidic. The flask is swirled to dissolve the residue and the mixture is transferred to a 125-mL separatory funnel. Dichloromethane (20 mL) is added to the round bottom flask and the flask is swirled to dislodge any residue adhering to the walls of the flask. The dichloromethane is poured into the separatory funnel and the mixture is shaken for approximately 30 seconds. The phases are allowed to separate and the dichloromethane layer is discarded. The dichloromethane partition is repeated with an additional portion of 20 mL, which is also used to rinse the round bottom flask before adding to the separatory funnel. The aqueous extract is then partitioned with diisopropyl ether (20 mL) and the aqueous layer is drained into the round bottom flask and set on the rotary evaporator, with the water bath at approximately 45 °C, to remove any traces of organic solvent. Only a few minutes are necessary to accomplish this.

7.4 Chem-Elat Clean-Up:

The acidic extract from 7.3 is basified by adding NaOH (10 mL of 10 N). The basified extract is poured into a 25-mL graduated cylinder and the volume is adjusted to 15-20 mL with distilled water. The solution is mixed and poured into the Chem-Elat column. After the solution has completely drained, the column is left to equilibrate for 5 minutes. The Chem-Elat column is eluted with 100 mL of diisopropyl ether (DIPE), collecting the eluate into a 250-mL round bottom flask. The diisopropyl ether is evaporated down to approximately 2 mL, using a rotary evaporator with a water bath set at approximately 40°C. The solution is transferred to a 5-mL volumetric flask with small portions of diisopropyl ether. The volume is adjusted to 5 mL with DIPE. An aliquot of this solution is transferred to an autosampler vial for analysis by GC/TSD or GC/Ion Trap MS.

Appendix IV (continued)

Xenos Method XAM-34

9

8.0 QUANTITATION GC/MS**8.1 Instrumentation:**

The Varian Saturn 3 Ion Trap Mass Spectrometer system consists of a model 3400 GC and a model 8200 autosampler. The MS was operated in the SECI mode.

8.2 Gas Chromatographic Conditions:

GC Column: DB-5MS 30 m x 25 mm ID, 0.25 μ m film
Carrier Gas Flow: approximately 50 cm³/min at 260 °C
Temperature:
 Inlet: Split/Splitless injector, 275 °C. Split Ratio 10:1
 Column: Initial: 70 °C, hold: 1.0 min.
 70-200 °C at 10 °C/min., hold: 0 min.
 200-300 °C at 50 °C/min., hold: 8.3 min
Transfer Line: 300 °C
Mass Spectrometer Interface: direct capillary
Injection Volume: 3.0 μ L
Injection Rate: 5 μ L/sec.
Retention Time: approximately 10 minutes
Limit of Quantitation: 0.025 ppm
Limit of Detection: 0.3 ng (S/N >20)

Representative chromatograms and calibration plot are included in Appendix II.

8.3 Mass Spectrometer Conditions:**8.3.1 Mass Spectrometer Acquisition Segments:**

A/M Amplitude:	4 V	Low Mass:	187 m/z
High Mass:	191 m/z	Scan Rate:	1000 ms
Segment Time:	11.4 min.	Peak Threshold:	0 counts
Filament Delay:	7.5 min.	Mass Defect:	50 m/z/100 u
Background Mass:	100 amu	Calibration Gas:	no
Scan Mode:	SECI	Ionization Control:	Automatic
Tune File:	PROPTUNE	Emission Current:	20 μ A
CI-Gas	Isobutane		



Appendix IV (continued)

Xmas Method XAM-34

10

8.4 GC/MS Calibration:

The GC/MS response (peak area) is determined for a series of calibration standards. The 189 m/z ion is used to generate extracted ion profiles for propanocarb. Detector response is non-linear and can be described by a logarithmic function of the form shown in Equation 1. For each analytical set, the calibration data was used to perform a power regression analysis. The natural logarithm of the amount injected (ng) was taken as the X-axis and the natural logarithm of the detector response (peak area) was taken as the Y-axis to give Equation 2.

$$\ln y = \ln xm + b \quad [\text{Eq. 1}]$$

where: y = peak area response for analyte in injected sample/standard
 m = slope of the regression line
 x = amount (ng) of analyte found in the sample/standard
 b = intercept of the regression line

$$\ln (\text{peak area}) = \ln (\text{ng in the sample/standard})m + b \quad [\text{Eq. 2}]$$

8.5 Sample Analysis:

For sample extract of unknown propanocarb content, the amount (ng) found may be calculated, from the observed peak area, using Equation 3.

$$\text{ng (x) in sample} = e^{\ln (\text{peak area}) - Mb} \quad [\text{Eq. 3}]$$

Both samples and standards must be analyzed under the same GC conditions and within the same analytical sequence.

9.0 QUANTITATION GC/TSD

9.1 Instrumentation:

The gas chromatographic system consisted of a Varian 3400 GC equipped with a thermionic specific detector and an 8200 autosampler.

Appendix IV (continued)

Xenon Method XAM-34

11

9.2 Gas Chromatographic Conditions:

The gas chromatographic system consists of a Varian 3400 GC equipped with a thermionic specific detector (TSD) and an 8200 autosampler.

GC Column:	DB-1, 15 m x 0.53 mm ID, 1.5 or 3.0 μ m film
Carrier Gas Flow:	4.5 mL/min
Temperatures,	
Inlet:	273 °C
Column:	Initial: 100 °C, hold: 0.2 min 100-140 °C at 5 °C/min., hold: 0 min. 140-230 °C at 30 °C/min., hold: 13.8 min
Detector:	300 °C
Injection Volume:	2.0 μ L, on-column
Injection Rate:	1.0 μ L/sec
Retention Time:	10.2 minutes
Limit of Quantitation:	0.025 ppm
Limit of Detection:	0.2 ng (S/N \geq 3)

Representative chromatograms and calibration plot are included in Appendix III.

9.3 TSD Calibration:

The GC/TSD response (peak area) is determined for a series of calibration standards. The detector response is linear and can be described by a function of the form shown in Equation 1. For each analytical set, the calibration data was used to perform a linear regression analysis. The amount injected (ng) was taken as the X-axis and the detector response (peak area) was taken as the Y-axis to give Equation 2.

$$y = xm + b \quad [\text{Eq. 1}]$$

where: y = peak area response for analyte in injected sample/standard

m = slope of the regression line

x = amount (ng) of analyte found in the sample/standard

b = intercept of the regression line

$$\text{peak area} = (\text{ng in the sample/standard})m + b \quad [\text{Eq. 2}]$$

Appendix IV (continued)

Xenos Method XAM-34

12

9.4 Sample Analysis:

The peak area response for Hoe 074189 is computed using the Varian Star Workstation. The amount of material is determined from the corresponding calibration plot. For sample extracts of unknown propamocarb content, the amount found (in ng) may be calculated, from the observed peak area, using Equation 3

$$\text{ng (x) in sample} = \frac{\text{peak area} \cdot k}{m} \quad [\text{Eq. 3}]$$

Both samples and standards must be analyzed under the same GC conditions and within the same analytical sequence.

10.0 CALCULATION OF RESIDUES

10.1 Calculation:

The amount of residues in a sample are expressed in parts per million. They are calculated using the following equation:

$$\text{ppm} = \frac{\text{Amt}}{B}$$

where: Amt = ng (x) analyte found from the standard curve

B = mg of sample injected

$$B = \frac{W_s (\text{g}) \cdot V_{inj} (\mu\text{L})}{V_f (\text{mL})}$$

where: W_s = weight of sample (g)

V_f = final volume of extract (mL)

V_{inj} = sample volume injected (μL)

$$\% \text{ Recovery} = \frac{\text{ppm found} \cdot 100}{\text{ppm added}}$$



Appendix IV (continued)

Xenos Method XAM-34

13

10.2 Example Calculation of Procedural Recovery-GC/MS:

The data for sample POT-F2-1, fortified at 0.50 ppm is given in GC data sheet, X9618V04, (Appendix II). The procedural recovery calculation for Hoe 074189 is shown below:

$$\text{Amt} = x = e^{(1.872 + 18.235)(1.344)} = e^{10.54} = 2,869 \text{ ng}$$

$$B = \frac{20 \text{ g} * 3.0 \text{ uL}}{10 \text{ mL}} = 6.0 \text{ mg}$$

ppm found = $\frac{\text{ng found from curve}}{\text{mg-equivalent injected}}$

$$\text{ppm found} = \frac{2,869 \text{ ng}}{6.0 \text{ mg}} = 0.478 \text{ ppm as Hoe 074189}$$

$$\text{Recovery} = \frac{\text{ppm found}}{\text{ppm fortification level}} = \frac{0.478 \text{ ppm} * 100}{0.50 \text{ ppm}} = 95.6\%$$

10.3 Example Calculation of Procedural Recovery-GC/TSD:

The data for sample POT-F1-3, fortified at 0.05 ppm is given in GC data sheet, X9618V01, (Appendix II). The procedural recovery calculation for Hoe 074189 is shown below:

$$\text{Amt} = x = \frac{4207 - 125 - 217.18}{10473.3} = 0.369 \text{ ng}$$

where: 125 is the peak area found in the control sample extract

$$B = \frac{20 \text{ g} * 2.0 \text{ uL}}{5 \text{ mL}} = 8.0 \text{ mg}$$

ppm found = $\frac{\text{ng found (from linear regression)}}{\text{mg-equivalent injected}}$

$$\text{ppm found} = \frac{0.369 \text{ ng}}{8 \text{ mg}} = 0.046 \text{ ppm as Hoe 074189}$$

$$\text{Recovery} = \frac{\text{ppm found}}{\text{ppm fortification level}} = \frac{0.046 \text{ ppm} * 100}{0.05 \text{ ppm}} = 92.3\%$$

Appendix IV (continued)

Xenon Method XAM-34

14

11.0 QUALITY CONTROL PROCEDURES**11.1 Laboratory Fortified Controls:**

To assure the quality of the analytical data, laboratory fortified controls are analyzed with each set of samples. These fortified controls should cover the expected range of residues in the samples and be at least 10% of the samples within an analytical set.

11.2 GC Analysis:

To verify the stability of the response, a standard plot is drawn for the levels of interest. Within an analytical set, analytical standards are typically injected after every sample, this serves as an ongoing quality control for detector sensitivity and analyte retention time. Calibration standards are analyzed at the beginning and end of an analytical set.

The lowest level analytical standard should correspond to 50 to 70% of the limit of quantitation. Residue results must not be determined by extrapolation of calibration data outside of the concentration range of the calibration standards (using $\pm 10\%$ tolerance, typically). Samples with residue levels greater than the calibrated range must be diluted and re-injected so that they do fall within the calibrated range.

The recovery of laboratory fortified controls should fall within a range of 70 - 120%. The analyte signal should be 2-3 times the background signal, and the intra-laboratory reproducibility as indicated by the relative standard deviation ($n \geq 3$) obtained from replicated analyses should fall within 20% (rel.) of the averaged result.

11.3 Sample Storage:

Field samples should be kept frozen until analyzed. After obtaining a representative analytical sample, the remaining material should be promptly refrozen and stored until authorization for disposal is received.

12.0 RESULTS

Potato samples were fortified at 0.050 and 0.500 ppm with Propamocarb. Five replicate samples were analyzed for each level and the recoveries determined. For GC/MS, the mean of the recoveries for the two levels of fortification were: $115 \pm 11.6\%$ and $97.2 \pm 3.3\%$ for 0.050 ppm (LOQ) and 0.50 ppm, respectively (Table I). The procedure was validated also using GC/TSD, the mean of the recoveries for the two levels of fortification were: $74.4 \pm 13.5\%$ and $94.4 \pm 5.8\%$ for 0.050 ppm (LOQ) and 0.50 ppm, respectively (Table II). These data validate an LOQ of 0.050 ppm for potato.



Appendix IV (continued)

Xenos Method XAM-34

15

Table I: Recoveries of Propamocarb from Fortified Potato - GC/MS

Sample No.	Fortification Level (ppm)	Found (ppm)	Recovery (%)
POT-F1-1	0.050	0.056	112
POT-F1-2	0.050	0.049	98.9
POT-F1-3	0.050	0.063	131
POT-F1-4	0.050	0.057	115
POT-F1-5	0.050	0.059	119
Mean Recovery ± std dev.			115 ± 11.6%
POT-F2-1	0.50	0.478	95.6
POT-F2-2	0.50	0.500	100
POT-F2-3	0.50	0.489	97.9
POT-F2-4	0.50	0.461	92.3
POT-F2-5	0.50	0.502	100
Mean Recovery ± std dev.			97.2 ± 3.3 %
Overall Mean Recovery (n = 10)			105 ± 12.4

Table II: Recoveries of Propamocarb from Fortified Potato - GC/TSD

Sample No.	Fortification Level (ppm)	Found (ppm)	Recovery (%)
POT-F1-1	0.050	0.037	73.2
POT-F1-2	0.050	0.039	78.6
POT-F1-3	0.050	0.046	92.3
POT-F1-4	0.050	0.027	54.7
POT-F1-5	0.050	0.037	73.2
Mean Recovery ± std dev.			74.4 ± 13.5
POT-F2-1	0.50	0.478	95.7
POT-F2-2	0.50	0.487	97.4
POT-F2-3	0.50	0.480	96.0
POT-F2-4	0.50	0.421	84.2
POT-F2-5	0.50	0.493	98.7
Mean Recovery ± std dev.			94.4 ± 5.8
Overall Mean Recovery (n = 10)			84.4 ± 14.4

Appendix IV (continued)

Xenos Method XAM-34

16

19.0 REFERENCES

1. Analytical Method for the Determination of Residues of Propamocarb x HCl in Soil, J. Moede; Report Number UPSR 54/91 - PA 66 752.5/14 Schering AG,
2. Dissipation of Propamocarb X HCl in Soil following application of Banol - USA 1990, A.Wrods-Rucker, Schering report no. UPSR 57/91 - PA 66 752.7/16, July 7, 1992)
3. Dissipation of Propamocarb HCl in Soil Following Application of Banol to Bare Plot, USA, 1993, Martin G. Cole; AgrEvo USA Registration reference, Propamocarb HCl/W132 A54951


Appendix IV (continued)

Krene Method XAM-34

17

Appendix I

Schematic Flow Diagram XAM-34



Appendix IV (continued)

Xmas Method XAM-34

18

Appendix II

GC/MS - Representative Ion Chromatograms, Typical Calibration Curve
and Chromatographic Data Sheet

- Figure 1: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 0.10 ng/ μ L, 3 μ L injected. Reference Data Sheet #: X9618V04.
- Figure 2: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 0.15 ng/ μ L, 3 μ L injected. Reference Data Sheet #: X9618V04
- Figure 3: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 0.25 ng/ μ L, 3 μ L injected. Reference Data Sheet #: X9618V04
- Figure 4: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 0.50 ng/ μ L, 3 μ L injected. Reference Data Sheet #: X9618V04
- Figure 5: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 0.75 ng/ μ L, 3 μ L injected. Reference Data Sheet #: X9618V04
- Figure 6: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 1.00 ng/ μ L, 3 μ L injected. Reference Data Sheet #: X9618V04
- Figure 7: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 1.50 ng/ μ L, 3 μ L injected. Reference Data Sheet #: X9618V04
- Figure 8: Ion Chromatogram (m/z 189) of control potato sample, 3 μ L of 5 mL injected.
POT-C1, Reference Data Sheet #: X9618V03
- Figure 9: Ion Chromatogram (m/z 189) of control potato sample, 3 μ L of 10 mL injected.
POT-C2, Reference Data Sheet #: X9618V04
- Figure 10: Ion Chromatogram (m/z 189) of fortified control potato sample with
Hoe 074189 at 0.030 ppm, 3 μ L of 5 mL injected.
POT-F1-2, Reference Data Sheet #: X9618V03
- Figure 11: Ion Chromatogram (m/z 189) of fortified control potato sample with
Hoe 074189 at 0.50 ppm, 3 μ L of 10 mL injected.
POT-F2-2, Reference Data Sheet #: X9618V04



Report No. AV96R001


Page 173

Appendix IV (continued)

Xmas Method XAM-34

19

Figure 1: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 0.10 ng/ μ L, 3 μ L injected. Reference Data Sheet #: X9618V04.



Integration Report Scan File: P803 Call File: P80F Entries: 1
Comment: X9618V04\PROPANOCARD\0.10 MC/L
Scanned via: Entry Number 1 IS Factor: 1.000 Multi: 1.000 Div: 1.000

Name of Compound	R/R Ratios	R Time	Scan#	Pk Height	Peak Area
PROPANOCARD	24.8	18.21	613	639	1.648

Xmas Report Number: XLEN97-37

Page 43




Report No. AV96R001
Page 174

Appendix IV (continued)

Xenos Method XAM-34

20

Figure 2: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 0.15 ng/ μ L, 3 μ L injected. Reference Data Sheet #: X9618V04



Integration Report Quan File: PE05 Cal File: PNAP Entries: 1
Comment: XEN9617\96ACBL\PROPANOCARB-0.15 NG/UL
Sorted via: Entry Number t IS Factor: 1.000 Mult: 1.000 Div: 1.000

Name of Compound	E/M Ratio	R Time	Scan#	Pk Height	Pk Area
PROPANOCARB	68.5	18.14	685	1.720	16.274

Xenos Report Number: XA97-37

Page 44



Report No. AV96R001


Page 175

Appendix IV (continued)

Xmass Method XAM-34

21

Figure 3: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 0.25 ng/ μ L.
3 μ L injected. Reference Data Sheet #: X9618V04



Integration Report Run File: PE07 Call File: PMP Entries: 1
Comment: XEN9617\SEAC01\PROPANOCARB\0.25 MC/L
Sorted via: Entry Number F IS Factor: 1.000 Mult: 1.000 Div: 1.000

Name of Compound	E/M Ratio	R Time	Scan#	Pk Height	Peak Area
PROPANOCARB	145.2	18.81	587	4.533	22.655

Xmass Report Number: X9618V04

Page 21




Appendix IV (continued)

Xemai Method XAM-34

23

Figure 4: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 0.50 ng/ μ L, 3 μ L injected. Reference Data Sheet #: X96189V04



Integration Report Quan File: P099 Cali File: PROP Entries: 1
Comment: XM19517/96ACB1-PROFAMOCARB/8.68 MC/UL
Sorted via: Entry Number f 12 Factor: 1.000 Multi: 1.000 Div: 1.000

Name of Compound	S/N Ratio	R Time	Scan#	Pk Height	Peak Area
PROPHYLACIDE	559.1	10.15	685	15,625	55,282




Report No AV96R001
Page 177

Appendix IV (continued)

Xmas Method XAM-34

23

Figure 5: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 0.75 ng/ μ L, 3 μ L injected. Reference Data Sheet #: X9618V04



Integration Report Quan File: PE91 Cal File: PHOP Entries: 1
Comment: X9618V04\96ACB1\PROPAMOCARD\0.75 NG/UL
Sorted via: Entry Number 1 IS Factor: 1.000 Multi: 1.000 Div: 1.000


Name of Compound	R/R Ratio	R Time	Scan#	Pk Height	Pk Area
PROPAMOCARD	799.9	18.49	606	24,519	82,080

Appendix IV (continued)

Xenos Method XAM-34

24

Figure 6: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 1.00 ng/ μ L,
 3 μ L injected, Reference Data Sheet #: X9618V04



Integration Report Run File: P093 Cali File: PROF Entries: 1
 Comment: XEN9617\96ACB\PROPANOCARB\1.0 NG\UL
 Sorted via: Entry Number ? IS Factor: 1.000 Mult: 1.000 Div: 1.000

Name of Compound	R/H Ratio	R Time	Scan#	Pk Height	Pk Area
PROPANOCARB	1305.3	10.89	605	36,668	118,053

Xenos Report Number: XEN9770


Page 10

Appendix IV (continued)

Xmas Method XAM-34

25

Figure 7: Ion Chromatogram (m/z 189) of Hoe 074189 standard at 1.50 ng/ μ L, 3 μ L injected. Reference Data Sheet #: X9618V04



Integration Report Scan File: PO05 Call File: PROD Entries: 1
Comment: X9617\96AC01\PROPHOCARD\1.5 NG\UL
Sorted via: Entry Number : IS Factor: 1,000 Mult: 1,000 Div: 1,000

Name of Compound	E/N Ratio	R Time	Scan#	Pk Height	Peak Area
PROPHOCARD	9537.0	18.00	605	74,234	299,578

Xmas Report Number: X9617

Page 71




Appendix IV (continued)

Xmas Method XAM-34

26

Figure 8: Ion Chromatogram (m/z 189) of control potato sample, 3 μ L of 5 mL injected.
POT-C1, Reference Data Sheet #: X9615V03



Integration Report Query File: Plots Call File: PROP Entries: 1
Comment: XEN9617/96ACN/PROPANOCARD/POT-C1
Selected via: Entry Number 1 IB Factor: 1.000 Relt: 1.000 Dis: 1.000

Name of Compound	Z/N Ratio	R Time	Scand	Pk Height	Peak Area
PROPANOCARD	20.7	18.15	605	643	4.347

Xmas Report Number: XEN97-3


Page 72

Appendix IV (continued)

Xenes Method XAM-34

27

Figure 9: Ion Chromatogram (m/z 189) of control potato sample, 3 μ L of 10 mL injected.
POT-C2, Reference Data Sheet #: X9618V04



Integration Report Run File: PS04 Call File: PROP Entries: 1
Comment: XEN9617/96ACSL\PROPNOCARB\POT-C2
Sorted via: Entry Number 1 IS Factor: 1.000 Mult: 1.000 Div: 1.000

Name of Compound	S/N Ratio	R Time	Scan#	Pk Height	Peak Area
PROPNOCARB	5.4	18.02	614	174	731

Xenes Report Number: XEN9737

Page 73




Appendix IV (continued)

Xenos Method XAM-34

28

Figure 10: Ion Chromatogram (m/z 189) of fortified control potato sample with Hoe 074189 at 0.05 ppm, 3 μ L of 5 mL injected.
POT-F1-2, Reference Data Sheet #: X9618V03



Integration Report Quan File: P115 Cali File: PROP Entries: 1
Comment: XEN9617\55ACB1\PROFMOCARD\POT-F1-2
Sorted via: Entry Number IS Factor: 1.000 Multi: 1.000 Div: 1.000

Name of Compound	R:N Ratio	R:Time	Scan#	Pk Height	Peak Area
PROFMOCARD	172.7	18.89	666	4.491	19.220

Xenos Report Number: XEN97-37

Page 74




Report No. AV96R001
Page 183

Appendix IV (continued)

Xmas Method XAM-34

29

Figure 11: Ion Chromatogram (m/z 189) of fortified control potato sample with Roc 074189 at 0.50 ppm, 3 μ L of 10 mL injected.
POT-F2-2, Reference Data Sheet #: X9618V04



Integration Report Quant File: PE08 Call File: PROP Entries: 1
Comment: XEN9617/96ACB1/PROPANOCARE/POT-F2-2
Sorted via: Entry Number ? IS Factor: 1.000 Multi: 1.000 Div: 1.000

Name of Compound	R/H Ratio	R Time	Scan#	Tk Height	Peak Area
PROPANOCARE	1723.8	10.88	605	36,318	125,207

Xmas Report Number: XEN961775




Report No. AV96R001
Page 184

Appendix IV (continued)

Xmas Method XAM-34

298

CALIBRATION CURVE FOR PROGAMOCARB UNDER GC/ION TRAP CONDITIONS
Datafile Name: X9615F04, R²=0.9922



Xmas Report Number: XEN97-37

Page 76



Report No. AV96R001

Page 185

Appendix IV (continued)

Xenos Method XAM-34

Xerox Report Number: XER97-37

Page 77

Appendix IV (continued)

Xenon Method XAM-34

30

Appendix III**GC/TSD - Representative Chromatograms, Typical Calibration Curve and Chromatographic Data Sheet**

- Figure 1:** Chromatogram of Hoe 074189 standard at 0.10 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01
- Figure 2:** Chromatogram of Hoe 074189 standard at 0.15 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01
- Figure 3:** Chromatogram of Hoe 074189 standard at 0.25 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01
- Figure 4:** Chromatogram of Hoe 074189 standard at 0.50 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01
- Figure 5:** Chromatogram of Hoe 074189 standard at 0.75 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01
- Figure 6:** Chromatogram of Hoe 074189 standard at 1.00 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01
- Figure 7:** Chromatogram of Hoe 074189 standard at 1.50 ng/ μ L, 2 μ L injected,
Reference Data Sheet #: X9618V01
- Figure 8:** Chromatogram of control potato sample, 2 μ L of 5 mL injected.
POT-C1, Reference Data Sheet #: X9618V01
- Figure 9:** Chromatogram of control potato sample, 2 μ L of 10 mL injected.
POT-C2, Reference Data Sheet #: X9618V02
- Figure 10:** Chromatogram of fortified control potato sample with Hoe 074189
at 0.05 ppm, 2 μ L of 5 mL injected.
POT-F1-2, Reference Data Sheet #: X9618V01
- Figure 11:** Chromatogram of fortified control potato sample with Hoe 074189
at 0.50 ppm, 2 μ L of 10 mL injected.
POT-F2-2, Reference Data Sheet #: X9618V02



Report No. AV96R001

Page 187

Appendix IV (continued)

Xmes Method XAM-34


21

Figure 1: Chromatogram of Hoe 074189 standard at 0.10 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01

```
Title      : Propamocarb in potato - XEN96-18, 96AC03
Run File   : C:\star\module1\X9618V01.RUN
Method File : C:\STAR\MODULES\PROFFOT.MTH
Sample ID   : 0.10 NG/UL

Injection Date: 31-OCT-96 1:11 AM    Calculation Date: 31-OCT-96 1:30 AM
Operator     : Xmes Laboratories        Detector Type: ADCB (1 Volt)
Workstation: MS-DOS 5                 Bus Address : 16
Instrument  : GC 3400                  Sample Rate : 10.00 Hz
Channel      : A = TSD                 Run Time   : 18.002 min

***** STAR Chromatography Workstation ***** Version 4.3 *****
Chart Speed = 1.00 cm/min Attenuation = 64  Zero Offset = 51
Start Time  = 0.000 min   End Time   = 17.200 min  Min / Tick = 1.00
```



Xmes Report Number: XEN97-07

Page 79




Report No. AV96R001
Page 188

Appendix IV (continued)

Xmas Method XMAS-34

32

Figure 2: Chromatogram of Hoe 074189 standard at 0.15 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01




Appendix IV (continued)

Xmas Method XAM-34

33

Figure 3: Chromatogram of Hoe 974189 standard at 0.25 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01






Appendix IV (continued)

Xenos Method XAM-34

34

Figure 4: Chromatogram of Hoe 074189 standard at 0.50 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01





Report No. AV96R001

Page 191

Appendix IV (continued)

Xmas Method XAM-34


35

Figure 5: Chromatogram of Hoe 07418 standard at 0.75 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01

```
Title      : Propiconazole in potato - XEN96-18, 96AC03
Run File   : c:\star\module1\X9618065.run
Method File : C:\STAR\MODULE1\PROPPOT.MTH
Sample ID  : 0.75 NG/UL

Injection Date: 31-OCT-96 8:08 AM    Calculation Date: 31-OCT-96 8:26 AM
Operator    : Xmas Laboratories        Detector Type: AODE (1 Volt)
Workstation: WS-0005                 Bus Address : 16
Instrument : GC 3400                  Sample Rate : 10.00 Hz
Channel    : A - TCD                 Run Time  : 18.002 min

***** Star Chromatography Workstation ***** Version 4.5 *****
Chart Speed = 1.00 cm/min Attenuation = 64    Zero Offset = 51
Start Time = 0.0000 min End Time = 17.200 min Min / Tick = 1.00
```



Xmas Report Number: XEN97-37

Page 35



Appendix IV (continued)

Xmas Method XAM-34


34

Figure 6: Chromatogram of Hoe 074189 standard at 1.00 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01

Title : Propamocarb in potato - XEN96-18, 96AC03
Run File : C:\STAR\module1\9618057.run
Method File : C:\STAR\MODULE18\PROPPOT.MTH
Sample ID : 1.0 NG/UL

Injection Date: 31-OCT-96 6:06 AM Calculation Date: 31-OCT-96 6:24 AM
Operator : Xmas Laboratories Detector Type: ADCB (1 Volt)
Workstation: HS-DOS 5 Bus Address : 18
Instrument : GC 340B Sample Rate : 10.00 Hz
Channel : A - TSD Run Time : 18.002 min

***** Star Chromatography Workstation ***** Version 4.5 *****
Chart Speed = 1.00 cm/min Attenuation = 64 Zero Offset = 51
Start Time = 0.000 min End Time = 17.200 min Min / Tick = 1.00



Xmas Report Number: XEN97-37

Page 14



Report No AV96R001

Page 193

Appendix IV (continued)


Xenos Method XAM-34

37

Figure 7: Chromatogram of Hos 07418 standard at 1.50 ng/ μ L, 2 μ L injected.
Reference Data Sheet #: X9618V01

Title : Propiconazole in potato - XEN96-18, 96AC03
Run File : c:\star\module1\X9618V01.run
Method File : C:\STAR\MODULE1\PROPPOT.MTH
Sample ID : 1.5 NG/G

Injection Date: 31-OCT-96 7:05 AM Calculation Date: 31-OCT-96 7:23 AM
Operator : Xenos Laboratories Detector Type: ADC8 (1 Volt)
Workstation: HS-POS 5 Bus Address : 18
Instrument : GC 3400 Sample Rate : 10.00 Hz
Channel : A = TSD Run Time : 18.002 min
***** Star Chromatography Workstation ***** Version 4.5 *****
Chart Speed = 1.00 cm/min Attenuation = 64 Zero Offset = 50
Start Time = 0.000 min End Time = 17.200 min Min / Tick = 1.00



Xenos Report Number: XEN97-37


Page 35

Appendix IV (continued)

Xenos Method XAM-34

38

Figure 8: Chromatogram of control potato sample, 2 μ L of 5 mL injected.
POT-C1, Reference Data Sheet #: X9618V01





Report No. AV96R001

Page 195


Appendix IV (continued)

Xmas Method XAM-34

39

Figure 9: Chromatogram of control potato sample, 2 μ L of 10 mL injected.
POT-C2, Reference Data Sheet #: X9618V02

Title : Propamocarb in potato - XEN96-18, 96AC03
Run File : c:\star\module1\9618044.run
Method File : C:\STAR\MODULE18\PROPPOT.MTH
Sample ID : POT-C2
Injection Date: 30-OCT-96 6:50 PM Calculation Date: 30-OCT-96 7:08 PM
Operator : Xmas Laboratories Detector Type: ADCB (1 Volt)
Workstation: MS-DOS 5 Bus Address : 18
Instrument : GC 3400 Sample Rate : 10.00 Hz
Channel : A - TSD Run Time : 18.002 min
***** star Chromatography Workstation ***** Version 4.5 *****
Chart Speed = 1.00 cm/min Attenuation = 64 Zero Offset = 54
Start Time = 0.000 min End Time = 17.200 min Min / Tick = 1.00



Xmas Report Number: X9617-37

Page 87


Appendix IV (continued)

Xmas Method XAM-34

40

Figure 10: Chromatogram of fortified control potato sample with Hec 074189
at 0.05 ppm, 2 μ L of 5 mL injected.
POT-F1-2, Reference Data Sheet #: X9618V01

Title : Propamocarb in potato - XEM96-18, 96AC03
Run File : c:\star\moduled1\9618062.run
Method File : C:\STAR\MODULE18\PROPOT.MTH
Sample ID : POT-F1-2
Injection Date: 31-OCT-96 3:39 AM Calculation Date: 31-OCT-96 3:38 AM
Operator : Xanos Laboratories Detector Type: ADCB (1 Volt)
Workstation: MS-DOS 5 Bus Address : 118
Instrument : GC 3400 Sample Rate : 10.00 Hz
Channel : A = TSD Run Time : 18.002 min
***** Star Chromatography Workstation ***** Version 4.5 *****
Chart Speed = 1.00 cm/min Attenuation = 64 Zero Offset = 58
Start Time = 0.000 min End Time = 17.200 min Min / Tick = 1.00



Xmas Report Number: X9607-07

Page 80


Appendix IV (continued)

Xmas Method XAM-34

41

Figure 11: Chromatogram of fortified control potato sample with Hoe 074189 at 0.50 ppm, 2 μ L of 10 mL injected.
POT-F2-2, Reference Data Sheet #: X9618VM2

Title : Propanocarb in potato - XEN96-18, 96Ac01
Run File : c:\star\module1\9618046.run
Method File : C:\STAR\MODULE1\PROPPOT.MTH
Sample ID : POT-F2-2
Injection Date: 30-OCT-96 8:47 PM Calculation Date : 30-OCT-96 9:06 PM
Operator : Xmas Laboratories Detector Type: A CB (Volt)
Workstation: MS-DOS 5 Bus Address : 1
Instrument : GC 3400 Sample Rate : 1.00
Channel : A = TSD Run Time : 1.002 min
***** Star Chromatography Workstation ***** arsi : 4.5 *****
Chart Speed = 1.00 cm/min Attenuation = 64 Zero Offset = 59
Start Time = 0.000 min End Time = 17.200 min Min / Tick = 1.00



Xmas Report Number: X96D7-37


Page 39

4

Appendix IV (continued)

Xmas Method XAM-34

41b



Xmas Report Number: XEN97-37

Page 98



Report No AV96R001

Page 199

Appendix IV (continued)

GAS CHROMATOGRAPHIC DATA SHEET

Xmas Project No: XXN96-18

Speaker Study No: Validation

Submitted: October 20, 1996


Received: October 26, 1996

Analyst: N. Bertoncello /26

Xmas Ref. No: X9611P91


Verifier:

GC Run No.	Sample ID	Sample Weight	Inj. Vol.	Inj. Vol.	Peak Area	Total Peak Area	Total Peak Area normalized	Peak Area (ppm)	Peak Area (ppm)	% Recovery
37	DIFP	0.100	10.00	1	2	8.0	20.74	ND	ND	ND
38	POT-C1	0.100	20.00	1	2	8.0	12.25	0.029	-0.001	100
39	POT-F1-1	0.100	20.00	1	2	8.0	25.24	0.300	0.300	100
40	POT-F1-2	0.200	20.00	1	2	8.0	34.07	0.393	0.371	95
41	POT-F1-3	0.500	20.00	1	2	8.0	53.89	0.590	0.590	100
42	POT-F1-4	0.750	20.00	1	2	8.0	74.70	0.787	0.787	100
43	POT-F1-5	1.000	20.00	1	2	8.0	95.53	0.984	0.984	100
44	POT-F1-6	1.500	20.00	1	2	8.0	126.36	1.252	1.252	100
45	POT-F1-7	0.500	20.00	1	2	8.0	42.07	0.489	0.489	100
46	POT-F1-8	0.750	20.00	1	2	8.0	59.97	0.667	0.667	100
47	POT-F1-9	1.000	20.00	1	2	8.0	76.83	0.835	0.835	100
48	POT-F1-10	1.500	20.00	1	2	8.0	103.13	1.000	1.000	100
49	POT-F1-11	2.000	20.00	1	2	8.0	134.00	1.293	1.293	100
50	POT-F1-12	3.000	20.00	1	2	8.0	211.92	1.800	1.800	100



Xmas Report Number: XXN97-37

Page 91

Appendix V Map of Trial Locations

EPA ADDENDUM

PP#6F04707

Propamocarb HCl on whole potatoes and tomatoes

1. ACB used a Hewlett-Packard 6890 *plus* Gas Chromatograph equipped with a nitrogen phosphorus detector (GC-NPD). The detector was interfaced to a TurboChrom (ver 6.1.1) data system. The analytical column was a fused silica DB-1, 15m, 0.53 mm i.d. (Megabore), 1.5 μ m film thickness (J&W Scientific). The column temperature was programmed according to the method. A Hewlett-Packard 7683 Automatic Liquid Sampler was used for sample injections. A 4mm splitless straight bore injection liner was used in the injection port. The injection port was maintained at 275 °C. The detector temperature was 300 °C. Helium was used as the carrier gas with a constant flow rate of ~5 mL/min @ 100 °C. The detector gas flows were set at the following: hydrogen @ 3 mL/min, air @ 60 mL/min, and nitrogen (make-up) @ 19 mL/min. A 2 μ L injection volume was used and the retention time under these conditions was approximately 8.8 min. The analysis time was 25 min.

2. The TurboVap®LV (Zymark Corporation, Hopkinton, MA) operated at 40 °C under a nitrogen stream was used in place of the vacuum rotary evaporators to concentrate samples during the final step only (the evaporation of 100 mL of DIPE to 2 mL). Initially, ACB used the TurboVap at 45 °C. to evaporated the methanol from the first extraction step rather than the Rotary Vacuum Evaporator set at 45 °C. This reduced the acidic aqueous to 5 mL (as recommended in the method) but took about 3 hours of time/sample. In a teleconference with the petitioner, it was suggested that 3 hours was too long and some loses might be occurring in that step. The recoveries were consistently lower when the TurboVap was used during the first step. Recoveries were improved when the rotary vacuum evaporator was used during the first step of the extraction. In our laboratory, it was difficult to achieve the 5 mL volume as recommended in the method. Our final volumes ranged from 10-12 mL. This created problems in the subsequent steps of the extraction.

3. ACB found, and the petitioner agreed, that the volume of 10N NaOH could be altered from the method (10 mL) and still yield good recoveries. Following the first step from the rotary vacuum evaporator, an additional 5 mL of 1N HCl is added to the extract to ensure a low pH, this acidic aqueous extract should be transferred to the 25 mL graduated cylinder before the base is added (the 20 mL volume cannot be exceeded). ACB found that 2-3 mL of base was enough to partition the propamocarb into the organic phase when the extract was poured onto the Chem Elute™ columns. If less than 2 ml of base was added to the acidic aqueous extract, low or no recoveries resulted. The Chem Elute™ columns are a gravity flow system (propamocarb is a very volatile compound and if a vacuum was applied to the Chem Elute™ columns,

EPA ADDENDUM

PP#6F04707

Propamocarb HCl on whole potatoes and tomatoes
Page 2

losses of analyte occurred). The 5 min equilibration timing should begin as soon as the meniscus from the 20 mL volume disappears onto the column.

4. The method suggests the use of calibration curves to determine sample concentrations. ACB ran calibration curves to demonstrate linearity over the dynamic range of the fortified samples. ACB determined sample concentrations from a ratio of sample responses to the average of standard responses that bracketed the samples during analysis. The bracketing standards were made to be the same concentration as the expected concentration of the fortified samples if 100% recovery were obtained. ACB observed some sample carry-over between injections. ACB found that it was important to use fresh washing solvent for both vial A and vial B in the autosampler tower at the beginning of a sample set. ACB suggests the use of blank injections of DIPE within a sample set to minimize carry-over. ACB used acetone for both washing solvents A and B, the method did not specify a solvent.

5. ACB determined that a good stopping point was after the first step of the extraction. The acidic aqueous extract could be stoppered and stored in the refrigerator overnight and the clean-up resumed the following day. Once the final extracts were prepared, they were assayed the same day (usually unattended overnight) with bracketing standards. The method did not give information regarding the stability of the final extract nor did they mention possible stopping points.

6. The fortification solutions and the calibration standard solutions are made from individual weighings of analytical reference standards. ACB recommends that both weighings occur under the same laboratory conditions (same day), since one set of solutions is measured against a different set of solutions.