

US EPA ARCHIVE DOCUMENT

3791



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

WASHINGTON, D.C. 20460

9/12/95
9/12/95 S. Ira C. Smith
Analytical Chemistry Section
Building 306, BARC-East
Beltsville, Maryland 20705

SEP 14 1995

OFFICE OF
PREVENTION, PESTICIDES AND
TOXIC SUBSTANCES

MEMORANDUM

SUBJECT: Glyphosate Trimesium in or on Corn, Soybeans, Citrus Fruit, Stone Fruit and the Nut Crop Group (except Almonds). (MRID#s: 432736-04 and -05. Barcode D216509. Chemical No 128501. CBTS# 15715.)

FROM: Douglas Swineford, Chemist ~~out~~
Analytical Chemistry Section

THRU: Harvey K. Hundley, Head
Analytical Chemistry Section

THRU: Donald A. Marlow, Chief ~~out~~
Analytical Chemistry Branch
Biological and Economic Analysis Division (7503W)

TO: F. D. Griffith, Acting Section Head
Tolerance Petition Section III
Chemistry Branch I, Tolerance Support
Health Effects Division (7509C)

Introduction

The Analytical Chemistry Laboratory (ACL) validated the method: "Determination of Residues of the Trimethylsulfonium Cation in Agricultural Crops by Gas Chromatography" (Study Number GLYP-93-AM-04/Document Number RR93-105B). The matrices validated were pecans, corn grain, corn fodder, soybean seed and oranges.

Method Summary

The samples were extracted with an aqueous phosphate buffer and filtered. The extract is treated with basic anion-exchange resin to neutralize acids followed by phenylisocyanate and barium hydroxide to remove amino acids. The samples were cleaned up with cation-exchange resin and dealkylated prior to quantitation by gas chromatography (GC) using a flame photometric detector (FPD).

Comments

1. Two experienced chemists can process a set of six samples in 12 hours, excluding the optional thermal cleanup (sec. 3.5) that may be necessary for some matrices. GC analysis time is approximately 5 minutes per injection.
2. The limit of detection from visual inspection of the chromatograms is estimated to be 0.01 ppm for the matrices tested.
3. No special safety hazards were noticed. Normal laboratory safety procedures were followed.
4. Standards were obtained from the registrant (Zeneca Ag Products, Richmond, CA 94840-0023). The EPA standards repository in RTP, North Carolina indicated they have a small amount of the standard available.
5. ACL has received various analytical methods this year for the validation of the trimethylsulfonium cation from Zeneca Ag Products for multiple raw agricultural commodities. Has the registrant considered the feasibility of combining the various approaches into a single method?
*Sept 12, 1991
Comments*
6. If the above comments are taken into consideration and incorporated with additional comments noted in the laboratory pre-review, the method would meet 40 CFR 158 and EPA's requirements as published in the Pesticide Assessment Guidelines, Subdivision "O" for Residue Chemistry, Part 171-4(b) as an enforcement method.
7. Recovery data, example chromatograms, pre-review checklist and pre-review comments are attached.

VALIDATION REPORT TRIMETHYLSULFONIUM CATION (TMS)

<u>Commodity</u>	<u>Chemical Added</u>	<u>PPM Added</u>	<u>PPM Found</u>	<u>%Recovery</u>
Corn Grain	Control	0	N.D.*	---
	Control	0	N.D.*	---
	TMS	0.05	0.0463	92.6
	TMS	0.05	0.0465	93.0
	TMS	0.10	0.0854	85.4
	TMS	0.10	0.0820	82.0

* N.D. = < 0.01 ppm

<u>Commodity</u>	<u>Chemical Added</u>	<u>PPM Added</u>	<u>PPM Found</u>	<u>%Recovery</u>
Corn Fodder	Control	0	N.D.*	---
	Control	0	N.D.*	---
	TMS	0.10	0.1030	103
	TMS	0.10	0.1021	102
	TMS	1.0	1.0119	101
	TMS	1.0	0.8960	89.6
	TMS	2.0	1.9367	96.8
	TMS	2.0	1.8282	91.4

N.D.* = < 0.01 ppm

VALIDATION REPORT TRIMETHYLSULFONIUM CATION (TMS) (continued)

<u>Commodity</u>	<u>Chemical Added</u>	<u>PPM Added</u>	<u>PPM Found</u>	<u>%Recovery</u>
Soybean Seed	Control	0	N.D.*	---
	Control	0	N.D.*	---
	TMS	0.05	0.0426	85.2
	TMS	0.05	0.0478	95.6
	TMS	0.50	0.4017	80.3
	TMS	0.50	0.4447	88.9
	TMS	1.0	0.9494	94.9
	TMS	1.0	0.9743	97.4

* N.D. = < 0.01 ppm

<u>Commodity</u>	<u>Chemical Added</u>	<u>PPM Added</u>	<u>PPM Found</u>	<u>%Recovery</u>
Oranges	Control	0	N.D.*	---
	Control	0	N.D.*	---
	TMS	0.05	0.0450	90.0
	TMS	0.05	0.0437	87.4
	TMS	0.50	0.5290	106
	TMS	0.50	0.5050	101
	TMS	1.0	1.1752	118
	TMS	1.0	1.1115	111

* N.D. = < 0.01 ppm

VALIDATION REPORT TRIMETHYLSULFONIUM CATION (TMS) (continued)

<u>Commodity</u>	<u>Chemical Added</u>	<u>PPM Added</u>	<u>PPM Found</u>	<u>%Recovery</u>
Pecans	Control	0	N.D.*	---
	Control	0	N.D.*	---
	TMS	0.05	0.0510	102
	TMS	0.05	0.0544	109
	TMS	0.10	0.1070	107
	TMS	0.10	0.1026	103

* N.D. = < 0.01 ppm

Modifications to method:

None.

Special precautions to be taken:

Handle all chemicals in a safe manner.

Source of analytical standard:

Zeneca Ag Products, Richmond, CA 94804-0023.

If derivatized standard, give source:

N/A

Instrument for quantitation:

GC/FPD

Instrument for confirmation:

N/A

If instrument parameters differ from those given in method, list parameters used:

N/A

Commercial sources of any special chemicals or apparatus:

N/A

Additional comments:

See report.

Chromatograms:

Copies attached.

TMV Pre-Review of Glyphosate- Trimesium

Reviewed By: Dallas Wright, Jr. *DW*

Date: June 15, 1995

Laboratory Assignment Number: B95-34,36,37,38,39

Analyte: Trimethylsulfonium Cation

Commodities: Pecans, Corn grain, *Corn Fodder*, Soybean Hay, Soybean Seed, and Oranges

Petitioner: Zeneca Ag Products

Method: "Determination of Residues of the Trimethylsulfonium Cation in Agricultural Crops by Gas Chromatography;" Dec. 30, 1993; RR 93-105B

Section 3: Analytical Procedure

3.2.2: In the note for this step there are suggestions for filtering the extract; however, there is no filtration step in the extraction for oily crops.

3.3.5: Is it necessary to continue the evacuation of the flask for 1-2 minutes as in step 3.3.4? Should the filter cake be removed before the second filtration? (0.95 hr. at they did 2.5 hr. regeneration)

3.5 Thermal Cleanup: Normally ACL safety regulations do not allow overnight heating or extractions. However, this step should be included in the validation. Is it possible to reduce the time to 8 or 9 hours? (1)

3.6 Dealkylation: 1) Since the toluene is added before the sample is heated, it is possible that the volume will change due to leaking during the heating. Since there is only 0.5 to 1 ml of toluene to start with, a small loss during heating could result in a large percentage change to the final volume used for calculations.

Section 5: Calculations: The instructions are misleading and seem to indicate that it would be appropriate to use the average response factor when using an FPD. The average response factor method for calculating cannot be used when the samples are quantitated on an FPD in sulfur mode, since the response is not linear. The method should clearly state that the analyst must use either a linear regression based on the square root of the response or closely match standard and sample responses.

Section 8.4 Lower Limit of Quantitation: 1) The LOQ for this method is 0.05 ppm for food crops and 0.1 ppm for feed items.

Section 9 Conclusions: The method has not been tested for samples less than the LOQ and greater than 10 times the LOQ. Some of the requested spiking levels are greater than 10 times the LOQ. The following table summarizes the data for each commodity requested, the LOQ, and the requested spiking levels:

Commodity	LOQ (PPM)	10x LOQ (PPM)	Data Provided (PPM)	Requested Levels (PPM)
Pecans	0.05	0.5	0.05; 0.1	0.05, 0.10
Corn Grain	0.05	0.5	0.05; 0.5	0.05, 0.10
Corn Fodder	0.1	1.0	0.1; 1.0	0.1; 1.0; 2.0
Soybean Seed	0.1	1.0	0.05; 0.1	0.05; 0.5; 1.0
Oranges	0.05	0.5	0.05; 0.09; 0.23	0.05; 0.5; 1.0

Recovery Data: Not all recovery data was within the range of 70-120%. The following list summarizes the data: Pecans: 62-112%; Corn grain: 81-118%; Corn Fodder: 85-115%; Soybean seed: 74-110%; Oranges: 104-125%

ILV Data: The Independent Laboratory Validation was run by Morse Laboratories, Inc. on corn grain (0.05 and 0.5 ppm), forage (0.1, and 1.0 ppm), and fodder (0.1, and 1.0 ppm). Acceptable recoveries were obtained on the first set of samples run.

ANALYTICAL CHEMISTRY BRANCH
SCREEN FOR RESIDUE METHODS FOR TMV

1. LABORATORY ASSIGNMENT NUMBER: B95-34; 3C-39
2. PP#: 9F03796, OF03866, 3F04238, 4F04343
3. TECHNICAL REVIEWER: Dallas Wright, Jr.
4. DATE: June 19, 1995
5. ANALYTES/LEVEL: See Review
6. COMMODITIES: Pecans, Corn Grain, Corn Fodder, Soybean Seed, Oranges
7. METHOD: Determination of Residues of the Trimethylsulfonium Cation in Agricultural crops by Gas Chromatography

The Analytical Chemistry Section has been asked to screen the residue chemistry methods submitted by the registrant in order to determine if they contain the essential requirements identified in the Residue Chemistry Guidelines. Full scientific review and laboratory evaluation of those methods will take place after the initial screen. The following items need to be resolved before the analytical method can be evaluated.

- | | <u>YES</u> | <u>NO</u> |
|--|------------|-----------|
| 1. Does the method use exotic equipment and/or supplies that are not commercially available in the U.S.? | <u> </u> | <u> </u> |
| 2. Does the method require any new equipment before the laboratory work begins? | <u> </u> | <u> </u> |
| 3. Are chromatograms included? <ol style="list-style-type: none"> a. Is (are) peak(s) of interest sufficiently resolved from other peaks? b. Has registrant included chromatograms of analyses at or below tolerance on all crop types for which tolerance is requested by HED? See Review c. Do the control samples have reasonably low levels of the analyte in relation to the proposed tolerance? d. Is the method sufficiently sensitive and specific to measure and identify the residues at levels specified by HED in the TMV request? | <u> </u> | <u> </u> |

YES NO

4. Has recovery data been provided to ACL for the residues that are specified in the TMV request? See Review ✓
5. Are recovery values between 70% and 120% at all levels and for all commodity types? See Review ✓
6. Are all procedures clearly written with no ambiguities so that the method can be run without communication with the registrant? See Review ✓
7. Does the method require correction for a sample of the untreated commodities or a blank? ✓
8. Does the method require the use of an internal or procedural standard to compensate for lost analyte during analysis? ✓
9. Are 2nd laboratory validation data provided with the method? ✓
10. Are there any deficiencies other than those covered above that would prevent ACS from conducting a method trial? ✓
11. Is this method suitable for validation testing? ✓

Any deficiencies/problems noted for any above items should be addressed in the full scientific review of this method to be attached as an addendum.

Dallen Wynn

Signature

6/19/95

Date

The following is to be completed by the analyst performing the TMV.

12. Can a set of 6 samples be run within 24 hours? ✓
13. a. Are standards available at RTP repository? ✓
- b. Are derivatized analytical reference standards available? N/A


B95-34-35

**Linearity of TMS
GLC/FPD
9/11/95**

ug/mL	PK. HT. 1	PK. HT.2	AVE PK. HT.	Sqr Avg Pk Ht	Sqr CALC. PK. HT.
0.1250	18082.0	17314.0	17698.0	133.0338	
0.2500	29935.0	29323.0	29629.0	172.1308	
0.5000	128467.0	146861.0	137664.0	371.0310	
1.0000	1223028.0	1223179.0	1223103.5	1105.9401	
0.0200					-69.6693
1.2000					1285.0712

Regression Output:

Constant	-92.6310293
Std Err of Y Est	104.8272386
R Squared	0.964211424
No. of Observations	4
Degrees of Freedom	2
X Coefficient(s)	1148.08523
Std Err of Coef.	156.4030464



RUN # 545 AUG 5, 1995 13:32:16
START

IF
IF
IF 1.163 ←

TIMETABLE STOP

RUN# 545 AUG 5, 1995 13:32:16

SAMPLE NAME: DMSS STD SAMPLE# 1
0.25 UG/ML

B95/34 THRU 39 (SOYBEAN SEED) ZUL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
<u>TMS</u> 1.163	49240	BU	.015	100.00000

TOTAL HEIGHT= 49240
MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT 2^ = 4
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 546 AUG 5, 1995 13:35:56
START

IF
IF
IF 1.215

TIMETABLE STOP

RUN# 546 AUG 5, 1995 13:35:56

SAMPLE NAME: CONTROL #1 SAMPLE# 2
FINAL VOL 1ML/5G ALIQUOT

B95/34 THRU 39 (SOYBEAN SEED) ZUL INJ

NO RUN PEAKS STORED

RUN PARAMETERS

13

RUN PARAMETERS

ZERO = 20
ATT 2^ = 4
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 546 AUG 5, 1995 13:35:56
START

IF
IF
IF 1.215

TIMETABLE STOP

RUN# 546 AUG 5, 1995 13:35:56

SAMPLE NAME: CONTROL #1 SAMPLE# 2
FINAL VOL 1ML/5G ALIQUOT

B95/34 THRU 39 (SOYBEAN SEED) 2UL INJ

NO RUN PEAKS STORED

RUN PARAMETERS

ZERO = 20
ATT 2^ = 4
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 547 AUG 5, 1995 13:39:36
START

IF
IF
IF 111852

TIMETABLE STOP

RUN# 547 AUG 5, 1995 13:39:36

SAMPLE NAME: CONTROL #1 SAMPLE# 2
FINAL VOL 1ML/5G ALIQUOT

B95/34 THRU 39 (SOYBEAN SEED) 2UL TNT

14

RUN PARAMETERS

ZERO = 20
ATT 2^A = 4
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 554 AUG 5, 1995 14:05:18
START

IF

IF

IF

1.155 ←

TIMETABLE STOP

RUN# 554 AUG 5, 1995 14:05:18

SAMPLE NAME: 0.05PPM SPK# SAMPLE# 6
FINAL VOL 1ML/5G ALIQUOT

B95/34 THRU 39 <SOYBEAN SEED> 2UL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.155	40259	PU	.018	100.0000

TOTAL HEIGHT= 40259
MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT 2^A = 4
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 555 AUG 5, 1995 14:09:01
START

IF

IF

IF

1.155

TIMETABLE STOP

RUN# 555 AUG 5, 1995 14:09:01

SAMPLE NAME: 0.05PPM SPK# SAMPLE# 6
FINAL VOL 1ML/5G ALIQUOT

TOTAL HEIGHT= 27393
MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT Z^ = 4
CHT SP = 1.0
RR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 560 AUG 5, 1995 14:27:22
START

IF
IF
IF

1.152 ←

TIMETABLE STOP

RUN# 560 AUG 5, 1995 14:27:22

SAMPLE NAME: 0.5PPM SPK#2 SAMPLE# 9
FINAL VOL 10ML/5G ALIQUOT

B95/34 THRU 39 (SOYBEAN SEED) 2UL INJ

HEIGHTX

RT	HEIGHT	TYPE	WIDTH	HEIGHTX
1.152	34198	PU	.017	100.00000

TOTAL HEIGHT= 34198
MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT Z^ = 4
CHT SP = 1.0
RR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 561 AUG 5, 1995 14:31:03
START

IF
IF
1.061 ← 1.152
IF

14

TOTAL HEIGHT= 32333
 MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
 ATT 2^ = 4
 CHT SP = 1.0
 AR REJ = 6000
 THRSH = 3
 PK WD = 0.02

RUN # 567 AUG 5, 1995 14:53:11
 START

IF

IF

IF

1.148



TIMETABLE STOP

RUN# 567 AUG 5, 1995 14:53:11

SAMPLE NAME: 1PPM SPK#2 SAMPLE# 12
 FINAL VOL 20ML/5G ALIQUOT

B95/34 THRU 39 <SOYBEAN SEED> 2UL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.148	32041	PU	.018	100.00000

TOTAL HEIGHT= 32041
 MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
 ATT 2^ = 4
 CHT SP = 1.0
 AR REJ = 6000
 THRSH = 3
 PK WD = 0.02

RUN # 568 AUG 5, 1995 14:56:48
 START

IF

IF

IF

1.148

TIMETABLE STOP

B95/34 THRU 39 < CORN FODDER > ZUL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.174	95550	UU	.015	100.00000

TOTAL HEIGHT= 95550

MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT 2^o = 4
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 371 JUL 27, 1995 16:34:23
START

IF

IF

IF

1.175

TIMETABLE STOP

RUN# 371 JUL 27, 1995 16:34:23

SAMPLE NAME: DMS STD SAMPLE# 7
0.25 UG/ML *zul inj*

B95/34 THRU 39 < CORN FODDER > ZUL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
<i>TMS</i>	1.175	84906	PB	.014 100.00000

TOTAL HEIGHT= 84906

MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT 2^o = 4
CHT SP = 1.0
-- -- --

18

RUN PARAMETERS

ZERO = 20
ATT Z^ = 4
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 361 JUL 27, 1995 15:39:45
START

IF

IF

11.1075 ←
IF

TIMETABLE STOP

RUN# 361 JUL 27, 1995 15:39:45

SAMPLE NAME: CONTROL #1 SAMPLE# 2
FINAL VOL 2ML/5G ALIQUOT

B95/34 THRU 39 (CORN FODDER) 2UL INJ

NO RUN PEAKS STORED

RUN PARAMETERS

ZERO = 20
ATT Z^ = 4
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 362 JUL 27, 1995 15:45:10
START

IF

IF

11.1085
IF

TIMETABLE STOP

RUN# 362 JUL 27, 1995 15:45:10

SAMPLE NAME: CONTROL #2 SAMPLE# 3

RUN PARAMETERS

ZERO = 20
ATT Z^ = 4
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 367 JUL 27, 1995 16:12:25
START

IF

IF

~~1.121~~

IF

1.172

TIMETABLE STOP

RUN# 367 JUL 27, 1995 16:12:25

SAMPLE NAME: 0.1 SPK#1 SAMPLE# 5
FINAL VOL 2ML/5G ALIQUOT

B95/34 THRU 39 <CORN FODDER> 2UL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.121	6857	PU	.015	7.12342
1.172	89403	UB	.014	92.87658

TOTAL HEIGHT= 96260

MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT Z^ = 4
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 368 JUL 27, 1995 16:17:57
START

IF

IF

IF

1.172

70

ATT 2^ = 1
 CHT SP = 1.0
 AR REJ = 6000
 THRSH = 3
 PK WD = 0.02

RUN # 373 JUL 27, 1995 16:45:17
 START

IF

IF

IF

1.176 ↙

TIMETABLE STOP

RUN# 373 JUL 27, 1995 16:45:17

SAMPLE NAME: 1PPM SPK#1 SAMPLE# 8
 FINAL VOL 20ML/5G ALIQUOT

B95/34 THRU 39 (CORN FODDER) 2UL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.176	87658	PB	.014	100.00000

TOTAL HEIGHT= 87658
 MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
 ATT 2^ = 4
 CHT SP = 1.0
 AR REJ = 6000
 THRSH = 3
 PK WD = 0.02

RUN # 374 JUL 27, 1995 16:50:45
 START

IF

IF

IF

1.174

TIMETABLE STOP

21

RUN PARAMETERS

ZERO = 20
ATT Z^ = 4
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 380 JUL 27, 1995 17:23:36
START

IF

IF

IF

1.180 ↙

TIMETABLE STOP

RUN# 380 JUL 27, 1995 17:23:36

SAMPLE NAME: 2PPM SPK#2 SAMPLE# 12
FINAL VOL 40ML/5G ALIQUOT

B95/34 THRU 39 (CORN FODDER) 2UL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.180	76587	PB	.013	100.00000

TOTAL HEIGHT= 76587
MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT Z^ = 4
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 381 JUL 27, 1995 17:29:03
START

IF

IF

IF

1.178 ↙

22

HED Records Center Series 361 Science Reviews - File R061972 - Page 191 of 409

RUN PARAMETERS
ZERO = 20
RTT 2^ = 3
CHT SP = 1.0
RR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 147 JUL 12, 1995 09:24:13
START

IF

IF

1.115

IF

1.165

← DMS (TMS)

TIMETABLE STOP

RUN# 147 JUL 12, 1995 09:24:13

SAMPLE NAME: TMS STD SAMPLE# 1
0.25 UG/ML *2 ul ing*

B95-34 THRU 39(CORN GRAIN)

HEIGHT%

TMS RT HEIGHT TYPE WIDTH HEIGHT%
1.165 28482 PB .015 100.0000

TOTAL HEIGHT= 28482

MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
RTT 2^ = 3
CHT SP = 1.0
RR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 148 JUL 12, 1995 09:29:43
START

IF

IF

1.170

IF


TIMETABLE STOP

23

RUN PARAMETERS

ZERO = 20
ATT 2^ = 3
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 148 JUL 12, 1995 09:29:43
START



TIMETABLE STOP

RUN# 148 JUL 12, 1995 09:29:43

SAMPLE NAME: CONTROL#1 SAMPLE# 2
2UL INJ/FINAL VOL 1ML/5G ALIQUOT

B95-34 THRU 39<CORN GRAIN>

NO RUN PEAKS STORED

RUN PARAMETERS

ZERO = 20
ATT 2^ = 3
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 149 JUL 12, 1995 09:35:12
START

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.164	25518	UU	.016	100.00000

TOTAL HEIGHT= 25518
MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
RTT 2^ = 3
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

RUN # 154 JUL 12, 1995 10:02:31
START

IF
IF
1.14867 1.164 ←
IF

TIMETABLE STOP

RUN# 154 JUL 12, 1995 10:02:31

SAMPLE NAME: 0.05PPM SPK1 .SAMPLE# 5
2UL INJ/ FINRL VOL 1ML/ 5G ALIQUOT

B95-34 THRU 39(CORN GRAIN)

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.164	23776	UU	.015	100.00000

TOTAL HEIGHT= 23776
MUL FACTOR=1.0000E+00

RUN PARAMETERS


ZERO = 20
RTT 2^ = 3
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK WD = 0.02

25

RUN PARAMETERS

ZERO = 20
ATT 2^ = 3
CHT SP = 1.0
AR REJ = 6000
THRSH = 3
PK W0 = 0.02

RUN # 161 JUL 12, 1995 10:41:03
START



RUN# 161 JUL 12, 1995 10:41:03

SAMPLE NAME: 0.1SPK1 SAMPLE# 8
2UL INJ/ FINAL VOL2ML/ 5G ALIQUOT,

B95-34 THRU 39(CORN GRAIN)

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.163	23813	PU	.016	100.00000

TOTAL HEIGHT= 23813
MUL FACTOR=1.0000E+00

RUN PARAMETERS


He

TOTAL HEIGHT= 158790
 MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
 ATT 2^ = 5
 CHT SP = 1.0
 RR REJ = 2000
 THRSH = 3
 PK WD = 0.06

RUN # 1217 AUG 29, 1995 10:54:38
 START



TIMETABLE STOP

Storing processed peaks to A:Q14398EF.PRO

RUN# 1217 AUG 29, 1995 10:54:38

SAMPLE NAME: 0.25UG/ML SAMPLE# 1
 DMS-STD

PERK FILE : A:Q14398EF.PRO

B95/34-39 (PECANS) ZUL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.013	8814	BB	.015	6.79301
<i>TMS</i>	118403	BB	.022	91.25402
1.252	2534	BB	.016	1.95297

TOTAL HEIGHT= 129751

MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
 ATT 2^ = 5
 CHT SP = 1.0
 RR REJ = 2000
 THRSH = 3
 PK WD = 0.06

RUN # 1218 AUG 29, 1995 10:59:49
 START

IF

27

RUN PARAMETERS

ZERO = 20
ATT 2^ = 5
CHT SP = 1.0
AR REJ = 2000
THRSH = 3
PK WD = 0.06

RUN # 1221 AUG 29, 1995 11:15:13
START

IF

1.088 ←
1.250
1.759

TIMETABLE STOP

Storing processed peaks to A:Q1439DC2.PRO

RUN# 1221 AUG 29, 1995 11:15:13

SAMPLE NAME: CONTROL#2 SAMPLE# 3
FINAL VOL 1ML/5G ALIQUOT

PEAK FILE : A:Q1439DC2.PRO

B95/34-39 (PECANS) 2UL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.088	10119	PB	.017	66.63373
1.250	2913	BV	.031	19.18214
1.759	2154	PU	.162	14.18412

TOTAL HEIGHT= 15186

MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT 2^ = 5
CHT SP = 1.0
AR REJ = 2000
THRSH = 3
PK WD = 0.06

RUN # 1222 AUG 29, 1995 11:20:24
START

IF

1.018 ←
1.250
1.674

1.139

TIMETABLE STOP

28

RUN PARAMETERS

ZERO = 20
ATT 2^ = 5
CHT SP = 1.0
AR REJ = 2000
THRSH = 3
PK WD = 0.06

RUN # 1232 AUG 29, 1995 12:34:18
START

IF

1.010

1.251

1.139

←

IT

TIMETABLE STOP

Storing processed peaks to R:Q143B04B.PRO

RUN# 1232 AUG 29, 1995 12:34:18

SAMPLE NAME: .0.05PPM SPK#/ SAMPLE# 5
FINAL VOL 1ML/5G ALIQUOT

PEAK FILE : R:Q143B04B.PRO

B95/34-39 <PECANS> ZUL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.010	12107	PP	.020	7.97368
1.088	9140	PP	.014	6.01961
<u>1.139</u>	<u>127017</u>	PB	.021	83.65354
1.251	3573	PB	.023	2.35318

TOTAL HEIGHT= 151837

MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT 2^ = 5
CHT SP = 1.0
AR REJ = 2000
THRSH = 3
PK WD = 0.06

RUN # 1233 AUG 29, 1995 12:39:39
START

IF

1.010

1.251

1.139


IT

29

RUN PARAMETERS

ZERO = 20
ATT 2^ = 5
CHT SP = 1.0
RR REJ = 2000
THRSH = 3
PK WD = 0.06

RUN # 1239 AUG 29, 1995 13:11:46
START



TIMETABLE STOP

Storing processed peaks to A:Q143B913.PRO

RUN# 1239 AUG 29, 1995 13:11:46

SAMPLE NAME: 0.1PPM SPK#1 SAMPLE# 8
FINAL VOL 2ML/5G ALIQUOT

PEAK FILE : A:Q143B913.PRO

B95/34-39 (PECANS) 2UL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.137	<u>136306</u>	PB	.022	97.70829
1.248	3197	BB	.013	2.29171

TOTAL HEIGHT = 139503

MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT 2^ = 5
CHT SP = 1.0
RR REJ = 2000
THRSH = 3
PK WD = 0.06

RUN # 1240 AUG 29, 1995 13:17:06
START



TIMETABLE STOP

30

1.165	110459	PB	.020	43.91850
1.653	124705	PB	.048	49.58272
2.207	1251	PP	.034	.49740

TOTAL HEIGHT= 251509
MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
RTT 2^ = 5
CHT SP = 1.0
AR REJ = 2000
THRSH = 3
PK WD = 0.04

RUN # 1075 AUG 19, 1995 13:40:50
START

IF

1.24517

1.163

1.651

2.203

TIMETABLE STOP

Storing processed peaks to A:Q13690E3.PRO

RUN# 1075 AUG 19, 1995 13:40:50

SAMPLE NAME: 0.25UG/ML SAMPLE# ?
DMS STD

PEAK FILE : A:Q13690E3.PRO

B95/34-39 (ORANGES) 2UL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.117	15068	PU	.017	6.15988
<i>TMS</i>	118768	UU	.021	48.55302
1.245	13158	UB	.021	5.37907
1.651	96244	PB	.027	39.34609
2.203	1377	PP	.052	.56293

TOTAL HEIGHT= 244615
MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
RTT 2^ = 5
CHT SP = 1.0

31

RUN # 1064 AUG 19, 1995 12:36:40
START

IF

0 .865
1.123 ←
1.251

1.660

2.214

TIMETABLE STOP

Storing processed peaks to R:Q1368109.PRO

RUN# 1064 AUG 19, 1995 12:36:40

SAMPLE NAME: CONTROL#1 SAMPLE# 2
FINAL VOLUME/5G ALIQUOT

PEAK FILE : R:Q1368109.PRO

B95/34-39 <ORANGES> 2UL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
.865	2353	BU	.142	1.65205
1.123	16199	UU	.016	11.37339
1.251	13339	UU	.021	9.36537
1.660	109064	PB	.052	76.57430
2.214	1474	UP	.039	1.03490

TOTAL HEIGHT= 142429

MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT 2^ = 5
CHT SP = 1.0
RR REJ = 2000
THRSH = 3
PK WD = 0.04

RUN # 1065 AUG 19, 1995 12:42:32
START

IF

0 .891
1.255

1.665

2.225

TIMETABLE STOP


32

TOTAL HEIGHT= 109984
 MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
 ATT Z^ = 5
 CHT SP = 1.0
 RR REJ = 2000
 THRSH = 3
 PK WD = 0.04

RUN # 1071 AUG 19, 1995 13:17:02
 START



TIMETABLE STOP

Storing processed peaks to A:Q1368B4F.PRO

RUN# 1071 AUG 19, 1995 13:17:02

SAMPLE NAME: 0.05PPM SPK# SAMPLE# 5
 FINAL VOL 1ML/5G ALIQUOT

PERK FILE : A:Q1368B4F.PRO

895/34-39 (ORANGES) 2UL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
.731	2239	UP	.135	.91650
1.168	103408	BU	.021	42.32862
1.250	8494	UB	.021	3.47690
1.656	128886	BB	.056	52.75770
2.209	1271	BU	.045	.52027

TOTAL HEIGHT= 244298

MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
 ATT Z^ = 5
 CHT SP = 1.0
 RR REJ = 2000
 THRSH = 3
 PK WD = 0.04

33

RUN PARAMETERS

ZERO = 20
ATT 2ⁿ = 5
CHT SP = 1.0
AR REJ = 2000
THRSH = 3
PK WD = 0.04

RUN # 1091 AUG 19, 1995 14:57:23
START

IF
0.537 ,
1.112 ←
1.239 ← 1.156 ←
TIMETABLE STOP 1.645

Storing processed peaks to A:Q136A2D4.PRO

RUN# 1091 AUG 19, 1995 14:57:23

SAMPLE NAME: 0.5PPM SPK#1* SAMPLE# 8
FINAL VOL 10ML/5G ALIQUOT

PEAK FILE : A:Q136A2D4.PRO

B95/34-39 <ORANGES> ZUL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
.537	1896	BP	.110	.79938
1.112	13316	PU	.017	5.61418
1.156	75239	UU	.023	31.72165
1.239	11545	UU	.023	4.06751
1.645	135109	BP	.122	56.99728

TOTAL HEIGHT= 237185
MUL FACTOR=1.0000E+00

RUN PARAMETERS

ZERO = 20
ATT 2ⁿ = 5
CHT SP = 1.0
AR REJ = 2000
THRSH = 3
PK WD = 0.04

RUN # 1092 AUG 19, 1995 15:01:23
START

IF
0.460

34

NO RUN PEAKS STORED

RUN PARAMETERS

ZERO = 20
ATT Z^ = 5
CHT SP = 1.0
AR REJ = 2000
THRSRH = 3
PK WD = 0.04

* EDIT SEQ

- 1 = ALS INFORMATION
- 2 = EQUILIBRATION TIME DELAY
- 3 = METHOD FILE SPECIFICATION
- 4 = SAMPLE INFORMATION TABLE

SECTION TO BE EDITED: 1

ALS INFORMATION

INET SAMPLER CONTROL [Y*/N]:

7673A SAMPLER:

LOOP ADDRESS: 8

FRONT INJECTOR

INJ/BOTTLE	2 -->	@
FIRST BOTTLE	7 -->	12 @
LAST BOTTLE	13 -->	BREAK

* SEQ START

RUN # 1099 AUG 19, 1995 15:27:39
START

IF

3.6 min
1.175055 1.098

TIMETABLE STOP

1.566

Storing processed peaks to A:Q136A9EC.PRO

RUN# 1099 AUG 19, 1995 15:27:39

SAMPLE NAME: 1PPM SPK#2 SAMPLE# 12
FINAL VOL 20ML/5G ALIQUOT

PEAK FILE : A:Q136A9EC.PRO

895/34-39 (ORANGES) 2UL INJ

HEIGHT%

RT	HEIGHT	TYPE	WIDTH	HEIGHT%
1.055	16919	PV	.016	7.69164
1.098	62111	UU	.023	28.23663
1.175	14877	UB	.019	6.39962
1.566	126859	PB	.056	57.67210

TOTAL VOLUME: 21.00000

35