US ERA ARCHIVE DOCUMENT



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

WASHINGTON, D.C. 20460

Analytical Chemistry Section Building 306, BARC-East Beltsville, Maryland 20705

JUN 1 0 1993

OFFICE OF PREVENTION, PESTICIDES AND TOXIC SUBSTANCES

MEMORANDUM

SUBJECT: Product Chemistry method validation of the analytical procedure for the determination of the Active Ingredient Sumilarv content in a Technical material sample.

Other names for Sumilarv are pyriproxyfen (BSI proposed) and S-31183. The chemical name is 4-phenoxyphenyl (RS)-2-(2-pyridyloxy)propyl ether, 2-[1-methyl-2-(4-phenoxyphemoxy)ethoxy]pyridine.

Registrant: Sumitomo Chemical ACS Project No. 1992-AI-07

FROM: Gregory P. Verdin, Physical Science Technician

Analytical Chemistry Section

THRU: Harvey K. Hundley, Head

Analytical Chemistry Section

THRU: Donald A. Marlow, Chief

Analtyical Chemistry Branch

TO: Alfred Smith, Chemist

Product Chemistry Review Branch Registration Support Branch Registration Division (H7505C)

As required by §62-3 and §64-1 of the <u>Pesticide Assessement Guidelines</u>, <u>Subdivison D</u>, Product Chemistry, regarding submission of analytical methods and samples for analysis, the Registrant, in compliance with the above requirements, has submitted the following:

- 1) "Method for the Determination of Sumilarv", from the Study:
 Analytical Methods to Verfiy Certified Limits of Sumilarv
 Technical Grade, Ken-ichi Takahashi, June 14,1988, Laboratory
 Project ID NO. NNA-90-0011, Sumitomo Chemical Co., Ltd.,
 Environmental Health Science Laboratory, Osaka, Japan



ACL reviewed and used the method submitted to assay the sample provided. A summary and any variations made will be noted under a Comments Section.

Five portions of the sample were analyzed by the method provided by the Registrant. The results are below.

<u>Material</u>	Assay Values(%)	Avg(%)	SD	<u>Precision</u>
Technical Sumilarv [™]	98.4, 98.6, 99.0, 99.1, 98.7	98.8	0.29	0.29

Comments:

The method utilizes HPLC analysis with a reverse-phase, C-18 column and an ultra-violet detector wavelength set at 254 nm. 4-Benzyldiphenyl is the internal standard. Active Ingredient quantitation is determined from response ratios between the compound of interest (active ingredient) and the internal standard. All data acquisition was by electronic integration.

ACL determined that the response ratios of the analytical standard to a constant weight of the internal standard are linear over a concentration range of 2.60 - 8.13 mg/ml, which the method recommended.All analyses above are within this range.

Acceptable results with good precision are achieved by ACL using a Zorbax-ODS C-18 column, which is equivalent to the Lichrosorb RP-18 column suggested in the method. An isocratic system of 2:1 acetonitrile/water mobile phase is used for analysis. The flow rate, injection volume, and other instrument parameters are adjusted to achieve reproducible retention times and separations.

Precision and results for the technical materials assayed by this method are within limits suggested by the <u>Pesticide Assessment</u> <u>Guidelines</u>, <u>Subdivision D</u>, Product Chemistry and the Protocol for the Active Ingredient Program.

This method appears to produce acceptable results for the assay of technical materials of Sumilarv $^{\text{m}}$.

ACL will dispose of the standards and samples used for this method validation in safe and appropriate manner thirty (30) days from the date of this report unless we receive other specific instructions from you.

Contact the analyst or Adrian W. Burns at (301)504-8133 if there are any questions concerning this method validation.

2