

US EPA ARCHIVE DOCUMENT

430444-01

STUDY TITLE

Revised Residue Analytical Method for Parent RH-7592 and
Its Lactone Metabolites RH-9129 and RH-9130 in Stonefruit

DATA REQUIREMENT

Guideline 171-4 (c)

AUTHOR

John J. Martin

REVISION by

Theodore F. Burnett

STUDY COMPLETED ON

December 8, 1993

PERFORMING LABORATORY

Rohm and Haas Company
727 Norristown Road
Spring House, PA 19477

LABORATORY PROJECT ID

Rohm and Haas Technical Report No. 34-90-47R

Revisions to this document are explained on Page 14A

Revision to Report No. 34-90-47 (MRID 418750-38)

Page 1 of 32

Including pages 2a, 11A, 14a

Rohm and Haas Company Report Number 34-90-47R
Revision Number 1 to Rohm and Haas Report Number 34-90-47
EPA MRID Number 418750-38

Statement of No Confidentiality Claims

No claim of confidentiality is made for any information contained in this document on the basis of its falling within the scope of FIFRA § 10(d)(1)(A), (B), or (C).


Company: Rohm and Haas Company
Company Agent: Richard D. Costlow, Ph.D., D.A.B.T.
Title: Product Registration Manager
Date: 8 December, 1993

SIGNATURE: _____


Richard D Costlow

GLP Compliance Statement

This revised report contains a tolerance enforcement method, and as per 40 CFR 160.3 method development is not required to be conducted in compliance with GLP. However the work was conducted in a laboratory facility that is in compliance with Good Laboratory Practice as defined by the United States Environmental Protection Agency.



T. F. Burnett Date
Revisionist/Study Director
Rohm and Haas Company


S. S. Stawinski Date
Applicant/Submitter
Rohm and Haas Company


S. S. Stawinski Date
Sponsor
Rohm and Haas Company

Quality Assurance Statement

This report revision has been reviewed by the Quality Assurance Unit of the Rohm and Haas Company Agricultural Research Division for compliance with relevant SOPs and has been verified as a true and accurate representation of the data collected.


R. S. Krick Date
QA Analyst
Rohm and Haas Company

Technical Report
No. 34-90-47R

Table of Contents

	Page
I. Summary	5
II. Introduction	8R
III. Method	8R
A. Chemicals/Supplies	8R
B. Preparation of Solutions	9
C. Equipment	9
D. Instrumentation	10
E. Analytical Procedure	10
1. Sample Processing	10
2. Extraction	10
3. Partition	11R
4. Cleanup	11R
a. Standardization of Silica Gel elution	11A
b. Silica Gel Column Chromatography	11A
c. Standardization of Florisil elution	11A
d. Florisil Column Chromatography	11A
5. Gas Chromatography	11A
a. Preparation of Standard Curves	12R
b. Sample Analysis	12R
6. Method of Calculation	12R
a. Fortification Recovery	12R
b. Component Residue Concentration	12R
IV. Results and Discussion	13
V. Project Information and Study Certification	14A
VI. Method Revisions	

Tables

I. Summary of Recovery Data	15
II. Recovery Data for RH-7592, RH-9130 and RH-9130 in Cherry	16

Figures

1. 1.0 µg/ml Calibration Standards	17
2. 0.5 µg/ml Calibration Standards	18
3. 0.1 µg/ml Calibration Standards	19
4. 0.05 µg/ml Calibration Standards	20
5. Typical Calibration Curve for RH-7592	21
6. Typical Calibration Curve for RH-9130	22
7. Typical Calibration Curve for RH-9129	23
8. Control Cherry Sample No. 90-0066-009	24
9. Control Cherry Sample No. 90-0066-009 Fortified at 0.4 ppm	25
10. Treated Cherry Sample No. 90-0066-008 (5 ml)	26
11. Treated Cherry Sample No. 90-0066-008 (20 ml)	27
12. Reagent B lank	28
13. Reagent Fortification	29

Technical Report
No. 34-90-47R

Appendix

A. Material Safety Data Sheet

PROPERTY OF ROHM AND HAAS COMPANY PHILADELPHIA

I. Summary

The RH-7592 residues (parent and the two metabolites) are extracted from stonefruit by blending with methanol. The methanol extract is transferred to a 500 ml separatory funnel and partitioned with 10% sodium chloride solution and methylene chloride. The methylene chloride is evaporated to dryness and the samples are cleaned up on silica gel and Florisil.

The concentrations of RH-7592 and the metabolites are determined by gas chromatography using a 0.53 mm ID capillary column (SPB-608) and a capillary thermionic specific detector optimized for nitrogen selectivity. A flow diagram of the method is shown on page 7.


The sensitivity of the method is 0.01 ppm for all compounds in all substrates.

II. Introduction


RH-7592 is an experimental triazole fungicide being developed for use on stonefruit, wheat, apples, almonds and other crops. In order to obtain commercial registration an analytical method is needed to determine the magnitude of the total toxic residues from the proposed uses of RH-7592 and to provide a means of obtaining data for the setting and enforcement of tolerances for residues in food and feed.

Field studies conducted with radiolabeled RH-7592 have shown that the residues of concern resulting from the use of RH-7592 on stonefruit are extractable by blending with methanol and consist mainly of parent RH-7592 and lesser amounts of the lactone, (RH-9129). Since there are two unsymmetrical substituted quaternary carbons in the lactone structure, this material can exist as two diastereomeric isomers which can have different physical and chemical properties. This method measures parent RH-7592 and the two lactone isomers RH-9129 and RH-9130. The structures of these compounds are shown on the following page.

Technical Report
No. 34-90-47R



α -(2-[4-chlorophenyl]ethyl)- α -phenyl 1-3-(1H-1,2,4-triazole)-1-propanenitrile



Flow Diagram

Blender Extraction with Methanol or Toluene/Methanol

Partition with 10% NaCl and Methylene Chloride

Silica Gel Column Chromatography

Florisil Column Chromatography

Gas Liquid Chromatography

Quantitation

Technical Report
No. 34-90-47R

III. METHOD

A. Chemicals/Supplies

Acetone, Pesticide Grade
Silica gel, 40-140 mesh (Activate at 200 °C
for 24 hours. Bottle and cap immediately
Store in a dessicator cabinet
Standardize before use, Section IV.)
Celite-545
Cotton, Sterile, Absorbent
Whatman #3 Filter Paper (7.0 cm, Qualitative)
Florisil 60-100 Mesh (Activate at 200°C
for 24 hours. Bottle and cap immediately
Store in a dessicator cabinet
Standardize before use. Section IV.)
Methanol, Pesticide Grade
Methylene Chloride, Pesticide Grade
Sodium Chloride, Certified A.C.S.
Sodium Sulfate, Anhydrous, Granular AR
Toluene, Pesticide Grade
RH-7592 Analytical Standard
RH-9129 Lactone Analytical Standard
RH-9130 Lactone Analytical Standard
Water Milli-Q

Fisher
Baker

Johns Mansville
Johnson & Johnson
Fisher
U.S. Silica Co.

Baker
Baker
Fisher
Rohm and Haas Co.
Baker
Rohm and Haas Co.
Rohm and Haas Co.
Rohm and Haas Co.
Millipore

B. Preparation of Solutions

1. Prepare a 9.1% sodium chloride solution by dissolving 400 g of sodium chloride in four liters of Milli-Q water.
2. Prepare 100/30, 100/10 and 100/5 (v/v) toluene/acetone solutions by adding 300, 100 and 50 ml of acetone to 1000 ml of toluene.
3. Prepare 100/3 (v/v) toluene/methanol solution by adding 30 ml of methanol to 1000 ml of toluene.
4. Prepare GLC standard solutions by carefully weighing on an analytical balance 0.100g, correcting for purity, of each of the analytical standards (RH-7592, RH-9129 and RH-9130) into individual 100 ml volumetric flasks. Add approximately 80 ml of the toluene/methanol (100/3) solution and sonicate until dissolution occurs. Bring to volume with toluene/methanol (100/3). These are the primary standards of 1000 µg/ml.

Technical Report No. 34-90-47R

Into a 100 ml volumetric flask carefully pipet 10 ml of each of the primary standards and bring to volume with toluene/methanol (100/3). This is the primary multi-component standard of 100 ug/ml. Make serial dilutions of the multi-component standard to 10, 5, 1, 0.5, 0.1, and 0.05 ug/ml for working standards.

5. Prepare fortification standard solutions by carefully weighing on an analytical balance 0.10 g of each analytical standard (RH-7592, RH-9129 and RH-9130) into individual 100 ml volumetric flasks. Dissolve in a small amount of methanol. Bring to volume with methanol. These are the primary standards of 1000ug/ml. Into a 100 ml volumetric flask carefully pipet 10 ml of each of the standards and bring to volume with methanol. This is the primary multi-component fortification standard of 100 ug/ml. Make serial dilutions of the multi-component standard to 25, 12.5, 2.5, 1.25 and 0.25 ug/ml for working standards.

C. Equipment

Blenders Explosion-Resistant, Waring EPI
Buchner Funnels, Porcelain, 83 mm ID
Chromatographic Columns, 19 mm ID x 30 cm
length with 250 ml integral reservoir at
the top, teflon stopcock (Code 5907-T-15)
Rotary Evaporator, Buchi Rotovap R with
dry ice trap
Round bottom flasks, 500, 300,
Separatory funnels, 500 ml
Standard laboratory equipment, balances,
beakers, etc.
Ultrasonic Cleaner
Volumetric flasks, 100 ml

A. H. Thomas
A. H. Thomas

Ace Glass, Inc.

Brinkman
Kimax
Pyrex

Branson
Kimble

D. Instrumentation

Varian 3500 Capillary Gas Chromatograph equipped with a Varian Model 8035 Autosampler, a 1040 Megabore Injector and a Capillary Thermionic Specific Detector. Data are obtained with an HP 300 Data Acquisition and Processing Station with Hewlett-Packard Extrachrom Software. Data are processed by Nelson Analytical Software.

Column: Fused silica capillary, SPB-608, 0.53 mm ID, 15 meters, 0.5 um df - Supelco

Temperatures:	Column	-	245°C
	Injector	-	265°C
	Detector	-	300°C

Technical Report
No. 34-90-47R

Flows: Air (zero grade) - 175 ml/min
Hydrogen (UPC) - 4.5 ml/min
Helium (UPC) - 18 ml/min, 10 ml/min purge

Bead Current: 3.3 amps (varies)

Under these conditions, the retention times are as follows:

RH-7592 - 4.18 minutes
RH-9130 - 5.24 minutes
RH-9129 - 5.82 minutes

E. Analytical Procedure

1. Sample Processing

The fruit samples are received fresh at the Rohm and Haas Research Laboratories. The samples are split and the stones removed and discarded. The fruit is well chopped in the Hobart Food Chopper with dry ice and thoroughly mixed. The dry ice is allowed to sublime overnight in a freezer (4°F).

Store the samples in the freezer (-10°C) until analysis.

2. Extraction

Allow the samples to thaw sufficiently to remove an aliquot for analysis. Weigh 25 g of the chopped fruit into a 1 pint blender jar. Add 5-10 g of Celite-545, 100ml of methanol, cap and blend for 5 minutes. Vacuum suction filter the blended sample through the 7.5 cm filter paper using a porcelain filter funnel and a 500 ml filter flask. Wash the blender jar with 25 ml of methanol and pour over the filter cake into the 500 ml flask. Transfer the filtrate to a 500 ml separatory funnel. Wash the filter flask with 10 ml of methanol and add it to the separatory funnel. Proceed to the partition step, Section 3.

3. Partition

To the 500 ml separatory funnel containing the methanol extract, add 150 ml of methylene chloride and 250 ml of 9.1% sodium chloride solution. Let the solution stand for one or two minutes to allow the release of pressure. Cap and shake for approximately 10 seconds. Invert the separatory funnel and release the pressure by opening the stopcock. Close the stopcock and shake vigorously for one minute. When the phases separate, draw down the lower methylene chloride phase into a 500 ml round bottom flask. Evaporate the methylene chloride on the rotary evaporator at 45-50°C at atmospheric pressure. Remove the final traces of methylene chloride under vacuum. Add 25 ml of toluene/acetone (100/10) and swirl to dissolve the residue. Proceed to the silica gel cleanup step.

Technical Report
No. 34-90-47R

4. Column Chromatography Cleanup

a. Standardization of Silica Gel elution

Place 10 µg/ml RH-7592, RH-9129, RH-9130 standard in a r.b. flask. Rotavap to dryness. Add 25 ml toluene/acetone(100/10) to the flask to dissolve the residue. Insert a small cotton plug into a 19mm ID column and dry pack the column with 15cc of activated Silica Gel. Top the column 1 inch of anhydrous sodium sulfate. Add the 25 ml of toluene/acetone (100/10) solution in the r.b. flask to the column and collect the eluent as Cut 1. Rinse the flask with 10 ml toluene/acetone (100/10), and add to the column when the previous addition just enters the top of the column and collect the eluent as Cut 2. Add 10 ml toluene/acetone (100/10) to the column in like fashion and collect the eluent as Cut 3. Add 8 x 25 ml toluene/acetone (100/30) in like fashion and collect eluents as Cuts 4, 5, 6, 7, 8, 9, 10 and 11. Each cut (1 to 11) is rotavap to dryness and dissolved in 2 ml toluene/methanol (100/3). Each cut is injected into the GC along with a set of standards. Further dilution might be necessary to give a response between the linear range of the standard curve. Cuts 4 to 9 should show the elution of RH-7592 and its metabolites. Quantitatively, recovery of the analytes should be about 85% or greater.

b. Silica Gel Column Chromatography

Insert a small cotton plug into a 19mm ID chromatographic column and dry pack the column with 15cc of the activated silica gel (measure in a 25cc graduate cylinder). Top the column with 1 inch of anhydrous sodium sulfate. Add the 25 ml of toluene/acetone (100/10) solution in the 500 ml round bottom flask from the partition step to the column and collect the eluent in a 150 ml Erlenmeyer flask. Rinse the 500 ml round bottom flask with 10 ml of toluene/acetone (100/10) and add to the column when the previous addition just enters the top of the column. Add 10 ml of toluene/acetone (100/10) to the column when the previous addition just enters the top of the column and collect the eluent in the 150 ml Erlenmeyer flask. Discard the combined eluent. Add 150 ml of toluene/acetone (100/30) to the column when the previous addition just enters the top of the column and collect the eluent in a 500 ml round bottom flask. Concentrate the eluent to dryness the rotary evaporator at 70°C under reduced pressure. Add 25 ml of toluene/acetone (100/5) to the flask and swirl to dissolve the residue and proceed to the Florisil column cleanup step.

Technical Report
No. 34-90-47R

c. Standardization of Florisil elution

Place 10 $\mu\text{g/ml}$ RH-7592, RH-9129, RH-9130 standard in a r.b. flask. Rotavap to dryness. Add 25 ml toluene/acetone (100/5) to the flask to dissolve the residue. Insert a small cotton plug into a 19mm ID column and dry pack the column with 15cc of activated Silica Gel. Top the column 1 inch of anhydrous sodium sulfate. Add the 25 ml of toluene/acetone (100/5) solution in the r.b. flask to the column and collect the eluent as Cut 1. Rinse the flask with 10 ml toluene/acetone (100/5), and add to the column when the previous addition just enters the top of the column and collect the eluent as Cut 2. Add 25 ml toluene/acetone (100/5) to the column in like fashion and collect the eluent as Cut 3. Add 8 x 25 ml toluene/acetone (100/30) in like fashion and collect eluents as Cuts 4, 5, 6, 7, 8, 9, 10 and 11. Each cut (1 to 11) is rotavap to dryness and dissolved in 2ml toluene/methanol (100/3). Each cut is injected into the GC along with a set of standards. Further dilution might be necessary to give a response between the linear range of the standard curve. Cuts 4 to 9 should show the elution of RH-7592 and its metabolites. Quantitatively, recovery of the analytes should be about 85% or greater.

d. Florisil Column Chromatography

Insert a small cotton plug into a 19mm ID chromatographic column and dry pack the column with 15cc of the activated Florisil (measure in a 25cc graduated cylinder). Top the column with 1 inch of anhydrous sodium sulfate. Add the 25 ml of the toluene/acetone (100/5) solution in the 500 ml round bottom flask to the column and collect the eluent in a 125 ml Erlenmeyer flask. Rinse the 500 ml round bottom flask with 10 ml of toluene/acetone (100/5) and add to the column when the previous addition just enters the top of the column. Add 25 ml of toluene/acetone (100/5) to the column in like fashion. Discard the eluent. Elute the column with 150 ml of toluene/acetone (100/30) and collect in a 300 ml round bottom flask. Concentrate the eluent to dryness on the rotary evaporator at 70°C under reduced pressure. Be sure the samples are free of residual toluene. Add the appropriate volume of toluene/methanol (100/3), Section 5.b., and proceed to the gas chromatography step.

Gas Chromatography

a. Preparation of Standard Curves

A 3 μl aliquot of each multi-component standard in the range of 1.0 to 0.05 $\mu\text{g/ml}$ is injected into the gas chromatograph. The resulting peak heights are measured and plotted vs concentration ($\mu\text{g/ml}$) to obtain three(3) standard calibration curves. Standard curves are prepared for each analysis day and are obtained by least squares fit of standard injection data.

Technical Report
No. 34-90-47R

b. Sample Analysis

A 3 μ l aliquot of the final GLC samples is injected into the gas chromatograph. All samples are run initially diluted to 5 ml and, if necessary, the sample is further diluted with toluene/methanol (100/3) to give a response within the standard curve range. This allows all compounds to be detected at the 0.01 ppm level. The peak heights are measured and the concentration of each component is determined from the standard curves. The limit of detection is set by the data system at one-half the value of the peak height of the lowest standard injected (0.5 X lowest peak height of the 0.05 μ g/ml standard).

6. Method of Calculation

The RH-7592, RH-9129 and RH-9130 residue concentration are determined as follows:

a. Fortification Recovery

For samples fortified with known amounts of RH-7592m RH-9129 and RH-9130 prior to extraction, measure the peak heights, determine the concentration (μ g/ml) from the standard curves and calculate the percent recovery from the equation 1.

$$\text{Eq. 1 } \frac{(\mu\text{g/ml Found}) \times \text{Final Sample Vol. (ml)} \times 100}{\text{Fortification } (\mu\text{g})} = \% \text{ Recovery}$$

b. Component Residue Concentration

The component residue concentration is determined as follows:

$$\text{Eq. 2 } \frac{\text{Final Sample Vol. (ml)} \times \text{Component Conc. } (\mu\text{g/ml}) \times 100}{\text{Sample Weight (g)}} = \text{ppm}$$

IV. RESULTS AND DISCUSSION

Recovery data are summarized in Table I. The values range from 78 to 100% for RH-7592, 81 to 99% for RH-9130 and 60 to 94% for RH-9129. Detailed recovery data are shown in Table II. The demonstrated sensitivity of the method is 0.01 ppm for all stonefruit analyzed (cherry, peach and plum). Typical chromatograms of standards, calibration curves, cherry control, fortified control, treated samples and reagents are shown in Figure 1 to 13.

It is important that each lot of silica gel and Florisil be evaluated for consistency by running a multi-component standard through the elution scheme, concentrating the 150 ml toluene/acetone (100/30) out to dryness and determining the recovery by GLC. This should be greater than 85%.

**Technical Report
No. 34-90-47R**

The SPB-608 capillary column should be conditioned overnight at 260°C before connecting it to the detector. It was also necessary to prime the column before analysis by making two injections of the 5 µg/ml standard solution followed by two injections of toluene/methanol (100/3) solvent. If this is not done, standards injected at the beginning of an analytical run give a lower detector response than standards injected at the end of a run.

The TSD bead should be conditioned and the bead current set following the manufacturers instructions. The bead current should then be adjusted to a setting that will give a detector response sufficient to reliably measure a standard of 0.05 µg/ml concentration. The bead current will need to be adjusted as the bead ages.

The extraction efficiency is very good as metabolism studies with field aged radiolabeled RH-7592 treated peaches have shown that blending with methanol extracts more than 80% of the total radioactive residue from the fruit.

A Material Safety Data Sheet for RH-7592 is included in Appendix A.

V. Project Information and Study Certification Statement

Sponsor and Testing Facility: Rohm and Haas Company
Research Laboratories
727 Norristown Road
Spring House, PA 19477

Dates: Work Initiated: October 27, 1986
Work Completed: August 10, 1990

Study Personnel: Study Director: John Martin
Technical Assistant: Theodore F. Burnett

Notebook References: JH-35 Rohm and Haas No. 048623
JH-36 Rohm and Haas No. 048674
JH-37 Rohm and Haas No. 51666
JH-38 Rohm and Haas No. 52305
JH-39 Rohm and Haas No. 53047
TB-4 Rohm and Haas No. 50808
TB-5 Rohm and Haas No. 52778
TB-6 Rohm and Haas No. 54730

Technical Report
No. 34-90-47R

Study Certificate Statement

This study was conducted in conformance with standards of Good Laboratory Practice as defined by the United States Environmental Protection Agency and is a true and accurate representation of the residue analytical method development.


John J. Martin


Theodore F. Burnett

Technical Report
No. 34-90-47R

VI. METHOD REVISIONS

This method was revised to include the following minor modifications:

1. Filter paper # --- Chemical/Supplies (Sec. 111. A)
2. Procedure for the standardization of the Silica Gel and Florisil column cleanup elution pattern --- Cleanup (Sec. 111.E.4a)
3. Removal of control background correction for fortification recovery --- Method of calculation (Sec. 111.E.6a)
4. Removal of % recovery correction for residue data --- Method of calculation (Sec. 111.E.6b)

Other changes in this report include:

- a. Title page -- revised
- b. GLP/QA Statement page was added
- c. Table of Contents was revised to reflect all changes and additions.

All modifications were given the letter "R" after the page number, and report number (eg.: 8R, TR 34-90-47R). Any pages that were added to this report were given the letter "A" after their page number (eg.: 2A).

Technical Report
No. 34-90-47R

Table I

Summary of Recovery Data

	<u>Average(%)</u>	<u>Max</u>	<u>Min</u>
RH-7592	89+/-6	100	78
RH-9129	90+/-6	99	81
RH-9130	82+/-10	94	60

Technical Report
No. 34-90-47R

Table II

Recovery Data for RH-7592, RH-9130, and RH-9129 in Cherries

<u>Sample</u>	<u>ppm Added</u>	<u>% Recovery</u>		
		<u>RH-7592</u>	<u>RH-9130</u>	<u>RH-9129</u>
90-0066-003	4.0	93	98	85
90-0066-003	2.0	90	95	87
90-0066-003	0.08	95	90	88
90-0066-003	0.01	93	88	86
90-0066-003	0.01	100	99	94
90-0066-006	0.10	87	85	74
90-0066-009	0.04	91	90	77
90-0066-009	0.40	84	93	90
90-0058-003	3.0	83	83	71
90-0058-003	3.0	78	81	60


TR 34-90-47R

017

Data file: J801901 Type: STANDARD
 Sample Name:
 Date: 1 Aug 1990 09:27 Method: K7592AM Operator: JJM
 Interface: 706 Cycle#: 1 Channel#: A

Instrument: VARIAN 3500(5469)
 Column: SPB-608 Column Length: 15 Meters
 Start Temp-Time (deg-min): 245 Ramp Hold (deg-min):
 Program Rate (deg/min): End Time-Temp (deg-min):
 Prog Slope (# or Linear): Inj Port Temp: 265
 Flowrate/Gas: 17.8 Split Ratios
 Det 1-Type & Temp: TSD/300 Det 2-Type & Temp:

Plot times: 0 to 20 minutes
 Plot range: 20 millivolts (-3.2 mv offset)



RH7592
 RH9130
 RH9129

Retention Time	Compound Name	PPM Injected	Area	Height
4.18	RH7592	1.00	1.630E+01	9.990E+02
5.24	RH9130	1.00	1.100E+01	5.670E+02
5.82	RH9129	1.00	1.120E+01	5.260E+02

Figure 1

3 ul injection of 1.0 ug/mL RH-7592, RH-9130
 and RH-9129 Standard Solution

TR 34-90-47R

018

Data file: J801902 Type: STANDARD

Sample Name: Cal. Curve: 08/01/90

Date: 1 Aug 1990 09:49 Method: H7592Am Operator: JJM

Interface: 706 Cycled: 2 Channel: A

Instrument: VARIAN 3500(5469)

Column: SPB-608 Column Length: 15 Meters

Start Temp-Time (deg-min): 245 Ramp Hold (deg-min):

Program Rate (deg/min): End Time-Temp (deg-min):


Prog Slope (S or Linear): Inj Port Temp: 265

Flowrate/Gas: 17.8 Split Ratio:

Det 1-Type & Temp: TSD/300 Det 2-Type & Temp:

Plot times: 0 to 20 minutes

Plot range: 20 millivolts (-3.2 mv offset)



RH7592
RH9130

Retention Time	Compound Name	PPH Injected	Area	Height
4.19	RH7592	0.500	8.410E+08	5.330E+02
5.25	RH9130	0.500	5.440E+08	2.850E+02
5.83	RH9129	0.500	5.500E+08	2.350E+02

Figure 2

3 ul injection of 0.5 ug/mL RH-7592, RH-9130
and RH-9129 Standard Solution

R 34-90-47R

019

Data file: J801903 Type: STANDARD

Sample Name: Cal. Curve: 08/01/90

Date: 1 Aug 1990 10:12 Method: H7592A Operator: JJM

Interface: 706 Cycled: 3 Channel: A

Instrument: VARIAN 3500(5469) Column Length: 15 Meters

Column: SPB-608

Start Temp-Time (deg-min): 245 Ramp Hold (deg-min):

Program Rate (deg/min): End Time-Temp (deg-min):


Prog Slope (S or Linear): Inj Port Temp: 265

Flowrate/Gas: 17.8 Split Ratio:

Det 1-Type & Temp: TSD/300 Det 2-Type & Temp:

Plot times: 0 to 20 minutes

Plot range: 20 millivolts (-3.2 mv offset)



Retention Time	Compound Name	PPM Injected	Area	Height
4.19	RH7592	0.100	1.750E+00	1.110E+02
5.25	RH9130	0.100	1.070E+00	5.620E+01
5.84	RH9129	0.100	1.080E+00	5.040E+01

Figure 3

3 ul injection of 0.1 ug/mL RH-7592, RH-9130 and RH-9129 Standard Solution

TR 34-90-47R

020

Data file: J801904

Type: STANDARD

Sample Name:

Cal. Curve: 08/01/90

Date: 1 Aug 1990 10:35

Method: H7592Am

Operator: JHM

Interface: 706

Cycled: 4

Channel: A

Instrument: VARIAN 3500(5469)

Column: SPB-608

Column Length: 15 Meters

Start Temp-Time (deg-min): 265

Ramp Hold (deg-min):

Program Rate (deg/min):

End Time-Temp (deg-min):

Prog Slope (S or Linear):

Inj Port Temp: 265

Flowrate/Gas: 17.8


Split Ratio:

Det 1-Type & Temp: TSD/300

Det 2-Type & Temp:

Plot times: 0 to 20 minutes

Plot range: 20 millivolts (-3.2 mv p/foot)




RH7592

RH9130

Retention Time	Compound Name	PPM Injected	Area	Height
4.19	RH7592	0.050	7.990E-01	5.110E+01
5.25	RH9130	0.050	4.670E-01	2.520E+01
5.84	RH9129	0.050	4.960E-01	2.340E+01

Figure 4


3 ul injection of 0.05 ug/mL RH-7592, RH-9130 and RH-9129 Standard Solution



Concentrations in report are calculated from equation:
 $HEIGHT = B0 + B1(UG/ML) + B2(UG/ML)^2$
obtained by least-squares fit of standard injection data.

Figure 5

Typical RH-7592 Calibration Curve




Concentrations in report are calculated from equation:
 $HEIGHT = B0 + B1(UG/ML) + B2(UG/ML)^2$
obtained by least-squares fit of standard injection data.

Figure 6

Typical RH-9130 Calibration Curve

HEIGHT RH9129 CALIBRATION CURVE

08/01/90

 $B0 = .32787$ $B1 = 441.28$ $B2 = 84.365$ 

Concentrations in report are calculated from equation:
 $HEIGHT = B0 + B1(UG/ML) + B2(UG/ML)^2$
obtained by least-squares fit of standard injection data.

Figure 7

Typical RH-9129 Calibration Curve

TR 34-90-47R

024


Date file: J801906
Method file: N7592Am
Type: SAMPLE

RAI number: 90-0066
Sample No: 009
Component: FRUIT

Sample Name: _____ Cal. Curve: 05/01/90
Date: 1 Aug 1990 11:21 Method: N7592Am Operator: JRM
Interface: 706 Cycles: 6 Channel: A

Instrument: VARIAN 3500(5469) Column Length: 15 Meters
Column: SPB-608 Ramp Hold (deg-min):
Start Temp-Time (deg-min): 245 End Time-Temp (deg-min):
Program Rate (deg/min): Inj Port Temp: 255
Prog Slope (S or Linear): Split Ratios:
Flowrate/Gas: 17.8 Det 1-Type & Temp: TSD/300 Det 2-Type & Temp:

Plot times: 0 to 20 minutes
Plot range: 20 millivolts (-3.2 mv offset)



RH7592
RH9130

Ret. Time	Compound Name	Peak Area	Peak Height	ug/ml Found	Volume (ml)	Samp Wt.	Recov. Fact.	PPM
4.17	RH7592	7.330E-02	.477E+01	0.00669	5.00	25.0	1.00	0.00134
5.25	RH9130	0.000E+00	.000E+00	0.00	5.00	25.0	1.00	0.00
5.84	RH9129	0.000E+00	.000E+00	0.00	5.00	25.0	1.00	0.00

Figure 8

Control Cherry Sample No. 90-0066-009

TR 34-9C-47R

025

Data file: J801907

BAR number: 90-0066

Method file: M7592Am

Sample No: 009

Type: FORTIFICATION

Component: FRUIT

Sample Name:

Cal. Curve: 08/01/90

Date: 1 Aug 1990 11:44 Method: M7592Am

Operator: JRM

Interface: 706

Cycled: 7

Channel: A

Instrument: VARIAN 3300(3469)

Column Length: 15 Meters

Column: SPB-608

Start Temp-Time (deg-min): 245

Ramp Hold (deg-min):

Program Rate (deg/min):

End Time-Temp (deg-min):

Prog Slope (0 or Linear):

Inj Port Temp: 265

Flowrate/Gas: 17.8


Split Ratio:

Det 1-Type & Temp: TSD/300

Det 2-Type & Temp:

Plot time: 0 to 20 minutes

Plot range: 20 millivolts (-3.2 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	ug/ml Found	Volume (ml)	Ctrl Corr.	ug Added	Pct Recov
4.19	RH7592	7.340E+00	.455E+03	0.449	20.0	.033	10.7	83.6
5.24	RH9130	4.750E+00	.244E+03	0.437	20.0	0.008	9.80	93.3
5.83	RH9129	4.590E+00	.213E+03	0.444	20.0	0.000	9.90	89.7

Figure 9
Control Cherry Sample No. 90-0066-009
Fortified with 0.4 ppm of RH-7592, RH-9130 and RH-9129

TR 34-90-47R

026


Data file: J801909
Method file: H7592Am
Type: SAMPLE

RAR number: 90-0066
Sample No: 008
Component: FRUIT

Sample Name: _____ Cal. Curve: 08/01/90
Date: 1 Aug 1990 12:29 Method: H7592Am Operator: JAM
Interface: 706 Cycled: 9 Channel: A

Instrument: VARIAN 3500(3449)
Column: SPB-608 Column Length: 15 Meters
Start Temp-Time (deg-min): 245 Ramp Hold (deg-min):
Program Rate (deg/min): End Time-Temp (deg-min):
Prog Slope (S or Linear): Inj Port Temp: 265
Flowrate/Gas: 17.8 Split Ratios:
Det 1-Type & Temp: TSD/300 Det 2-Type & Temp:

Plot times: 0 to 20 minutes
Plot range: 20 millivolts (-3.2 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	ug/ml Found	Volume (ml)	Samp Wt.	Recov. Fact.	PPM
4.18	RH7592	2.340E+01	.133E+04	1.35	5.00	25.0	1.00	0.270
5.28	RH9130	8.990E+02	.497E+01	0.0111	5.00	25.0	1.00	0.00222
5.83	RH9129	2.790E+01	.130E+02	0.0304	5.00	25.0	1.00	0.00408

Figure 10
Treated Cherry Sample No. 90-0066-008 (5ml)

TR 34-90-47R

027


Data file: J801908
 Method file: H7592Am
 Type: SAMPLE

RAR number: 90-0066
 Sample No: 008
 Component: FRUIT

Sample Name: _____ Cal. Curve: 08/01/90
 Date: 1 Aug 1990 12:07 Method: H7592Am Operator: JLR
 Interface: 706 Cycle#: 8 Channel#: A

Instrument: VARIAN 3500(5469)
 Column: SPB-608 Column Length: 15 Meters
 Start Temp-Time (deg-min): 245 Ramp Hold (deg-min):
 Program Rate (deg/min): End Time-Temp (deg-min):
 Prog Slope (S or Linear): Inj Port Temp: 265
 Flowrate/Gas: 17.8 Split Ratio:
 Det 1-Type & Temp: TSD/300 Det 2-Type & Temp:

Plot times: 0 to 20 minutes
 Plot range: 20 millivolts (-3.2 mv offset)



RH7592
 RH9130

Ret. Time	Compound Name	Peak Area	Peak Height	ug/ml Found	Volume (ml)	Samp Wt.	Recov. Fact.	PPM
4.19	RH7592	6.920E+00	.427E+03	0.421	20.0	25.0	1.00	0.337
5.63	RH9129	8.580E-02	.496E+01	0.0105	20.0	25.0	1.00	0.0084
5.25	RH9130	0.000E+00	.000E+00	0.00	20.0	25.0	1.00	0.00

Figure 11

Treated Cherry Sample No. 90-C066-008 (20ml)

TR 34-90-47R

028


Data file: J8019011
Method file: H7592Am
Types: SAMPLE

RAR number: REAGENT
Sample No: 001
Component:

Sample Name:
Date: 1 Aug 1990 13:15 Method: H7592Am Operator: JMM
Interface: 706 Cycles: 11 Channel: A

Instrument: VARIAN 3500(3469)
Column: SPB-608 Column Length: 15 Meters
Start Temp-Time (deg-min): 245 Ramp Hold (deg-min):
Program Rate (deg/min): End Time-Temp (deg-min):
Prog Slope (S or Linear): Inj Port Temp: 265
Flowrate/Max: 17.8 Split Ratio:
Det 1-Type & Temp: TSD/308 Det 2-Type & Temp:

Plot times: 0 to 20 minutes
Plot ranges: 20 millivolts (-3.2 mv offset)



RH7592
RH8138

Ret. Time	Compound Name	Peak Area	Peak Height	ug/ml Found	Volume (ml)	Samp Wt.	Recov. Fact.	PPM
4.19	RH7592	0.000E+00	0.000E+00	0.00	5.00	25.0	1.00	0.00
5.25	RH8138	0.000E+00	0.000E+00	0.00	5.00	25.0	1.00	0.00
5.84	RH8129	0.000E+00	0.000E+00	0.00	5.00	25.0	1.00	0.00

Figure 12
Reagent Blank

TR 34-90-47R

029

Date File: J8019012

RAR number: REAGENT

Method File: H7592Am

Sample No: 001

Type: FORTIFICATION

Component:

Sample Name:

Cal. Curve: 08/01/90

Date: 1 Aug 1990 13:38 Method: H7592Am

Operator: JAM

Interface: 706

Cycled: 12

Channel: A

Instrument: VARIAN 3500(5469)

Column Length: 15 Meters

Column: SPB-608

Start Temp-Time (deg-min): 245

Ramp Hold (deg-min):

Program Rate (deg/min):

End Time-Temp (deg-min):

Prog Slope (# or Linear):

Inj Port Temp: 265

Flowrate/Min: 17.8


Split Ratios:

Det 1-Type & Temp: TSD/300

Det 2-Type & Temp:

Plot times: 0 to 20 minutes

Plot range: 20 millivolts (-3.2 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	ug/ml Found	Volume (ml)	Conc. Corr.	ug Added	Pct Recov
4.79	RH7592	7.580E+00	.470E+03	0.464	20.0	0.000	10.7	86.7
5.25	RH9130	4.460E+00	.228E+03	0.428	20.0	0.000	9.89	87.3
5.83	RH9129	4.680E+00	.218E+03	0.434	20.0	0.000	9.98	91.7

Figure 13

Reagent Fortification, 0.4 ppm

TR 34-90-47R

030

Appendix A

Material Safety Data Sheet for RH-7592

PROPERTY OF ROHM AND HAAS COMPANY PHILADELPHIA

TR 34-60-478

ROHM AND HAAS COMPANY

CORPORATE PRODUCT INTEGRITY DEPARTMENT
INDEPENDENCE WALL WEST
PHILADELPHIA, PA 19103EMERGENCY TELEPHONE
215-582-7000 ROHM AND HAAS
800-424-9300 CHEMTRECHAZARD RATING
1-EXTREME
2-HIGH
3-MODERATE
4-SLIGHT
5-INCONSIDERABLE
6-INCONSIDERABLE
7-SEE SECTION IV

CLOS

LIST 5

MATERIAL SAFETY DATA SHEET

NOT OSHA HAZARDOUS

MATERIAL RH-7592 TECHNICAL MATERIAL	CODE NONE	REV 892353-3	DOT HAZARD CLASS NONREGULATED
	DATE ISSUED 02/02/88		

FORMULA NA	CHEMICAL NAME OR SYNONYMS Proprietary triazole fungicide
---------------	---

I - COMPOSITIONAL INFORMATION

	CAS Reg. No.	APPROX WT %	TWA/TLV		
Triazole fungicide	TR SECRET #1	95-100	REL	OSHA	ACGIH
Related reaction products	NE	0-5	NE	NE	NE
			NE	NE	NE
			NE=none established		

II - PHYSICAL PROPERTY INFORMATION

APPEARANCE - COLOR - PH White solid; slight	VISCOSITY NA	
MELTING OR FREEZING POINT 116-118C / 241-244F	BOILING POINT No data	VAPOR PRESSURE mm Hg 12200 / 392F
SOLUBILITY IN WATER Insoluble	PERCENT VOLATILE BY WEIGHT 0	SPECIFIC GRAVITY WATER @ 1.20
		EVAPORATION RATE BUTYL ACETATE NA

III - FIRE AND EXPLOSION HAZARD INFORMATION

FLASH POINT NA	AUTO IGNITION TEMPERATURE NA	LOWER EXPLOSION LIMIT % NA	UPPER EXPLOSION LIMIT % NA
-------------------	---------------------------------	-------------------------------	-------------------------------

EXTINGUISHING MEDIA
☐ FOAM ☐ "ALCONOL" FOAM ☒ CO₂ ☒ DRY CHEMICAL ☐ WATER ☐ OTHER

SPECIAL FIRE FIGHTING PROCEDURES

Wear respirator (pressure-demand, self-contained breathing apparatus, MSHA/NIOSH-approved equivalent) and full protective gear. Run-off water should be contained.

UNUSUAL FIRE AND EXPLOSION HAZARDS

Toxic fumes are evolved when material is exposed to fire. Personnel should remain upwind and avoid exposure to smoke, because pesticides particulate may become airborne.

IV - HEALTH HAZARD INFORMATION

ROHM AND HAAS RECOMMENDED WORK PLACE EXPOSURE LIMITS

TWA—See Section I.

EFFECTS OF OVEREXPOSURE

Skin Contact: Irritating to skin upon repeated or prolonged exposure.

Eye Contact: Possibly irritating to eyes.

Ingestion: Possibly harmful if swallowed.

EMERGENCY AND FIRST AID PROCEDURES

Inhalation: Move subject to fresh air.

Eye and Skin Contact: Flush eyes with large amounts of water for at least 15 minutes. Consult a physician if irritation persists. Wash affected skin areas with soap and water.

Ingestion: If swallowed dilute by giving 2 glasses of water to drink and call a physician. Never give anything by mouth to an unconscious person.

TR 34-90-47R

V - REACTIVITY INFORMATION

032

STABILITY <input checked="" type="checkbox"/> STABLE <input type="checkbox"/> UNSTABLE		CONDITIONS TO AVOID NA
HAZARDOUS DECOMPOSITION PRODUCTS None known		
HAZARDOUS POLYMERIZATION <input type="checkbox"/> MAY OCCUR <input checked="" type="checkbox"/> WILL NOT OCCUR		CONDITIONS TO AVOID NA
INCOMPATIBILITY MATERIALS TO AVOID <input type="checkbox"/> WATER <input checked="" type="checkbox"/> OTHER Oxidizing materials, acids.		

VI - SPILL OR LEAK PROCEDURE INFORMATION

STEPS TO BE TAKEN IN CASE MATERIAL IS RELEASED OR SPILLED
Scoop or shovel into containers for disposal or recovery. Keep dusting to a minimum. Flush contaminated area with a large amount of water to a chemical sewer. Remove contaminated clothing and wash affected skin areas with soap and water. Wash clothing before reuse. Keep spills and cleaning run-offs out of the municipal sewers and open bodies of water.

WASTE DISPOSAL METHODS Pesticide, spray mixture, and rinsate that cannot be used according to label directions must be disposed of by incineration at a permitted facility according to state and local regulations.

VII - SPECIAL PROTECTION INFORMATION

VENTILATION TYPE Mechanical local exhaust at point of contaminant release.	
RESPIRATORY PROTECTION None required if good ventilation is maintained. Wear suitable MSHA/NIOSH approved respirator or equivalent where high mist concentrations are encountered.	
PROTECTIVE GLOVES Impervious	EYE PROTECTION Chemical splash goggles (ANSI Z87.1)
OTHER PROTECTIVE EQUIPMENT Eyewash facility, safety shower, protective clothing	

VIII - STORAGE AND HANDLING INFORMATION

STORAGE TEMPERATURE MAX 4 MIN	INDOOR YES	HEATED NO	REFRIGERATED NO	OUTDOOR NO
Indoor storage should be in a cool, dry, well ventilated area. Store away from food or feed.				

IX - TOXICITY INFORMATION

Acute oral LD50 (rat): >5 g/kg (male)
Acute dermal LD50 (rat) >5 g/kg
Skin rabbit: practically non-irritating, Eye rabbit: inconsequentially irritating.

X - MISCELLANEOUS INFORMATION

CONTAINER DISPOSAL Triple rinse (or equivalent) or puncture and dispose of in a sanitary landfill or by other state or federally approved procedures, such as, incineration.

NA - NOT APPLICABLE C - CHEMICAL VALUE	KEY 892353-3	DATE OF ISSUE 02/02/89	SUPERSEDES 04/11/85
THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED RELIABLE. HOWEVER NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA OR THE RESULTS TO BE OBTAINED THEREFROM.		ROHM AND HAAS COMPANY ASSUMES NO RESPONSIBILITY FOR PERSONAL INJURY OR PROPERTY DAMAGE TO WHATEVER DEGREE OR TO WHATEVER CAUSE BY THE USER'S, EMPLOYEE'S OR AGENT'S MISUSE OF THIS PRODUCT.	