

US EPA ARCHIVE DOCUMENT

*Valente*  
*BAS 199*

439235-01

BASF CORPORATION  
Agricultural Products  
P.O. Box 13528  
Research Triangle Park, NC 27709-3528

Study Title:

Validation of BASF Method 343/1 for the Determination  
of BAS 9075 I (Hexythiazox) in Dried Hops

Data Requirement:

EPA Guideline Numbers 171-4C

Author:

Dr. N. Tilting

Study Completed On:

December 1995

Performing Laboratory:

BASF Aktiengesellschaft  
Ecology and Environmental Analytics  
Environmental Chemistry  
Postfach 120  
D-67114 Limburgerhof

Registration Document No. BASF:

95/11149

This Report Consists of 49 Pages

95/11149 0001

**PR 86-5 DATA CONFIDENTIALITY CLAIM**

No claim of confidentiality is made for any information contained in this study on the basis of its falling within the scope of FIFRA 10 (d) (1) (A), (B) or (C).

Company: BASF Corporation, Agricultural Products

Company Agent: Rodney C. Akers Date: Feb 6, 1996

Title: Registration Scientist Signature: Rodney C. Akers


35/11149 0002

## **PR 86-5 GOOD LABORATORY PRACTICE CERTIFICATION**

This study does not meet the requirements for 40 CFR 160, Good laboratory Practice Standards and differs in the following way:

1. This study was performed in accordance with the GLP-Regulation of Appendix 1 to § 19a, Section 1, Chemikaliengesetz of 25-July-1994 (Official Bulletin/Federal Republic of Germany, I 1994, p. 1703).

Study Director: See page 6



Sponsor/ Submitter


Rodney Akers

Registration Scientist

BASF CORPORATION

P.O. Box 13528

Research Triangle Park, NC 27709-3528



Date

95/11149 0003

**Reg.Doc.# BASF 95/11149**

**Study Title**

Validation of BASF method 343/1 for the determination of BAS 9075 I (hexythiazox) in dried hops

**Data requirement**

EPA Guidelines Subdivision O, 171-4c

**Author**

Dr. N. Tilting

**Study completed in**

December 1995

**Testing facility**

BASF Aktiengesellschaft  
Ecology and Environmental Analytics  
Environmental Chemistry,  
Postfach 120  
D-67114 Limburgerhof

**Study Code**

30796

**Number of Pages**

46

95/11149 0004

# RheinlandPfalz



Ministerium für Arbeit, Soziales und Gesundheit  
Postfach 50 11 19421 Mainz

Ministerium für Arbeit,  
Soziales und Gesundheit

Rauhhoferstraße 9, 55116 Mainz  
Telefon (06131) 16-0 2405  
Telefax (06131) 162452  
Bemerkens/Beobachter: Herr Kail  
Aktenzeichen: 6324-80 283-0  
Mainz 20. Dezember 1994

## GLP - BESCHEINIGUNG

### Bescheinigung

Hiermit wird bestätigt, daß die Prüfeinrichtung

APS/U  
der Firma BASF AG  
- Landwirtschaftliche  
Versuchsstation -

57114 Limburgerhof  
Bundesrepublik Deutschland

am 01. und 14. September 1993  
und vom 15.-16. Dezember 1993  
von mir für die Überwachung zuständigen Behörde über die Einhaltung der Grundsätze der Guten Laborpraxis inspiziert worden ist.

Es wird hiermit bestätigt, daß die Prüfungen in dieser Prüfeinrichtung nach den Grundsätzen der Guten Laborpraxis durchgeführt werden.

In Antrag

*W. Fresenius*  
Dr. Fresenius



### Certificate

It is hereby certified that  
the test facility

APS/U  
der Firma BASF AG  
- Landwirtschaftliche  
Versuchsstation -

57114 Limburgerhof

Federal Republic of Germany.

at 01. and 14. September 1993  
and from 15.-16. December 1993  
was inspected by the  
competent authority regarding  
compliance with the principles  
of Good Laboratory  
Practice.

It is hereby certified that  
studies in this test facility  
is conducted in compliance  
with the principles of Good  
Laboratory Practice.

### Anlage zur GLP-Bescheinigung für die Prüfeinrichtung APS/U

Die Prüfeinrichtung führt in Übereinstimmung mit dem GLP-Grundzügen Prüfungen in folgenden Prüfkategorien durch:

1. Prüfungen auf physikalisch-chemische Eigenschaften und Gehaltsbestimmungen
4. Umwelttoxikologische Prüfungen zu Auswirkungen auf aquatische und terrestrische Organismen
5. Prüfungen zum Verhalten im Boden, Wasser und in der Luft; Bioakkumulation; Metabolismus
6. Prüfungen auf Rückstände
7. Prüfungen zu Auswirkungen auf Mikroben und natürliche Ökosysteme
8. Analytische Untersuchungen an biologischen Materialien

95/11149 0005

**GOOD LABORATORY PRACTICE**

**STATEMENT OF COMPLIANCE**

The study was conducted in compliance with Good Laboratory Practice Regulations:  
Appendix 1 to § 19 a, Section 1, Chemikaliengesetz of 25<sup>th</sup> July 1994 (Official  
Bulletin/Federal Republic of Germany, I 1994, p. 1703).

Study Director:

Date: 23.12.95

Signature: T. Lütting

95/11149-0006

B A S F Aktiengesellschaft

BASF Landw. Versuchsstation - Postfach 120 - 67114 Limburgerhof

**STATEMENT OF THE QUALITY ASSURANCE UNIT**

Study Code: 30796

The quality assurance unit of the testing facility inspected the study plan, study and audited the final report. Findings were reported to the study director and to management.

| Date of inspections | Report to study director<br>and to management |
|---------------------|---|
| 11-Jan-1995         |   |
| 02-Feb-1995         |   |
| 18-Dec-1995         |   |

67114 Limburgerhof

Scheper 19.12.95  
Signature QAU

95/11149 0007


### Quality Assurance Statement

I hereby attest to the authenticity of the study and guarantee that the data are correct and accurate to the best of our knowledge and that the study was performed by the procedures described.

Study Director and Author:

Dr. Norman Tilting  
APS/UT  
BASF Aktiengesellschaft  
Landwirtschaftliche Versuchsstation  
Postfach 120  
D-67114 Limburgerhof  
Federal Republic of Germany  
Tel.: 0621/60-27398

Date: 20. XII. 95

Signature: 

**SUMMARY**

The residue analytical method 343 for the gas chromatographic determination of hexythiazox was modified to improve performance of dry hops analysis. A GPC clean-up step was added to remove extracted plant material. It is possible to quantify residues of hexythiazox and its major metabolites. The analytes were extracted with methanol and cleaned-up by partitioning and column chromatography on florisil/charcoal. On injection into the hot injection port of the gas chromatograph the metabolite PT 1-3 was formed ( $T = 300^\circ\text{C}$ ) in almost quantitative yield which was quantified against a PT 1-3 calibration standard.

The limit of quantitation of the method was 0.5 mg/kg for dried hops. The mean recovery at the limit of quantitation was 102.8 % with a coefficient of variation of 19%.

Samples fortified with hexythiazox metabolites PT 1-3, PT 1-8 and PT 1-9 at the 2 mg/kg level were also analysed according to the modified method. Mean recoveries ranged from 80.7 % to 103.1% with coefficients of variation from 2.2 to 5.1 ( $n = 4$ ).

**TABLE OF CONTENTS**

|  | Page   |
|--|--------|
| COVER PAGE                                   | 1      |
| GLP CERTIFICATE                              | 2      |
| STATEMENT OF COMPLIANCE                      | 3      |
| STATEMENT OF THE QUALITY ASSURANCE UNIT      | 4      |
| QUALITY ASSURANCE STATEMENT                  | 5      |
| SUMMARY                                      | 6      |
| TABLE OF CONTENTS                            | 7      |
| LIST OF TABLES                               | 9      |
| LIST OF FIGURES                              | 9      |
| GENERAL                                      | 10     |
| <br>1. TEST SUBSTANCES                       | <br>11 |
| 1.1 Test Substances Used for Fortification   | 11     |
| 1.2 Test and Reference Substance (PT 1-3)    | 12     |
| 1.3 Principle of the Method                  | 13     |
| <br>2. MATERIALS AND METHODS                 | <br>13 |
| 2.1 Equipment                                | 13     |
| 2.2 Reagents and Standards                   | 14     |
| <br>3. ANALYTICAL PROCEDURE                  | <br>15 |
| 3.1 Preparation of Samples                   | 15     |
| 3.2 Clean-up Procedure                       | 15     |
| 3.2.1 Extraction                             | 15     |
| 3.2.2 Partition                              | 16     |
| 3.2.3 Preparation of GPC-column              | 16     |
| 3.2.4 GPC column clean-up                    | 16     |
| 3.2.5 Florisil Column clean-up               | 16     |
| 3.3 Instrumentation and Operating Conditions | 17     |
| 3.4 Calibration Procedures                   | 17     |
| 3.5 Interferences                            | 17     |
| 3.6 Confirmatory Technique                   | 17     |
| 3.7 Potential Problems                       | 18     |
| <br>4. METHOD OF CALCULATION                 | <br>18 |
| 4.1 Calibration                              | 18     |
| 4.2 Analyte in Samples                       | 18     |
| 4.2.1 Principle                              | 18     |
| 4.2.2 Calculation of Residues                | 18     |
| 4.2.3 Calculation of Recoveries              | 19     |
| <br>5. RECOVERIES, ACCURACY, AND PRECISION   | <br>19 |
| 5.1 Limit of Quantification                  | 19     |

|                 |                               |           |
|-----------------|-------------------------------|-----------|
| <b>6.</b>       | <b>ADDITIONAL INFORMATION</b> | <b>20</b> |
| 6.1             | Ruggedness Testing            | 20        |
| 6.2             | Limitations                   | 20        |
| <b>7.</b>       | <b>REFERENCES</b>             | <b>20</b> |
| <b>Appendix</b> |                               | <b>21</b> |

**LIST OF TABLES**

|          |   |    |
|----------|---|----|
| Table 1: | Recoveries, accuracy and precision of the method          | 19 |
| Table 2: | Analytical raw data of analytical series ME95006          | 22 |
| Table 3: | Analytical raw data of analytical series ME95009          | 23 |
| Table 4: | Analytical raw data of analytical series ME95015          | 24 |
| Table 5: | Analytical raw data of analytical series ME95026 (PT 1-8) | 25 |
| Table 6: | Analytical raw data of analytical series ME95026 (PT 1-9) | 26 |

**LIST OF FIGURES**

|            |   |    |
|------------|---|----|
| Figure 1:  | Flow Chart of Analytical Procedure                          | 21 |
| Figure 2:  | PT 1-3 Standard 0.25 µg/ml                                  | 27 |
| Figure 3:  | PT 1-3 Standard 0.5 µg/ml                                   | 27 |
| Figure 4:  | PT 1-3 Standard 1.0 µg/ml                                   | 28 |
| Figure 5:  | Example of reagent blank (no.: 9500062)                     | 28 |
| Figure 6:  | PT 1-3 Example of calibration curve                         | 29 |
| Figure 7:  | Dried hops control sample no.: 9500054                      | 30 |
| Figure 8:  | Dried hops control sample no.: 9500055                      | 30 |
| Figure 9:  | Dried hops control sample no.: 9500056                      | 31 |
| Figure 10: | Dried hops control sample no.: 9500057                      | 31 |
| Figure 11: | Dried hops spiked with 0.5 mg/kg hexythiazox (no.: 9500058) | 32 |
| Figure 12: | Dried hops spiked with 0.5 mg/kg hexythiazox (no.: 9500059) | 32 |
| Figure 13: | Dried hops spiked with 0.5 mg/kg hexythiazox (no.: 9500060) | 33 |
| Figure 14: | Dried hops spiked with 0.5 mg/kg hexythiazox (no.: 9500061) | 33 |
| Figure 15: | Dried hops spiked with 2.0 mg/kg hexythiazox (no.: 9500063) | 34 |
| Figure 16: | Dried hops spiked with 2.0 mg/kg hexythiazox (no.: 9500064) | 34 |
| Figure 17: | Dried hops spiked with 2.0 mg/kg hexythiazox (no.: 9500065) | 35 |
| Figure 18: | Dried hops spiked with 2.0 mg/kg hexythiazox (no.: 9500066) | 35 |
| Figure 19: | Dried hops spiked with 2.0 mg/kg PT 1-3 (no.: 9500356)      | 36 |
| Figure 20: | Dried hops spiked with 2.0 mg/kg PT 1-3 (no.: 9500357)      | 36 |
| Figure 21: | Dried hops spiked with 2.0 mg/kg PT 1-3 (no.: 9500358)      | 37 |
| Figure 22: | Dried hops spiked with 2.0 mg/kg PT 1-3 (no.: 9500359)      | 37 |
| Figure 23: | Dried hops spiked with 2.0 mg/kg PT 1-8 (no.: 9500360)      | 38 |
| Figure 24: | Dried hops spiked with 2.0 mg/kg PT 1-8 (no.: 9500361)      | 38 |
| Figure 25: | Dried hops spiked with 2.0 mg/kg PT 1-8 (no.: 9500362)      | 39 |
| Figure 26: | Dried hops spiked with 2.0 mg/kg PT 1-8 (no.: 9500363)      | 39 |
| Figure 27: | Dried hops spiked with 2.0 mg/kg PT 1-9 (no.: 9500364)      | 40 |
| Figure 28: | Dried hops spiked with 2.0 mg/kg PT 1-9 (no.: 9500365)      | 40 |
| Figure 29: | Dried hops spiked with 2.0 mg/kg PT 1-9 (no.: 9500366)      | 41 |
| Figure 30: | Dried hops spiked with 2.0 mg/kg PT 1-9 (no.: 9500367)      | 41 |
| Figure 31: | Example of GC-MS run of PT 1-3 with NH <sub>3</sub> Cl      | 42 |
| Figure 32: | Analytical certificates of hexythiazox                      | 43 |
| Figure 33: | Analytical certificates of PT 1-3                           | 44 |
| Figure 34: | Analytical certificates of PT 1-8                           | 45 |
| Figure 35: | Analytical certificates of PT 1-9                           | 46 |

**GENERAL****Performing Laboratory**

BASF Aktiengesellschaft  
Ecology and Environmental Analytics  
Plant Metabolism and Residue Studies  
Postfach 120  
D-67114 Limburgerhof

**Guidelines covered**

EPA Guidelines Subdivision O, 171-4c.

GLP Consensus Document: The Application of the GLP Principles to Field Studies, Env. Monograph No. 50, OECD, Paris 1992

For Good Laboratory Practice: Appendix 1 to § 19a, Section 1, Chemikaliengesetz of 25 July 1994 (Official Bulletin, Federal Republic of Germany, I 1994, p 1703)

**Location of raw data**

Archive belonging to APS/U, building Li 444

**Filing of reports, raw data, samples and specimen**

Study plan, raw data, original report or certified duplicates, samples and specimen, will be stored in an archive at BASF Aktiengesellschaft, Landwirtschaftliche Versuchsstation, D-67117 Limburgerhof, at least for the period of time specified in the GLP regulations.

**Schedule**

|                              |               |
|------------------------------|---------------|
| Protocol signed:             | January 1995  |
| Start of experimental part:  | January 1995  |
| Completion of experim. part: | June 1995     |
| Final report:                | December 1995 |

**Technical staff**

Manfred Metzger

**Test system**


Dried hops

## 1. TEST SUBSTANCES

### 1.1 Test Substances Used for Fortification

#### *Hexythiazox*


Common name: Hexythiazox  
Code name: Reg. No. 282502  
Chemical name: (4RS,5RS)-5-(4-chlorophenyl)-N-cyclohexyl-4-methyl-2-oxo-1,3-thiazolidine-3-carboxamide  
Structural formula:



Empirical formula C<sub>17</sub>H<sub>21</sub>ClN<sub>2</sub>O<sub>2</sub>S  
Molecular weight: 352.89  
Solubility: Soluble in n-hexane, methanol, dichloromethane

#### *PT 1-8*

Code name: PT 1-8, Reg. No. 248650  
Chemical name: (4RS,5RS)-5-(4-chlorophenyl)-N-4-hydroxycyclohexyl-4-methyl-2-oxo-1,3-thiazolidine-3-carboxamide  
Structural formula:



Empirical formula: C<sub>17</sub>H<sub>21</sub>ClN<sub>2</sub>O<sub>3</sub>S  
Molecular weight: 368.99

**PT 1-9**


Code name :

PT 1-9, Reg. No. 248667

Chemical name:

5-(4-chlorophenyl)-N-4-oxo-cyclohexyl-4-methyl-2-oxo-1,3-thiazolidine-3-carboxamide

Structural formula:



Empirical formula:

 $C_{17}H_{19}ClN_2O_3S$ 

Molecular weight:

366.87

**1.2 Test and Reference Substance (PT 1-3)**


Code name:

PT 1-3, Reg. No. 248649

Chemical name:

trans-5-(4-chlorophenyl)-4-methyl-2-oxothiazolidone

Structural formula:



Molecular formula:

 $C_{10}H_{10}ClNOS$ 

Molecular weight:

227.7

Solubility:

Soluble in acetone, dichloromethane, methanol  
slightly soluble in n-hexane

### 1.3 Principle of the Method

Samples are extracted with methanol. After n-hexane partition, GPC and Florisil column clean-up, the residue of hexythiazox is determined by GC using NPD-detection. On injection, hexythiazox is pyrolyzed to PT 1-3, which is used as reference compound for the determination.

## 2. MATERIALS AND METHODS

Equipment and reagents in the following lists are examples and can be replaced by equivalent ones. (Suppliers, batch no., etc describe products that were in use at BASF's laboratories during method development.)

### 2.1 Equipment

|   |  |
|---|--|
| Glass bottle (wide neck)                          | 500 ml   |
| Round bottom flasks with standard ground joint 29 | 100, 250 and 500 ml  |
| Büchner funnels                                   | 90 mm i.d.   |
| Separatory funnels                                | 500 ml and 1000 ml   |
| Funnels   | glass, 90 mm i.d., 30 mm   |
| Chromatography columns                            | glass, 20 mm i.d.  |
| Tapered flasks with standard ground joint 14.5    | 10 ml  |
| Volumetric pipettes                               | 10 ml, 50 ml, 100 ml   |
| Volumetric flasks                                 |  |
| Glass syringe 10 ml luer-fitting                  | 90 mm i.d., Art. No. 300 009,<br>Schleicher & Schüll, Dassel, FRG<br>e.g. Rotavapor, Büchi, Flawil, CH |
| Round filters                                     |  |
| Rotary evaporator with water bath                 | e.g. MZZC, Vacuubrand GmbH & Co,<br>Wertheim, FRG  |
| Vacuum controller and pump                        |  |
| Instrumentation for GPC                           | Analytical Bio-Chemistry Laboratories  |
| GPC Autoprep 1002A                                | Columbia M.O. 65205, USA   |
| Column  | Analytical Bio-Chemistry   |
| Stationary phase                                  | Bio Bead, S-X3, Bio Rad Laboratories,<br>Richmond, CA 9404, USA 8                                      |
| Filter: Millex-HV 0.45 µm                         | Millipore, SLHV025NB, Millipore Corp., Bedford<br>M.A. 01730, USA                                      |

## 2.2 Reagents and Standards

It is recommended to use pesticide grade solvents.

|                               |  |
|-------------------------------|--|
| Acetone dist.                 |  |
| Cyclohexane                   |  |
| Dichloromethane dist.         |  |
| Ethylacetate                  |  |
| n-Hexane                      |  |
| Methanol dist.                |  |
| Celite                        | Art. No. 16391, Serva, Heidelberg, FRG       |
| Florisil, 0.15 - 0.25 mm      | Art. No. 81534, Macherey & Nagel, Düren, FRG |
| Sodium sulfate                |  |
| Cotton wool                   |  |
| Activated charcoal (powdered) | Merck, Darmstadt, FRG                        |

### Solutions and Solvent Mixtures

Dichloromethane + acetone 83 + 17 (v/v)

Cyclohexane + ethylacetate 1 + 1 (v/v)

### Standard substance for fortification

Hexythiazox (structural formula see 1.1)

|          |         |         |
|----------|---------|---------|
| Reg. No. | 282 502 | 282 502 |
| Purity   | > 99.9% | 98.9%   |
| Lot No.  | L43-251 | 2EI 72D |

PT 1-8 (structural formula see 1.1)

|          |         |
|----------|---------|
| Reg. No. | 248 650 |
| Purity   | 99.9%   |
| Lot No.  | L43-265 |

PT 1-9 (structural formula see 1.1)

|          |         |
|----------|---------|
| Reg. No. | 248 667 |
| Purity   | 99.9%   |
| Lot No.  | L43-267 |

Standards supplied by Dr. Roschger, BASF Aktiengesellschaft, APS/UD, 67114 Limburgerhof, Phone 0621/60-27103.

Preparation of standard solutions for fortification (example other dilution procedures can also be used):

Dissolve 10.0 mg hexythiazox in 100 ml methanol to yield a solution of 100 µg/ml. Dilute 1 ml of this solution with 99 ml of methanol to get a solution of 1 µg/ml. Use appropriate dilution for spiking, depending on sample weight and spiking level. Spike samples for fortification purposes prior to the extraction step (3.2.1).

Standard substance for Calibration

PT 1-3 (structural formula see 1.2)

|          |         |
|----------|---------|
| Reg. No. | 248 649 |
| Purity   | 99.9 %  |
| Lot No.  | L43-255 |

Standards supplied by Dr. Roschger, BASF Aktiengesellschaft, APS/UD, 67114 Limburgerhof, Phone 0621/60-27103.

Preparation of standard solutions for GC calibration (recommended procedure)

Dissolve 10 mg of PT 1-3 in 10 ml of acetone to yield a 1 mg/ml solution. Dilute 1 ml of this 1 mg/ml solution to 100 ml to yield a 10 µg/ml solution.

To prepare calibration standards in the range of 0.25 µg/ml to 2.0 µg/ml proceed as follows (other concentrations can also be used):

Dilute 2.0 ml of the 10 µg/ml solution to 10 ml = 2.0 µg/ml

Dilute 1 ml of the 10 µg/ml solution to 10 ml = 1.0 µg/ml

Dilute 0.5 ml of the 10 µg/ml solution to 10 ml = 0.5 µg/ml

Dilute 1 ml of the 2.5 µg/ml solution to 10 ml = 0.25 µg/ml

✓ All dilutions with hexane + acetone (9 + 1 v/v)

← should be 2.0 µg/ml

The solutions must be stored in a refrigerator when not in use. At 4°C in the dark Hexythiazox and PT 1-3 solutions are stable for 30 days (ref. 2 and ref. 3); PT 1-8 and PT 1-9 solutions were prepared fresh before spiking.

### 3. ANALYTICAL PROCEDURE

#### Note:

The symbols in brackets are identical with the ones used in the formula for the calculation of the analytical results (see 4.2.2). The analytical procedure is described for one sample.

#### 3.1 Preparation of Samples

Dried hops must be powdered before taking the laboratory sample.

The samples have to be stored at appr. -20°C until analysis.

#### 3.2 Clean-up Procedure

##### 3.2.1 Extraction

Sample weight: hops (dried, powdered): 5 g

Weigh sample into a wide mouth glass bottle, add 200 ml methanol and macerate with an Ultra Turrax for 3 minutes. Place a filter paper covered with a 0.5 cm layer

of Celite into a Büchner funnel and filter the suspension. Wash the filter cake with 50 ml methanol.

### 3.2.2 Partition

Transfer the extract from 3.2.1 into a 1000 ml separatory funnel. Add 200 ml of sodium chloride solution (200 g/l) and extract twice with 50 ml n-hexane each. *After* ~~Filtrate~~ the n-hexane phase through a glass funnel plugged with cotton wool and sodium sulfate. Repeat the extraction with dichloromethane (2 times 50 ml), filter over sodium sulfate, wash the sodium sulfate with 50 ml dichloromethane and combine the organic phases. Evaporate the organic phase to dryness, discard the water phase.

### 3.2.3 Preparation of GPC-column

A mixture of cyclohexane + ethylacetate is used as mobile phase in GPC. To prepare the column, weigh 40 g of SX3 Bio-Bead into a flask and add the solvent mixture that the beads are completely covered by solvent. Allow to swell over night. Pour the gel into the column and adjust the column adapters so that they touch the gel surface. Connect column to GPC system. Remove air from the system according to operating instructions from the manufacturer. Readjust column adaptors. As the retention volume of the analyte is dependent on column performance test it elution parameters (collect time) for the compound of interest (see 6). *↑  
is  
the*

### 3.2.4 GPC column clean-up

Dissolve extract from 3.2.2 in 5 ml ethylacetate. Sonicate if necessary. Add 5 ml cyclohexane. To remove undissolved particles, place a small piece of cotton wool into a funnel (i.d. 3 cm) and filter the solution into a 10 ml volumetric flask. Rinse cotton wool with some drops of cyclohexane-ethylacetate mixture and make up to the mark.

Fill filtrate into a 10 ml glass syringe, place a 0,45 µm filter disk onto the luer lock fitting and fill the injection loop of the GPC instrument. Make sure not to inject air into the injection loop. Start chromatographic run according to the instructions of the instrument manufacturer. The instrument injects only 5 ml (50%) of the sample. With a GPC instrument prepared according to 3.2.3, the first 125 ml (25-min) must be dumped, the active ingredient elutes within the next 60 ml (12 min). The column should be washed for 10 min (50 ml) before injection of the next sample. Evaporate eluent to dryness. *✓*

### 3.2.5 Florisil column clean-up

Place a cotton wool plug at the lower end of a chromatographic column and fill the column with dichloromethane + acetone 83 + 17 (v/v). Pour 5 g Florisil into the column. After settling of the Florisil lower the solvent meniscus to the Florisil surface.

Mix 15 g of Florisil with 1 g of charcoal and suspend in 35 ml of the dichloromethane acetone mixture. Transfer the slurry to the column and allow to settle. Lower the solvent meniscus to the Florisil surface. Cover the Florisil surface with a small amount of cotton wool. Quantitatively transfer the residue from 3.2.4 onto the column using 3 times 5 ml dichloromethane. Lower the solvent meniscus to the florisil surface and elute with 100 ml dichloromethane + acetone 83 + 17 (v/v). Collect the eluate and evaporate to dryness. Dissolve the residue in 1 ml or other appropriate volume (=  $V_E$ ) hexane + acetone (9 + 1 v/v).

### 3.3 Instrumentation and Operating Conditions

|                      |   |
|----------------------|---|
| Gas chromatograph    | Perkin Elmer 8420   |
| Detector             | (NPD)   |
| Column               | L = 30 m, i.d. = 0.32 mm, film = 0.25 µm  |
| Stationary phase     | DB 1701   |
| Injection port temp. | 300 °C  |
| Detector temperature | 300 °C  |
| Oven program         | Start: 160 °C, 0 minutes<br>rate: 6 °/ minute<br>final temperature 250 °C/ 10 minutes |
| Injection volume     | 1 µl  |
| Carrier gas          | Helium 1000 mbar  |
| Retention time       |   |
| PT-1-3               | 14 - 17 min   |

← NPD

### 3.4 Calibration Procedures

Calculation of results is based on peak height measurements using calibration curves. To obtain this standard curve inject 1 µl of standard solutions containing 0.25 to 2.0 µg/ml of PT 1-3 into the gas chromatograph.

For routine analysis at least every third injection should be a standard. Plot peak heights versus amount (ng) of injected standard. Hexythiazox and its metabolites are derivatized in the injection port to PT 1-3 so that calibration against PT 1-3 includes the derivatization efficiency.

### 3.5 Interferences

Interferences from the matrix, other pesticides, solvents or labware are not known to date.

### 3.6 Confirmatory Technique

To confirm the identity of the PT 1-3 peak in the chromatogram, GC-MS can be used. An example of a GC- MS run with NH<sub>3</sub> chemical ionization can be found in

figure 31. (non-GLP, not produced under this study code).

### 3.7 Potential Problems

There can be variations in the properties of the Florisil from different lots. Elution volume for column clean-up should be checked (see 6.).

It is recommended to condition the GC column by several injections of standard substance or control extract before starting the analytical series. The derivatization reaction in the injection port of the GC must also be checked, a dependence on injector geometry and filling with glass wool could be observed.

## 4. METHOD OF CALCULATION

### 4.1 Calibration

For establishing the calibration curves see 3.4.

### 4.2 Analyte in Samples

#### 4.2.1 Principle

Calculation of results is based on peak height measurements. Determine the amount of substance measured in the injected volume from the calibration curves and use the equation described in 4.2.2. Calculation can also be achieved by a suitable computer program.

#### 4.2.2 Calculation of Residues

The residues R in mg/kg are calculated as follows:

$$R = \frac{V_E \cdot W_A \cdot D}{G \cdot V_I \cdot 0.5}$$

- V<sub>E</sub> = Final volume (ml) after all dilution steps  
G = Weight (g) of sample extracted  
W<sub>A</sub> = Amount of analyte (ng) read from calibration curve  
V<sub>I</sub> = Volume (μl) injected into GC  
D = Derivatization factor  
0.5 = Aliquotation factor from GPC injection
- To convert PT 1-3 to hexythiazox the derivatization factor is 1.548.  
For the conversion of PT 1-3 to PT 1-8 the factor is 1.62 and for the conversion PT 1-3 to PT 1-9 it is 1.611.
- 353  
228
- 367  
228
- 361  
228

#### 4.2.3 Calculation of Recoveries

Determine the recovery factor F from fortification experiments

$$\% \text{ Recovery} = \frac{(\text{mg/kg in fortified sample} - \text{mg/kg in control}) \cdot 100}{\text{mg/kg analyte added}}$$

$$F = \frac{100\%}{\% \text{ Recovery}}$$

To correct results for recoveries, multiply the result from 4.2.2 with the recovery factor. The formula for corrected residues is:

$$\text{Residue corrected for 100\% recovery} = \frac{V_E \cdot W_A \cdot D}{G \cdot V_I \cdot 0.5} \cdot F$$

### 5. RECOVERIES, ACCURACY, AND PRECISION

#### 5.1 Limit of Quantitation

The limit of quantitation which is defined to be the lowest concentration for which accuracy and precision of the method was tested by recovery trials, is 0.5 mg/kg for dry hops. For hexythiazox two analytical series of two controls and two samples at 3 different spiking levels were run at different times. For the metabolites series of one control and four spiked samples were analysed. The results of the trials was listed in table 1. Details of the analytical series were listed in the raw data spreadsheets of tables 2 to 7.

Table 1: Recoveries, accuracy and precision of the method

| Compound added | Fort.level (mg/kg) | Recoveries (%)           | Mean (%) | Std. Dev. | Coeff. of Variation |
|----------------|--------------------|--------------------------|----------|-----------|---------------------|
| Hexythiazox    | 0.5                | 96.7; 83.1; 101.3; 130.0 | 102.8    | 19.7      | 19.2                |
|                | 2.0                | 99.8; 90.1; 97.0; 92.6   | 94.9     | 4.3       | 4.6                 |
|                | 5.0                | 99.4; 93.2; 77.6; 84.1   | 88.6     | 9.6       | 10.9                |
| PT 1-3         | 2.0                | 98.6; 108.5; 106.6; 98.5 | 103.1    | 5.3       | 5.1                 |
| PT 1-8*        | 2.0                | 81.9; 81.2; 81.6; 78.1   | 80.7     | 1.8       | 2.2                 |
| PT 1-9*        | 2.0                | 90.8; 85.7; 94.7; 86.5   | 89.4     | 4.2       | 4.7                 |

\* new GPC column

35/11148 0022

## 6. ADDITIONAL INFORMATION

During analytical work, the GPC column had to be refilled. The new column had a slightly different elution behaviour. The collect time had to be readjusted. The first 150 ml of the eluate (30 min run time) were dumped, the next 75 ml (15 min run time) were collected and the wash phase remained at 10 minutes.

Due to the varying activity of florisil the elution volume of the column described in 3.2.5 was changed to 150 ml (DCM+acetone 83 + 17).

### 6.1 Ruggedness Testing


Not performed.

### 6.2 Limitations

Not known to date.

## 7. REFERENCES

1. Determination of Hexythiazox and its Metabolites as PT 1-3, Institut Fresenius Method 0012, Januar 1988, BASF Fremdmethode 925 (Mai 1988)
2. Stability of Hexythiazox in Methanol, Acetonitrile and Acetonitrile/Water (1+1, v/v), BASF Report No. 3803, Dr. Ch. Mackenroth, April 1994, Reg. Doc. No. 94/10250
3. Stability of PT 1-3 in Acetone and n-Hexane, Dr. N. Tilting, December 1994, Reg. Doc. No. 94/11071
4. Determination of Hexythiazox and its Metabolites in Soil, Apples, Strawberries, Hops, Beer, Yeast and Brewing Wastes by GC with NPD (Method 343), BASF Reg. No. 3944, Dr. N. Tilting, April 1994, Reg. Doc. No. 94/10226

**Appendix.****Figure 1:** Flow Chart of Analytical Procedure

95/1114S 0024

Study code 30796

Study Code: 30796  
Method validation with hexythiazox  
Queuefile name: ME95006  
Matrix:Dried hops

Lab Journal no.00842

| Sample number | Type | Sample weight (g) | Final dilution(mL) | Sample injected (mg) | Peak height | Std. Conc. (pg/mL) | Substance injected (ng) | Residue (mg/kg) | Spiking level (ng) | Recovery % | Date of extraction | Date of measurement. |
|---------------|------|-------------------|--------------------|----------------------|-------------|--------------------|-------------------------|-----------------|--------------------|------------|--------------------|----------------------|
| 9500054       | C*   | 5                 | 1                  | 2.5                  | 159         | 0.25               | 0.07996                 | < 0.5           |                    | 31.1.95    | 31.1.95            |                      |
| ME038         | S    |                   |                    |                      | 581         | 0.50               | 0.23956                 |                 |                    |            | 31.1.95            |                      |
| ME041         | S    | 5                 | 1                  | 2.5                  | 1208        | 0.50               | 0.47668                 |                 |                    |            | 31.1.95            |                      |
| 9500055       | C    | 5                 | 1                  | 2.5                  | 148         | 1.00               | 0.07580                 | < 0.5           |                    |            | 31.1.95            |                      |
| ME04          | S    |                   |                    |                      | 2640        | 1.00               | 1.01826                 |                 |                    |            | 31.1.95            |                      |
| 9500058       | F    | 5                 | 1                  | 2.5                  | 2321        | 0.25               | 0.89761                 |                 | 2.5                | 101.3      | 31.1.95            |                      |
| ME039         | S    |                   |                    |                      | 712         | 0.25               | 0.28910                 |                 |                    |            | 31.1.95            |                      |
| 9500059       | F    | 5                 | 1                  | 2.5                  | 2934        | 0.50               | 1.12944                 |                 | 2.5                | 130        | 31.1.95            |                      |
| ME042         | S    |                   |                    |                      | 1331        | 0.50               | 0.52301                 |                 |                    |            | 31.1.95            |                      |
| 9500063       | F    | 5                 | 4                  | 0.625                | 2072        | 1.00               | 0.80344                 |                 | 10                 | 97         | 31.1.95            |                      |
| ME045         | S    |                   |                    |                      | 2510        | 1.00               | 0.96909                 |                 |                    |            | 31.1.95            |                      |
| 9500084       | F    | 5                 | 4                  | 0.625                | 1977        | 0.25               | 0.76751                 |                 | 10                 | 92.6       | 31.1.95            |                      |
| ME040         | S    |                   |                    |                      | 612         | 0.50               | 0.25128                 |                 |                    |            | 31.1.95            |                      |
| 9500057       | F    | 5                 | 10                 | 0.25                 | 1626        | 0.25               | 0.63477                 |                 | 25                 | 77.6       | 31.1.95            |                      |
| ME043         | S    |                   |                    |                      | 1180        | 0.50               | 0.46869                 |                 |                    |            | 31.1.95            |                      |
| 9500088       | F    | 5                 | 10                 | 0.25                 | 1765        | 1.00               | 0.68734                 |                 | 25                 | 84.1       | 31.1.95            |                      |
| ME046         | S    |                   |                    |                      | 2536        | 1.00               | 1.01674                 |                 |                    |            | 31.1.95            |                      |
| 9500062       | R    |                   |                    | 1                    | 313         | 0.25               | 0.1382003               |                 |                    |            | 31.1.95            |                      |

C = Control R = Reagent blank F = Fortified sample S = Standard

Injection volume: 1.0  $\mu$ L Aliquotation 50%

Calibration: linear Intercept: 0.1983e-1 Coeff. of correlation: 0.99938

SD+-0.2788e-1 \* This blank was subtracted from recoveries

Deriv. factor: 1.648  
Slope: 0.3782e-3

Table 3: Analytical raw data of analytical series ME95009

| Sample number | Type | Sample weight (g) | Final dilution(mL) | Sample injected (mg) | Peak height | Std Conc. (µg/mL) | Substance injected (ng) | Residue (mg/g) | Spiking level (µg) | Recovery % | Date of extraction | Date of measurement |
|---------------|------|-------------------|--------------------|----------------------|-------------|-------------------|-------------------------|----------------|--------------------|------------|--------------------|---------------------|
| 9500056       | C    | 5                 | 2                  | 1.25                 | 0           | 0.00000           | <0.5                    |                |                    | 132.95     | 16.2.95            |                     |
| ME050         | S    |                   |                    |                      |             | 0.50343           |                         |                |                    |            | 16.2.95            |                     |
| 9500057       | C    | 5                 | 2                  | 1.25                 | 0           | 0.00000           | <0.5                    |                |                    | 132.95     | 16.2.95            |                     |
| ME053         | S    |                   |                    |                      |             | 1.00634           |                         |                |                    |            | 16.2.95            |                     |
| 9500060       | F    | 5                 | 2                  | 1.25                 | 500         | 0.39024           |                         | 2.5            | 96.7               | 132.95     | 16.2.95            |                     |
| ME048         | S    |                   |                    |                      |             | 0.26651           |                         |                |                    |            | 16.2.95            |                     |
| 9500061       | F    | 5                 | 1                  | 2.5                  | 927         | 0.67125           |                         | 2.5            | 83.1               | 132.95     | 16.2.95            |                     |
| ME057         | S    |                   |                    |                      |             | 1.03536           |                         |                |                    |            | 16.2.95            |                     |
| ME061         | S    |                   |                    |                      |             | 0.26847           |                         |                |                    |            | 16.2.95            |                     |
| ME051         | S    |                   |                    |                      |             | 0.46347           |                         |                |                    |            | 16.2.95            |                     |
| 9500065       | F    | 5                 | 4                  | 0.625                | 899         | 0.86625           |                         | 10             | 99.8               | 132.95     | 16.2.95            |                     |
| ME054         | S    |                   |                    |                      |             | 1.03448           |                         |                |                    |            | 16.2.95            |                     |
| ME052         | S    |                   |                    |                      |             | 0.27634           |                         |                |                    |            | 16.2.95            |                     |
| 9500066       | F    | 5                 | 4                  | 0.625                | 899         | 0.72755           |                         | 10             | 90.1               | 132.95     | 16.2.95            |                     |
| ME065         | S    |                   |                    |                      |             | 0.48796           |                         |                |                    |            | 16.2.95            |                     |
| ME058         | S    |                   |                    |                      |             | 0.93392           |                         |                |                    |            | 16.2.95            |                     |
| ME049         | S    |                   |                    |                      |             | 0.24299           |                         |                |                    |            | 16.2.95            |                     |
| ME053         | S    |                   |                    |                      |             | 0.27483           |                         |                |                    |            | 16.2.95            |                     |
| 9500069       | F    | 5                 | 10                 | 0.25                 | 249         | 0.25              |                         | 25             | 99.4               | 132.95     | 16.2.95            |                     |
| ME052         | S    |                   |                    |                      |             | 0.80246           |                         |                |                    |            | 16.2.95            |                     |
| ME056         | S    |                   |                    |                      |             | 0.49379           |                         |                |                    |            | 16.2.95            |                     |
| 9500070       | F    | 5                 | 10                 | 0.25                 | 476         | 0.50              |                         | 0.48704        |                    |            |                    |                     |
| ME055         | S    |                   |                    |                      |             | 0.75230           |                         | 25             | 93.2               | 132.95     | 17.2.95            |                     |
| ME059         | S    |                   |                    |                      |             | 1.03878           |                         |                |                    |            | 17.2.95            |                     |
|               |      |                   |                    |                      |             | 0.96258           |                         |                |                    |            |                    |                     |

C = Control F = Fortified sample S = Standard

Injection volume: 1µL

Calibration: linear

Coeff. of correlation: 0.9966/0.9973/0.9935

Aliquotation: 50%

Deriv. factor: 1.548

Intercept: 0.6118e-10/0.2013e-10/0.3484e-1

SD+: 0.2246e-10/-0.317e-10/0.3073e-1

Slope: 0.6581e-30/0.8744e-30/0.9846e-3

5 / 11149 0026

Table 4: Analytical raw data of analytical series ME95015

| Sample number | Type | Sample weight (g) | Final dilution(rat) | Sample injected (mg) | Peak height | Std. Conc. [µg/ml] | Substance injected (ng) | Residue (mg/kg) | Spiking level (µg) | Recovery % | Date of extraction | Date of measurement. |
|---------------|------|-------------------|---------------------|----------------------|-------------|--------------------|-------------------------|-----------------|--------------------|------------|--------------------|----------------------|
| 9500053       | C    | 5                 | 1                   | 2.5                  | 317         | 0.14171            | < 0.5                   |                 |                    | 7.395      |                    | 103.95               |
| ME142         | S    |                   |                     |                      | 1184        | 0.50               | 0.47788                 |                 |                    |            |                    | 103.95               |
| ME140         | S    |                   |                     |                      | 650         | 0.25               | 0.27083                 |                 |                    |            |                    | 103.95               |
| ME141         | S    |                   |                     |                      | 624         | 0.25               | 0.26075                 |                 |                    |            |                    | 103.95               |
| 9500356       | F    | 5                 | 5                   | 0.5                  | 2567        | 1.01412            | 10                      |                 |                    | 98.6       | 7.395              | 103.95               |
| ME143         | S    |                   |                     |                      | 1291        | 0.50               | 0.51987                 |                 |                    |            |                    | 103.95               |
| ME145         | S    |                   |                     |                      | 2582        | 1.00               | 1.01994                 |                 |                    |            |                    | 103.95               |
| 9500357       | F    | 5                 | 6                   | 0.417                | 2345        | 0.92805            | 10                      |                 |                    | 108.5      | 7.395              | 103.95               |
| ME158         | S    |                   |                     |                      | 697         | 0.25               | 0.28905                 |                 |                    |            |                    | 103.95               |
| ME160         | S    |                   |                     |                      | 1133        | 0.50               | 0.45811                 |                 |                    |            |                    | 103.95               |
| 9500358       | F    | 5                 | 6                   | 0.417                | 2303        | 0.91176            | 10                      |                 |                    | 106.6      | 7.395              | 103.95               |
| ME162         | S    |                   |                     |                      | 2602        | 1.00               | 1.02769                 |                 |                    |            |                    | 103.95               |
| ME159         | S    |                   |                     |                      | 596         | 0.25               | 0.24989                 |                 |                    |            |                    | 103.95               |
| 9500359       | F    | 5                 | 6                   | 0.417                | 2130        | 0.84468            | 10                      |                 |                    | 98.5       | 7.395              | 103.95               |
| ME161         | S    |                   |                     |                      | 1128        | 0.50               | 0.45617                 |                 |                    |            |                    | 103.95               |
| ME163         | S    |                   |                     |                      | 2454        | 1.00               | 0.97031                 |                 |                    |            |                    | 103.95               |

C = Control      T = Treated      F = Fortified sample      S = Standard

Injection volume: 1.0 µl

Calibration: linear

Coef. of correlation: 0.9911

Daily factor: 1

Slope: 0.3877e-3

Intercept: 0.1880e-1

SD4: 0.3102e-1

5 / 1140 0027

**Table 5:** Analytical raw data of analytical series ME95026 (PT 1-8)

| Sample number | Type | Sample weight (g) | Final dilution(ml) | Sample injected (mg) | Peak height | Std Conc. (µg/ml) | Substance injected (ng) | Residue (mg/kg) | Spiking level (µg) | Recovery % | Date of extraction | Date of measurement. |
|---------------|------|-------------------|--------------------|----------------------|-------------|-------------------|-------------------------|-----------------|--------------------|------------|--------------------|----------------------|
| 9500354       | C    | 5                 | 1                  | 2.5                  | 215         | 0.06696           | <0.5                    |                 |                    | 29.595     | 31.5.95            |                      |
| ME46          | S    |                   |                    |                      | 629         | 0.25              | 0.22602                 |                 |                    |            | 31.5.95            |                      |
| ME48          | S    | 5                 | 3                  | 0.833                | 1449        | 0.50              | 0.54107                 |                 |                    |            | 31.5.95            |                      |
| 9500360       | F    | 5                 | 3                  | 0.833                | 2293        | 1.00              | 0.86534                 |                 | 10                 | 81.9       | 29.595             |                      |
| ME50          | S    |                   |                    |                      | 2599        |                   | 0.98290                 |                 |                    |            | 31.5.95            |                      |
| 9500361       | F    | 5                 | 3                  | 0.833                | 2272        |                   | 0.85727                 |                 | 10                 | 81.2       | 29.595             |                      |
| ME47          | S    |                   |                    |                      | 598         | 0.25              | 0.24468                 |                 |                    |            | 31.5.95            |                      |
| 9500362       | F    | 5                 | 3                  | 0.833                | 2070        |                   | 0.86728                 |                 | 10                 | 81.6       | 29.595             |                      |
| ME49          | S    |                   |                    |                      | 1221        | 0.50              | 0.50819                 |                 |                    |            | 31.5.95            |                      |
| 9500363       | F    | 5                 | 3                  | 0.833                | 1985        |                   | 0.88133                 |                 | 10                 | 78.1       | 29.595             |                      |
| ME151         | S    |                   |                    |                      | 2377        | 1.00              | 0.98713                 |                 |                    |            | 31.5.95            |                      |

C = Control   F = Fortified sample   S = Standard

Injection volume: 1.0 µl

Aliquotation 50%

Deriv. factor: 1.62

Intercept: -0.15649-1/0.82529-2

Slope: 0.38429-30.4230e-3

SD+: 0.50546-10.10186-1

Coef. of correlation: 0.99120.9996

05/11/14 0028

Table 6: Analytical raw data of analytical series ME95026 (PT 1-9)

Study Code: 30976  
Method validation with PT 1-9



Chromfile name: ME95026  
Matrix: Dried hops

| Sample number | Type | Sample weight [g] | Final dilution [ml] | Sample injected [mg] | Peak height | Std. Conc. [µg/ml] | Substance injected [ng] | Residue (mg/kg) | Spiking level [µg] | Recovery % | Date of extraction | Date of measurement |
|---------------|------|-------------------|---------------------|----------------------|-------------|--------------------|-------------------------|-----------------|--------------------|------------|--------------------|---------------------|
| 9500355       | C    | 5                 | 1                   | 2.5                  | 165         | 0.25               | 0.10618                 | < 0.5           |                    | 29.5.95    | 31.5.95            |                     |
| ME164         | S    |                   |                     |                      | 500         | 0.50               | 0.23952                 |                 |                    |            | 31.5.95            |                     |
| ME168         | S    |                   |                     |                      | 1135        | 0.50               | 0.49225                 |                 |                    |            | 31.5.95            |                     |
| 9500364       | F    | 5                 | 3                   | 0.833                | 2348        |                    | 0.97505                 |                 | 10                 | 90.8       | 29.5.95            |                     |
| ME172         | S    |                   |                     |                      | 2377        | 1.00               | 0.98659                 |                 |                    |            | 1.6.95             |                     |
| 9500365       | F    | 5                 | 3                   | 0.833                | 2215        |                    | 0.92211                 |                 | 10                 | 85.7       | 29.5.95            |                     |
| ME165         | S    |                   |                     |                      | 594         | 0.25               | 0.27693                 |                 |                    |            | 1.6.95             |                     |
| 9500366       | F    | 5                 | 3                   | 0.833                | 2449        |                    | 1.01525                 |                 | 10                 | 94.7       | 29.5.95            |                     |
| ME169         | S    |                   |                     |                      | 1120        | 0.50               | 0.48628                 |                 |                    |            | 1.6.95             |                     |
| 9500367       | F    | 5                 | 3                   | 0.833                | 2235        |                    | 0.93007                 |                 | 10                 | 86.5       | 29.5.95            |                     |
| ME173         | S    |                   |                     |                      | 2457        | 1.00               | 1.01843                 |                 |                    |            | 1.6.95             |                     |

C = Control   F = Fortified sample   S = Standard

Aliquotation 50%  
Injection volume: 1.0 µl  
Calibration: linear  
Coeff. of correlation: 0.9973

Dev. factor: 1.611  
Slope: 0.3980e-3  
SD<sub>y</sub>: 0.2002e-1

**Figure 2:** PT 1-3 Standard 0.25 µg/ml**Figure 3:** PT 1-3 Standard 0.5 µg/ml

85/11149-0030



**Figure 4:** PT 1-3 Standard 1.0 µg/ml**Figure 5:** Example of reagent blank (no.: 9500062)

Figure 6: PT 1-3 Example of calibration curve




Figure 7: Dried hops control sample no.: 9500054







Figure 8: Dried-hops control sample no.: 9500055



35/11149.0033

**Figure 9:** Dried hops control sample no.: 9500056**Figure 10:** Dried hops control sample no.: 9500057

95/11149 0034

Figure 11: Dried hops spiked with 0.5 mg/kg hexythiazox (no.: 9500058)







Figure 12: Dried hops spiked with 0.5 mg/kg hexythiazox (no.: 9500059)




95/1114S 0035


**Figure 13:** Dried hops spiked with 0.5 mg/kg hexythiazox (no.: 9500060)**Figure 14:** Dried hops spiked with 0.5 mg/kg hexythiazox (no.: 9500061)

95/11149 0036

**Figure 15:** Dried hops spiked with 2.0 mg/kg hexythiazox (no.: 9500063)



**Figure 16:** Dried hops spiked with 2.0 mg/kg hexythiazox (no.: 9500064)



95/11149 0037

Study code 30796

Figure 17: Dried hops spiked with 2.0 mg/kg hexythiazox (no.: 9500065)





Figure 18: Dried hops spiked with 2.0 mg/kg hexythiazox (no.: 9500066)



95/11148 0038

Figure 19: Dried hops spiked with 2.0 mg/kg PT 1-3 (no.: 9500356)









Figure 20: Dried hops spiked with 2.0 mg/kg PT 1-3 (no.: 9500357)



95/11149 0039


**Figure 21:** Dried hops spiked with 2.0 mg/kg PT 1-3 (no.: 9500358)**Figure 22:** Dried hops spiked with 2.0 mg/kg PT 1-3 (no.: 9500359)

95/11140 0040


**Figure 23:** Dried hops spiked with 2.0 mg/kg PT 1-8 (no.: 9500360)**Figure 24:** Dried hops spiked with 2.0 mg/kg PT 1-8 (no.: 9500361)

95/11146-0041



**Figure 25:** Dried hops spiked with 2.0 mg/kg PT 1-8 (no.: 9500362)



**Figure 26:** Dried hops spiked with 2.0 mg/kg PT 1-8 (no.: 9500363)



05/11148. 0042

**Figure 27:** Dried hops spiked with 2.0 mg/kg PT 1-9 (no.: 9500364)**Figure 28:** Dried hops spiked with 2.0 mg/kg PT 1-9 (no.: 9500365)

35/11149 0043

Figure 29: Dried hops spiked with 2.0 mg/kg PT 1-9 (no.: 9500366)






Figure 30: Dried hops spiked with 2.0 mg/kg PT 1-9 (no.: 9500367)



95/11143 0044

Figure 31: Example of GC-MS run of PT 1-3 with NH<sub>3</sub> Cl

25/11149.0045

Figure 32: Analytical certificates of hexythiazox

**BASF****BASF Aktiengesellschaft**

Landeschematische Versuchsstation der BASF • Postfach 120 • 6711 Ludwigshafen

Landeschematische Versuchsstation der BASF • Postfach 120 • 6711 Ludwigshafen

**Purity Statement**

External Reference Standard -

APS/UP  
Product Chemistry

Reg. No.: 282502


Lot-No.: L43-251

Common name: Hexythiazox

Supplier: Nippon Soda Co., Ltd.  
Chemical name: (4RS, 5RS)-5-(4-chlorophenyl)-N-cyclohexyl-1-4-methyl-1-2-oxo-1,3-thiazolidine-3-carboxamideEmpirical formula: C<sub>12</sub>H<sub>21</sub>CIN<sub>2</sub>O<sub>2</sub>S

Molecular weight: 352.89

Structure:



Purity stated: &gt;99.9%, homogeneous

Date: 04.11.1993

Purity stated: 98.9%, homogeneous

Date: 04.11.1993

Methods used for purity check:

HPLC:  
Stationary Phase: Zorbax BP ODS, 250 x 4.6mm  
Mobile Phase: CH<sub>3</sub>CN + H<sub>2</sub>O = 75 + 25  
Wavelength: 240nm  
Flowrate: 1.5ml/minHandling advice:  
Keep in refrigerator.  
Before opening cooled sample vial must be warmed to room temperature in desiccator.Stability:  
Expected to be stable at least for 2.5 years at room temperature (PCP1404).

Purity Statement valid until: 04/96

Date: 24 SEP 98

Signature:

**Purity Statement**APS/UP  
Product Chemistry

Lot-No.: 2E172D


Reg. No.: 282502

Study Code: PCP00526

Common name: Hexythiazox  
Chemical name: (4RS, 5RS)-5-(4-chlorophenyl)-N-cyclohexyl-1-4-methyl-1-2-oxo-1,3-thiazolidine-3-carboxamideEmpirical formula: C<sub>12</sub>H<sub>21</sub>CIN<sub>2</sub>O<sub>2</sub>S

Molecular weight: 352.89

Structure:



Purity stated: 98.9%

Date: 07.12.1998

Purity stated: 98.9%

Date: 20.10.1994 (PCP03305)

Method of Purification: Recrystallization from diisopropylether

Methods used for purity check:

HPLC

Storage advice:  
Keep in refrigerator.  
Before opening cooled sample vial must be warmed to room temperature.Stability:  
Expected to be stable at least for 6 years at room temperature.

Purity Statement valid until: 11/2000

Date: 12 MAI 1995

Signature:

05/141400046

Figure 33: Analytical certificates of PT 1-3

BASF Aktiengesellschaft

**BASF**

Landwirtschaftliche Versuchsstation der BASF - Postfach 120 - 67114 Limburgerhof

### Purity Statement

#### - External Reference Substance -

APS/UP  
Product Chemistry

Reg.No.: 248649


Lot-No.: L43-255

Supplier: Nippon Soda Co., Ltd.

Chemical name: --

Empirical formula: C<sub>10</sub>H<sub>10</sub>ClNOS  
Molecular weight: 227.71**Structure:**

PT-1-3



Purity stated: &gt;99.9%, homogeneous

Date: 27.12.1993

## Methods used for purity check:

**HPLC:**

Stationary Phase: CAPCELL PAK C18, 250 x 4.6mm

Mobile Phase: CH<sub>3</sub>OH + H<sub>2</sub>O = 7 + 3

Wavelength: 223.235nm

Flowrate: 1.0ml/min

Result: &gt;99.9%

**Storage advice:**

Keep in freezer.

Before opening cooled sample vial must be warmed to room temperature.

**Stability:**Expected to be stable at least for 6 years stored in freezer  
(PCP03304).

Purity Statement valid until: 01/2000

Date: 12. MAI 1995

Signature:

R. Reijer

95/11149 0047

Figure 34: Analytical certificates of PT 1-8

BASF Aktiengesellschaft

**BASF**

Landwirtschaftliche Versuchsstation der BASF - Postfach 120 - 67114 Limburgerhof

**Purity Statement**

- External Reference Substance -

APS/UP  
Product Chemistry

Reg.No.: 248650

Lot-No.: L43-365

Supplier: Nippon Soda Co., Ltd.


Common name: --/-

Chemical name: --/-

Empirical formula: C<sub>17</sub>H<sub>21</sub>ClN<sub>2</sub>O<sub>3</sub>S  
Molecular weight: 368.89

Structure:

Code: PT-1-8-t



Purity stated: 99.9%, homogeneous

Date: 27.12.1993

Methods used for purity check:

HPLC:

Stationary Phase: CAPCELL PAK C18, 250 x 4.6mm

Mobile Phase: CH<sub>3</sub>OH + H<sub>2</sub>O = 75 + 25

Wavelength: 225.240nm

Flowrate: 1.0ml/min

Result: 99.9%

The analysis is conducted in compliance with GLP.

Storage advice:

Keep in freezer.

Stability: --/-

Expiry Date: --/-

Date: 09 MAI 1995

Signature:

R. [Signature] 35/11148 0048

Figure 35: Analytical certificates of PT 1-9

BASF Aktiengesellschaft

Landwirtschaftliche Versuchsstation der BASF - Postfach 120 - 67114 Limburgerhof



## Purity Statement

- External Reference Substance -

APS/UP  
Product Chemistry

Reg. No.: 248667

Lot-No.: L43-267

Supplier: Nippon Soda Co., Ltd.


Common name: -/-

Chemical name: -/-

Empirical formula: C<sub>17</sub>H<sub>19</sub>ClN<sub>2</sub>O<sub>3</sub>S  
Molecular weight: 366.87

Structure:

Code: PT-i-9



Purity stated: 99.9%, homogeneous

Date: 27.12.1993

Methods used for purity check:

HPLC:

Stationary Phase: CAPCELL PAK C18, 250 x 4.6mm

Mobile Phase: CH<sub>3</sub>OH + H<sub>2</sub>O = 7 + 3

Wavelength: 225.240nm

Flowrate: 1.0ml/min

Result: 99.9%

The analysis is conducted in compliance with GLP.

Storage advice:

Keep in freezer.

Stability: -/-

Expiry Date: -/-

Date: 14 MAI 1995

Signature:

R. Reichenbach

95/11143 0049