

US EPA ARCHIVE DOCUMENT

Application of Bicyclic Sesquiterpanes in Forensic Identification of Light Petroleum Products Spilled on Water



Chun Yang, Zhendi Wang, Bruce Hollebone, Carl E. Brown
Mike Landriault, Zeyu Yang

Emergencies Science and Technology Section
Science and Technology Branch, Environment Canada
335 River Road, Ottawa, Ontario, Canada K1A 0H3
E-mail: chun.yang@ec.gc.ca




What is Bicyclic Sesquiterpanes?


Decalin



Sesquiterpanes



- Polymethyl-substituted decalins
- Containing three isoprene subunits
- Containing two hexyl-ring (or cyclohexyl-).



What is Bicyclic Sesquiterpanes?


- The bicyclic biomarker sesquiterpanes with the drimane skeleton (C_{14} - C_{16}) are ubiquitous components of crude oils and ancient sediments.
- Most sesquiterpanes probably originate not only from *high plants* but also from *algae or bacteria*
- The concentration of C_{14} sesquiterpanes is higher at the immature stage, while those of C_{15} drimane and C_{16} homodrimane are relatively lower.

What is Bicyclic Sesquiterpanes?

- Early studies mainly focused on geological application of sesquiterpane compounds.
- The abundance of major compounds in the m/z 123 chromatograms has been employed to differentiate the organic matter input from various sedimentary environments.
- The naturally occurring bicyclic sesquiterpanes are stable in biodegradation, therefore, have potential applications in oil-source correlation and differentiation.



Why Sesquiterpanes?

- Biomarkers play a very important role in characterization, correlation, differentiation, and source identification in environmental forensic investigations of oil spills.
- The commonly used biomarkers include pentacyclic triterpanes (e.g., hopanes), regular and rearranged steranes, and mono- and tri-aromatic steranes.





Why Sesquiterpanes?

- The high MW biomarkers occur within the residual range. For lighter petroleum products like jet fuels and diesel, the refining processes remove most of these compounds from the original oil feedstock.
- The smaller bicyclic sesquiterpanes are concentrated in these products.



GC/MS Determination of bicyclic sesquiterpanes







Sesquiterpanes in Crude Oils and Petroleum Products



Sesquiterpanes in Crude Oils

Sesquiterpanes
(m/z 123)

Biomarker terpanes
(m/z 191)









Sesquiterpanes in Crude Oils

Oil samples	Total (μg/g oil)
Mississippi Canyon	917
Maya	969
Mars TLP	976
West Delta Canyon	1,010
Platform Elly	1,129
Orinoco bitumen	1,302
Arabian Heavy	1,904
Cold Lake bitumen	2,124
Alaska North Slope	2,179
Cook Inlet	2,211
Prudhoe Bay	2,551
Federated	3,670
Troll	5,119
South Louisiana	5,155

Sesquiterpanes in Petroleum Products

Sesquiterpanes

Biomarker terpanes



Sesquiterpanes in Petroleum Products

Oil samples	Total (µg/g oil)
Winter gasoline	ND
Kerosene	0.7
Aviation gasoline	30.6
Jet-A	6,770
Diesel-Pioneer	7,628
Diesel-Shell	8,703
Diesel-Stinson	8,447
Diesel-Ottawa	6,028
Korea diesel #1	7,658
Korea diesel #2	5,819
Korea diesel #3	6,129
Fuel #4	4,221
IFO-180	2,033
Fuel No. 5	1,449
Bunker C	255
Valvoline-10W-30	10.1
Pennzoil synthetic-10W-30	15.3
Extreme pressure gear	15.6

Summaries: Sesquiterpanes in Oils


- The presence of bicyclic sesquiterpanes in crudes and petroleum products are apparent.
- The oils from different sources demonstrate differences in both the absolute concentrations and relative distribution patterns of sesquiterpanes.
- The pentacyclic triterpanes and steranes are generally absent in lighter petroleum products, but bicyclic sesquiterpanes are present in diesels and certain jet fuels with significant abundances.

Sesquiterpanes in Distillation Fractions


Alaska north slope crude oil	F1	F2	F3	F4
Mass (%)	20.6	15.5	31.3	29.5
Boiling point range (°C)	IBP-180	180-287	287-481	>481
Main carbon range	<C ₁₂	C ₉ – C ₁₈	C ₁₄ – C ₃₄	>C ₃₄
Hydrocarbon groups (%)				
Total petroleum hydrocarbons	11.1	22.5	43.4	23.0
Total saturated hydrocarbons	13.3	26.8	41.2	18.8
Total aromatic hydrocarbons	4.0	9.3	50.4	36.2
n-Alkanes	20.2	34.8	41.8	3.2
Biomarker terpanes	0.0	0.3	53.4	46.3
Biomarker steranes	0.0	0.2	90.4	9.5
C ₃₀ bb-Hopane	0.0	0.1	51.3	48.6

03/13/2008

Biomarker Terpanes



Sesquiterpanes











Sesquiterpanes in Distillation Fractions

Sesquiterpanes (%)	F1	F2	F3	F4
BS1	1.0	89.4	9.6	0.0
BS2	1.3	86.6	12.1	0.0
BS3	0.8	71.3	28.0	0.0
BS4	0.6	62.5	36.9	0.0
BS5	0.4	56.1	43.4	0.0
BS6	0.4	52.6	47.0	0.0
BS7	0.3	48.4	51.3	0.0
BS8	0.2	43.5	56.3	0.0
BS9	0.3	40.4	59.3	0.0
BS10	0.4	37.5	62.1	0.0

Summaries: Sesquiterpanes in Distillation Fractions

- Biomarker terpanes and steranes were found almost entirely in heavier fractions 3 and 4 due to their high boiling points
- Sesquiterpanes were largely found in fractions 2 and 3.
- For fraction 2, all ten target sesquiterpanes were concentrated up to approximately five-fold compared to the native crude oils
- Sesquiterpanes were in low concentration in the lightest fraction 1 and rarely found in the heaviest fraction 4.

Weathering Effects on Sesquiterpane Distributions



Summaries: Weathering Effects on Sesquiterpane Distributions

- In general, **light to medium** weathering has little effects on distribution pattern and diagnostic ratio values of sesquiterpanes.
- Smaller C₁₄ sesquiterpanes (BS1 and BS2) start to be lost as oil subjected to high percentages of weathering .
- Most of sesquiterpane ratios are robust for correlation and differentiation of oil samples in that have only lightly to moderately weathered
- The **most heavily weathered** might result in certain changes for such ratios as BS1/BS5, BS3/BS10 and BS5/BS10.

Application of Bicyclic Sesquiterpanes in Forensic Identification of Oil Spill

Case Study:

Oil Spill in Northbay, Ontario

Oil Spill in Northbay, Ontario


- In 2008, Wagg's Petroleum Equipment Ltd. was hired by the Department of National Defence (DND) and Canadian Forces Exchange System (Canex) at Canadian Force Base (CFB) at North Bay to clean up a hydrocarbon spill at one of the sites.
- The site is a fueling site for gasoline only and this site has had 3 spills; one in 1998, one in 2002 and one in 2008
- Preliminary laboratory test indicates the spilled oil is not pure gasoline, could be a mixture with other type of fuel.

Oil Spill in Northbay, Ontario

- Spill site: Northbay, Ontario
- Two water samples:
- 2008/04/22/-1290,

2008/04/22/-1291


- Oil layer on water surface was clearly observed.



Oil Fingerprinting Analysis





- The water samples were separated into two parts: oil phase and water phase respectively;
- Oil layer: Approximately 0.5 mL of oil was removed from the top oil layer of two water samples. About 16.0 mg of oil was used for fingerprinting analysis.
- Water layer: The water samples were weighed and then transferred to separatory funnels and extracted (liquid/liquid extraction) with dichloromethane (DCM) three times (50/50/50 mL).

Oil Fingerprinting Analysis




Hydrocarbon Group Analysis Results

ESTD code		2008/04/22-1290Oil	2008/04/22-1290L	2008/04/22-1291Oil	2008/04/22-1291L	Diesel
Customer sample info.		Oil Layer	Water	Oil Layer	Water	Shell [2003]
TPH (mg/g)		768	803	805	705	957
TSH/TPH (%)		75.8	80.2	84.4	77.9	89.8
TAH/TPH (%)		24.2	19.8	15.6	22.1	13.6
Resolved/TPH (%)		57.0	53.5	53.9	52.0	22.3
Σn-Alkanes (mg/g)		72.4	77.4	76.5	67.9	128
ΣEPA PAHs (mg/g)		0.24	0.30	0.25	0.27	0.61
ΣAPAHs (mg/g)		13.4	16.6	13.6	14.9	19.9
ΣBiomarkers (mg/g)		0.052	0.067	0.060	0.060	0.014




N-Alkane Distribution


2008/04/22-1290Oil



2008/04/22-1291Oil




BTEX and PAHs Analysis




BTEX and PAHs Analysis



2008/04/22-1290Oil



2008/04/22-1291Oil



Biomarkers Analysis





Answer

- The spilled oil is a mixture of a Gasoline and a Diesel !!!

Questions

- How much *Gasoline* ?
- How much *Diesel* ?

Bicyclic Sesquiterpanes Analysis



Sesquiterpanes Analysis

Oil samples	Total (µg/g)
2008/04/22-1290 Oil	3,476
2008/04/22-1290 Water	3,973
2008/04/22-1291 Oil	4,022
2008/04/22-1291 Water	3,772
Mean	3811
Gasoline	ND
Kerosene	0.7
Diesel-Pioneer	7,628
Diesel-Shell	8,703
Diesel-Stinson	8,447
Diesel-Ottawa	6,028
Mean	7702

Answer

- The spilled oil is a mixture of a Gasoline and a Diesel
- Gasoline : Diesel = ~50:50 !!!

Simulation of Spilled Oils


Gasoline: from Petro Canada gas station,
Ottawa

Diesel: from Shell gas station, Ottawa


Gasoline : Diesel = 100:0, 90:10, ..., 0:100 (v:v).

Simulation of Spilled Oils


1290 F3




Diesel




Gasoline




Gasoline:Diesel=20:80 (v/v)



Gasoline:Diesel=50:50 (v/v)



Gasoline:Diesel=70:30 (v/v)



Conclusions

- Bicyclic sesquiterpanes are widely present in crude oils and diesel range petroleum products with significant abundances.
- Sesquiterpanes offer potential criteria for both oil identification and oil-source correlation in cases where the traditional biomarkers are absent due to refining of petroleum product or environmental weathering.

Conclusions

- The oil samples are gasoline/diesel type fuel.
- The oil samples are probably mixtures of gasoline and diesel (No. 2) at an approximate ratio of 50:50 (v:v).

Thank You!