


US EPA ARCHIVE DOCUMENT

Pathogen Association with Drinking Water Distribution Systems

T. H. (Helen) Nguyen and W.T. Liu, University of Illinois at Urbana-Champaign

R. M. Hozalski, Univ. of Minnesota



Project Objectives


- To identify the factors that influence the persistence of a model pathogen, *Legionella*, in multi-species biofilms grown on different pipe materials (Illinois).
- To characterize the microbial communities and the surface properties of the biofilms grown on different pipe materials and chlorine level (both Illinois and Minnesota).
- To further investigate novel chemical and enzymatic treatments to weaken the biofilm to promote detachment and improved cleaning (Minnesota).

Methodology (Illinois)


Biofilms are grown in CDC reactors using groundwater, tap water and dechlorinated tap water




Biofilm analysis






E. coli or
Legionella
deposition





Microbial community analysis using Terminal Restricted Fragment Length Polymorphism (TRFLP) showed relatively stable biofilm after 16 weeks




gfp-tagged environmental *E. coli* is used in deposition experiment.



Legionella pneumophila (ATCC33152) tagged with GFP plasmid pBG307 (Chen et al., Plasmid, 2006)



Parallel plate flow cells are used to observe *E. coli* or *Legionella* deposition on PVC and biofilms.




Sherwood number calculation for *E. coli* or *Legionella* adhesion to biofilm surface

Deposition rate coefficient (k_D)(m/s)


$$k_D = \frac{J}{C_0}$$

$$\text{Sherwood number } Sh = \frac{J \times a}{C_0 \times D}$$




Normalized
by viewing
area (296 X
222 μm)

a = bacteria diameter (m)
 D = diffusion coefficient
(m^2/s)




Sherwood
number

$$Sh = \frac{J \times a}{C_0 \times D}; \quad \text{confocal } J = \frac{N}{A \times t}$$



- Env. *E. coli* on PVC
- Env. *E. coli* on 2 week GW biofilms
- Env. *E. coli* on 8 week GW biofilms
- Env. *E. coli* on 14 week GW biofilms
- Env. *E. coli* on 20 week GW biofilms



- Env. *E. coli* on 2 week GW biofilms
- Env. *E. coli* on 8 week GW biofilms
- Env. *E. coli* on 14 week GW biofilms

Clean PVC

2w biofilm



8w biofilm

14w biofilm

20w biofilm

- Higher *E. coli* deposition on older biofilms compared to on PVC.
- *E. coli* deposition on old biofilms did not increase with ionic strength

The trend in *Legionella* deposition on biofilms is similar to the *E. coli* trend



Clean PVC

2w biofilm


4w biofilm

8w biofilm


12w biofilm

14w biofilm

Optical coherence tomography (OCT) captures micrometer-resolution, three-dimensional images from within optical scattering media.




Biofilm surface roughness is calculated from analysis of optical coherence tomography (OCT) images.



Mean thickness $42 \mu\text{m}$

Quantitative biofilm characterization



Mean Biofilm Thickness

$$\bar{z} = \frac{1}{n} \sum_{i=1}^N z_i$$


Absolute roughness

$$R_a = \frac{1}{n} \sum_{i=1}^N \left(\left| z_i - \bar{z} \right| \right)$$


Relative roughness coefficient

$$R'_a = \frac{1}{n} \sum_{i=1}^N \left(\frac{\left| z_i - \bar{z} \right|}{\bar{z}} \right)$$

Legionella deposition correlate with biofilm roughness





- Adhesion experiments with microscope
- Deposition experiments with CLSM





- Adhesion experiments with microscope
- Deposition experiments with CLSM

OCT and SEM images show no difference in untreated and chloramine treated biofilms



E. Coli attaches more on rougher biofilm



Conclusion and future work for deposition study


- Higher *Legionella* and *E. coli* deposition on biofilms compared to on PVC surface was observed.
- *Legionella* and *E. coli* deposition is controlled by biofilm roughness not the water ionic concentration.
- 20 week of exposure to chloramine did not change the biofilm roughness. Similar deposition was observed on biofilms with and without exposure to chloramine.

Conclusion and future work for deposition study

- Deposition will be conducted with biofilms grown on concrete coupons and coupons coated with calcite.
- Deposition will be conducted with biofilms treated with biofilm weakening agents.

UMn Tasks

- Investigate effect of pipe material (and chlorine level) on biofilm development, community composition, and detachment/removal
 - CDC Reactors
 - 24 Well Plates
- Test effect of treatments on biofilm weakening and removal
 - 96 Well Plates
 - Metal chelators: polyphosphate, citrate



A: Tap water

B: Pump (3.9 ml/min)

C: Dechlorination reactor

**D: Reservoir (Sodium Thiosulfate
1.17 g/L)**

E: Pump (0.1 ml/min)


F: Pump (3.9 ml/min)

**G1, G2, G3: CDC Biofilm Reactor
(125 rpm, dechlorinated tap water)**




**H: CDC Biofilm Reactor
(125 rpm, tap water)**




CDC Reactors



CDC Reactors – *Particle Counter*



Assimilable Organic Carbon(AOC)



$$AOC (\mu g/L) = [AOC_{P17} + AOC_{NOX}]$$

$$AOC (\mu g/L) = \left[\frac{\text{max density P17 } (\frac{cfu}{ml})}{\text{yield P17 } (\frac{cfu}{\mu g})} + \frac{\text{max density NOX } (\frac{cfu}{ml})}{\text{yield NOX } (\frac{cfu}{\mu g})} \right] \times 1000 (\frac{ml}{L})$$

Expediting the AOC test



PARTICLE COUNTER



FLOW CYTOMETRY




ATP ASSAY







PLATE COUNTS

Results (Pooled data)



Conceptual Models of Biofilm Cohesion



Courtesy of Phil Stewart at CBE

96 Well Plate Testing

✓ Preparatio



Culture for 24hrs

→
Innoculation



Culture for 24hrs




Washing using Milli-Q (use
pipet and washed two times


✓ Treatment and staining




Treatment for 2-24
hrs



Washing two times
using tap water



Staining for 15
min using crystal
violet



Washing three times
using tap water




Add 33% acetic
acid




Measure absorbance
at OD 600

96 well plates - *S. epidermidis*



24 Well Plate Testing


- Following 6 coupon materials were used
 - Cement
 - Glass
 - PVC
 - HDPE
 - Polycarbonate
 - Polypropylene
- Following 2 bacterial species were used :
 - *S. epidermidis*
 - *E. faecalis*



Experimental Plan

	1	2	3	4	5	6	
A	CE	GL	PV	HD	PC	PP	Cells
B	CE	GL	PV	HD	PC	PP	Controls
C	CE	GL	PV	HD	PC	PP	
D	CE	GL	PV	HD	PC	PP	

Glossary	
CE	Cement
GL	Glass
PV	PVC
HD	HDPE
PC	Polycarbonate
PP	Polypropylene



Summary

- Higher planktonic growth/activity in dechlorinated reactors
- Flow cytometry a good tool for expediting the current AOC protocol
- Effects of pipe surfaces and treatments on biofilm development/removal being explored using well-plate method
- Roughness matters (cement coupons)

Future Work

- CDC biofilm reactors. Test coupons for:
 - Biomass (total protein, DNA)
 - Thickness, roughness, coverage (CLSM)
 - Community analysis
 - Effect of treatment (weakening, detachment)
- 96 and 24 Well plate tests
 - Exploring polyphosphate and citrate treatments
 - Cultures: *S. epidermidis*, CDC reactor effluent
- Submit paper comparing flow cytometry, particle counting, and HPC for bacterial enumeration in the AOC test