US ERA ARCHIVE DOCUMENT




The Fate, Transport, Transformation and Toxicity of Manufactured Nanomaterials in Drinking Water

Arizona State University Investigators

Paul Westerhoff (p.westerhoff@su.edu) – Associate Professor Yongsheng Chen (yongsheng.chen@asu.edu) – Research Professor John Crittenden (j.crittenden@asu.edu) – Professor Yang Zhang & Kiril Hristovski – PhD students Civil and Environmental Engineering

David Capco (dcapco@asu.edu) – Professor Brian Koeneman – PhD Student Cellular and Molecular Biosciences Faculty, School of Life Sciences





Project Objectives

Goal: To understand the fate and significance of nanomaterials in drinking water

The objectives of this project are:

- 1) to characterize the fundamental properties of nanomaterials in aquatic environments
- 2) to examine the interactions between nanomaterials and pollutants and pathogens
- 3) to evaluate the removal efficiency of nanomaterials by drinking water unit processes
- 4) to test the toxicity of nanomaterials in drinking water using cell culture model system of the epithelium.



Project Timetable

Tasks	Scheduled Time		
Literature preparation	January 2005		
Nanomaterial Detection Methods	Oct 2004 – Mar 2005		
Characterization of Nanomaterials in Water	Jan 2005 – Dec 2005		
Adsorption of Dissolved Pollutants onto Nanoparticles	July 2006 – June 2007		
Aggregation and Coagulation of Nanomaterials	Jan 2005 – Mar 2007		
Nanoparticle Adsorption and Disinfectant Shielding of Virus	June 2006 – Mar 2007		
Nanoparticle Toxicity Screening for Drinking Water	Jan 2005 – Dec 2006		
Final Report	July 2007 - Oct 2007		



I. Nanoparticle Characterization

Interesting note: Several papers investigate environmental applications of nanoparticles, but use 0.45 um filtration to remove the "nano" particles.

We have found for metal oxide nanoparticles:

- Nanoparticles placed in water are aggregated
- Aggregation due to electrostatics of dry powders, manufacturing process, and/or aggregation in solution
- Sonication temporarily dis-aggregates some nanoparticles
- Surfactants and/or solvents promote some dis-aggregation
- Nanoparticles purchased in solutions are less, but still, aggregated
- Solution = Producing nanoparticles in laboratory rather than commercial sources

Challenge for all nanoparticle research: commercial nanoparticles in water are NOT nanoparticles (< 100 nm in at least one dimension)




Nanomaterials in Study

Nanoparticle	Vendor Reported Mean Particle Size	Density	Source	
Titanium dioxide	< 40 nm	3.9g/mL	5wt% Dispersion in water	
Aluminum oxide	< 20 nm	3.97g/mL	5wt% Dispersion in water	
Zinc Oxide	50 ~ 70 nm	5.61g/mL	Nanopowder	
Titanium dioxide	15 nm	3.9g/mL	99.7% Nanopowder	
Iron (III) Oxide	5 ~ 25 nm	5.24g/mL	Nanopowder	
Nickel Oxide	10 ~ 20 nm	6.67g/mL	99.8% Nanopowder	
Silica	10 nm	2.6g/mL	99.5% Nanopowder	

- All nanoparticles purchased as powders and liquid suspensions
- All values in above table are reported by Vendor.
- Other nanomaterials currently in use:
 - Several other commercial metal oxide nanoparticles
 - Carbon nanotubes and fullerenes
 - TiO2 nanotubes and nanoparticles fabricated in ASU laboratory
 - Gold and/or cadmium quantum dots




Size of nanoparticles in water



- 10mg/L and 5mg/L nanoparticles in Nanopure water sonicated for 15 min
- Filter paper with 3 µm pore size and 110nm diameter
- 100ml and 50ml suspensions for filtration
- Concentrations of particles analyzed by digestion/AAS.
- Conclusion:
 - Significant mass of "nanoparticles" are aggregated and > 3um
 - DLS particle size instrument only measures particles < 3 um



At high nanoparticle concentration (1g/L) more aggregation occurs



C(NP) = 1 g/L

10 mM NaHCO₃

 $pH = 8.4 \pm 0.2$

Analysis by:

Digestion/GF-AAS


Centrifugation:

G > 1300



10 mg/L TiO₂ in Water

(ZetaPALS dynamic light scattering; < 3um only)




Conclusion: Nanoparticles do not appear to be discrete NP in water




Size of nanoparticles in water


- Particle Size of nanoparticles were analyzed by Dynamic Light Scattering.
- 10mg/L nanoparticles in Nanopure water
- Sonicated for 15 min at 200 W/L and 20 kHz.
- Instrument range: 2nm
 3μm.
- Mean particle size is greater than reported by manufacturer.





SEM Analysis of Commercial TiO₂ NPs





TiO₂ NPs from **5wt% Dispersion in water**

TiO₂ NPs from a **dry powder**

Nanoparticles do not appear to be discrete NPs



Dispersion of aggregated NPs in water


- Sodium hexameta phosphate,
 Sodium dodecyl sulfate,
 Isopropanol, Acetone,
 Butonanone, Methanol, and
 Ethylene Glycol selected as a dispersant solvent.
- All suspensions sonicated for 15 mins
- ➤ With all above dispersants, mean particle sizes of TiO2 nanoparticles in Nanopure water still were much more than size of discrete nanoparticles.
- ➤ It is very difficult to disaggregate these NPs and obtain homogenous discrete nanoparticles.

Dispersant	Mean Particle Size
SDS (Sodium dodecyl sulfate)	530~570nm
(NaPO3)6 (Sodium hexameta phosphate)	560~580nm
Ethylene Glycol (10~100%)	300~550nm
Isopropanol (10%)	550~1050nm
Acetone (10%)	560~580nm
Butonanone (10%)	580~640nm
Methanol (10%)	550~570nm
Nanopure water (no dispersant)	540~600nm

Mean particle size of TiO2 by Dynamic Light Scattering in Nanopure water w/o dispersants


II. Removal of Nanoparticles in Simulated Water Treatment Systems








Further aggregation of nanoparticles in water




Mean particle size of 10mg/L nanoparticle suspensions with time. a.TiO2 b. NiO c. Fe2O3 d.ZnO



Electric Double Layer Compresssion





Removal of nanoparticles in nano pure water without coagulant and with 100mM MgCl2

Mean nanoparticle size in nano pure water without coagulant and with 100mM MgCl2


Test Conditions:

- 10 mg/L NPs in Nanopure water buffered by 10mM NaHCO3 (pH=8.1+ 0.2)
- 100 mM MgCl₂ for EDL compression and NP destabilization
- NPs in supernatant measured by digestion/AAS after sedimentation period Conclusions:
- Sedimentation removes aggregated NPs
- EDL compression leads to more aggregation for some NPs
- > 40% of NPs remain after sedimentation






TiO2 (5wt% Dispersion)






Mean particle size during coagulation experiment TiO2 (5wt% dispersion) in Nanopure water





III. Adsorption of Pollutants onto Nanoparticles

- Arsenic adsorption evaluated (see graphic)
- Adsorption of bacteriaphage (MS2 & PRD1):
 - Adsorption occurs onto positively charged NPs
 - Charge based upon zeta potential measurements



 $AS_{Initial} = 1 mg/L$

 $C(NP) = 1 g/L pH = 8.0 \pm 0.4$




IV. Toxicity of Nanoparticles

Objectives:

- Evaluate transport and necrotic effects of nanoparticles across epithelial layers (esophageal & intestinal)
- Uniform cells are established on semipermeable support. Continuous cell structure leads to conductivity gradient across biofilm. Disruption of conductivity gradient inferred as detrimental impact to cells or biofilm structure

Human intestinal tissue model (Caco-2BBe) – from ATCC:

- Cells were transferred to semi-porous membranes and allowed to anchor and form tight junctions (three days).
- Transepithelial electrical resistance (TEER) was utilized to monitor the density and junctional complex of the cell monolayer.
- Cells were maintained in DMEM containing 10% fetal calf serum, penicillin/streptomycin/fungizone, and transferrin at 37°C, humidified air containing 10% CO2. Subsequently, the medium was changed every day after seeding onto the membrane.



Filter Inserts



Voltohmmeter



Testing Procedures

Rejected Methodology

- 5% CO2
- Media absent of transferrin
- Membranes composed of polycarbonate
- Membranes composed of polyester
- 12mm membrane diameter
- 24mm membrane diameter
- 3.0 µm pore size in membrane
- Change media every three days
- Change media every two days

Adopted Methodology

- 10% CO2
- Media containing 10µg/mL transferrin
- Pre-coated collagen membranes
- 6.5mm membrane diameter
- 0.4µm pore size in membrane
- Change media daily



Optimization of Culture Conditions

Condition	Relative Growth						
	0.05	1	2	2	3	4	5
CO2 Level	5%						
	10%						
Transferrin (μg/mL)	0		10				
Filter diameter (mm)	24			12	2		6.5
Filter pore size (μm)	3	0.4					
Filter material	Polycarb	Poly	ester				
Filter casing	Untreated					Collag	ated-PTFE
Collogen coated	By hand		Prec	oated			
Frequency for media change-out (days)	3	2		1			





Optimization of Culture Conditions

- CO₂ Level: no difference with 5% or 10%
- Transferrin addition: 10 μg/L significantly better than no addition (costly but required)
- Filter diameter: 6.5 mm best (6.5 > 12 > 24 mm)
- Filter material: Polyester better than polycarbonate
- Filter coating:
 - Collagen coated better than untreated
 - Commercially precoated better than hand precoated
- Frequency for media changeout: 1 day better than 2 or 3 days




Caco-2 cells
9 days of growth in culture (10x)
Confluent cells are present







Scanning Electron Microscopy Nanoparticles may flatten down microvilla

Control 10ppm TiO₂









Confocal Scanning Laser Microscope


Red = gamma catenin



Green = Nuclei



Green =
Possible TiO₂
Particle
(no stain for nuclei)





Plan for Year 2

- I. Nanoparticle Characterization
 - Measure surface areas of commercial NPs
 - Synthesize in-lab true NP and nanotubes
 - Develop procedure to disaggregate commercial NPs
 - Measure DLS on sample from WTP effluent
- II. Nanoparticle removal during simulated water treatment
 - Evaluate role of size, surface charge density, hydrophobicity, shape
 - Removal after coag/sed and paper filter
- III. Adsorption of pollutants onto nanoparticles
 - Continue MS2/PRD1 work
 - Evaluate adsorption of 1hydrophobic and 1 hydrophilic SOC
- IV. Nanoparticle Toxicity
 - Test the toxicity of various nanoparticles
 - Determine effects of chronic and acute exposure to nanoparticles
 - Examine viability (live vs. dead) of the cells after treatment
 - Design esophageal and stomach models