

US EPA ARCHIVE DOCUMENT

3.13.81

BD-1061
TKR-365

Testing completed on 10/10/77.
EPA Acc. No. 236764

000365

000365

Mutagenicity Testing of Paraquat

Inveresk Research International Edinburgh, Scotland
Report No. 877
IRI Project No. 409042, dated October, 1977

Author : DB McGregor (IRI)

Report addressed to Dr. E. Longstaff, ICI Ltd.
Macclesfield, Cheshire

BEST AVAILABLE COPY

Notes taken from microfiche, received from Mr. Robert Taylor (PN; RD) on
3/13/81

1/20/85

(2)

000365

SUMMARY

"Paraquat dichloride was accepted from ICI Ltd. for mutagenicity testing with Salmonella typhimurium strains TA 1535, TA 1537, TA 1538, TA 98 and TA 100. The tests were conducted on agar plates both in the presence and absence of a preparation from the livers of male rats treated with Anachlor 1254 and the cofactors required for mixed function enzyme reactions.

Doses of paraquat dichloride used were between 1 µg and 1 mg/plate.

Clear toxic effects were observed at 1 mg/plate, while at 333 µg/plate there appeared to be inhibition of mutant colony growth but normal growth of non-mutant microcolonies. At lower dose levels, no mutagenic effect was observed.

It was concluded that paraquat dichloride was without mutagenic potential detectable by these tests."

REST AVAILABLE COPY

3

2

①

requires metabolic activation by oxidative reactions.

000365

MATERIALS AND METHODS

Test mat.: paraquat dichloride from ICI, 99.9% pure; white powder was dissolved in 0.05M phosphate buffer, pH 7.4.

2-Aminooxazirine (^{carcinogen} ~~known mutagen~~) was dissolved in dimethyl sulphoxide.

Animals Male rats, 200-300 g, were injected i.p. with Aroclor 1254 in corn oil, 500 mg Aroclor 1254 / kg body wt., 5 days before they were killed. Food was withdrawn 15 hr before they were killed.

Bacteria

All 5 strains had mutations in the histidine operon (had to have histidine for growth). The 3 mutations were: his G 46 in TA 1535, TA 100; his C 3076 in TA 1537; his D 3052 in TA 1538 & TA 98. The first type of mutation is reversed by a variety of mutagens that cause base-pair substitutions. The second type mutation is reversed by 9-aminocridine, Icl-191 & epoxides of polycyclic hydrocarbons. The 3rd type is reverted by aromatic amines & derivatives.

9000 g supernatant

BEST AVAILABLE COPY

Homogenized livers in 3 volumes their vol of 0.15M KCl (cold). Spin at 9000g for 10 min at 0° to get supernatant + pellet (mostly whole cells, nuclei & mito.)

Testing

Bacteria are grown in nutrient broth. A co-factor solution (0.05M PO₄ buffer, pH 7.4 + NADP-Na salt + glucose-6-PO₄ - diK-salt, + MgCl₂ · 6H₂O + KCl) is mixed with 9000 g supernatant in the ratio of 9:1 (S-9 mix). Then S-9 mix, bacteria, solvent or test solution & diluted agar are mixed. (Agar had histidine & biotin added to it). This mixture was poured on agar plates containing 2% glucose. The plates are incubated for 2 days at 37° & the colonies were counted. The plates

(4)

000365

were also examined for precipitates & microscopically, for microcolony growth.

RESULTS

Positive control system was highly sensitive to 2-Aminanthracene ($0.5 \mu\text{g}/\text{plate}$)

Paraquat dichloride No indications of mutagenic activity either with or without the mixed function oxidase system. Bacteria were almost completely killed at $1 \text{mg} \text{ per. dichloride}/\text{plate}$. At $333 \mu\text{g}/\text{plate}$ and, to a lesser extent, $100 \mu\text{g}/\text{plate}$, there was inhibition of mutant colony growth.

Conclusion "Paraquat dichloride was not found to be mutagenic in these tests, which achieved a satisfactory sensitivity. The substance was, however, toxic at a dose of $1 \text{mg}/\text{plate}$." Levels used: $10, 33, 100, 333, 1000 \mu\text{g}/\text{plate}$ of para. dichloride.

(Apparently, number of mutant colonies / plate is counted & compared with control, containing buffer in place of the test material. 2-Aminanthracene greatly increased colony growth). "At least a doubling of the control values was looked for at some concentration of the test substance if a mutagenic effect was to be suspected".

NOEL : $33 \mu\text{g}/\text{plate}$ (paraquat dichloride)

LEL: $100 \mu\text{g}/\text{plate}$. (paraquat dichloride)

BEST AVAILABLE COPY

Care - Minimum