US ERA ARCHIVE DOCUMENT

BAS 510 F BAS 510F and M510F01 PMRA a.i. code (CCH) Multiresidue Method OPPTS 860.1360 DACO7.2.4 PC Code: MRID: 128008

MRID: 45405107 Submission # 2001-1027, 1036, 1043



# UNITED STATES ENVIRONMENTAL PROTECTION AGENCY WASHINGTON, D.C. 20460

OFFICE OF PREVENTION, PESTICIDES AND TOXIC SUBSTANCES

# MEMORANDUM

Date:

July 2, 2003

Reviewers:

William T. Drew, Chemist

Reviewer

RAB2/HED (7509C)

Henri P. Bietlot, Chemist

Peer reviewer

FREAS, HED, PMRA

Richard A. Loranger

Branch Senior Scientist RAB2/HED (7509C) Ann Auy Section Head

FREAS, HED, PMRA

DP Barcode: D278386

Petition:

1F06313

Citation:

45405107 Fomenko, J. (2001) PAM I Multiresidue Testing for BAS 510 F and its

Hydroxy Metabolite (M 510 F01): Lab Project Number: A054.006:

2001/5000875: 46647. Unpublished study prepared by Maxim Technologies, Inc.

89 pages.

Sponsor:

**BASF** Corporation

## Background

The information contained herein was compiled by Dynamac Corporation (20440 Century Boulevard, Suite 100, Germantown MD 20874), contractor, under the supervision of RAB2/HED. This DER has undergone secondary review by RAB2, and reflects current HED and Office of Pesticide Programs (OPP) policies. This DER has also been peer-reviewed by PMRA/Canada.

## **Executive Summary**

Residues of BAS 510 F and its metabolite M510F01 were not adequately recovered using the multiresidue methods. BAS 510 F and M510F01 were not evaluated by Protocol A because

BAS 510 F	•
<b>BAS 510F</b>	and M510F01
PMRA a.i.	code (CCH)

Multiresidue Method **OPPTS 860.1360** DAC07.2.4

PC Code: MRID: Submission # 2001-1027, 1036, 1043

128008 45405107

neither analyte incorporates an N-methylcarbamate structure, and BAS 510 F was not evaluated by Protocol B since it is not an acid or a phenol. Methylated M510F01 was recovered at ~55% from Protocol B Florisil, and M510F01 was recovered at -95% from Protocol B GPC. However, the response of methylated M510F01 at levels necessary to quantify recovery under this protocol was inconsistent and non-linear, therefore no additional work was performed. BAS 510 F and M510F01 had good responses with GC/ECD on a DB-1 column at standard and elevated temperature GC conditions under Protocol C. BAS 510 F and M510F01 were not recovered at ≥30% using Protocols D, E, and F. These data will be forwarded to the FDA for further evaluation.

# **GLP Compliance**

Signed and dated GLP, quality assurance, and data confidentiality statements were provided. The petitioner stated that reference standards of chlorpyrifos, ethion, p,p'-DDT, phosalone, and permethrin for the testing of Protocol C were purchased commercially and were not characterized according to GLP standards. The standards used were believed to be of highest quality; this deviation was not expected to impact the study results or their interpretation.

#### 1. Test Substances

Table 1.1.1. List of Analytes Tested.				
Common Name:	Nicobifen, proposed (parent compound)	Hydroxy metabolite		
IUPAC Name:	2-Chloro-N-(4'-chlorobiphenyl-2-yl)- nicotinamide	2-Chloro-N-(4'-chloro-5-hydroxybiphenyl-2-yl)-nicotinamide		
CAS Name:	3-Pyridinecarboxamide, 2-chloro-N- (4'chloro[1,1'-biphenyl]-2-yl)-	Not available		
CAS Number:	188425-85-6	Not available		
Company Name:	BAS 510 F	M510F01		
Other Synonyms:	BASF Registry No. 300355	BASF Registry No. 398794		

#### 2. Results

Table 2.1. Results of Multiresidue Methods Testing with BAS 510 F and its metabolite M510F01.				
PAM I Protocol	Analyte	Results	Comments	
A	BAS 510 F	Not evaluated because neither analyte incorporates an N-methylcarbamate structure.		
	M510F01			
В	BAS 510 F	Not evaluated because it is not an acid or a phenol.		

BAS 510 F BAS 510F and M510F01 PMRA a.i. code (CCH)

Multiresidue Method **OPPTS 860.1360** DAC07.2.4

PC Code: MRID:

128008 45405107 Submission # 2001-1027, 1036, 1043

Table 2.1. Results of Multiresidue Methods Testing with BAS 510 F and its metabolite M510F01.				
PAM I Protocol	Analyte	Results	Comments	
	M510F01 (and methylated M510F01)	Methylated M510F01 recovered ~55% from Protocol B Florisil; M510F01 recovered ~95% from Protocol B GPC.	The response of methylated M510F01 at levels necessary to quantify recovery under this protocol was inconsistent and non-linear, therefore no additional work was performed.	
С		Good response with GC/ECD on a DB- 1 column at standard and elevated	Only the DB-1 column in combination with ECD was used because responses	
M510F01 (and methylated M510F01) temperature GC conditions to the conditions of the	temperature GC conditions.	were good for both analytes using this system.		
D		Not recovered at ≥30% through Protocol D Florisil using either	No additional work was performed with BAS 510 F or M510F01.	
	M510F01	methylene chloride or mixed ether elution systems.		
E	BAS 510 F	Not recovered at ≥30% through Protocol E Florisil using either	No additional work was performed with BAS 510 F or M510F01.	
	M510F01	methylene chloride or mixed ether elution systems.		
F	BAS 510 F	Not recovered at ≥30% through Florisil	No additional work was performed with BAS 510 F or M510F01.	
	M510F01	using either methylene chloride or mixed ether elution systems.		

# 3. Discussion

Residues of BAS 510 F and its metabolite M510F01 were not adequately recovered using the multiresidue methods. BAS 510 F and M510F01 were not evaluated by Protocol A because both analytes do not contain an N-methylcarbamate structure, and BAS 510 F was not evaluated by Protocol B since it is not an acid or a phenol. Methylated M510F01 was recovered at ~55% from Protocol B Florisil, and M510F01 was recovered at ~95% from Protocol B GPC. However, the response of methylated M510F01 at levels necessary to quantify recovery under this protocol was inconsistent and non-linear, therefore no additional work was performed. BAS 510 F and M510F01 had good responses with GC/ECD on a DB-1 column at standard and elevated temperature GC conditions under Protocol C. BAS 510 F and M510F01 were not recovered at ≥30% using Protocols D, E, and F. These data will be forwarded to the FDA for further evaluation.

## 4. Deficiencies

None.

## 5. References

BAS 510 F BAS 510F and M510F01 PMRA a.i. code (CCH) Multiresidue Method OPPTS 860.1360 DACO7.2.4 PC Code: 128008 MRID: 45405107 Submission # 2001-1027, 1036, 1043

None.