

US EPA ARCHIVE DOCUMENT



OECD Quantitative Method for Evaluating Mycobactericidal Activity of Microbicides

Review of the 2013 Collaborative Study

Rebecca Pines, Biologist
Microbiology Laboratory Branch
19 February 2014

Disclaimer:

This information is distributed solely for the purpose of pre-dissemination peer review under applicable information quality guidelines. It has not been formally disseminated by EPA. It does not represent and should not be construed to represent any Agency Determination or Policy.

Background

- In 2013, 6 labs completed a method performance evaluation of the OECD method using *Mycobacterium terrae*
- Labs followed EPA SOP-MB-25-01, based on the OECD test guidelines dated 08/05/11
- A log reduction (LR) of 4 has been proposed as the quantitative performance standard for tuberculocides
 - Target range for the untreated control counts: 4.5-5.5 logs

Terminology

- Control LD: \log_{10} of the total CFUs on a control carrier
- TestLD*: mean of the 4 Control LDs enumerated per test day
- LR: difference between the *TestLD* and the mean of the \log_{10} CFUs for the 3 treated carriers enumerated on the same test day
- S_r : repeatability standard deviation
- S_R : reproducibility standard deviation

Test Conditions and Treatments

- ◉ Test microbe: *Mycobacterium terrae*
 - ATCC # 15755, strain W-45
- ◉ Standard test parameters
 - 5 minute contact time
 - 375 ppm hard water as the diluent
 - 3-part soil load added to the inoculum
- ◉ Treatments
 - Sodium hypochlorite (reference standard)
 - Quaternary ammonium compounds
 - Phenol
 - Citric acid
 - Hydrogen peroxide/peroxyacetic acid


Collaborative Study Plan

- ◉ Step 1: Conducting the OECD method using one reference standard at two levels of presumed efficacy
 - 500 ppm and 2,000 ppm NaOCl
- ◉ Step 2: Conducting the neutralization confirmation assay
 - Each lab evaluated one test chemical
- ◉ Step 3: Conducting the main OECD method performance evaluation
 - 1 concentrated product, 3 ready to use products

Step 1: Reference Standard


- ⦿ 2 efficacy levels for sodium hypochlorite (at presumed high and low efficacy levels)
 - High efficacy level = 2,000 ppm
 - Low efficacy level = 500 ppm

Reference Standard – TestLDs



- Each point is the *TestLD* = mean of the 4 Control LDs from a single test.
- Solid horizontal line indicates the mean *TestLD* of 5.11.
- Dashed lines at 4.5 and 5.5, the limits of validity for *TestLD* specified in the protocol.

Reference Standard – TestLRs




- Each point is a log reduction for a single test.
- Solid horizontal lines indicate the mean LR for each treatment.

Reference Standard – Summary by Lab

Lab	<i>TestLD</i>		LR: Low Efficacy		LR: High Efficacy		Responsiveness	
	Mean	S _r	Mean	S _r	Mean	S _r	Mean	S _r
1	5.602	0.092	1.115	0.119	3.796	0.607	2.681	0.488
2	4.677	0.069	1.098	0.418	3.628	0.565	2.530	0.226
3	5.454	0.055	1.083	0.120	4.667	0.138	3.584	0.101
4	4.844	0.548	0.964	0.326	2.562	0.141	1.598	0.368
5	4.869	0.089	0.628	0.149	2.435	0.347	1.807	0.211
6	5.251	0.130	0.968	0.120	2.626	0.426	1.658	0.306

Responsiveness for each lab is:
 (mean LR for all high efficacy tests – mean LR for all low efficacy tests)

Responsiveness




- Each point is the responsiveness for a single test (LR for high efficacy treatment – LR for low efficacy treatment).
- Solid horizontal line indicates mean responsiveness .

Reference Standard – Conclusions


- ⦿ S_r and S_R values for *TestLD* were small.
- ⦿ For both efficacy levels, S_r and S_R values for LR were small.
- ⦿ The mean LR at the high efficacy level was statistically significantly larger than the mean LR at the low efficacy level ($p\text{-value} \leq 0.001$), confirming responsiveness.
- ⦿ The greater variability associated with the high efficacy treatment may be due to the moderate mean LR (3.28) achieved.

Method Performance – Test LDs



- Each point is the *TestLD* = mean of the 4 Control LDs from a single test.
- Solid horizontal line indicates the mean *TestLD* of 5.09.
- Dashed lines at 4.5 and 5.5, the limits of validity for *TestLD* specified in the protocol.

Method Performance – LRs



- Each point is a log reduction for a single test.
- Solid horizontal lines indicate the mean LR for each treatment.

Method Performance – Summary by Lab


Lab	<i>TestLD</i>		LR-A		LR-B		LR-C		LR-D	
	Mean	S _r	Mean	S _r	Mean	S _r	Mean	S _r	Mean	S _r
1	4.908	0.248	4.187	0.914	5.209	0.248	3.210	0.596	3.858	1.379
2	5.108	0.262	2.194	0.398	5.409	0.262	0.997	0.168	1.734	0.114
3	5.404	0.058	2.327	0.228	4.729	0.825	1.540	0.252	1.818	0.077
4	4.981	0.266	1.248	0.134	5.282	0.266	1.533	0.116	1.834	0.311
5	4.978	0.132	0.975	0.389	5.279	0.132	0.710	0.176	3.038	0.412
6	5.092	0.079	2.421	0.477	5.393	0.079	1.225	0.147	1.598	0.153

Method Performance – Conclusions

- ⦿ S_r and S_R for *TestLD* were small.
- ⦿ For all treatments, S_r and S_R for LR were small.

Next Steps

- ⦿ Increase the control carrier count range to 5.0-6.0 CFU/carrier.
- ⦿ Determine an optical density range in which to harvest the test culture to achieve carrier counts in the appropriate range.
- ⦿ Optimize the growth conditions for the *M. terrae* culture
 - Grow culture with agitation (150 rpm)
- ⦿ Optimize harvesting procedure for *M. terrae*
 - Homogenize the culture prior to carrier inoculation



Acknowledgements

◉ Collaborating laboratories

- EPA, Microbiology Laboratory Branch
- Lonza, Inc.
- ATS Labs
- Clorox Services Company
- MICROBIOTEST
- Stepan Company
- Metrex
- Eurofins

◉ Statistical Analysis

- Al Parker, Montana State University

Questions?