

US EPA ARCHIVE DOCUMENT

Quantitative Petri Plate Method for Vegetative Bacteria: Time Recording
 Sheet for Carrier Transfers
 OPP Microbiology Laboratory

TEST INFORMATION	
EPA Reg. No.	
Product Name	
Sample No.	
Test Organism	
Test Date	
Contact Time	

Exposure Start Time	
Exposure End Time	

Initials/date	Petri Plate No.	Wipe Interval	Carrier Wipe Start Time*	Carrier Wipe End Time (last carrier neutralized by addition of 20mL neutralizer)**
			Timer	Timer
	1			
	2			
	3			
	4			
	5			

*Contact time begins following the entire wiping procedure. Wiping process is conducted within 6±2 seconds.

**At end of contact time neutralizer is added within ±5 seconds for products with contact time of 1-10 minutes. For products with contact times ≤ 1 minute the neutralizer is added within ±2 seconds after contact time has elapsed.