US ERA ARCHIVE DOCUMENT

# Text Searchable File

PMRA Submission #	•	xsulam on the Algae, Selenastrum capricornutum  EPA MRID #: 45834805
Data Requirement:	PMRA DATA COD EPA DP Barcode OECD Data Point EPA MRID EPA Guideline	E {
Test material: Common name: Chemical name:	Penoxsulam XDE-638 IUPAC: Benzenesulfonamide,2-(2 dimethoxy[1,2,4]triazolo[pyrimid CAS name: Not reported CAS No.: Not reported Synonyms: Not reported	
Primary Reviewer: Staff Scientist, Dyna		Signature: Vicac bepon Date: 12/29/03
QC Reviewer: Dana Staff Scientist, Dyna		Signature: Larva works to Date: 12/29/03
Primary Reviewer: {EPA/OECD/PMRA	Bill Erigkson  A) J. GOODYE/	Date: {
Secondary Reviewe {EPA/OECD/PMRA	r(s):{} A}	Date: {}
Company Code { Active Code { EPA PC Code 19	[For PMRA] [For PMRA] [9031   119 03	

CITATION: H.D. Kirk, M.M Gilles, E.L. McClymont, and L.G. McFadden. 2000. Effects of XDE-638 on the Growth of the freshwater green alga, *Selenastrum capricornutum* PRINTZ. Unpublished study performed by Toxicology & Environmental Research and Consulting, The Dow Chemical Company, Midland, Michigan. Laboratory Project Identification No. 991177. Study submitted by Dow AgroSciences, LLC, Indianapolis, Indiana. Experimental start date was July 19, 1999 and the experimental termination date was July 23, 1999. The final report issued January 7, 2000.

**Date Evaluation Completed:** {dd-mmm-yyyy}



#### **EXECUTIVE SUMMARY:**

In a 96-hour acute toxicity study, cultures of *Selenastrum capricornutum* were exposed to Penoxsulam, as XDE-638, under static conditions. The nominal concentrations were 0 (negative and solvent controls), 4.69, 9.38, 18.8, 37.5, 75, 150, and 300  $\mu$ g a.i./. The mean measured concentrations were 4.62, 11.3, 14.6, 34.9, 74.3, 122, and 233  $\mu$ g a.i./L. The 96-hour cell density percent inhibitions were 9.7, 10.1, 16.2, 20.6, 34.6, 44.2, 58.3, and 75.5% for the 4.62, 11.3, 14.6, 34.9, 74.3, 122, and 233  $\mu$ g a.i./L treatment groups, respectively. The solvent control had 9.7% inhibition compared to the negative (dilution water) control. Cell density was significantly reduced at treatment levels equal to and greater than 11.3  $\mu$ g a.i./L. The cell density EC<sub>50</sub> was 92.0  $\mu$ g a.i./L, the EC<sub>05</sub> was 6.5  $\mu$ g a.i./L, and the NOAEC was 4.62  $\mu$ g a.i./L.

The study is scientifically sound and satisfies the U.S. EPA Guideline Subdivision J, §123-2 for an aquatic nonvascular plant study with *Selenastrum capricornutum*. This study is classified as Core.

#### **Results Synopsis**

Test Organism: Selenastrum capricornutum

Test Type: Static

#### **Cell Density:**

NOAEC:  $4.62 \mu g \ a.i./L$ 

EC $_{05}$ : 6.5 µg a.i./L 95% C.I.: 3.2-13 µg a.i./L EC $_{50}$ : 92 µg a.i./L 95% C.I.: 74-110 µg a.i./L

Slope: 1.42±0.150

#### Growth rate:

NOAEC/EC<sub>05</sub>: Not reported

EC<sub>50</sub>: 10,615 μg a.i./L 95% C.I.: Not reported

## Area Under the Growth Curve (Biomass, percent inhibition):

NOAEC/EC<sub>05</sub>: Not reported

EC<sub>50</sub>: 106 μg a.i./L 95% C.I.: 40.2-277 μg a.i./L

Endpoint(s) Affected: cell density

#### I. MATERIALS AND METHODS

GUIDELINE FOLLOWED: The test was based on the following guideline: U.S. EPA-FIFRA Pesticide

Assessment Guidelines, Subdivision J, Hazard Evaluation: Nontarget Plants Guideline 123-2, Growth and Reproduction of Aquatic Plants Tier 2. The

following deviation from U.S. EPA Guideline, §123-2 was noted:

1. The values of pH at test initiation and termination were not specified, but a range was reported.

This deviation did not affect the acceptability or the validity of the study.

**COMPLIANCE:** 

Signed and dated GLP, Quality Assurance and No Data Confidentiality statements

were provided. The study was conducted according to U.S. EPA, OECD, and EC

standards.

A. MATERIALS:

1. Test Material

Penoxsulam, XDE-638

**Description:** 

Pink, solid powder

Lot No./Batch No.:

ND05167938

**Purity:** 

97.5%

#### Stability of Compound

Under Test Conditions: The mean measured concentrations of XDE-638 were 78.0-142% of nominal at hour 0 and 71.3-116% of nominal at hour 96 (Table 3, pp. 26-27). The 96 hour recoveries were 70.6-142% of hour 0 recoveries.

(OECD requires water solubility, stability in water and light, pKa, Pow, vapor pressure of test compound)

Storage conditions of test chemicals: Not reported

#### 2. Test organism:

Name: Selenastrum capricornutum

EPA requires a nonvascular species: For tier I testing, only one species, S. capricornutum, to be tested; for tier II testing, S. costatum, A. flos-aquae, S. capricornutum, and a freshwater diatom is tested

OECD suggests the following species are considered suitable: S. capricornutum, S. subspicatus, and C. vulgaris. If other species are used, the strain should be reported

Strain: 1648

Source: Originally from Starr Algal Collection at the University of Texas, Austin, Texas. Current in-house

laboratory cultures.

Age of inoculum: 5 days old

Method of cultivation: Algal Assay Medium (Appendix A, p. 41).

## **B. STUDY DESIGN:**

a) Range-finding Study: A 96-hour range-finding study with XDE-638 was conducted in order to estimate the nominal test concentrations for the definitive study. The range-finder test concentrations ranged from 9.4 to 300  $\mu$ g a.i./L. The 96-hour EC<sub>50</sub> value and NOAEC were reported as 104  $\mu$ g a.i./L and 19.6  $\mu$ g a.i./L, respectively.

## b) Definitive Study

**Table 1. Experimental Parameters** 

		Remarks
Parameter	Details	Criteria
Acclimation period: culturing media and conditions: (same	Continuous Algal Assay Medium	Inoculum used in test was taken from stock culture and transferred to fresh medium 5 days before testing.
as test or not)	(Appendix A, p. 41); same as test, except for chelant used in cultures.	EPA recommends two week acclimation period.
health: (any toxicity observed)	Not reported	OECD recommends an amount of algae suitable for the inoculation of test cultures and incubated under the conditions of the test and used when still exponentially growing, normally after an incubation period of about 3 days. When the algal cultures contain deformed or abnormal cells, they must be discarded.
Test system static/static renewal: renewal rate for static renewal:	Static	
Incubation facility	Incubator	
Duration of the test	96 hours	
		EPA requires: 96 - 120 hours
		OECD: 72 hours
Test vessel material: (glass/polystyrene) size: fill volume:	Borosilicate Erlenmeyer flasks 250 mL 100 mL	OECD recommends 250 ml conical flasks are suitable when the volume of the test solution is 100 ml or use a culturing apparatus.

		Remarks
Parameter	Details	Criteria
Details of growth medium name: pH at test initiation: pH at test termination: Chelator used: Carbon source: Salinity (for marine algae):	Algal Assay Medium 6.9-7.6 (during entire test) Not reported No NaHCO <sub>3</sub> N/A	The values of pH at test initiation and termination were not reported.  OECD recommends the medium pH after equilibration with air is ~8 with less than .001 mmol/l of chelator if used.  EPA recommends 20X-AAP medium.
If non-standard nutrient medium was used, detailed composition provided (Yes/No)	N/A	
Dilution water source: type: pH: salinity (for marine algae): water pretreatment (if any): Total Organic Carbon: particulate matter: metals: pesticides: chlorine:	Deionized water Not reported 7.0-7.5 N/A None Not reported Not reported Not reported Not detected Not reported	EPA pH: Skeletonema costatum= ~8.0 Others = ~7.5 from beginning to end of the test. EPA salinity: 30- 35 ppt. EPA is against the use of dechlorinated water.  OECD: pH is measured at beginning of the test and at 72 hours, it should not normally deviate by more than one unit during the test.
Indicate how the test material is added to the medium (added directly or used stock solution)	Stock solutions	
Aeration or agitation	Agitation, 100 rpm	EPA recommends agitation only for Selenastrum at 100 cycles per min and Skeletonema at ~60 cycles per min. Aeration is not recommended.
Initial cells density	Approximately 10,000 cells/mL (actual range: 3,691-14,765 cells/mL)	EPA requires an initial number of 3,000 - 10,000 cells/mL. For Selenastrum capricornutum, cell counts on day 2 are not required.  OECD recommends that the initial cell concentration be approximately

		Remarks
Parameter	Details	Criteria
		10,000 cells/ml for S. capricornutum and S. subspicatus. When other species are used the biomass should be comparable.
Number of replicates control: solvent control: treated ones:	3 3 3	Three replicates with plants, one replicate without plants.  EPA requires a negative and/or solvent control with 3 or more replicates per doses. Navicula sp. tests should be conducted with four replicates.  OECD preferably three replicates at each test concentration and ideally twice that number of controls. When a vehicle is used to solubilize the test substance, additional controls containing the vehicle at the highest concentration used in the test cultures should be included in the test.
Test concentrations nominal:  measured:	0 (negative and solvent controls), 4.69, 9.38, 18.8, 37.5, 75, 150, and 300 μg a.i./L  2-40 μg a.i./L (LOQ, controls), 4.62, 11.3, 14.6, 34.9, 74.3, 122, and 233 μg a.i./L	The day 0 LOQ was 40 µg a.i./L and the day 4 LOQ was 2 µg a.i./L. EPA requires at least 5 test concentrations, with each at least 60% of the next higher one.  OECD recommends at least five
	122, and 233 μg a.i./L	concentrations arranged in a geometric series, with the lowest concentration tested should have no observed effect on the growth of the algae. The highest concentration tested should inhibit growth by at least 50% relatively to the control and, preferably, stop growth completely.
Solvent (type, percentage, if used)	Acetone, 100 μL/L	
Method and interval of analytical verification	HPLC; 0 and 96 hours	
Test conditions temperature:	25.3-25.4°C	EPA temperature: Skeletonema:

		Remarks
Parameter	Details	Criteria
photoperiod: light intensity and quality:	Continuous 3230-5400 lux	20°C, Others: 24-25°C; EPA photoperiod: S. costatum 14 hr light/ 10 hr dark, Others: Continuous; EPA light: Anabaena: 2.0 Klux (±15%), Others: 4 - 5 Klux (±15%)  OECD recommended the temperature in the range of 21 to25°C maintained at ± 2°C and continuous uniform illumination provided at approximately 8000 Lux measured with a spherical collector.
Reference chemical {if used} name: concentrations:	N/A	
Other parameters, if any	None	

## 2. Observations:

**Table 2: Observation parameters** 

Parameters	Details	Remarks/Criteria
Parameters measured including the growth inhibition/other toxicity symptoms	Cell count (area under the growth curve and growth rates were calculated).	EPA recommends the growth of the algae expressed as the cell count per mL, biomass per volume, or degree of growth as determined by spectrophotometric means.
Measurement technique for cell density and other end points	Electron particle counting using a Coulter Multisizer.	EPA recommends the measurement technique of cell counts or chlorophyll a  OECD recommends the electronic particle counter, microscope with counting chamber, fluorimeter, spectrophotometer, and colorimeter.

Parameters	Details	Remarks/Criteria
		(note: in order to provide useful measurements at low cell concentrations when using a spectrophotometer, it may be necessary to use cuvettes with a light path of at least 4 cm).
Observation intervals	Every 24 hours	Observations were not conducted every 24 hours.  EPA and OECD: every 24 hours.
Other observations, if any	None	
Indicate whether there was exponential growth in the control	Yes, dilution water group cell densities at test termination was 140.7X greater than the dilution water control group cell densities at test initiation.	EPA requires control cell count at termination to be ≥2X initial count or by a factor of at least 16 during the test.  OECD: cell concentration in control cultures should have increased by a factor of at least 16 within three days.
Were raw data included?	Yes	

#### II. RESULTS and DISCUSSION:

#### A. INHIBITORY EFFECTS:

The 96-hour cell density percent inhibitions were 9.7, 10.1, 16.2, 20.6, 34.6, 44.2, 58.3, and 75.5% for the 4.62, 11.3, 14.6, 34.9, 74.3, 122, and 233 mg a.i./L treatment groups, respectively. The solvent control had 9.7% inhibition compared to the negative (dilution water) control.

Table 3: Effect of Penoxsulam, XDE-638, on Algae (Selenastrum capricornutum)

Treatment mean	Initial cell	Mean Cell density (cells/mL) at		
measured and nominal concentrations a	density (cells/mL) 24 hour		24 hours 96 hours	
(μg a.i./L)			cell count	% inhibition
Dilution water control	7534	32337	1060157	
Solvent control	5675	33998	957742	9.7
4.62 (4.69)	5180	32034	952561	10.1
11.3 (9.38)	6488	31339	888315	16.2
14.6 (18.8)	5586	29747	842261	20.6
34.9 (37.5)	6338	31522	693774	34.6
74.3 (75.0)	6843	28578	591618	44.2
122 (150)	10296	23720	442233	58.3
233 (300)	5075	20475	259843	75.5
Reference chemical	N/A	N/A	N/A	N/A

<sup>&</sup>lt;sup>a</sup> The nominal test concentrations are presented in parentheses.

Table 4: Effect of Penoxsulam, XDE-638, on the Algae Selenastrum capricornutum

Mean Measured and Nominal Treatment Concentrations <sup>a</sup> (µg a.i./L)	Initial cell density (cells/mL)	Mean Growth Rate per day	% inhibition (Mean Growth Rate per day)	Mean Area Under Growth Curve	% inhibition (Mean Area Under Growth Curve)
Dilution water control	7534	NR	NR	NR	NR
Solvent control	5675	NR	NR	NR	NR
4.62 (4.69)	5180	NR	NR	NR	NR
11.3 (9.38)	6488	NR	NR	NR	NR
14.6 (18.8)	5586	NR	NR	NR	NR
34.9 (37.5)	6338	NR	NR	NR	NR _
74.3 (75.0)	6843	NR	NR	NR	NR

Mean Measured and Nominal Treatment Concentrations <sup>a</sup> (µg a.i./L)	Initial cell density (cells/mL)	Mean Growth Rate per day	% inhibition (Mean Growth Rate per day)	Mean Area Under Growth Curve	% inhibition (Mean Area Under Growth Curve)
122 (150)	10296	NR	NR	NR	NR
233 (300)	5075	NR	NR	NR	NR
Reference chemical	N/A	N/A	N/A	N/A	N/A

<sup>&</sup>lt;sup>a</sup> The nominal test concentrations are presented in parentheses.

NR-Not Reported

Table 5: Statistical endpoint values.

Statistical Endpoint	Biomass (percent inhibition)	Growth rate	Cell density
NOAEC or EC <sub>05</sub> (μg a.i./L)	NR	NR	11.3
EC <sub>50</sub> (μg a.i./L)	106	10,615	86.4
IC <sub>50</sub> or EC <sub>50</sub> (μg a.i./L) (95% C.I.)	40.2-277	NR	33.6-223
IC <sub>25</sub> /EC <sub>25</sub> (μg a.i./L) (95% C.I.)	NR	NR	20.7 (8.08-52.8)
Reference chemical, if used NOAEC IC <sub>50</sub> /EC <sub>50</sub>	N/A	N/A	N/A

NR-Not Reported

N/A = Not applicable.

#### **B. REPORTED STATISTICS:**

Statistical Method: The  $EC_{25}$  and  $EC_{50}$  values were calculated using least squares linear regression for algal cell counts. The NOAEC was determined using analysis of variance and the Dunnett's t-test. All statistical calculations were performed using the mean measured concentrations.

### **Cell Density:**

NOAEC/EC<sub>05</sub>: 11.3 µg a.i./L

EC<sub>50</sub>: 86.4 μg a.i./L 95% C.I.: 33.6-638 μg a.i./L

#### Growth rate:

NOAEC/EC<sub>05</sub>: Not reported

EC<sub>50</sub>: 10,615 μg a.i./L 95% C.I.: Not reported

## Area Under the Growth Curve (Biomass, percent inhibition):

NOAEC/EC<sub>05</sub>: Not reported

EC<sub>50</sub>: 106 µg a.i./L 95% C.I.: 40.2-277 µg a.i./L

Endpoint(s) Affected: cell density

#### C. VERIFICATION OF STATISTICAL RESULTS:

Statistical Method: Cell density data satisfied the assumptions of ANOVA. The negative and solvent controls were compared using a Student's t-test and, upon finding no difference, the two were pooled for comparison to treatment. The NOAEC was determined using this test, followed by William's multiple comparison test via TOXSTAT statistical software. The EC<sub>x</sub> values were determined using the Probit method via Nuthatch statistical software. Results could not be verified for biomass or growth rate, as replicate data were not provided. The reviewer used the mean measured concentrations to calculate toxicity values.

#### Cell Density:

NOAEC: 4.62 µg a.i./L

EC<sub>05</sub>: 6.5 μg a.i./L EC<sub>50</sub>: 92 μg a.i./L 95% C.I.: 3.2-13 μg a.i./L 95% C.I.: 74-110 μg a.i./L

Slope:  $1.42 \pm 0.150$ 

757

#### D. STUDY DEFICIENCIES:

The deviation was minor, so it did not affect the acceptability or validity of the study.

## **E. REVIEWER'S COMMENTS:**

The reviewer's toxicity values for cell density differed from those of the study authors. Because the reviewer's  $EC_{50}$  estimate was associated with a narrower 95% confidence interval, it is reported in the Executive Summary and Conclusions sections (as is the reviewer's NOAEC estimate).

F. CONCLUSIONS: The study is scientifically sound and satisfies the guidelines for an aquatic nonvascular plant study with *Selenastrum capricornutum* [§123-2]. This study is classified as Core. Cell density was significantly reduced at treatment levels equal to and greater than 11.3  $\mu g$  a.i./L. The cell density EC<sub>50</sub> was 92.0  $\mu g$  a.i./L and the NOAEC was 4.62  $\mu g$  a.i./L.

#### Cell Density:

NOAEC: 4.62 µg a.i./L

EC<sub>05</sub>: 6.5 μg a.i./L EC<sub>50</sub>: 92 μg a.i./L 95% C.I.: 3.2-13 μg a.i./L 95% C.I.: 74-110 μg a.i./L

Slope:  $1.42 \pm 0.150$ 

Growth rate:

NOAEC/EC<sub>05</sub>: Not reported

EC<sub>50</sub>: 10,615 μg a.i./L 95% C.I.: Not reported

## Area Under the Growth Curve (Biomass, percent inhibition):

NOAEC/EC<sub>05</sub>: Not reported

EC<sub>50</sub>: 106  $\mu$ g a.i./L 95% C.I.: 40.2-277  $\mu$ g a.i./L

Endpoint(s) Affected: cell density

## III. REFERENCES:

- Holst, R.W. and T.C. Ellwanger, 1982, Pesticide Assessment Guidelines Subdivision J Hazard Evaluation: Non-target Plants, EPA 540/9-82-020, Washington, D.C.
- Holst, R.W., 1986, Hazard Evaluation Division: Standard Evaluation Procedure Non-Target Plants: Growth and Reproduction of Aquatic Plants Tiers 1 and 2. EPA 540/9-86-134, Washington, D.C.
- Organization of Economic Cooperation and Development (OECD). OECD Guideline for Testing of Chemicals. Algal Growth, Inhibition Test. Number 201. Adopted 7 June, 1984.
- Environmental Protection Agency-FIFRA GLPs. Title 40 CFR, 160-Federal Insecticide, Fungicide and Rodenticide Act (FIFRA); Good Laboratory Practice Standards, Final Rule.
- OECD Series on Principles of Good Laboratory Practice Compliance and Monitoring, Number 1. OECD Principles on Good Laboratory Practice (as revised in 1997) ENV/MC/CHEM(98)17.
- EC Directive 99/11/EC of 8 March 1999 (OJ No. L 77/8-21, 23/3/1999).
- Smith, A.J., "Certificate of Analysis for Test/Reference/Control/Substances Analytical Report FA & PC Number 993090. 20 May, 1999.
- A.J. Smith, Purity Report for XDE-638, FA &PC 993090, May 20, 1999.
- Miller, W.E., Green, J.C. and Shiroyama, T. (1978). The *Selenastrum capricornutum* Printz Algal Assay Bottle Test. EPA-600/9-78-018.5.
- Kirk, H.D. and Gilles, M.M., "Phase Two Aquatic Screening of the Experimental Sulfonamide XR-638 with the Freshwater Green Alga, *Selenastrum capricornutum* Printz." Report 980324. The Dow Chemical Company. 27 August, 1998.
- Neter, J., Wasserman, W. and Kutner, M.H. (1983). Applied Linear Regression Models. Richard D. Irwin Inc., Homewood, Illinois.
- Winer, B.J. (1971). Statistical Principles on Experimental Design. 2<sup>nd</sup> Ed., McGraw Hill, Co. New York, New York.

## APPENDIX I. OUTPUT OF REVIEWER'S STATISTICAL VERIFICATION: cell density

File: 4805cd

Transform: NO TRANSFORMATION

ANOVA TABLE

SOURCE DF SS MS F

Between 7 1682495151742.000 240356450248.750

34.287

Within (Error) 19 133193629408.000 7010191021.477

Total

26 1815688781150.000

Critical F value = 2.54 (0.05,7,19)

Since F > Critical F REJECT Ho:All groups equal

cell density

File: 4805cd

Transform: NO TRANSFORMATION

BONFERRONI T-TEST - TABLE 1 OF 2

Ho:Control<Treatment

TRANSFORMED MEAN CALCULATED IN

GROUP	IDENTIFICATION	MEAN

ORIGINAL UNITS T STAT SIG

1	GRPS 1&2	POOLED 100	04449.833	1004449.833	
2	4.62	952560.667	952560.6	67 0.876	
3	11.3	888315.000	888315.0	00 1.962	
4	14.6	842261.333	842261.3	33 2.739	*
5	34.9	693773.667	693773.6	67 5.248	*
6	74.3	591618.333	591618.3	33 6.973	*
7	122	442232.667	442232.6	67 9.496	*
8	233	259843.333	259843.3	33 12.577	7 *

Bonferroni T table value = 2.70 (1 Tailed Value, P=0.05, df=19,7)

cell density

File: 4805cd

Transform: NO TRANSFORMATION

BONFERRONI T-TEST - TABLE 2 OF 2 Ho:Control<Treatment

NUM OF Minimum Sig Diff % of DIFFERENCE GROUP IDENTIFICATION REPS (IN ORIG. UNITS) CONTROL FROM CONTROL

1 GRPS 1&2 POOLED 6

2 4.62 3 159672.776 15.9 51889.167

11.3 3 159672.776 15.9 116134.833 3

4	14.6	3	159672.776	15.9	162188.500
5	34.9	3	159672.776	15.9	310676.167
6	74.3	3	159672.776	15.9	412831.500
7	122	3	159672.776	15.9	562217.167
8	233	3	159672.776	15.9	744606.500
6 7	74.3 122	3	159672.776 159672.776	15.9 15.9	412831.500 562217.167

cell density

File: 4805cd

Transform: NO TRANSFORMATION

WILLIAMS TEST (Isotonic regression model) TABLE 1 OF 2

CDOID	·	ODIODIAL	TD ANGE	CONTED	ICOTONIATED
GROUP		ORIGINAL	TRANSF	ORMED	ISOTONIZED
II	DENTIFICATION	N ME	CAN M	EAN	MEAN
1	GRPS 1&2 POOL	ED 6 10044	449.833 100	04449.833	1004449.833
2	4.62 3 95	2560.667	952560.667	952560.6	67
3	11.3 3 88	8315.000	888315.000	888315.0	00
4	14.6 3 84	2261.333	842261.333	842261.3	33
5	34.9 3 69	3773.667	693773.667	693773.6	67
6	74.3 3 59	1618.333	591618.333	591618.3	33
7	122 3 44	2232.667	442232.667	442232.6	67
8	233 3 25	9843.333	259843.333	259843.3	33

cell density

File: 4805cd

Transform: NO TRANSFORMATION

WILLIAMS TEST (Isotonic regression model) TABLE 2 OF 2

ISOTONIZED CALC. SIG TABLE DEGREES OF IDENTIFICATION MEAN WILLIAMS P=.05 WILLIAMS FREEDOM

GRPS 1&2 POOLED1004449.833

4.62 952560.667 0.876 1.73 k= 1, v=19

11.3 888315.000 1.962 \* 1.81 k= 2, v=19

14.6 842261.333 2.739 \* 1.84 k= 3, v=19

14.6 842261.333 2.739 \* 1.84 k= 3, v=19
34.9 693773.667 5.248 \* 1.85 k= 4, v=19
74.3 591618.333 6.973 \* 1.86 k= 5, v=19
122 442232.667 9.496 \* 1.87 k= 6, v=19
233 259843.333 12.577 \* 1.87 k= 7, v=19

s = 83726.884

Note: df used for table values are approximate when v > 20.

**Estimates of EC%** 

Παγε 15 οφ 16

Parameter	er Estimate		95% Bounds		Std.Err.	Lower Bound
	Low	er (	Jpper	/E	stimate	
EC5	6.5	3.2	13.	0.15	0.49	
EC10	12.	6.5	21.	0.12	0.56	
EC25	31.	21.	46.	0.082	0.68	
EC50	92.	74.	1.1E+0	0.0	<b>46</b> 0.	80

Slope = 1.42 Std.Err. = 0.150

Goodness of fit: p = 0.74 based on DF= 19.

4805CD : cell density

Observed vs. Predicted Treatment Group Means Dose #Reps. Obs. Pred. Obs. Pred. %Change Mean Mean -Pred. %Control 0.00 6.00 1.00e+06 9.91e+05 1.31e+04 100. 0.00 4.62 3.00 9.53e+05 9.60e+05 -7.09e+03 96.8 3.20 11.3 3.00 8.88e+05 8.95e+05 -6.87e+03 90.3 9.70 14.6 3.00 8.42e+05 8.65e+05 -2.31e+04 87.3 12.7 34.9 3.00 6.94e+05 7.20e+05 -2.61e+04 72.6 27.4 74.3 3.00 5.92e+05 5.48e+05 4.32e+04 55.3 44.7 122. 3.00 4.42e+05 4.28e+05 1.47e+04 43.1 56.9 233. 3.00 2.60e+05 2.81e+05 -2.09e+04 28.3 71.7