

US EPA ARCHIVE DOCUMENT



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

Office of Research and Development
Washington, D.C. 20460



**ENVIRONMENTAL TECHNOLOGY VERIFICATION PROGRAM
VERIFICATION STATEMENT**

TECHNOLOGY TYPE:	POLYCHLORINATED BIPHENYL (PCB) FIELD ANALYTICAL TECHNIQUES
APPLICATION:	MEASUREMENT OF PCBs IN SOILS AND SOLVENT EXTRACTS
TECHNOLOGY NAME:	EnviroGard PCB TEST KIT
COMPANY: ADDRESS:	STRATEGIC DIAGNOSTICS INC. 111 PENCADER DRIVE NEWARK, DE 19702-3322
PHONE:	(302) 456-6789

The U.S. Environmental Protection Agency (EPA) has created a program to facilitate the deployment of innovative technologies through performance verification and information dissemination. The goal of the Environmental Technology Verification (ETV) Program is to further environmental protection by substantially accelerating the acceptance and use of improved and more cost-effective technologies. The ETV Program is intended to assist and inform those involved in the design, distribution, permitting, and purchase of environmental technologies. This document summarizes the results of a demonstration of the Strategic Diagnostics Inc. (SDI) EnviroGard PCB test kit.

PROGRAM OPERATION

EPA, in partnership with recognized testing organizations, objectively and systematically evaluates the performance of innovative technologies. Together, with the full participation of the technology developer, they develop plans, conduct tests, collect and analyze data, and report findings. The evaluations are conducted according to a rigorous demonstration plan and established protocols for quality assurance. EPA's National Exposure Research Laboratory, which conducts demonstrations of field characterization and monitoring technologies, with the support of the U.S. Department of Energy's (DOE) Environmental Management (EM) program, selected Oak Ridge National Laboratory (ORNL) as the testing organization for the performance verification of polychlorinated biphenyls (PCBs) field analytical techniques.

DEMONSTRATION DESCRIPTION

In July 1997, the performance of six PCB field analytical techniques was determined under field conditions. Each technology was independently evaluated by comparing field analysis results with those obtained using approved reference methods. Performance evaluation (PE) samples were also used to assess independently the accuracy and comparability of each technology.

The demonstration was designed to detect and measure PCBs in soil and solvent extracts. The demonstration was conducted at ORNL in Oak Ridge, Tennessee, from July 22 through July 29, 1997. The study was conducted under two environmental conditions. The first site was outdoors, with naturally fluctuating temperatures and relative humidity conditions. The second site was inside a controlled environmental chamber, with generally cooler temperatures and lower relative humidities. Multiple soil types, collected from sites in Ohio, Kentucky, and Tennessee, were analyzed in this

US EPA ARCHIVE DOCUMENT

study. Solutions of PCBs were also analyzed to simulate extracted surface wipe samples. The results of the soil and extract analyses conducted under field conditions by the technology were compared with results from analyses of homogenous replicate samples conducted by conventional EPA SW-846 methodology in an approved reference laboratory. Details of the demonstration, including a data summary and discussion of results, may be found in the report entitled *Environmental Technology Verification Report: Immunoassay Kit, Strategic Diagnostics Inc., EnviroGard PCB Test Kit*, EPA/600/R-98/113.

TECHNOLOGY DESCRIPTION

The EnviroGard PCB test kit is a competitive binding enzyme immunoassay that performs rapid, interval testing for PCBs in soils and solutions at specified action levels of 1, 5, 10, and 50 parts per million (ppm). These results are reported in intervals (e.g., < 1 ppm, 1 to 5 ppm, etc.) rather than in specific quantities (e.g., 6.7 ppm). The test kit is standardized using Aroclor 1248, but it can also detect Aroclors 1016, 1242, 1254, and 1260. The following items are needed to run a test: the EnviroGard PCB test kit, the Soil Extraction Bottle kit, the equipment contained in the Soil Field Lab, methanol, and water. The test procedure entails collecting a 5-g soil sample and extracting the PCBs from it using methanol. To initiate the PCB test, PCB-enzyme conjugate is added to the antibody-coated test tubes. The soil extract sample is then added to the test tube. After a 15-min incubation period, the tubes are rinsed and a color developing solution is added. Color development is inversely related to the PCB concentration (e.g. the darker the color, the less analyte PCB is present in the sample). PCBs are detected using a photometer that measures the absorbance of each tube.

VERIFICATION OF PERFORMANCE

The following performance characteristics of the EnviroGard PCB test kit were observed:

Throughput: Throughput was 18 samples/hour under outdoor conditions and 9 to 10 samples/hour under chamber conditions. This rate included sample preparation and analysis.

Ease of Use: Three operators analyzed samples during the demonstration, but the technology can be run by a single trained operator. Minimal training (2 to 4 h) is required to operate the EnviroGard kit, provided the user has a fundamental understanding of basic chemical and field analytical techniques.

Completeness: The EnviroGard kit generated results for all 232 PCB samples for a completeness of 100%.

Blank results: All of the blank soil samples were reported as the lowest reporting interval, which included zero; therefore, the percentage of false positive results was 0%. One false positive result (13%) was reported for the extract samples. The EnviroGard kit reported no false negative results.

Precision: The overall precision—based on the percentage of combined sample sets where all four replicates were reported as the same interval—was 38% for the PE soils, 47% for the environmental soils, and 67% for the extracts.

Accuracy: Accuracy was assessed using PE soil and extract samples. Accuracy, defined as the percentage of EnviroGard results that agreed with the accepted concentration, was 51% for PE soils and 58% for extracts. In general, the percentage of samples that was biased high was much greater (47% for PE soils and 38% for extracts) than the percentage biased low (1% for PE soils and 4% for extracts).

Comparability: Comparability, like accuracy, was defined as the percentage of results that agreed with, was above (i.e., biased high), or was below (i.e., biased low) the reference laboratory result. The percentage of samples that agreed with the reference laboratory results was 53% for all soils (PE and environmental) and 63% for extracts. The percentage of samples that was biased high was again much greater (45% for soils and 38% for extracts) than the percentage that was

biased low (2% for soils and 0% for extracts).

Regulatory Decision-making: One objective of this demonstration was to assess the technology's ability to perform at regulatory decision-making levels for PCBs, specifically 50 ppm for soils and 100 $\mu\text{g}/100\text{cm}^2$ for surface wipes. For PE and environmental soil samples in the range of 40 to 60 ppm, 39% of the EnviroGard results agreed with the reference laboratory. In contrast, 59% were biased high and 2% were biased low. For extract samples representing surface wipe sample concentrations of 100 $\mu\text{g}/100\text{cm}^2$ and 1000 $\mu\text{g}/100\text{cm}^2$ (assuming a 100 cm^2 wipe sample), 63% of the EnviroGard results agreed with the extract spike concentration. In comparison, the percentage of extract samples that was biased high was 38%, and the percentage of samples that was biased low was 4%.

Data quality levels: The performance of the EnviroGard PCB test kit was characterized as biased and imprecise about 50% of the time, because nearly half of the data were biased relative to the accepted concentration values (in terms of accuracy) and had replicate results that were not reported as the same interval (in terms of precision). It should be noted that there was an increased likelihood of results being biased high as a result of the conservatism that the manufacturer has incorporated into the calculation of results.

The results of the demonstration show that the EnviroGard PCB test kit can provide useful, cost-effective data for environmental problem-solving and decision-making. Undoubtedly, it will be employed in a variety of applications, ranging from serving as a complement to data generated in a fixed analytical laboratory to generating data that will stand alone in the decision-making process. As with any technology selection, the user must determine if this technology is appropriate for the application and the project data quality objectives. For more information on this and other verified technologies, visit the ETV web site at <http://www.epa.gov/etv>.

Gary J. Foley, Ph.D.
Director
National Exposure Research Laboratory
Office of Research and Development

NOTICE: EPA verifications are based on an evaluation of technology performance under specific, predetermined criteria and the appropriate quality assurance procedures. EPA makes no expressed or implied warranties as to the performance of the technology and does not certify that a technology will always, under circumstances other than those tested, operate at the levels verified. The end user is solely responsible for complying with any and all applicable Federal, State, and Local requirements.

